

1

2

3

4

5

6

7 **Dialysis Machine HL7 Implementation**

8 **Guidelines**

9

10

11

12

13 Copyright Licenses – HL7®

14 HL7 licenses its standards and select IP free of charge. If you did not acquire a free license from HL7 for this
15 document, you are not authorized to access or make any use of it. To obtain a free license, please visit
16 <http://www.HL7.org/implement/standards/index.cfm>.

17 If you are the individual that obtained the license for this HL7 Standard, specification or other freely licensed work (in
18 each and every instance "Specified Material"), the following describes the permitted uses of the Material.

19 **A. HL7 INDIVIDUAL, STUDENT AND HEALTH PROFESSIONAL MEMBERS**, who register and agree to the terms of
20 HL7's license, are authorized, without additional charge, to read, and to use Specified Material to develop and sell
21 products and services that implement, but do not directly incorporate, the Specified Material in whole or in part without
22 paying license fees to HL7.

23 **INDIVIDUAL, STUDENT AND HEALTH PROFESSIONAL MEMBERS** wishing to incorporate additional items of
24 Special Material in whole or part, into products and services, or to enjoy additional authorizations granted to HL7
25 ORGANIZATIONAL MEMBERS as noted below, must become ORGANIZATIONAL MEMBERS of HL7.

26 **B. HL7 ORGANIZATION MEMBERS**, who register and agree to the terms of HL7's License, are authorized, without
27 additional charge, on a perpetual (except as provided for in the full license terms governing the Material), non-
28 exclusive and worldwide basis, the right to (a) download, copy (for internal purposes only) and share this Material
29 with your employees and consultants for study purposes, and (b) utilize the Material for the purpose of developing,
30 making, having made, using, marketing, importing, offering to sell or license, and selling or licensing, and to
31 otherwise distribute, Compliant Products, in all cases subject to the conditions set forth in this Agreement and any
32 relevant patent and other intellectual property rights of third parties (which may include members of HL7). No other
33 license, sublicense, or other rights of any kind are granted under this Agreement.

34 **C. NON-MEMBERS**, who register and agree to the terms of HL7's IP policy for Specified Material, are authorized,
35 without additional charge, to read and use the Specified Material for evaluating whether to implement, or in
36 implementing, the Specified Material, and to use Specified Material to develop and sell products and services that
37 implement, but do not directly incorporate, the Specified Material in whole or in part.

38 **NON-MEMBERS** wishing to incorporate additional items of Specified Material in whole or part, into products and
39 services, or to enjoy the additional authorizations granted to HL7 ORGANIZATIONAL MEMBERS, as noted above,
40 must become ORGANIZATIONAL MEMBERS of HL7.

41 Please see <http://www.HL7.org/legal/ippolicy.cfm> for the full license terms governing the Material.
42

43 Copyright Licenses - IHE

44 IHE International hereby grants to each Member Organization, and to any other user of these documents, an
45 irrevocable, worldwide, perpetual, royalty-free, nontransferable, nonexclusive, non-sublicensable license under its
46 copyrights in any IHE profiles and Technical Framework documents, as well as any additional copyrighted materials
47 that will be owned by IHE International and will be made available for use by Member Organizations, to reproduce
48 and distribute (in any and all print, electronic or other means of reproduction, storage or transmission) such IHE
49 Technical Documents.

50 The licenses covered by this Copyright License are only to those copyrights owned or controlled by IHE International
51 itself. If parts of the Technical Framework are included in products that also include materials owned or controlled by
52 other parties, licenses to use those products are beyond the scope of this IHE document and would have to be obtained
53 from that other party.
54

55 Copyright of Base Standards

56 IHE technical documents refer to and make use of a number of standards developed and published by several
57 standards development organizations. All rights for their respective base standards are reserved by these
58 organizations. This agreement does not supersede any copyright provisions applicable to such base standards.

59 Health Level Seven® has granted permission to IHE to reproduce tables from the HL7 standard. The HL7 tables in this
60 document are copyrighted by Health Level Seven. All rights reserved. Material drawn from these documents is credited
61 where used.
62

63 Trademark

64 IHE® and the IHE logo are trademarks of the Healthcare Information Management Systems Society in the United
65 States and trademarks of IHE Europe in the European Community. They may only be used with the written consent
66 of the IHE International Board Operations Committee, which may be given to a Member Organization in broad terms
67 for any use that is consistent with the IHE mission and operating principles.
68

69 Disclaimer Regarding Patent Rights

70 Attention is called to the possibility that implementation of the specifications in this document may require use of
71 subject matter covered by patent rights. By publication of this document, no position is taken with respect to the
72 existence or validity of any patent rights in connection
73

Revision History

Version	Change Description	Date (mm/yy)
1.0	Initial Release	TBD

74

75

76

77

Table of Contents

78			
79	1	Introduction & Scope	6
80	1.1	Purpose.....	6
81	1.2	Audience.....	6
82	1.3	Scope	6
83	1.4	Assumptions	6
84	1.5	Conventions.....	6
85	1.6	Referenced Standards	7
86	1.7	Glossary	7
87	2	Patient Identification.....	8
88	3	Time Synchronization	8
89	4	Reporting Treatment Information.....	8
90	4.1	Overview	8
91	4.2	Device Request.....	9
92	4.2.1	Message Structure.....	9
93	4.2.2	Example 1 – Minimal Message	10
94	4.2.3	Example 2 – Full Message.....	11
95	4.3	EMR Response.....	16
96	4.3.1	Message Structure.....	16
97	4.3.2	Example 1 – Accepted Response	16
98	4.3.3	Example 2 – Rejected Response	16
99	5	Reporting Alarm Information	17
100	5.1	Overview	17
101	5.2	Device Request.....	18
102	5.2.1	Message Structure.....	18
103	5.2.2	Example 1 - Alarm Initiation	25
104	5.2.3	Example 2 - Alarm Cessation	25
105	5.3	EMR Response.....	25
106	5.3.1	Message Structure.....	25
107	5.3.1	Example 1 – Accepted Response	26
108	6	HL7 Data Elements	26
109	6.1	Segments.....	26
110	6.1.1	MSA – Message Acknowledgement Segment	26
111	6.1.2	MSH – Message Header Segment.....	27
112	6.1.3	OBR Segment – Observation Request	29
113	6.1.4	OBX Segments – Dialysis Observation Results	30
114	6.1.5	PID – Patient Identification Segment	33
115	6.1.6	PV1 Segment.....	34
116	6.2	HL7 Tables	36
117	7	Dialysis Data Elements	39
118	7.1	Tables	39
119	7.2	Message Element Terminology	43
120	7.3	Alarms/Alerts	57
121			
122			

1 Introduction & Scope

1.1 Purpose

The Dialysis Machine Implementation Guide is a messaging specification intended to standardize the electronic reporting of patient, clinical, and machine data to electronic health record (EHR) systems and/or electronic medical record (EMR) systems found in most care settings today. The goal of this specification is to provide an applicable set of data definitions and structures that can be adopted as a worldwide industry standard, thereby obviating the need for proprietary solutions for dialysis healthcare interoperability.

1.2 Audience

The standards defined in this guide are intended for use by data exchange systems found in Renal Care groups listed below:

- Renal Care Providers
- Dialysis System Manufacturers
- Other groups or agencies that employ a dialysis machine to EHR/EMR interface

1.3 Scope

This specification focuses exclusively on the electronic reporting of device treatment data from a dialysis machine to an EHR/EMR for hemodialysis therapies based on HL7 v2.6 messaging. Specifically, the patient info, treatment system readings, thresholds, and alarms which may be generated during machine usage can be communicated.

This specification does NOT cover peritoneal dialysis or communication from an EHR/EMR to the dialysis machine.

1.4 Assumptions

This document assumes that users have a need to exchange treatment data between dialysis devices and EHRs. It also assumes that the user is familiar with the HL7 v2.6 standard as well as the IHE PCD specifications.

1.5 Conventions

This specification is primarily based on the Integrating the Healthcare Enterprise (IHE), Patient Care Device (PCD) technical framework, which is based on a profile as defined in the Health Level 7 (HL7)'s v2.6 Normative Standard. A message profile is an unambiguous specification of an HL7 message type intended for a particular use case, as defined in Section 2.B of the HL7 v2.6 Standard. Additionally, the message structures and terminology are based on the ISO/IEEE 11073 Standard for Medical Device Communication.

An HL7 message profile defines both the *dynamic* aspects of information interchange (i.e., the systems that participate in such interchanges and the real-world events that trigger messaging) as well as the *static* aspects of messaging (i.e., the structure and contents of the electronic messages that are exchanged).

Optionality (aka. Usage) is defined as one of the values in the following table.

Usage Values

Value	Description
B	Left for backwards compatibility.
C	Conditional
CE	Conditional but maybe empty.
O	Optional
R	Required
RE	Required but maybe Empty. The field or data type component description must stipulate when the field or data type component may be empty.

Usage Values

Value	Description
W	Withdrawn
X	Not used

The IHE Patient Care Device Technical Framework uses an information model and a nomenclature from the IEEE 11073. The information model is defined in ISO/IEEE 11073- 10201 Health Informatics – Point-of-care medical device communication – Part 10201: Domain Information Model. The nomenclature is defined in ISO/IEEE 11073-10101 Health Informatics – Point -of-care medical device communication – Part 10101: Nomenclature. Familiarity with these standards is necessary for implementers of the Device Observation Reporter and Device Observation Consumer Actors.

HL7 V2.6 Chapter 7 Observation Reporting defines the general HL7 syntax and coding requirements related to observation reporting, used for PCD data communications in the PCD TF. Familiarity with HL7 Chapter 7 is necessary for implementers of the PCD TF transactions.

This Technical Framework specifies conventions that are used to represent the information model hierarchy for medical devices embodied in the IEEE 11073 Domain Information Model within the syntactic and semantic conventions of HL7 v. 2.6

1.6 Referenced Standards

Ref	Title
[1]	IHE Patient Care Device, Technical Framework Volume 1 - Profiles, Rev 9.0 – Final Text, December 12, 2019
[2]	IHE Patient Care Device, Technical Framework Volume 2 - Transactions, Rev 9.0 – Final Text, December 12, 2019
[3]	IHE Patient Care Device, Technical Framework, Volume 3 – Semantic Content, Rev 9.0 – Final Text, December 12, 2019
[4]	IHE IT Infrastructure (ITI) Technical Framework Volume 1 - Integration Profiles, Rev 14.0, July 21, 2017
[5]	HL7 v2.6 Normative Standard, October 2007
[6]	ISO/IEEE 11073-10101R Health informatics – Point-of-care, Medical Device Communication – Part 10101: Nomenclature, 2015-09
[7]	ISO/IEEE 11073-10201 Health informatics – Point-of-care, Medical Device Communication – Part 10201: Domain information model, First Edition, 2004-12-15
[8]	ISO/IEEE 11073-10101R Health informatics – Point-of-care, Medical Device Communication – Part 10101b: Nomenclature
[9]	ISO/IEEE 11073-10301 Health informatics – Point-of-care, Medical Device Communication – Part 10201: Domain information model, First Edition, 2004-12-15
[10]	ISO/IEEE 11073-20101 Health informatics – Point-of-care, Medical Device Communication – Part 20101: Application Profiles – Base Standard, First Edition, 2004-12-15

1.7 Glossary

DOC: Device Observation Consumer

DOR: Device Observation Report

IHE: Integrating the Healthcare Enterprise – An initiative by healthcare professionals and industry to improve the way computer systems in healthcare share information.

MDS: Medical Device System

MDC: Medical Device Communications

MLLP: Minimal Lower Layer Protocol – is used for transferring HL7 messages over Ethernet. It defines delimiters which identify the beginning and ends of the HL7 message.

178 **PCD:** Patient Care Device – An IHE domain formed to address the integration of medical devices
179 into the healthcare enterprise.

180 **VMD:** Virtual Medical Device

181 **UOM:** Unit of Measure

182 **CARD:** Cardinality

184 2 Patient Identification

185 Patient Identification is perhaps the most essential infrastructural component of any interoperability and
186 communication process, particularly when PCD data is exported to the enterprise. It is the key element in medical
187 device, communication, data analysis, reporting and record keeping. Automation of the entry of patient identification
188 to patient care device has the potential for improving throughput, reducing errors, increasing safety and device and
189 drug effectiveness, and efficiency. It is strongly recommended that implementations use IHE Patient Demographics
190 Query (PDQ) Integration Profile compliant transactions for acquisition of Patient Identification credentials. These
191 transactions include: ITI-21, ITI-30 and ITI-31. Other mechanisms such as bar code or RFID are also valid alternatives or
192 complements.

193 If the Dialysis Machine has no way of obtaining a Patient Identifier, it will report a string value that is
194 Model/Serial_Number.

195 3 Time Synchronization

196 The data reported by a dialysis machine must be reconcilable with data from other parts of the health care enterprise.
197 In order for this to occur, time must be synchronized between the machines involved in the transactions. To this end,
198 dialysis machines compliant with this standard will use the IHE Consistent Time (CT) Protocol. The Consistent Time
199 Profile specifies the use of the Network Time Protocol (NTP) defined in RFC1305.


200 4 Reporting Treatment Information

201 4.1 Overview

202 The sending of Treatment Information to the EMR System is accomplished using the Patient Care Device (PCD)
203 domain in the Device Enterprise Communication (DEC) IHE profile. In doing so, the Dialysis Machine acts as the
204 Device Observation Report (DOR) and the EMR is the Device Observation Consumer (DOC).

205

Figure 1 - Communicate PCD Data Interaction Diagram



206

207 In general, the Dialysis Machine sends periodic reports at an interval of between several times per minute (high
 208 acuity) and a maximum interval of 24 hours (chronic, home health). The minimum and maximum intervals are
 209 configured at implementation. The Dialysis Machine may also send aperiodic reports for "event type" information.

210 This specification includes 4 different observation types based on differing real-world triggers. The elements for
 211 each observation are indicated in Appendix B – Dialysis Data Elements. The types are:

- 212 • Status – Sent initially, terminally, and when the treatment status changes
- 213 • Parameter – Sent initially and when any parameter changes
- 214 • Identifier – Sent with all messages
- 215 • Blood Pressure – Sent when communicating a measurement

216

217 A Status Report will contain all mandatory elements, the applicable conditional elements, and the optional
 218 elements supported by the dialysis machine.

219

220 4.2 Device Request

221 4.2.1 Message Structure

222 Segments shown in gray are not used in this standard.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
[[SFT]]	Software Segment	X	[0..0]

Segment	Meaning	Use	Card
[UAC]	User Authentication Credential	X	[0..0]
{	--- PATIENT_RESULT begin		
[--- PATIENT begin		
PID	Patient Identification	R	[1..1]
[PD1]	Additional Demographics	X	[0..0]
{{PRT}}	Participation	X	[0..0]
{{NTE}}	Notes and Comments	X	[0..0]
{{NK1}}	Next of Kin/Associated Parties	O	[0..3]
[--- VISIT begin		
PV1	Patient Visit	R	[1..1]
[PV2]	Patient Visit – Additional Info	X	[0..0]
{{PRT}}	Participation	X	[0..0]
]	--- VISIT end		
]	--- PATIENT end		
{	---ORDER_OBSERVATION begin		
[ORC]	Order Common	X	[0..0]
OBR	Observation Request	R	[1..1]
{{NTE}}	Notes and Comments	O	[0..1]
{{PRT}}	Participation	X	[0..0]
{{	--- TIMING_QTY begin		
TQ1	Timing/Quantity	X	[0..0]
{{TQ2}}	Timing/Quantity Order Sequence	X	[0..0]
}}	--- TIMING_QTY end		
[CTD]	Contact Data	X	[0..0]
{{	--- OBSERVATION begin		
OBX	Observation Result	R	[1..N]
{{PRT}}	Participation	X	[0..0]
{{NTE}}	Notes and comments	X	[0..0]
}}	--- OBSERVATION end		
{{FT1}}	Financial Transaction	X	[0..0]
{{CTI}}	Clinical Trial Identification	X	[0..0]
{{	--- SPECIMEN begin		
SPM	Specimen	X	[0..0]
{{OBX}}	Observation related to Specimen	X	[0..0]
}}	--- SPECIMEN end		
}	--- ORDER_OBSERVATION end		
}	--- PATIENT_RESULT end		
[DSC]	Continuation Pointer	X	[0..0]

223

224

4.2.2 Example 1 – Minimal Message

225

This message is a periodic report of a dialysis treatment status. It contains the minimum number of reported items.

226

227

Note, this standard reports Boolean event values in the PCD-01 transaction. This is done so that a complete record of the machine state is captured in the PCD-01. The state of the events is not considered clinically actionable. Clinically actionable events are communicated in the PCD-04 transaction. In the message below,

228

229

these events are highlighted to bring them to the readers attention and the foster a discussion.

```

230
231 MSH|^~\&|Dialysis
232 Machine|||20191003092006+0000||ORU^R01^ORU_R01|20191003092005|P|2.6||AL|NE||||IHE_PCD_001^
233 IHEPCD^1.3.6.1.4.12559.11.1.1.129^ISO

234 PID|||Scrubber 2000/SC678932^^^"^U||^"^^^"^^U

235 OBR|1|||0^MDC_DEV_DIALYSIS^MDC|||20191003092005+0000

236 OBX|1|ST|0^MDC_ATTR_ID_MODEL^MDC|1.1.1.1|NxStage System One|||||F
237 OBX|2|ST|7682^MDC_ID_MODEL_MANUFACTURER^MDC|1.1.1.2|NxStage|||||F
238 OBX|3|ST|7681^MDC_ID_MODEL_NUMBER^MDC|1.1.1.3|System One|||||F
239 OBX|4|ST|7684^MDC_ID_PROD_SPEC_SERIAL^MDC|1.1.1.4|1000478|||||F
240 OBX|5|ST|7667^MDC_ID_PROD_SPEC_SW^MDC|1.1.1.5|1.2.3.4|||||F
241 OBX|6|DTM|0^MDCX_HDIALY_MACH_TIME^MDC|1.1.1.1|20191003092005+0000|||||F
242 OBX|7|ST|0^MDCX_HDIALY_MACH_MODE_OF_OPERATION^MDC|1.1.1.3|Treatment|||||F
243 OBX|8|ST|0^MDCX_HDIALY_MACH_BLD_PUMP_ON^MDC|1.1.1.6|T|||||F
244 OBX|9|ST|0^MDCX_HDIALY_MACH_TX_FLUID_BYPASS^MDC|1.1.1.7|F|||||F
245 OBX|10|ST|0^MDCX_HDIALY_MACH_TX_MODALITY^MDC|1.1.1.8|HDF|||||F
246 OBX|11|NM|0^MDCX_HDIALY_MACH_THERAPY_TIME^MDC|1.1.1.9|180|min^minutes^UCUM|||||F
247 OBX|12|NM|0^MDCX_HDIALY_MACH_TIME_REMAIN^MDC|1.1.1.10|600|min^minutes^UCUM|||||F
248 OBX|13|NM|0^MDCX_HDIALY_ANTICOAG_PUMP_NUMBER^MDC|1.1.2.7|1|||||F
249 OBX|14|NM|0^MDCX_HDIALY_BLD_PRESS_ART^MDC|1.1.3.2|-75|mm[Hg]^Millimeters of Mercury^UCUM|< -
250 200|||||F

251 OBX|15|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_AVG_SETTING^MDC|1.1.3.4|250|ml/min^milliliter
252 per minute^UCUM|||||F

253 OBX|16|ST|0^MDCX_HDIALY_BLD_PUMP_MODE^MDC|1.1.3.5|Double_Needle|||||F

254 OBX|17|ST|0^MDCX_EVT_HDIALY_BLD_PUMP_STOP^MDC|1.1.3.6|F|||||F

255 OBX|18|NM|0^MDCX_HDIALY_BLD_PUMP_PRESS_VEN^MDC|1.1.3.15|200|mm[Hg]^Millimeters of
256 Mercury^UCUM|20-400|||||F

257 OBX|19|NM|0^MDCX_HDIALY_DIALYSATE_COND^MDC|1.1.4.3|13.81|mS/cm^millsiemens per
258 centimeter^UCUM|||||F

259 OBX|20|ST|0^MDCX_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|F|||||F

260 OBX|21|NM|0^MDCX_HDIALY_FILTER_TRANSMEMBRANE_PRESS^MDC|1.1.5.2|35|mm[Hg]^Millimeters of
261 Mercury^UCUM|||||F

262 OBX|22|ST|0^MDCX_EVT_HDIALY_SAFETY_SYSTEM_GENERAL^MDC|1.1.7.4|F|||||F

263 OBX|23|ST|0^MDCX_EVT_HDIALY_SAFETY_SELF_TEST^MDC|1.1.7.5|F|||||F

264 OBX|24|ST|0^MDCX_EVT_HDIALY_SAFETY_VEN_AIR_DETECT^MDC|1.1.7.7|F|||||F

265 OBX|25|ST|0^MDCX_HDIALY_THERAPY_COMPLETE_METHOD^MDC|1.1.8.18|UF Removed|||||F
266 OBX|26|NM|0^MDCX_HDIALY_UF_TARGET_VOL_TO_REMOVE^MDC|1.1.9.4|2000|ml^milliliter^UCUM|||||F
267 OBX|27|NM|0^MDCX_HDIALY_UF_ACTUAL_REMOVED_VOL^MDC|1.1.9.5|555|ml^milliliter^UCUM|||||F
268 OBX|28|ST|0^MDCX_HDIALY_UF_MODE^MDC|1.1.9.8|CONSTANT|||||F
269 OBX|29|NM|0^MDCX_HDIALY_UF_RATE^MDC|1.1.9.9|100|ml/h^milliliter per hour^UCUM|||||F
270 OBX|30|ST|0^MDCX_EVT_HDIALY_UF_RATE_RANGE^MDC|1.1.9.10|F|||||F
271

```

4.2.3 Example 2 – Full Message

This message reports all mandatory and optional fields for a pre/post dilution HDF therapy.

```

273
274 MSH|^~\&|Dialysis
275 Machine|||20191003092025+0000||ORU^R01^ORU_R01|20191003092024|P|2.6||AL|NE||||IHE_PCD_001^
276 IHEPCD^1.3.6.1.4.12559.11.1.1.129^ISO

277 PID|||Scrubber 2000/SC678932^^^"^U||^"^^^"^^U

```

278 OBR|1||0^MDC_DEV_DIALYSIS^MDC||20191003092024+0000

279 OBX|1|ST|0^MDC_ATTR_ID_MODEL^MDC|1.1.1.1|NxStage System One|||||F

280 OBX|2|ST|7682^MDC_ID_MODEL_MANUFACTURER^MDC|1.1.1.2|NxStage|||||F

281 OBX|3|ST|7681^MDC_ID_MODEL_NUMBER^MDC|1.1.1.3|System One|||||F

282 OBX|4|ST|7684^MDC_ID_PROD_SPEC_SERIAL^MDC|1.1.1.4|1000478|||||F

283 OBX|5|ST|7667^MDC_ID_PROD_SPEC_SW^MDC|1.1.1.5|1.2.3.4|||||F

284 OBX|6|ST|0^MDC_ATTR_ID_UDI^MDC|1.1.1.6|+M535NX10003A0/\$\$+735241/16D20180305J|||||F

285 OBX|7|DTM|0^MDCX_HDIALY_MACH_TIME^MDC|1.1.1.1|20191003092024+0000|||||F

286 OBX|8|ST|0^MDCX_HDIALY_MACH_MODE_DESCRIPTION^MDC|1.1.1.2|NxStage System One|||||F

287 OBX|9|ST|0^MDCX_HDIALY_MACH_MODE_OF_OPERATION^MDC|1.1.1.3|Treatment|||||F

288 OBX|10|NM|0^MDCX_HDIALY_MACH_PM_TIME_REMAIN^MDC|1.1.1.4|12|h^hours^UCUM|||||F

289 OBX|11|NM|0^MDCX_HDIALY_MACH_MAINT_TX_REMAIN^MDC|1.1.1.5|45|||||F

290 OBX|12|ST|0^MDCX_HDIALY_MACH_BLD_PUMP_ON^MDC|1.1.1.6|T|||||F

291 OBX|13|ST|0^MDCX_HDIALY_MACH_TX_FLUID_BYPASS^MDC|1.1.1.7|F|||||F

292 OBX|14|ST|0^MDCX_HDIALY_MACH_TX_MODALITY^MDC|1.1.1.8|HDF|||||F

293 OBX|15|NM|0^MDCX_HDIALY_MACH_THERAPY_TIME^MDC|1.1.1.9|180|min^minutes^UCUM|||||F

294 OBX|16|NM|0^MDCX_HDIALY_MACH_TIME_REMAIN^MDC|1.1.1.10|600|min^minutes^UCUM|||||F

295 OBX|17|ST|0^MDCX_EVT_HDIALY_ANTICOAG_STOP^MDC|1.1.2.4|F|||||F

296 OBX|18|ST|0^MDCX_HDIALY_ANTICOAG_TYPE^MDC|1.1.2.5|heparin|||||F

297 OBX|19|NM|0^MDCX_HDIALY_ANTICOAG_PUMP_NUMBER^MDC|1.1.2.7|1|||||F

298 OBX|20|ST|0^MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_EMPTY^MDC|1.1.2.8|F|||||F

299 OBX|21|ST|0^MDCX_HDIALY_ANTICOAG_SYRINGE_BRAND^MDC|1.1.2.9|Fishman|||||F

300 OBX|22|NM|0^MDCX_HDIALY_ANTICOAG_SYRINGE_VOL^MDC|1.1.2.10|60|ml^milliliter^UCUM|||||F

301 OBX|23|ST|0^MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_SIZE^MDC|1.1.2.11|F|||||F

302 OBX|24|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE^MDC|1.1.3.1|250|ml/min^milliliter per
303 minute^UCUM|||||F

304 OBX|25|NM|0^MDCX_HDIALY_BLD_PRESS_ART^MDC|1.1.3.2|-75|mm[Hg]^Millimeters of Mercury^UCUM|< -
305 200|||||F

306 OBX|26|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_AVG^MDC|1.1.3.3|250|ml/min^milliliter per
307 minute^UCUM|||||F

308 OBX|27|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_AVG_SETTING^MDC|1.1.3.4|250|ml/min^milliliter
309 per minute^UCUM|||||F

310 OBX|28|ST|0^MDCX_HDIALY_BLD_PUMP_MODE^MDC|1.1.3.5|Double Needle|||||F

311 OBX|29|ST|0^MDCX_EVT_HDIALY_BLD_PUMP_STOP^MDC|1.1.3.6|F|||||F

312 OBX|30|ST|0^MDCX_HDIALY_BLD_PUMP_TUBING_SIZE^MDC|1.1.3.7|8 mm|||||F

313 OBX|31|NM|0^MDCX_HDIALY_BLD_PUMP_PAT_BODY_TEMP_BY_BLOOD^MDC|1.1.3.8|39.1|Cel^degree
314 Celsius^UCUM|||||F

315 OBX|32|NM|0^MDCX_HDIALY_BLD_PUMP_CHANGE_IN_ENERGY^MDC|1.1.3.9|1|kJ/h^Kilojoule per
316 hour^UCUM|||||F

317 OBX|33|NM|0^MDCX_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP^MDC|1.1.3.10|100|mm[Hg]^Millimeters of
318 Mercury^UCUM|||||F

319 OBX|34|NM|0^MDCX_HDIALY_BLD_PUMP_PRIMING_VOL^MDC|1.1.3.11|191|ml^milliliter^UCUM|||||F

320 OBX|35|NM|0^MDCX_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS^MDC|1.1.3.12|200|mm[Hg]^Millimeters of
321 Mercury^UCUM|||||F

322 OBX|36|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_PROCESSED_TOTAL^MDC|1.1.3.14|120|L^liter^UCUM|||||F

323 OBX|37|NM|0^MDCX_HDIALY_BLD_PUMP_PRESS_VEN^MDC|1.1.3.15|200|mm[Hg]^Millimeters of
324 Mercury^UCUM|20-400|||||F

325 OBX|38|NM|0^MDCX_HDIALY_BLD_PUMP_BLOOD_TEMP_VEN^MDC|1.1.3.16|39.1|Cel^degree Celsius^UCUM|||||F

326 OBX|39|NM|0^MDCX_HDIALY_BICARB_COND^MDC|1.1.4.1|13.81|mS/cm^millisiemens per
327 centimeter^UCUM|||||F

328 OBX|40|NM|0^MDCX_HDIALY_BICARB_COND_SETTING^MDC|1.1.4.2|13.81|mS/cm^millisiemens per
329 centimeter^UCUM||||F

330 OBX|41|NM|0^MDCX_HDIALY_DIALYSATE_COND^MDC|1.1.4.3|13.81|mS/cm^millisiemens per
331 centimeter^UCUM||||F

332 OBX|42|NM|0^MDCX_HDIALY_DIALYSATE_COND_SETTING^MDC|1.1.4.4|13.81|mS/cm^millisiemens per
333 centimeter^UCUM||||F

334 OBX|43|NM|0^MDCX_HDIALY_DIALYSATE_FLOW_RATE^MDC|1.1.4.5|99|ml/min^milliliter per
335 minute^UCUM||||F

336 OBX|44|NM|0^MDCX_HDIALY_DIALYSATE_FLOW_RATE_SETTING^MDC|1.1.4.6|100|ml/min^milliliter per
337 minute^UCUM||||F

338 OBX|45|ST|0^MDCX_HDIALY_DIALYSATE_FLOW_MODE^MDC|1.1.4.7|CONSTANT|||||F

339 OBX|46|NM|0^MDCX_HDIALY_DIALYSATE_AMMONIA^MDC|1.1.4.10|0.00|[ppm]^parts per million^UCUM|||||F

340 OBX|47|NM|0^MDCX_HDIALY_DIALYSATE_FLOW_RATE_AVG^MDC|1.1.4.11|100|ml/min^milliliter per
341 minute^UCUM||||F

342 OBX|49|NM|0^MDCX_HDIALY_CONC_HCO3_SETTING^MDC|1.1.4.13|32.0|mmol/L^MilliMolesPerLiter^UCUM|||||F

343 OBX|50|ST|0^MDCX_HDIALY_CONC_HCO3_MODE^MDC|1.1.4.14|PROFILE|||||F

344 OBX|51|ST|0^MDCX_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|F|||||F

345 OBX|52|ST|0^MDCX_HDIALY_DIALYSATE_NAME^MDC|1.1.4.16|FMC smartbag 111.5||Bbraun Duosol 35||||F

346 OBX|53|NM|0^MDCX_HDIALY_DIALYSATE_CONC_ACETATE^MDC|1.1.4.17|3.00|mmol/L^MilliMolesPerLiter^UCUM|
347 ||||F

348 OBX|54|NM|0^MDCX_HDIALY_DIALYSATE_CONC_HCO3^MDC|1.1.4.18|35.0|mmol/L^MilliMolesPerLiter^UCUM|35
349 mEq/L||||F

350 OBX|55|NM|0^MDCX_HDIALY_DIALYSATE_CONC_CHLORIDE^MDC|1.1.4.19|108.00|mmol/L^MilliMolesPerLiter^UCUM|109
351 mEq/L||||F|||20191003085024+0000

352 OBX|56|NM|0^MDCX_HDIALY_DIALYSATE_CONC_MG^MDC|1.1.4.20|0.5|mmol/L^MilliMolesPerLiter^UCUM|1.0
353 mEq/L||||F|||20191003085024+0000

354 OBX|57|NM|0^MDCX_HDIALY_DIALYSATE_CONC_CA^MDC|1.1.4.21|1.50|mmol/L^MilliMolesPerLiter^UCUM|3
355 mEq/L||||F|||20191003085024+0000

356 OBX|58|NM|0^MDCX_HDIALY_DIALYSATE_CONC_CITRATE^MDC|1.1.4.22|0.0|mmol/L^MilliMolesPerLiter^UCUM|
357 |||F

358 OBX|59|NM|0^MDCX_HDIALY_DIALYSATE_CONC_GLU^MDC|1.1.4.23|1.00|mmol/L^MilliMolesPerLiter^UCUM|1
359 g/L||||F

360 OBX|60|NM|0^MDCX_HDIALY_DIALYSATE_CONC_K^MDC|1.1.4.24|1.00|mmol/L^MilliMolesPerLiter^UCUM|0
361 mEq/L||||F|||20191003085024+0000

362 OBX|61|NM|0^MDCX_HDIALY_DIALYSATE_CONC_NA^MDC|1.1.4.25|138|mmol/L^MilliMolesPerLiter^UCUM|140
363 mEq/L||||F|||20191003085024+0000

364 OBX|62|NM|0^MDCX_HDIALY_DIALYSATE_CONC_NA_SETTING^MDC|1.1.4.26|140|mmol/L^MilliMolesPerLiter^UCU
365 M||||F

366 OBX|63|ST|0^MDCX_HDIALY_DIALYSATE_CONC_NA_MODE^MDC|1.1.4.27|CONSTANT|||||F

367 OBX|64|NM|0^MDCX_HDIALY_DIALYSATE_CONC_PH^MDC|1.1.4.28|7.0|[pH]^pH^UCUM|||||F

368 OBX|65|NM|0^MDCX_HDIALY_DIALYSATE_VOL_DELIV^MDC|1.1.4.29|24.34|L^liter^UCUM|||||F

369 OBX|66|ST|0^MDCX_HDIALY_FILTER_NAME^MDC|1.1.5.1|NxStage CAR125|||||F

370 OBX|67|NM|0^MDCX_HDIALY_FILTER_TRANSMEMBRANE_PRESS^MDC|1.1.5.2|35|mm[Hg]^Millimeters of
371 Mercury^UCUM|||||F

372 OBX|68|ST|0^MDCX_HDIALY_RF_DILUTION_LOCATION^MDC|1.1.6.1|PREF-POSTF|||||F

373 OBX|69|ST|0^MDCX_HDIALY_RF_POST_FILTER_NAME^MDC|1.1.6.2|multiBIC||| |||F

374 OBX|70|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_ACETATE^MDC|1.1.6.3|0|mmol/L^MilliMolesPerLiter^UCUM
375 |||||F

376 OBX|71|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_HCO3^MDC|1.1.6.4|35.0|mmol/L^MilliMolesPerLiter^UCUM
377 |35 mEq/L||||F

378 OBX|72|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_CHLORIDE^MDC|1.1.6.5|111.00|mmol/L^MilliMolesPerLite
379 r^UCUM|109 mEq/L||||F|||20191003085024+0000

380 OBX|73|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_MG^MDC|1.1.6.6|0.5|mmol/L^MilliMolesPerLiter^UCUM|1.
381 0 mEq/L||||F|||20191003085024+0000

382 OBX|74|NM|0^MDCX_HDIALY_DIALYSATE_CONC_CA^MDC|1.1.6.7|1.50|mmol/L^MilliMolesPerLiter^UCUM|3
383 mEq/L||||F|||20191003085024+0000

384 OBX|75|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_CITRATE^MDC|1.1.6.8|0.0|mmol/L^MilliMolesPerLiter^UC
385 UM||||F

386 OBX|76|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_GLU^MDC|1.1.6.9|5.55|mmol/L^MilliMolesPerLiter^UCUM|
387 1 g/L||||F

388 OBX|77|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_K^MDC|1.1.6.10|2.00|mmol/L^MilliMolesPerLiter^UCUM|0
389 mEq/L||||F|||20191003085024+0000

390 OBX|78|NM|0^MDCX_HDIALY_RF_POST_FILTER_CONC_NA^MDC|1.1.6.11|140|mmol/L^MilliMolesPerLiter^UCUM|1
391 40 mEq/L||||F|||20191003085024+0000

392 OBX|79|NM|0^MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE^MDC|1.1.6.12|60|ml/min^milliliter per
393 minute^UCUM||||F

394 OBX|80|NM|0^MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE_SETTING^MDC|1.1.6.13|60|ml/min^milliliter per
395 minute^UCUM||||F

396 OBX|81|NM|0^MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE_AVG^MDC|1.1.6.14|60|ml/min^milliliter per
397 minute^UCUM||||F

398 OBX|82|NM|0^MDCX_HDIALY_RF_POST_FILTER_TEMP^MDC|1.1.6.15|39.9|Cel^degree Celsius^UCUM||||F

399 OBX|83|NM|0^MDCX_HDIALY_RF_POST_FILTER_TEMP_SETTING^MDC|1.1.6.16|41.0|Cel^degree
400 Celsius^UCUM||||F

401 OBX|84|NM|0^MDCX_HDIALY_RF_POST_FILTER_VOL^MDC|1.1.6.17|6.00|L^liter^UCUM||||F

402 OBX|85|NM|0^MDCX_HDIALY_RF_POST_FILTER_VOL_SETTING^MDC|1.1.6.18|12.00|L^liter^UCUM||||F

403 OBX|86|ST|0^MDCX_HDIALY_RF_PRE_FILTER_NAME^MDC|1.1.6.19|multiBIC|||F

404 OBX|87|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_ACETATE^MDC|1.1.6.20|0|mmol/L^MilliMolesPerLiter^UCUM
405 ||||F

406 OBX|88|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_HCO3^MDC|1.1.6.21|35.0|mmol/L^MilliMolesPerLiter^UCUM
407 |35 mEq/L||||F

408 OBX|89|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_CHLORIDE^MDC|1.1.6.22|111.00|mmol/L^MilliMolesPerLite
409 r^UCUM|109 mEq/L||||F|||20191003085024+0000

410 OBX|90|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_MG^MDC|1.1.6.23|0.5|mmol/L^MilliMolesPerLiter^UCUM|1.
411 0 mEq/L||||F|||20191003085024+0000

412 OBX|91|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_CA^MDC|1.1.6.24|1.50|mmol/L^MilliMolesPerLiter^UCUM|3
413 mEq/L||||F|||20191003085024+0000

414 OBX|92|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_CITRATE^MDC|1.1.6.25|0.0|mmol/L^MilliMolesPerLiter^UC
415 UM||||F

416 OBX|93|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_GLU^MDC|1.1.6.26|5.55|mmol/L^MilliMolesPerLiter^UCUM|
417 1 g/L||||F

418 OBX|94|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_K^MDC|1.1.6.27|2.00|mmol/L^MilliMolesPerLiter^UCUM|0
419 mEq/L||||F|||20191003085024+0000

420 OBX|95|NM|0^MDCX_HDIALY_RF_PRE_FILTER_CONC_NA^MDC|1.1.6.28|140|mmol/L^MilliMolesPerLiter^UCUM|14
421 0 mEq/L||||F|||20191003085024+0000

422 OBX|96|NM|0^MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE^MDC|1.1.6.29|60|ml/min^milliliter per
423 minute^UCUM||||F

424 OBX|97|NM|0^MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE_SETTING^MDC|1.1.6.30|60|ml/min^milliliter per
425 minute^UCUM||||F

426 OBX|98|NM|0^MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE_AVG^MDC|1.1.6.31|60|ml/min^milliliter per
427 minute^UCUM||||F

428 OBX|99|NM|0^MDCX_HDIALY_RF_PRE_FILTER_TEMP^MDC|1.1.6.32|39.9|Cel^degree Celsius^UCUM||||F

429 OBX|100|NM|0^MDCX_HDIALY_RF_PRE_FILTER_TEMP_SETTING^MDC|1.1.6.33|41.0|Cel^degree
430 Celsius^UCUM||||F

431 OBX|101|NM|0^MDCX_HDIALY_RF_PRE_FILTER_VOL^MDC|1.1.6.34|1.23|L^liter^UCUM||||F

432 OBX|102|NM|0^MDCX_HDIALY_RF_PRE_FILTER_VOL_SETTING^MDC|1.1.6.35|3.00|L^liter^UCUM||||F

433 OBX|103|NM|0^MDCX_HDIALY_RF_CONV_CLEARANCE^MDC|1.1.6.36|31.6|ml/min^milliliter per
434 minute^UCUM||||F

435 OBX|104|ST|0^MDCX_EVT_HDIALY_RF_EXCESS_DELIV^MDC|1.1.6.37|F||||F

436 OBX|105|ST|0^MDCX_EVT_HDIALY_RF_INSUFF_DELIV^MDC|1.1.6.38|F||||F

437 OBX|106|NM|0^MDCX_HDIALY_RF_PRE_POST_FLOW_RATIO_SETTING^MDC|1.1.6.39|3.00| |||||F

438 OBX|107|NM|0^MDCX_HDIALY_RF_BOLUS_RATE^MDC|1.1.6.40|0|ml/min^milliliter per minute^UCUM| |||||F

439 OBX|108|NM|0^MDCX_HDIALY_RF_BOLUS_VOL^MDC|1.1.6.41|0.00|ml^milliliter^UCUM| |||||F

440 OBX|109|NM|0^MDCX_HDIALY_RF_BOLUS_VOL_SETTING^MDC|1.1.6.42|3.00|ml^milliliter^UCUM| |||||F

441 OBX|110|ST|0^MDCX_HDIALY_RF_FLOW_MODE^MDC|1.1.6.43|CONSTANT| |||||F

442 OBX|111|ST|0^MDCX_HDIALY_RF_SOURCE^MDC|1.1.6.44|BAG| |||||F

443 OBX|112|ST|0^MDCX_EVT_HDIALY_SAFETY_DIALYSATE_AIR_DETECT^MDC|1.1.7.2|F| |||||F

444 OBX|113|ST|0^MDCX_EVT_HDIALY_SAFETY_DIALYSATE_COMPOSITION^MDC|1.1.7.3|F| |||||F

445 OBX|114|ST|0^MDCX_EVT_HDIALY_SAFETY_SYSTEM_GENERAL^MDC|1.1.7.4|F| |||||F

446 OBX|115|ST|0^MDCX_EVT_HDIALY_SAFETY_SELF_TEST^MDC|1.1.7.5|F| |||||F

447 OBX|116|ST|0^MDCX_EVT_HDIALY_SAFETY_VEN_ACCESS^MDC|1.1.7.6|F| |||||F

448 OBX|117|ST|0^MDCX_EVT_HDIALY_SAFETY_VEN_AIR_DETECT^MDC|1.1.7.7|F| |||||F

449 OBX|118|ST|0^MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ALERT^MDC|1.1.7.8|F| |||||F

450 OBX|119|ST|0^MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ERROR^MDC|1.1.7.9|F| |||||F

451 OBX|120|NM|0^MDCX_HDIALY_THERAPY_MASS_TRF_AREA_COEFF^MDC|1.1.8.1|700| |||||F

452 OBX|121|NM|0^MDCX_HDIALY_THERAPY_ACCESS_FLOW^MDC|1.1.8.2|250|ml/min^milliliter per
453 minute^UCUM| |||||F

454 OBX|122|NM|0^MDCX_HDIALY_THERAPY_BLOOD_VOL^MDC|1.1.8.3|20.0|{%vol}^VolumePercent^UCUM| |||||F

455 OBX|123|NM|0^MDCX_HDIALY_THERAPY_RATIO_EKT_OVER_V_DELIVERED^MDC|1.1.8.4|1.1|%^Percent^UCUM
456 |||||F

457 OBX|124|NM|0^MDCX_HDIALY_THERAPY_KT_DELIVERED^MDC|1.1.8.5|42.0|L^liter^UCUM| |||||F

458 OBX|125|NM|0^MDCX_HDIALY_THERAPY_RATIO_SPKT_OVER_V_DELIVERED^MDC|1.1.8.6|1.1|%^Percent^UCUM
459 |||||F

460 OBX|126|ST|0^MDCX_EVT_HDIALY_THERAPY_PAT_TX_GENERAL^MDC|1.1.8.7|F| |||||F

461 OBX|127|NM|0^MDCX_HDIALY_THERAPY_RATIO_KT_OVER_V_GOAL^MDC|1.1.8.10|1.21|%^Percent^UCUM | |||||F

462 OBX|128|NM|0^MDCX_HDIALY_THERAPY_MEAN_UREA_CLEARANCE^MDC|1.1.8.11|70| |||||F

463 OBX|129|NM|0^MDCX_HDIALY_THERAPY_BODY_START_WT^MDC|1.1.8.12|75.9|kg^kilograms^UCUM| |||||F

464 OBX|130|NM|0^MDCX_HDIALY_THERAPY_PCT_RECIRC^MDC|1.1.8.13|20|%^Percent^UCUM| |||||F

465 OBX|131|NM|0^MDCX_HDIALY_THERAPY_PLASMA_NA_CONC^MDC|1.1.8.14|140|mmol/L^MilliMolesPerLiter^UCUM|
466 |||||F

467 OBX|132|NM|0^MDCX_HDIALY_THERAPY_RATIO_SPKT_OVER_V_PROJECTED^MDC|1.1.8.16|1.1|%^Percent^UCUM
468 |||||F

469 OBX|133|ST|0^MDCX_HDIALY_THERAPY_COMPLETE_METHOD^MDC|1.1.8.18|UF Removed| |||||F

470 OBX|134|ST|0^MDCX_EVT_HDIALY_THERAPY_TX_END_TIME^MDC|1.1.8.19|F| |||||F

471 OBX|135|NM|0^MDCX_HDIALY_THERAPY_UREA_CLEARANCE^MDC|1.1.8.20|196| |||||F

472 OBX|136|NM|0^MDCX_HDIALY_THERAPY_BODY_END_WT^MDC|1.1.8.21|75.9|kg^kilograms^UCUM| |||||F

473 OBX|137|ST|0^MDCX_EVT_HDIALY_UF_LO^MDC|1.1.9.1|F| |||||F

474 OBX|138|ST|0^MDCX_EVT_HDIALY_UF_NEG^MDC|1.1.9.2|F| |||||F

475 OBX|139|NM|0^MDCX_HDIALY_UF_TIME_TO_TARGET^MDC|1.1.9.3|55|min^minutes^UCUM| |||||F

476 OBX|140|NM|0^MDCX_HDIALY_UF_TARGET_VOL_TO_REMOVE^MDC|1.1.9.4|2000|ml^milliliter^UCUM| |||||F

477 OBX|141|NM|0^MDCX_HDIALY_UF_ACTUAL_REMOVED_VOL^MDC|1.1.9.5|555|ml^milliliter^UCUM| |||||F

478 OBX|142|NM|0^MDCX_HDIALY_UF_RATE_MAX_NON_ALM_SETTING^MDC|1.1.9.6|100|ml/h^milliliter per
479 hour^UCUM| |||||F

480 OBX|143|ST|0^MDCX_EVT_HDIALY_UF_GOAL_MET^MDC|1.1.9.7|F| |||||F

481 OBX|144|ST|0^MDCX_HDIALY_UF_MODE^MDC|1.1.9.8|CONSTANT| |||||F

482 OBX|145|NM|0^MDCX_HDIALY_UF_RATE^MDC|1.1.9.9|100|ml/h^milliliter per hour^UCUM| |||||F

483 OBX|146|ST|0^MDCX_EVT_HDIALY_UF_RATE_RANGE^MDC|1.1.9.10|F| |||||F

484 OBX|147|NM|0^MDC_ATTR_TIME_PD_MSMT^MDC|1.2.1.1|30|min^minutes^UCUM| |||||F|||20191003085024+0000

485 OBX|148|ST|0^MDCX_ATTR_BODY_POSN^MDC|1.2.1.2|SUPINE|||||F|||20191003085024+0000
 486 OBX|149|NM|150022^MDC_PRESS_BLD_NONINV_DIA^MDC|1.2.1.3|80|mm[Hg]^Millimeters of
 487 Mercury^UCUM|||||F|||20191003085024+0000
 488 OBX|150|NM|149546^MDC_PULS_RATE_NON_INV^MDC|1.2.1.4|70|{Beats}/min^Beats Per
 489 Minuts^UCUM|||||F|||20191003085024+0000
 490 OBX|151|NM|150023^MDC_PRESS_BLD_NONINV_MEAN^MDC|1.2.1.5|100|mm[Hg]^Millimeters of
 491 Mercury^UCUM|||||F|||20191003085024+0000
 492 OBX|152|NM|150021^MDC_PRESS_BLD_NONINV_SYS^MDC|1.2.1.6|120|mm[Hg]^Millimeters of
 493 Mercury^UCUM|||||F|||20191003085024+0000
 494 OBX|153|NM|150456^MDC_PULS_OXIM_SAT_O2^MDC|1.3.1.1|98|%^Percent^UCUM|||||F
 495 OBX|154|NM|0^MDC_PULS_OXIM_PULS_RATE^MDC|1.3.1.2|67|{Beats}/min^Beats Per Minuts^UCUM|||||F
 496 OBX|155|ST|0^MDCX_EVT_PULS_OXIM_ERROR^MDC|1.3.1.3|F|||||F
 497 OBX|156|ST|0^MDCX_HDIALY_PLASMA_VOL_MARKER^MDC|1.4.1.1|NONE|||||F
 498 OBX|157|ST|0^MDCX_HDIALY_PLASMA_VOL_PROFILE^MDC|1.4.1.2|A|||||F
 499 OBX|158|NM|0^MDC_CONC_HCT^MDC|1.4.1.3|0.45|{%vol}^VolumePercent^UCUM|||||F
 500 OBX|159|NM|0^MDC_CONC_HB^MDC|1.4.1.4|13.6|g/dL^GramsPerDeciLiter^UCUM|||||F
 501 OBX|160|NM|0^MDCX_HDIALY_REL_BLOOD_VOL^MDC|1.4.1.5|35|%^Percent^UCUM|||||F
 502 OBX|161|NM|150456^MDC_SAT_O2^MDC|1.4.1.6|98|%^Percent^UCUM|||||F
 503

4.3 EMR Response

The EMR responds with a simple Acknowledgement Message (ACK_R01_ACK).

4.3.1 Message Structure

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
[(SFT)]	Software Segment	X	[0..0]
[UAC]	User Authentication Credential	X	[0..0]
MSA	Message Acknowledgement	R	[1..1]
ERR	Error	RE	[0..100]

4.3.2 Example 1 – Accepted Response

The following message is an example of an accepted PCD-01 message.

```
MSH|^~\&|EMR|||20191003092025+0000||ACK^R01^ACK|XX3657|P|2.6
MSA|AA|20191003092024
```

4.3.3 Example 2 – Rejected Response

The following message is an example of a rejected PCD-01 message.


```
MSH|^~\&|EMR|||20191003092025+0000||ACK^R01^ACK|XX3657|P|2.6
MSA|AR|20191003092024
ERR||PID^1^11^^9|103|E<cr>
```


5 Reporting Alarm Information

5.1 Overview

The sending of Alarm Information to the EMR System is accomplished using the PCD-04 Transaction of the Patient Care Device (PCD) domain in the Device Enterprise Communication (DEC) IHE profile. In doing so, the Dialysis Machine acts as the Alert Reporter (AR) and the EMR is the Alert Manager (AM) and/or the Alert Consumer (ACON).

Figure 2 - AR sends Report Alert to AM and/or ACON as an HL7 ORU message.



The Dialysis Machine will report the alarm when it starts, clears, and periodically while it is active. The dialysis machine will send the appropriate parameters and measurements with each alarm. For example, when an Arterial Pressure High Alarm is reported, the message will contain the Arterial pressure Measurement and Arterial Pressure Upper Limit.

The ORU^R40 message represents a unitary alert, which is to be acknowledged as a whole by an ORA message. Multiple alerts requiring separate acknowledgement will be sent as individual messages.

The ORDER_OBSERVATION Segment Group which has OBR-49 value A (Alert provider when abnormal) conveys the alert observation(s). One or more OBX segments in this Segment Group will typically have OBX-8 Interpretation Codes value of LL, HH, or AA. At least one OBR segment shall have OBR-49 value A. Other ORDER_OBSERVATION Segment Groups within the message shall be considered supporting information for the alert observation(s).

5.2 Device Request

The PCD-04 Report Alert message is used to communicate ACM data from an Alert Reporter (AR) to Alert Manager (AM) and/or Alert Consumer (ACON). In addition, the Dialysis Machine sends a PCD-04 Message to the Device Observation Consumer (DOC) so that the information can be logged with the PCD-01 Message.

Common HL7 segments are defined in Appendix A – HL7 Data Elements. There are sections discussing considerations specific to PCD-04 where applicable.

While there can be multiple OBR segments per PCD-04 transaction (in support of inclusion of alert common containment and evidentiary data) there is at most one alert per PCD-04 transaction.

5.2.1 Message Structure

Message Type - ORU^R40^ORU_R40

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
[[SFT]]	Software Segment	X	[0..0]
{	--- ALERT_begin		
[--- PATIENT begin		
PID	Patient Identification	R	[1..1]
[--- LOCATION begin		
PV1	Alert Location	O	[0..1]
]	--- LOCATION end		
]	--- PATIENT end		
{	--- ALERT_IDENTIFICATION begin		
[ORC]	Alert Order Common	X	[0..0]
{OBR}	Alert Identification	R	[1..*]
[{	--- ALERT_OBSERVATION begin		
{OBX}	Alert observations relative to OBR	R	[1..*]
{[NTE]}	Notes and Comments	X	[0..0]
}}	--- ALERT OBSERVATION end		
}	--- ALERT_IDENTIFICATION end		
}	--- ALERT end		

A single Report Alert [PCD-04] transaction contains at most one alert for a given patient. The PCD-04 message contains a single OBR segment that contains five OBX segments.

- OBX 1 identifies the type of alarm
- OBX 2 identifies the event source. For numeric values this include the object value and limits.
- OBX 3 identifies the event phase
- OBX 4 identifies the alarm state
- OBX 5 identifies the alarm activity state

557

The first OBX (ID=1) identifies the type of alert.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

558

OBX-1 Set ID

559

This field contains the sequence number. For this segment it will always be 1.

560

OBX-2 Value Type

561

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be “CWE”.

562

OBX-3 Observation Identifier

563

This field will be one of the following.

Field Value	Usage
196670^MDC_EVT_LO^MDC	This value is used when the alarm is generated because a measure value is below a numeric limit.
196648^MDC_EVT_HI^MDC	This value is used when the alarm is generated because a measure value is above a numeric limit.
196616^MDC_EVT_ALARM^MDC	This value is used when the alarm is not generated due to a numeric value.

564

OBX-4 Observation Sub-ID

565

This value will always be “1.0.0.0.1”.

566

OBX-5 Observation Value

This field identifies the alarm being generate. For numeric values, this field is the value in the Source Column in Table 2. For non-numeric alarms, this field is the Event Identifier Column in Table 2.

For example, "0^MDCX_HDIALY_ALARM_PRESS_ART^MDC"

"0^ MDCX_EVT_HDIALY_ALARM_BLD_PUMP_STOP^MDC"

OBX-8 Interpretation Codes

This field is used to provide the Alarm Priority and Alarm Type, and optionally the Abnormality Type.

Examples,

High Priority, Technical, Upper Limit Alarm H~ST~PH

Non-Numeric High Priority, Technical Alarm ST~PH

OBX-11 Observation Result Status

For this segment the value will always be "F".

The second OBX (ID=2) identifies the parameter and the reference range.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

OBX-1 Set ID

For this segment the value will always be 2.

OBX-2 Value Type

This field defines the data type of OBX-5, Observation Value. For numeric alarms, this value will be "NM". For non-numeric alarms, this value will be "CWE".

OBX-3 Observation Identifier

For numeric alarms, this field identifies the object being monitored. This object is identified in the Source Column in Table 2.

For non-numeric alarms, this field will always be "68480^MDC_ATTR_ALERT_SOURCE^MDC".

OBX-4 Observation Sub-ID

For numeric alarms, this field identifies the object being measured. This value is the Observation Sub ID of the object identified in the Source Column in Table 2.

For non-numeric alarms, this field identifies the VMD of the device generating the alarm. This value is the Observation Sub-ID of the object identified in the Source Column in Table 2.

In both cases, the observation sub-ID has ".2" appended to signify facet 2.

For example, the Therapy Outcomes VMD would be reported as "1.1.8.2".

OBX-5 Observation Value

For numeric alarms, this field contains the value of the object that resulted in the alarm.

For non-numeric alarms, this field contains the identifier of the VMD generating the alarm. This is identified in the Source Column in Table 2.

Example, 0000000^MDCX_HDIALY_THERAPY_OUTCOMES_CHAN^MDC.

OBX-6 Observation Units

For numeric alarms, this field contains the units of measure for the object.

For non-numeric alarms, this field is blank.

OBX-7 Reference Range

This field contains the alarm limits for numeric alarms. For non-numeric alarms this field is empty. Ranges can be specified in one of three ways.

	Format	Example
Upper Limit Only	> Upper	> 400
Lower Limit Only	< Lower	< -200
Upper and Lower Limit	Lower – Upper	20 - 400

OBX-11 Observation Result Status

For this segment the value will always be "F".

The third OBX (ID=3) identifies the event phase.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

OBX-1 Set ID

This field contains the sequence number. For this segment it will always be 3.

OBX-2 Value Type

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "ST".

OBX-3 Observation Identifier

This field will always be "68481^MDC_ATTR_EVENT_PHASE^MDC".

OBX-4 Observation Sub-ID

This value will always be "1.0.0.0.3".

OBX-5 Observation Value

This value will always be one of the following.

- start start of an interval event/alert – an end is expected
- continue continuation of an ongoing interval event/alert
- end end of an interval event/alert

OBX-11 Observation Result Status

For this segment the value will always be "F".

631 The third OBX (ID=4) identifies the alarm start.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

632 **OBX-1 Set ID**

633 This field contains the sequence number. For this segment it will always be 4.

634 **OBX-2 Value Type**

635 This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "ST".

636 **OBX-3 Observation Identifier**

637 This field will always be "68482^MDC_ATTR_ALARM_STATE^MDC".

638 **OBX-4 Observation Sub-ID**

639 This value will always be "1.0.0.0.4".

640 **OBX-5 Observation Value**

641 This value will one of the following:

- 642 • off
- 643 • inactive
- 644 • active
- 645 • latched

OBX-11 Observation Result Status

For this segment the value will always be “F”.

The third OBX (ID=5) identifies the alarm activity state.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

OBX-1 Set ID

This field contains the sequence number. For this segment it will always be 5.

OBX-2 Value Type

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be “ST”.

OBX-3 Observation Identifier

This field will always be “68483^MDC_ATTR_ALARM_INACTIVIATION_STATE^MDC”.

OBX-4 Observation Sub-ID

This value will always be “1.0.0.0.5”.

OBX-5 Observation Value

This value will one of the following:

- enabled
- audio-paused

- audio-off
- alarm-paused
- alarm-off
- alert-acknowledged (IEC 60601-1-8 Amendment 1; CP126)

OBX-11 Observation Result Status

For this segment the value will always be “F”.

5.2.2 Example 1 - Alarm Initiation

The following message is an example of a Venous Pressure Low Alarm.

```
MSH|^~\&|Dialysis Machine|||20191003092025+0000||
ORU^R40^ORU_R40|20191003092024|P|2.6||AL|NE|||IHE_PCD_001^IHEPCD^1.3.6.1.4.12559.11.1.1.12
9^ISO
PID|||Scrubber 2000/SC678932^^^"U||^"^^^U
OBR|1|||0^MDC_DEV_DIALYSIS^MDC|||20191003092024+0000
OBX|1|CWE|196670^MDC_EVT_LO^MDC|1.0.0.0.1|0^MDCX_HDIALY_ALARM_PRESS_VEN^MDC|PH~SP|||F
OBX|2|NM|0^MDCX_HDIALY_BLD_PUMP_PRESS_VEN^MDC|1.1.3.15.2|15|mm[Hg]^Millimeters of Mercury^UCUM
|20-400|||F
OBX|3|ST|68481^MDC_ATTR_EVENT_PHASE^MDC|1.0.0.0.3|start|||F
OBX|4|ST|68482^MDC_ATTR_ALARM_STATE^MDC|1.0.0.0.4|active|||F
OBX|5|ST|68483^MDC_ATTR_ALARM_INACTIVIATION_STATE^MDC|1.0.0.0.5|enabled|||F
```

5.2.3 Example 2 - Alarm Cessation

The following message shows the cessation of a Blood Pump Stopped alarm.

```
MSH|^~\&|Dialysis Machine|||20191003092025+0000||
ORU^R40^ORU_R40|20191003092024|P|2.6||AL|NE|||IHE_PCD_001^IHEPCD^1.3.6.1.4.12559.11.1.1.12
9^ISO
PID|||Scrubber 2000/SC678932^^^"U||^"^^^U
OBR|1|||0^MDC_DEV_DIALYSIS^MDC|||20191003092024+0000
OBX|1|CWE|196616^MDC_EVT_ALARM^MDC|1.0.0.0.1|0^MDCX_EVT_HDIALY_ALARM_BLD_PUMP_STOP^MDC|||F
OBX|2|CWE|68480^MDC_ATTR_ALERT_SOURCE^MDC|1.1.3.0.2|||F
OBX|3|ST|68481^MDC_ATTR_EVENT_PHASE^MDC|1.0.0.0.3|end|||F
OBX|4|ST|68482^MDC_ATTR_ALARM_STATE^MDC|1.0.0.0.4|inactive|||F
OBX|5|ST|68483^MDC_ATTR_ALARM_INACTIVIATION_STATE^MDC|1.0.0.0.5|enabled|||F
```

5.3 EMR Response

The EMR responds with a Report Alert Acknowledgement (ORA^R41^ORA_R41).

5.3.1 Message Structure

Message Type - ORA^R41^ORA_R41

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]

[[SFT]]	Software Segment	X	[0..0]
[UAC]	User Authentication Credential	X	[0..0]
MSA	Message Acknowledgement	R	[1..1]
ERR	Error	RE	[0..100]
[[PRT]]	Participation (Acknowledging User)	O	[0..1]

5.3.1 Example 1 – Accepted Response

The following message is an example of an accepted PCD-04 message.

```
MSH|^~\&|EMR|||20191003092025+0000||ORA^R41^ORA_R41|XX3657|P|2.6
MSA|AA|20191003092024
```

6 HL7 Data Elements

The information contained in this section is a subset of the full HL7 Specification. Only segments used by this standard are shown. Restrictions that are specific to this specification are identified for each segment.

6.1 Segments

6.1.1 MSA – Message Acknowledgement Segment

The MSA segment contains INFORMATION sent while acknowledging another message. This message is sent from the EMR System to the Dialysis Machine. Although other fields may be present, the Dialysis Machine processes MSA-1 Acknowledgement Code and MSA-2 Message Control ID.

HL7 Attribute Table - MSA Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ID	R	0008	Acknowledgment Code
2	ST	R		Message Control ID
3		W		Text Message
4	NM	O		Expected Sequence Number
5		W		Delayed Acknowledgment Type
6		W		Error Condition
7	NM	O		Message Waiting Number
8	ID	O	0520	Message Waiting Priority

MSA-1 Acknowledgment Code (ID)

HL7 Definition: This field contains an acknowledgment code, see message processing rules. Refer to HL7 Table 0008 - Acknowledgment Code for valid values.

In this table, the value “CA” corresponds to “Enhanced mode: Accept acknowledgment: Commit Accept.” This value means that the receiving system has committed the message to safe storage in a manner that releases the sending system from the need to resend the message. “CR” means that a message was rejected for a message uniquely identified by the information in MSH-9, MSH-12, MSH-11, MSH-21. “CE” is a message processing error or failure condition.

This Specification: The Dialysis Machine accepts both original and enhanced mode responses. The response of the Dialysis Machine to an Error or Rejected message is device and manufacturer specific and is not covered by this guide.

MSA-2 Message Control ID (ST)

HL7 Definition: This field contains the message control ID of the message sent by the sending system. It allows the sending system to associate this response with the message for which it is intended.

This Specification: Note that the combination of the values in *MSA-2 Message Control ID* and *MSH-6 Receiving Facility* should be unique to the recipient of the acknowledgement message (i.e., the Order Fulfiller).

Note on Element Length: The element length for MSA-2 has been extended to 50 characters from the HL7-prescribed length of 20 characters. This element has been extended to accommodate the extended element length of MSH-10 Message Control ID in the ORU message types.

6.1.2 MSH – Message Header Segment

The MSH segment defines the intent, source, destination, and some specifics of the syntax of a message.

HL7 Attribute Table - MSH Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	R		Field Separator
2	ST	R		Encoding Characters
3	HD	X	0361	Sending Application
4	HD	X	0362	Sending Facility
5	HD	X	0361	Receiving Application
6	HD	X	0362	Receiving Facility
7	DTM	R		Date/Time of Message
8	ST	X		Security
9	MSG	R		Message Type
10	ST	R		Message Control ID
11	PT	R		Processing ID
12	VID	R		Version ID
13	NM	X		Sequence Number
14	ST	X		Continuation Pointer
15	ID	X	0155	Accept Acknowledgment Type
16	ID	X	0155	Application Acknowledgment Type
17	ID	X	0399	Country Code
18	ID	X	0211	Character Set
19	CWE	X		Principal Language Of Message
20	ID	X	0356	Alternate Character Set Handling Scheme
21	EI	R		Message Profile Identifier
22	XON	X		Sending Responsible Organization
23	XON	X		Receiving Responsible Organization
24	HD	X		Sending Network Address
25	HD	X		Receiving Network Address

MSH-1 Field separator (ST)

HL7 Definition: This field contains the separator between the segment ID and the first real field. As such it serves as the separator and defines the character to be used as a separator for the rest of the message. Recommended value is | (ASCII 124).

This Specification: This value will be | (ASCII 124).

MSH-2 Encoding characters (ST)

HL7 Definition: This field contains four characters in the following order: the component separator, repetition separator, escape character, and subcomponent separator. Recommended values are ^~\& (ASCII 94, 126, 92, and 38, respectively).

This Specification: This value will be recommended values are ^~\& (ASCII 94, 126, 92, and 38, respectively).

MSH-7 Date/time of message (TS)

HL7 Definition: This field contains the date/time that the sending system created the message. If the time zone is specified, it will be used throughout the message as the default time zone.

MSH-9 Message type (MSG)

HL7 Definition: This element contains the message type, trigger event, and the message structure ID for the message.

CM Components: <message type (ID)> ^ <trigger event (ID)> ^ <message structure (ID)>

The allowed components of this element are listed in several tables maintained by HL7 (*HL7 Table 0076 - Message type, HL7 Table 0003 - Event type, and HL7 Table 0354 - Message structure*).

This Specification: In the MT-ACK-1 message type, this element shall be hard coded to ACK^R01^ACK

MSH-10 Message control ID (ST)

HL7 Definition: This element contains a number or other identifier that uniquely identifies the message.

This Specification: The sending system shall assign an identifier for the message that is unique within the namespace of the sending facility and/or application.

Note that the value of this element is *not* a reference to the message that is being acknowledged. The Message control ID of the acknowledged message appears in MSA-2 Message Control ID.

Note on Element Length: The element length for MSH-10 has been extended to 50 characters from the HL7-prescribed length of 20 characters. This extension allows sending systems to use globally unique identifiers (such as GUIDs) for Message IDs, an increasingly common practice.

MSH-11 Processing ID (PT)

HL7 Definition: This field is used to decide whether to process the message as defined in HL7 Application (level 7) Processing rules.

This Specification: This value will be 'P'.

MSH-12 Version ID (VID)

HL7 Definition: This field identifies the version.

This Specification: This value will be "2.6".

MSH-21 Message Profile Identifier (EI)

HL7 Definition: Sites may use this element to assert adherence to a Conformance Statement published by HL7 or by a site. Conformance Statements contain detailed explanations of grammar, syntax, and usage for a particular message or set of messages. This element is treated like all EI data type information, if the underscores or dashes are used as encoding characters in an interface, those characters shall be escaped to be sent in the above element.

This Specification: In the MT-ACK-1 message type, MSH-21 Conformance statement ID should be: DIALY_MT-ACK-1_R1

786
787
788
789
790

6.1.3 OBR Segment – Observation Request

The ORU message may include discrete OBX segments for individual observations reported. An OBR Segment will be used for each set of such OBX segments to establish the equipment context for the observations. All observation dates and times reported here should match OBX segments that report the same information.

HL7 Attribute Table - OBR Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBR
2	EI	X		Placer Order Number
3	EI	X		Filler Order Number
4	CWE	R	9999	Universal Service Identifier
5		X		Priority
6		X		Requested Date/Time
7	DTM	R		Observation Date/Time #
8	DTM	X		Observation End Date/Time #
9	CQ	X		Collection Volume *
10	XCN	X		Collector Identifier *
11	ID	X	0065	Specimen Action Code *
12	CWE	X	9999	Danger Code
13	ST	X		Relevant Clinical Information
14		X		Specimen Received Date/Time *
15		X		Specimen Source
16	XCN	X		Ordering Provider
17	XTN	X		Order Callback Phone Number
18	ST	X		Placer Field 1
19	ST	X		Placer Field 2
20	ST	X		Filler Field 1 +
21	ST	X		Filler Field 2 +
22	DTM	X		Results Rpt/Status Chng – Date/Time +
23	MOC	X		Charge to Practice +
24	ID	X	0074	Diagnostic Serv Sect ID
25	ID	C	0123	Result Status +
26	PRL	X		Parent Result +
27		X		Quantity/Timing
28	XCN	X		Result Copies To
29	EIP	X		Parent
30	ID	X	0124	Transportation Mode
31	CWE	X	9999	Reason for Study
32	NDL	X		Principal Result Interpreter +
33	NDL	X		Assistant Result Interpreter +
34	NDL	X		Technician +
35	NDL	X		Transcriptionist +
36	DTM	X		Scheduled Date/Time +
37	NM	X		Number of Sample Containers *
38	CWE	X	9999	Transport Logistics of Collected Sample *
39	CWE	X	9999	Collector's Comment *
40	CWE	X	9999	Transport Arrangement Responsibility
41	ID	X	0224	Transport Arranged
42	ID	X	0225	Escort Required
43	CWE	X	9999	Planned Patient Transport Comment

HL7 Attribute Table - OBR Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
44	CNE	X	0088	Procedure Code
45	CNE	X	0340	Procedure Code Modifier
46	CWE	X	0411	Placer Supplemental Service Information
47	CWE	X	0411	Filler Supplemental Service Information
48	CWE	X	0476	Medically Necessary Duplicate Procedure Reason
49	CWE	X	0507	Result Handling
50	CWE	X		Parent Universal Service Identifier
51	EI	X		Observation Group ID
52	EI	X		Parent Observation Group ID
53	CX	X		Alternate Placer Order Number

791

OBR-4 Universal Service ID

792

HL7 Definition: This field contains the identifier code for the requested observation/test/battery.

793

This specification: Identifier and Text can identify unique OBR segments that partition observations. The values for this field will be taken from the 11073_10103 MDC_IDC_SESS_TYPE enumerator MDC_IDC_ENUM_SESS_TYPE.

794

795

796

797

OBR-7 Observation Date/Time

798

HL7 Definition: This field is the clinically relevant date/time of the observation. In the case of observations taken directly from a subject, it is the actual date and time the observation was obtained. In the case of a specimen-associated study, this field shall represent the date and time the specimen was collected or obtained. (This is a results-only field except when the placer or a third party has already drawn the specimen.) This field is conditionally required. When the OBR is transmitted as part of a report message, the field **must** be filled in. If it is transmitted as part of a request **and** a sample has been sent along as part of the request, this field must be filled in because this specimen time is the physiologically relevant date/time of the observation.

799

800

801

802

803

804

805

This Specification: This value will be the same as the message time, MSH-7.

806

OBR-25 Result Status

807

HL7 Definition: This field contains the status of results for this order. This conditional field is required whenever the OBR is contained in a report message. See HL7 Table 0085.

808

809

This Specification: This value will always be Final ('F').

810

811

6.1.4 OBX Segments – Dialysis Observation Results

812

Discrete OBX segments for individual observations will be encoded into separate OBX segments as individual observations or measurements. These OBX segments will be preceded by an appropriate OBR segment to set the context for observations dealing with the dialysis device, characteristics, or treatment measurements.

813

814

815

816

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R	9999	Observation Identifier
4	ST	R		Observation Sub-ID

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
5	varies	R		Observation Value
6	CWE	R	9999	Units
7	ST	CE		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	CE		Date/Time of the Observation
15	CWE	X	9999	Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X	9999	Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

817

OBX-1 Set ID

818

HL7 Specification: This field contains the sequence number.

819

OBX-2 Value Type

820

HL7 Specification: This field defines the data type of OBX-5, Observation Value.

821

This specification: Only value identified in the table below will be used. The value is based upon the value type as identified in the Dialysis Data Elements Table.

822

823

824

IEEE to HL7 Data Type Matching

Applicable IEEE 11073 MDC_IDC types	HL7 v2 data type
String	ST
Enumerated	CWE or CNE
Date Time	DTM
Numeric	NM

825

OBX-3 Observation Identifier

826

HL7 Specification: This field contains a unique identifier for the observation.

827

This specification: The format is that of the Coded Element (CWE). Example: "0^MDCX_HDIALY_DIALYSATE_FLOW_MODE^MDC".

828

829

OBX-4 Observation Sub-ID

830

HL7 Specification: This field is used to distinguish between multiple OBX segments with the same observation ID organized under one OBR.

831

832

This specification: This value is expressed in a dotted notation consisting of 4 parts, MDS, VMD, CHAN, and Parameter. These values are identified in Dialysis Data Elements Table.

OBX-5 Observation Value

HL7 Specification: This field contains the value observed by the observation producer.

This specification: This is the actual value is expressed a string value.

OBX-6 Unit

HL7 Specification: This field contains the units of measurement for the value.

This specification: Will be coded with the MDC_IDC Nomenclature (based on UCUM) Unit for associated observation. Example, "ml/min^milliliter per minute^UCUM"

OBX-7 Reference Range

HL7 Specification: Reference range for the value.

- a) lower limit-upper limit (when both lower and upper limits are defined, e.g., for Venous Pressure, 20-400)
- b) > lower limit (if no upper limit, e.g., >10)
- c) < upper limit (if no lower limit, e.g., <15)

This specification: For Dialysis Machines' this field is used to provide the alarm limits for a numeric parameter.

OBX-8 Interpretation Codes

HL7 Specification: One or more codes specifying a categorical assessment of the observation value (OBX-5), such as "Normal", "Abnormal", "Positive", "Negative", "Resistant", "Susceptible", etc.

This specification: This field is used to provide zero or more codes separated by the tilde ("~"). A dialysis machine provides two mandatory codes, Alarm priority and Alarm Type, and one optional code, Abnormality Type.

The alarm priority is dialysis machine specific and may vary by manufacturer.

Alarm Priority	Code
High	PH
Medium	PM
Low	PL

The Alarm Type is specific in the Alert Type Column in the Dialysate Data Elements Table.

Alarm Type	Code
Physiological	SP
Technical	ST
Advisory	SA

Abnormality Type	Code
Below low normal	L
Above high normal	H

OBX-11 Observation Result Status

HL7 Specification: This field contains the observation result status, *HL7 Table 0085*.

This specification: This field will always be Final ('F').

OBX-14 Date/Time of Observation

HL7 Specification: This field is required in two circumstances. The first is when the observations reported beneath one report header (OBR) have different dates/times. This could occur in the case of queries, timed test sequences, or clearance studies where one measurement within a battery may have a different time than another measurement.

This specification: Dialysis Machines' use this field to report the time and date of spot check values such as non-invasive blood pressure.

6.1.5 PID – Patient Identification Segment

The PID segment is used by all applications as the primary means of communicating patient identification information. This segment contains permanent patient identifying and demographic information that, for the most part, is not likely to change frequently.

HL7 Attribute Table - PID Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID - PID
2		X		Patient ID
3	CX	R		Patient Identifier List
4		X		Alternate Patient ID - PID
5	XPN	R	0200	Patient Name
6	XPN	X		Mother's Maiden Name
7	DTM	X		Date/Time of Birth
8	CWE	X	0001	Administrative Sex
9		X		Patient Alias
10	CWE	X	0005	Race
11	XAD	X		Patient Address
12		X		County Code
13	XTN	X		Phone Number - Home
14	XTN	X		Phone Number - Business
15	CWE	X	0296	Primary Language
16	CWE	X	0002	Marital Status
17	CWE	X	0006	Religion
18	CX	X	0061	Patient Account Number
19		X		SSN Number - Patient
20		X		Driver's License Number - Patient
21	CX	X	0061	Mother's Identifier
22	CWE	X	0189	Ethnic Group
23	ST	X		Birth Place
24	ID	X	0136	Multiple Birth Indicator
25	NM	X		Birth Order
26	CWE	X	0171	Citizenship
27	CWE	X	0172	Veterans Military Status
28		X		Nationality
29	DTM	X		Patient Death Date and Time
30	ID	X	0136	Patient Death Indicator
31	ID	X	0136	Identity Unknown Indicator
32	CWE	X	0445	Identity Reliability Code
33	DTM	X		Last Update Date/Time
34	HD	X		Last Update Facility
35	CWE	X	0446	Species Code

HL7 Attribute Table - PID Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
36	CWE	X	0447	Breed Code
37	ST	X		Strain
38	CWE	X	0429	Production Class Code
39	CWE	X	0171	Tribal Citizenship
40	XTN	X		Patient Telecommunication Information

873

874

PID-3.1 Patient Identifier List

875

876

877

878

879

880

881

HL7 Table 0203

Code	Description	Notes	Use
U	Model and Serial Number of Device MDC_ID_MODEL_NUMBER and MDC_ID_PROD_SPEC_SERIAL	Model and Serial number will be concatenated together and will be unique within an Assigning Authority. The format of the ID will be following: "model:xxx/serial:yyy" Example: model:XZY987/serial:abc123	R
SS	Patient Social Security Number	Social Security number will be included if known.	RE

882

883

6.1.6 PV1 Segment

884

885

886

887

The PV1 segment is used by Registration/Patient Administration applications to communicate information on an account or visit-specific basis. Because it is required in the PCD-01 Message, it is included in this standard. The Dialysis Machine Provides the bare minimum of data to be compliant.

HL7 Attribute Table - PV1 Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID - PV1
2	CWE	R	0004	Patient Class
3	PL	X		Assigned Patient Location
4	CWE	X	0007	Admission Type
5	CX	X		Preadmit Number
6	PL	X		Prior Patient Location
7	XCN	X	0010	Attending Doctor
8	XCN	X	0010	Referring Doctor
9	XCN	X		Consulting Doctor
10	CWE	X	0069	Hospital Service
11	PL	X		Temporary Location
12	CWE	X	0087	Preadmit Test Indicator
13	CWE	X	0092	Re-admission Indicator
14	CWE	X	0023	Admit Source
15	CWE	X	0009	Ambulatory Status

HL7 Attribute Table - PV1 Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
16	CWE	X	0099	VIP Indicator
17	XCN	X	0010	Admitting Doctor
18	CWE	X	0018	Patient Type
19	CX	RE		Visit Number
20	FC	X	0064	Financial Class
21	CWE	X	0032	Charge Price Indicator
22	CWE	X	0045	Courtesy Code
23	CWE	X	0046	Credit Rating
24	CWE	X	0044	Contract Code
25	DT	X		Contract Effective Date
26	NM	X		Contract Amount
27	NM	X		Contract Period
28	CWE	X	0073	Interest Code
29	CWE	X	0110	Transfer to Bad Debt Code
30	DT	X		Transfer to Bad Debt Date
31	CWE	X	0021	Bad Debt Agency Code
32	NM	X		Bad Debt Transfer Amount
33	NM	X		Bad Debt Recovery Amount
34	CWE	X	0111	Delete Account Indicator
35	DT	X		Delete Account Date
36	CWE	X	0112	Discharge Disposition
37	DLD	X	0113	Discharged to Location
38	CWE	X	0114	Diet Type
39	CWE	X	0115	Servicing Facility
40		X		Bed Status
41	CWE	X	0117	Account Status
42	PL	X		Pending Location
43	PL	X		Prior Temporary Location
44	DTM	X		Admit Date/Time
45	DTM	X		Discharge Date/Time
46	NM	X		Current Patient Balance
47	NM	X		Total Charges
48	NM	X		Total Adjustments
49	NM	X		Total Payments
50	CX	X	0203	Alternate Visit ID
51	CWE	X	0326	Visit Indicator
52		X		Other Healthcare Provider
53	ST	X		Service Episode Description
54	CX	X		Service Episode Identifier

888

889

PV1-2 Patient Class (IS)

890

HL7 Definition: This field is used by systems to categorize patients by site. See HL7 Table 0004.

891

This specification: Dialysis machines will always use the Unknown Patient Class.

892

PV1-19 Visit Number (CX)

893

HL7 Definition: This field contains the unique number assigned to each patient visit.

894

This specification: This field will be empty.

6.2 HL7 Tables

The values contained in the tables in this section may be a subset of the values supported by the HL7 specification. Only values supported or used by Dialysis Machines are shown.

HL7 Table 0008 - Acknowledgment Code

Value	Description
AA	Original mode: Application Accept - Enhanced mode: Application acknowledgment: Accept
AE	Original mode: Application Error - Enhanced mode: Application acknowledgment: Error
AR	Original mode: Application Reject - Enhanced mode: Application acknowledgment: Reject
CA	Enhanced mode: Accept acknowledgment: Commit Accept
CE	Enhanced mode: Accept acknowledgment: Commit Error
CR	Enhanced mode: Accept acknowledgment: Commit Reject

HL7 Table 0038 - Order status

Value	Description
A	Some, but not all, results available
CA	Order was canceled
CM	Order is completed
DC	Order was discontinued
ER	Error, order not found
HD	Order is on hold
IP	In process, unspecified
RP	Order has been replaced
SC	In process, scheduled

HL7 Table 0085 - Observation result status codes interpretation

Value	Description
F	Final results; Can only be changed with a corrected result.

HL7 Table 0123 - Result status for OBR segment

Value	Description
O	Order received; specimen not yet received
I	No results available; specimen received, procedure incomplete
S	No results available; procedure scheduled, but not done
A	Some, but not all, results available
P	Preliminary: A verified early result is available, final results not yet obtained

Value	Description
C	Correction to results
R	Results stored; not yet verified
F	Final results; results stored and verified. Can only be changed with a corrected result.
X	No results available; Order canceled.
Y	No order on record for this test. (Used only on queries)
Z	No record of this patient. (Used only on queries)

906

907

HL7 Table 0125 - Value type

Value	Description
DT	Date
NM	Numeric
ST	String Data.
TM	Time

908

909

HL7 Table 0155 - Accept/application acknowledgment conditions

Value	Description
AL	Always

910

911

912

913

User Defined Table 0396 – Coding Systems

Value	Description
MDC	ISO/IEEE
UCUM	Unified Code for Units of Measure (UCUM)

914

915

HL7 table 0357 – Message Error Condition Codes

Value	Description	Comment
0	Message accepted	Success. Optional, as the AA conveys success. Used for systems that shall always return a status code.
100	Segment sequence error	Error: The message segments were not in the proper order, or required segments are missing.
101	Required field missing	Error: A required field is missing from a segment
102	Data type error	Error: The field contained data of the wrong data type, e.g. an NM field contained "FOO".

Value	Description	Comment
103	Table value not found	Error: A field of data type ID or IS was compared against the corresponding table, and no match was found.
200	Unsupported message type	Rejection: The Message Type is not supported.
201	Unsupported event code	Rejection: The Event Code is not supported.
202	Unsupported processing id	Rejection: The Processing ID is not supported.
203	Unsupported version id	Rejection: The Version ID is not supported.
204	Unknown key identifier	Rejection: The ID of the patient, order, etc., was not found. Used for transactions <i>other than</i> additions, e.g. transfer of a non-existent patient.
205	Duplicate key identifier	Rejection: The ID of the patient, order, etc., already exists. Used in response to addition
206	Application record locked	Rejection: The transaction could not be performed at the application storage level, e.g., database locked.
207	Application internal error	Rejection: A catchall for internal errors not explicitly covered by other codes.

916

917

HL7 table 0516 – Error Severity

Value	Description	Comment
W	Warning	Transaction successful, but there may issues
I	Information	Transaction was successful but includes information e.g., inform patient
E	Error	Transaction was unsuccessful

918

919

7 Dialysis Data Elements

7.1 Tables

This section defines tables of values used by the message elements.

_TBL_01 – Mode of Operation

Value	Description
PRETX	Pre-Treatment
TX	Treatment
POSTTX	Post-Treatment
DIS	Disinfection
IDL	Idle
SVC	Service

_TBL_02 – Treatment Modality

Value	Description
HD	Hemodialysis
HDF	Hemodiafiltration
HF	Hemofiltration
SLED	Sustained low efficiency dialysis.
IUF	Isolated Ultrafiltration
HP	Hemoperfusion

_TBL_03 – Anticoagulation Mode

Value	Description
BOL	Bolus
CON	Continuous
BOLCON	Bolus and Continuous
PRO	Profile
BOLPRO	Bolus and Profile
NONE	None

_TBL_04 – Patient Position

Value	Description
SITTING	Sitting
STANDING	Standing
SUPINE	Supine

_TBL_05 – Blood Pump Mode

Value	Description
2N	Double Needle
1N1P	Single Needle/Single Pump
1N2P	Single Needle/Double Pump

_TBL_06 – Dialysis Flow Mode

Value	Description
CONST	Constant
AUTO	Auto
PRO	Profile
STBY	Standby

_TBL_07 – Replacement Fluid Delivery Mode

Value	Description
CONST	Constant
PRO	Profile
NONE	None, Not Applicable
BOLUS	Bolus
AUTO	Automatic

_TBL_08 – RF Delivery Location

Value	Description
PREF	Pre-Filter
POSTF	Post-Filter
PREF-POSTF	Pre-Filter/Post-Filter
PREP-POSTF	Pre-Pump/Post-Filter

_TBL_09 – Fluid Source

Value	Description
BAG	Bag
ONLINE	Online

_TBL_10 – Plasma Volume Marker

Value	Description
NONE	None – No selection made at time of reading/reporting
SYMP	Symptomatic – Noting a symptom for continue monitoring.
INTER	Intervention – Action was taken or a machine setting was changed.

_TBL_11 – Plasma Volume Profile

Value	Description
FLAT	Change less than or equal to 3% per hour
GRAD	Change greater than 3% and less than or equal to 6.5% per hour
STEEP	Change greater than 6.5% per hour.

_TBL_12 – Completion Method

Value	Description
CT	Clock Time
AT	Actual Treatment Time
UF	UF Removed
KTV	Kt/V

_TBL_13 – Ultrafiltration Mode

Value	Description
CIUF	Constant Isolated UF
PIUF	Profile Isolated UF
CUF	Constant UF
PUF	Profile UF
NONE	None


_TBL_14 – Bicarb Delivery Mode

Value	Description
CONST	Constant
PRO	Profile
NA	Not Applicable

_TBL_15 – Sodium Delivery Mode

Value	Description
CONST	Constant
PRO	Profile

7.2 Message Element Terminology



In the table below when usage is conditional (C), a number is given. The number refers to one of the following conditions.

- 1) if Anticoagulation Mode is not None
- 2) if Anticoagulation Mode includes Bolus
- 3) if Anticoagulation Mode includes Continuous or Profile

- 4) One for each anticoagulant pump
- 5) If Blood Pump Mode is not double needle
- 6) if mode is HD or HDF
- 7) If Mode is HF or HDF and Dilution is Post or Pre-Post
- 8) If Mode is HDF or HF
- 9) if Replacement Fluid Mode is Bolus

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_DEV_HDIALY_MACHINE_MDS					1								
MDC_ATTR_ID_MODEL		Dialysis Device	Dialysis Device Model and Manufacturer	String containing device model and manufacturer's name	1.1.1.1	All	Identifier	All	String	Alphanumeric	N/A	1..1	M
MDC_ID_MODEL_MANUFACTURER		Dialysis Device	Dialysis Device Manufacturer	String containing device manufacturer's name (or use MDC_ATTR_ID_MODEL)	1.1.1.2	All	Identifier	All	String	Alphanumeric	N/A	1..1	M
MDC_ID_MODEL_NUMBER		Dialysis Device	Dialysis Device Model	String containing device manufacturer's model identifier for the device (or use MDC_ATTR_ID_MODEL)	1.1.1.3	All	Identifier	All	String	Alphanumeric	N/A	1..1	M
MDC_ID_PROD_SPEC_SERIAL		Dialysis Device	Dialysis Device Serial Number	String containing the device's serial number	1.1.1.4	All	Identifier	All	String	Alphanumeric	N/A	1..1	M
MDC_ID_PROD_SPEC_SW		Dialysis Device	Device Software Version	String defining the version of the software on the machine	1.1.1.5	All	Identifier	All	String	Alphanumeric	N/A	1..1	M
PRT-10 PRT-16-20 MDC_ATTR_ID_UDI		Dialysis Device	Dialysis Device UDI	String containing the device's FDA required Universal Device Identifier [Note: MDC_ATTR_ID_UDI contains three components: UdiAuthority, UdiIssuer, UdiLabel to support non-FDA entities, but currently there is no approved way of encoding this using HL7 V2; PRT-10 and PRT-16-20 is used instead; UdiLabel = PRT-10.]		All	Identifier	All	String	Alphanumeric	N/A	1..1	O
MDCX_DEV_HDIALY_VMD					1.1								
MDCX_DEV_HDIALY_MACH_CONFIG_CHAN					1.1.1								
MDCX_HDIALY_MACH_TIME		Dialysis Device	Dialysis Device Time	Date and time as recorded on the dialysis device's internal clock with offset to UTC (based on location and DST)	1.1.1.1	All	Identifier	All	Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	M
MDCX_HDIALY_MACH_MODE_DESCRIPTION		Dialysis Device	Machine Mode Description	Manufacturer-specific description of the machine mode of operation	1.1.1.2	All	Parameter	Episodic	String	Alphanumeric	N/A	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_MACH_MODE_OF_OPERATION		Dialysis Device	Machine Mode Of Operation	Process the machine is currently performing	1.1.1.3	All	Parameter	Episodic	_TBL_01	Alphanumeric	N/A	1..1	M
MDCX_HDIALY_MACH_PM_TIME_REMAIN		Dialysis Device	Preventive Maintenance Time Remaining	Hours of operation remaining until next preventative maintenance is required	1.1.1.4	All	Status	Episodic	Numeric	XXXX	h	0..1	O
MDCX_HDIALY_MACH_MAINT_TX_REMAIN		Dialysis Device	Preventive Maintenance Treatments Remaining	Number of treatments remaining until next preventative maintenance is required	1.1.1.5	All	Status	Episodic	Numeric	XXXX	N/A	0..1	O
MDCX_HDIALY_MACH_BLD_PUMP_ON		Dialysis Device	Treatment Blood Pump On	Status indicating if blood pump is running	1.1.1.6	Intradialytic	Parameter	Episodic	Bool	T / F	N/A	1..1	M
MDCX_HDIALY_MACH_TX_FLUID_BYPASS		Dialysis Device	Treatment Fluid Bypass	Fluid pathway state to indicate if fluid is NOT flowing through dialyzer	1.1.1.7	Intradialytic	Parameter	Episodic	Bool	T / F	N/A	1..1	M
MDCX_HDIALY_MACH_TX_MODALITY		Dialysis Device	Treatment Modality	The way a patient receives dialysis	1.1.1.8	All	Parameter	Episodic	_TBL_02	Alphanumeric	N/A	1..1	M
MDCX_HDIALY_MACH_THE_RAPY_TIME		Dialysis Device State	Elapsed Treatment Time	Accumulated time that the patient has been actively treated	1.1.1.9	Intradialytic	Status	Periodic	Numeric	XXXX	min	1..1	M
MDCX_HDIALY_MACH_TIME_REMAIN		Dialysis Device State	Remaining Treatment Time	Number of minutes left for current dialysis session	1.1.1.10	Intradialytic	Status	Periodic	Numeric	XXXX	min	1..1	M
MDCX_DEV_HDIALY_ANTI COAG_PUMP_CHAN					1.1.2								
MDCX_HDIALY_ANTI COAG_ACCUM_DELIV		Anticoagulant Pump	Accumulated Delivered Anticoagulant	Total amount of anticoagulation drug administered to patient so far	1.1.2.1	Intradialytic	Status	Both	Numeric	XXX.X	mL	0..1	C1
MDCX_HDIALY_ANTI COAG_BOLUS_VOL		Anticoagulant Pump	Anticoagulant Bolus Volume	User defined amount of anticoagulation drug to be administered in this single large dose	1.1.2.2	Intradialytic	Parameter	Episodic	Numeric	XX.X	mL	0..1	C2
MDCX_HDIALY_ANTI COAG_INFUS_RATE		Anticoagulant Pump	Anticoagulant Infusion Rate	User defined amount of anticoagulation drug to be administered per unit time	1.1.2.3	Intradialytic	Parameter	Episodic	Numeric	XX.X	mL/h	0..1	C3
MDCX_EVT_HDIALY_ANTI COAG_STOP	tech	Anticoagulant Pump	Anticoagulant Pump Stop	Notification that the anticoagulant pump has stopped	1.1.2.4	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_ANTI COAG_TYPE		Anticoagulant Pump	Anticoagulant Type	Name of the anticoagulant	1.1.2.5	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	0..1	O
MDCX_HDIALY_ANTI COAG_MODE		Anticoagulant Pump	Anticoagulation Mode	Method of anticoagulant administration during dialysis treatment	1.1.2.6	Intradialytic	Parameter	Episodic	_TBL_03	Alphanumeric	N/A	1..1	C4
MDCX_HDIALY_ANTI COAG_PUMP_NUMBER		Anticoagulant Pump	Number of Anticoagulant Pumps	The number of anticoagulant pumps being used.	1.1.2.7	Intradialytic	Parameter	Episodic	Numeric	X	N/A	1..1	M

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_EMPTY	tech	Anticoagulant Pump	Syringe Empty	Notification that the syringe plunger has reached the end of its functional range	1.1.2.8	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_ANTICOAG_SYRINGE_BRAND		Anticoagulant Pump	Syringe Name	Name of the syringe	1.1.2.9	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	0..1	O
MDCX_HDIALY_ANTICOAG_SYRINGE_VOL		Anticoagulant Pump	Syringe Volume	Volume of the syringe with anticoagulant.	1.1.2.10	Intradialytic	Parameter	Episodic	Numeric	XXX	mL	0..1	O
MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_SIZE	tech	Anticoagulant Pump	Wrong Syringe Size	Notification that the size of the syringe is not the size defined for anticoagulant administration	1.1.2.11	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_DEV_HDIALY_BLOOD_PUMP_CHAN					1.1.3								
MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE		Blood Pump	Actual Blood Flow Rate	Adjusted blood flow rate based on the blood flow rate setting and pressure drop caused by blood line, needle and vascular access	1.1.3.1	Intradialytic	Status	Periodic	numeric	XXX	mL/min	0..1	O
MDCX_HDIALY_BLD_PRESS_ART	phys tech high low thr	Blood Pump	Arterial Pressure	Pressure of arterial access line pre blood pump	1.1.3.2	Intradialytic	Status	Periodic	Numeric	±XXX	mm[Hg]	1..1	M
MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_AVG		Blood Pump	Average Blood Flow Rate	Average of the actual blood flow rate over the course of treatment.	1.1.3.3	Intradialytic	Status	Periodic	Numeric	XXX	mL/min	0..1	O
MDCX_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_AVG_SETTING		Blood Pump	Blood Flow Rate Setting	The rate at which the user programmed the blood flow	1.1.3.4	Intradialytic	Parameter	Episodic	Numeric	XXX	mL/min	1..1	M
MDCX_HDIALY_BLD_PUMP_MODE		Blood Pump	Blood Pump Mode	Therapy method in which blood is retrieved and returned to the patient.	1.1.3.5	Intradialytic	Parameter	Episodic	_TBL_05	Alphanumeric	N/A	1..1	M
MDCX_EVT_HDIALY_BLD_PUMP_STOP	tech	Blood Pump	Blood Pump Stop	Notification that the blood pump has stopped	1.1.3.6	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_HDIALY_BLD_PUMP_TUBING_SIZE		Blood Pump	Blood Tubing Size	Inner diameter of pump segment of tubing used in extracorporeal circulation	1.1.3.7	Intradialytic	Parameter	Episodic	String	alphanumeric	N/A	0..1	O
MDCX_HDIALY_BLD_PUMP_PAT_BODY_TEMP_BY_BLOOD		Blood Pump	Body Temperature	Patient body temperature, directly measured or calculated based on temperature of blood	1.1.3.8	Intradialytic	Status	Periodic	Numeric	XX.X	Cel	0..1	O
MDCX_HDIALY_BLD_PUMP_CHANGE_IN_ENERGY		Blood Pump	Change In Energy	Measurement of the energy (temperate) flux to and from the patient	1.1.3.9	Intradialytic	Status	Periodic	Numeric	±XXX.X	kJ/h	0..1	O
MDCX_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP	tech high low thr	Blood Pump	Post Pump Arterial Pressure	Pressure of arterial access line post blood pump	1.1.3.10	Intradialytic	Status	Periodic	Numeric	±XXX	mm[Hg]	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_BLD_PUMP_PRIMING_VOL		Blood Pump	Priming Volume	Volume of fluid used to prime the extracorporeal circuit	1.1.3.11	Intradialytic	Parameter	Episodic	Numeric	XXX	mL	0..1	O
MDCX_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS	tech high low thr	Blood Pump	Single Needle Switching Pressure	In single needle mode, the pressure used to switch phases at high pressure	1.1.3.12	Intradialytic	Parameter	Episodic	Numeric	± XXX	mm[Hg]	0..1	O
MDCX_HDIALY_BLD_PUMP_STROKE_VOL		Blood Pump	Stroke Volume	The volume of blood passing through the dialyzer during each single needle cycle.	1.1.3.13	Intradialytic	Status	Periodic	Numeric	XXX	mL	0..1	C5
MDCX_HDIALY_BLD_PUMP_BLOOD_PROCESSED_TOTAL		Blood Pump	Total Blood Processed	Accumulated volume of blood circulated by extracorporeal circulation.	1.1.3.14	Intradialytic	Status	Periodic	Numeric	XXXX.XX	L	0..1	O
MDCX_HDIALY_BLD_PUMP_PRESS_VEN	phys tech high low thr	Blood Pump	Venous Pressure	Pressure of the venous access line	1.1.3.15	Intradialytic	Status	Periodic	Numeric	±XXX	mm[Hg]	1..1	M
MDCX_HDIALY_BLD_PUMP_BLOOD_TEMP_VEN		Blood Pump	Venous Temperature	Temperature of the blood measured in the venous access line	1.1.3.16	Intradialytic	Status	Periodic	Numeric	XX.X	Cel	0..1	O
MDCX_DEV_HDIALY_FLUID_CHAN					1.1.4								
MDCX_HDIALY_BICARB_COND		Dialysis Fluid	Actual Bicarbonate Conductivity	Measured conductivity of bicarbonate	1.1.4.1	Intradialytic	Status	Periodic	Numeric	XX.XX	mS/cm	0..1	O
MDCX_HDIALY_BICARB_COND_SETTING		Dialysis Fluid	Bicarbonate Conductivity Setting	Bicarbonate conductivity	1.1.4.2	Intradialytic	Parameter	Episodic	Numeric	XX.XX	mS/cm	0..1	O
MDCX_HDIALY_DIALYSATE_COND	tech high low	Dialysis Fluid	Actual Conductivity	Measured conductivity of dialysate	1.1.4.3	Intradialytic	Status	Periodic	Numeric	XX.XX	mS/cm	1..1	M
MDCX_HDIALY_DIALYSATE_COND_SETTING		Dialysis Fluid	Dialysate Conductivity Setting	Acid and Bicarbonate conductivity combined	1.1.4.4	Intradialytic	Parameter	Episodic	Numeric	XX.XX	mS/cm	0..1	O
MDCX_HDIALY_DIALYSATE_FLOW_RATE	tech low	Dialysis Fluid	Actual Dialysate Flow Rate	Rate of dialysate flow at any given point during a treatment	1.1.4.5	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	1..1	C6
MDCX_HDIALY_DIALYSATE_FLOW_RATE_SETTING		Dialysis Fluid	Dialysate Flow Rate Setting	Rate at which the user programmed the dialysate flow	1.1.4.6	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL/min	0..1	C6
MDCX_HDIALY_DIALYSATE_FLOW_MODE		Dialysis Fluid	Dialysate Flow Mode	The mode of controlling dialysate flow rate	1.1.4.7	All	Parameter	Episodic	_TBL_06	alphanumeric	N/A	1..1	C6
MDCX_HDIALY_DIALYSATE_TEMP	tech high low	Dialysis Fluid	Actual Dialysate Temperature	Temperature of dialysate	1.1.4.8	Intradialytic	Status	Periodic	Numeric	XX.X	Cel	1..1	O
MDCX_HDIALY_DIALYSATE_TEMP_SETTING		Dialysis Fluid	Dialysate Temperature Setting	Temperature at which dialysate is to be delivered to patient	1.1.4.9	Intradialytic	Parameter	Episodic	Numeric	XX.X	Cel	1..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_DIALYSATE_AMMONIA		Dialysis Fluid	Ammonia Concentration	Measured ammonia concentration of the dialysate	1.1.4.10	Intradialytic	Status	Periodic	Numeric	XX.XX	[ppm]	0..1	O
MDCX_HDIALY_DIALYSATE_FLOW_RATE_AVG		Dialysis Fluid	Average Dialysate Flow Rate	Average rate of dialysate flow during course of treatment	1.1.4.11	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	0..1	O
MDCX_HDIALY_CONC_HCO3_SETTING		Dialysis Fluid	Bicarbonate Concentration Setting	User programmed Bicarbonate concentration in dialysate.	1.1.4.13	Intradialytic	Parameter	Episodic	Numeric	XXX.X	mmol/L	1..1	O
MDCX_HDIALY_CONC_HCO3_MODE		Dialysis Fluid	Bicarbonate Mode	Process for delivering adjusted bicarbonate concentration in dialysate	1.1.4.14	Intradialytic	Parameter	Episodic	_TBL_14	Alphanumeric	N/A	1..1	O
MDCX_EVT_HDIALY_BLOOD_LEAK	tech	Dialysis Fluid	Blood Leak	Notification that blood has been detected in the dialysate fluid	1.1.4.15	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_HDIALY_DIALYSATE_NAME		Dialysis Fluid	Dialysate Name	String containing a descriptive name of the dialysate	1.1.4.16	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	1..1	O
MDCX_HDIALY_DIALYSATE_CONC_ACETATE		Dialysis Fluid	Dialysate Acetate Concentration	Acetate concentration in dialysate	1.1.4.17	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_HCO3		Dialysis Fluid	Dialysate Bicarbonate Concentration	Bicarbonate concentration in dialysate	1.1.4.18	Intradialytic	Parameter	Episodic	Numeric	XXX.X	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_CHLORIDE		Dialysis Fluid	Dialysate Chloride Concentration	Chloride concentration in dialysate	1.1.4.19	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_MG		Dialysis Fluid	Dialysate Magnesium Concentration	Magnesium concentration in dialysate	1.1.4.20	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_CA		Dialysis Fluid	Dialysate Calcium Concentration	Calcium concentration in dialysate	1.1.4.21	Intradialytic	Parameter	Episodic	Numeric	X.XX	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_CITRATE		Dialysis Fluid	Dialysate Citrate Concentration	Citrate concentration in dialysate	1.1.4.22	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_GLU		Dialysis Fluid	Dialysate Glucose Concentration	Glucose concentration in dialysate	1.1.4.23	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_K		Dialysis Fluid	Dialysate Potassium Concentration	Potassium concentration in dialysate	1.1.4.24	Intradialytic	Parameter	Episodic	Numeric	X.X	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_NA		Dialysis Fluid	Dialysate Sodium Concentration	Sodium concentration in dialysate	1.1.4.25	Intradialytic	Parameter	Episodic	Numeric	XXXX	mmol/L	0..1	O
MDCX_HDIALY_DIALYSATE_CONC_NA_SETTING		Dialysis Fluid	Dialysate Sodium concentration Setting	User programmed sodium concentration in dialysate	1.1.4.26	Intradialytic	Parameter	Episodic	Numeric	XXXX	mmol/L	1..1	O
MDCX_HDIALY_DIALYSATE_CONC_NA_MODE		Dialysis Fluid	Dialysate Sodium Mode	Process for delivering adjusted sodium concentration in dialysate	1.1.4.27	Intradialytic	Parameter	Episodic	_TBL_15	Alphanumeric	N/A	1..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_DIALYSATE_CONC_PH		Dialysis Fluid	pH	Measured pH of the dialysate	1.1.4.28	Intradialytic	Status	Periodic	Numeric	XX.X	[pH]	0..1	O
MDCX_HDIALY_DIALYSATE_VOL_DELIV		Dialysis Fluid	Total Dialysate Volume	Total amount of dialysate used during a treatment	1.1.4.29	Intradialytic	Status	Periodic	Numeric	XXX.XX	L	0..1	O
MDCX_DEV_HDIALY_FILTER_CHAN					1.1.5								
MDCX_HDIALY_FILTER_NAME		Filter (Dialyzer)	Dialyzer Name	Dialyzer manufacturer and model	1.1.5.1	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	0..1	O
MDCX_HDIALY_FILTER_TRANSMEMBRANE_PRESS	tech high low	Filter (Dialyzer)	Actual Transmembrane Pressure	Difference in pressure between blood compartment and dialysate compartment of dialyzer	1.1.5.2	Intradialytic	Status	Periodic	Numeric	±XXX	mm[Hg]	1..1	M
MDCX_DEV_HDIALY_CONVECTIVE_CHAN					1.1.6								
MDCX_HDIALY_RF_DILUTION_LOCATION		Replacement Fluid	Replacement Fluid Dilution Location	Delivery location of the replacement fluid	1.1.6.1	Intradialytic	Parameter	Episodic	_TBL_08	Alphanumeric	N/A	1..1	C8
MDCX_HDIALY_RF_POST_FILTER_NAME		Replacement Fluid	Post Filter RF Name	String containing a descriptive name of the Post Filter RF	1.1.6.2	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	1..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_ACETATE		Replacement Fluid	Post Filter RF Acetate Concentration	Acetate concentration in the Post Filter RF	1.1.6.3	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_HCO3		Replacement Fluid	Post Filter RF Bicarbonate Concentration	Bicarbonate concentration in the Post Filter RF	1.1.6.4	Intradialytic	Parameter	Episodic	Numeric	XXX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_CHLORIDE		Replacement Fluid	Post Filter RF Chloride Concentration	Chloride concentration in the Post Filter RF	1.1.6.5	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_MG		Replacement Fluid	Post Filter RF Magnesium Concentration	Magnesium concentration in the Post Filter RF	1.1.6.6	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_CA		Replacement Fluid	Post Filter RF Calcium Concentration	Calcium concentration in the Post Filter RF	1.1.6.7	Intradialytic	Parameter	Episodic	Numeric	X.XX	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_CITRATE		Replacement Fluid	Post Filter RF Citrate Concentration	Citrate concentration in the Post Filter RF	1.1.6.8	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_GLU		Replacement Fluid	Post Filter RF Glucose Concentration	Glucose concentration in the Post Filter RF	1.1.6.9	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_K		Replacement Fluid	Post Filter RF Potassium Concentration	Potassium concentration in the Post Filter RF	1.1.6.10	Intradialytic	Parameter	Episodic	Numeric	X.X	mmol/L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_CONC_NA		Replacement Fluid	Post Filter RF Sodium Concentration	Sodium concentration in the Post Filter RF	1.1.6.11	Intradialytic	Parameter	Episodic	Numeric	XXXX	mmol/L	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE		Replacement Fluid	Actual Post-Filter Replacement Fluid Flow Rate	Rate of post -filter replacement fluid flow at any given point during a treatment	1.1.6.12	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	1..1	C7
MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE_SETTING		Replacement Fluid	Post-Filter Replacement Fluid Flow Rate Setting	Rate at which the user programmed the post -filter replacement fluid flow	1.1.6.13	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL/min	0..1	C7
MDCX_HDIALY_RF_POST_FILTER_FLOW_RATE_AVG		Replacement Fluid	Average Post-Filter Replacement Fluid Flow Rate	Average rate of post -filter replacement fluid flow during course of treatment.	1.1.6.14	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	0..1	O
MDCX_HDIALY_RF_POST_FILTER_TEMP		Replacement Fluid	Actual Post-Filter Replacement Fluid Temperature	Temperature of post -filter replacement fluid	1.1.6.15	Intradialytic	Status	Periodic	Numeric	XX.X	Cel	1..1	C7
MDCX_HDIALY_RF_POST_FILTER_TEMP_SETTING		Replacement Fluid	Post-Filter Replacement Fluid Temperature Setting	Temperature at which post -filter replacement fluid is to be delivered to patient	1.1.6.16	Intradialytic	Parameter	Episodic	Numeric	XX.X	Cel	1..1	C7
MDCX_HDIALY_RF_POST_FILTER_VOL		Replacement Fluid	Total Post-Filter Replacement Fluid Volume	Total amount of post -filter replacement fluid used during a treatment	1.1.6.17	Intradialytic	Status	Periodic	Numeric	XXX.XX	L	0..1	O
MDCX_HDIALY_RF_POST_FILTER_VOL_SETTING		Replacement Fluid	Post-Filter Replacement Fluid Volume setting	Programmed total amount of post -filter replacement fluid volume	1.1.6.18	Intradialytic	Parameter	Episodic	Numeric	XXX.XX	L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_NAME		Replacement Fluid	Pre-Filter RF Name	String containing a descriptive name of the Pre-Filter RF	1.1.6.19	Intradialytic	Parameter	Episodic	String	Alphanumeric	N/A	1..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_ACETATE		Replacement Fluid	Pre-Filter RF Acetate Concentration	Acetate concentration in the Pre-Filter RF	1.1.6.20	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_HCO3		Replacement Fluid	Pre-Filter RF Bicarbonate Concentration	Bicarbonate concentration in the Pre-Filter RF	1.1.6.21	Intradialytic	Parameter	Episodic	Numeric	XXX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_CHLORIDE		Replacement Fluid	Pre-Filter RF Chloride Concentration	Chloride concentration in the Pre-Filter RF	1.1.6.22	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_MG		Replacement Fluid	Pre-Filter RF Magnesium Concentration	Magnesium concentration in the Pre-Filter RF	1.1.6.23	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_CA		Replacement Fluid	Pre-Filter RF Calcium Concentration	Calcium concentration in the Pre-Filter RF	1.1.6.24	Intradialytic	Parameter	Episodic	Numeric	X.XX	mmol/L	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_RF_PRE_FILTER_CONC_CITRATE		Replacement Fluid	Pre-Filter RF Citrate Concentration	Citrate concentration in the Pre-Filter RF	1.1.6.25	Intradialytic	Parameter	Episodic	Numeric	XX.X	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_GLU		Replacement Fluid	Pre-Filter RF Glucose Concentration	Glucose concentration in the Pre-Filter RF	1.1.6.26	Intradialytic	Parameter	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_K		Replacement Fluid	Pre-Filter RF Potassium Concentration	Potassium concentration in the Pre-Filter RF	1.1.6.27	Intradialytic	Parameter	Episodic	Numeric	X.X	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_CONC_NA		Replacement Fluid	Pre-Filter RF Sodium Concentration	Sodium concentration in the Pre-Filter RF	1.1.6.28	Intradialytic	Parameter	Episodic	Numeric	XXXX	mmol/L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE		Replacement Fluid	Actual Pre-Filter Replacement Fluid Flow Rate	Rate of pre-filter replacement fluid flow at any given point during a treatment	1.1.6.29	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	1..1	C7
MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE_SETTING		Replacement Fluid	Pre-Filter Replacement Fluid Flow Rate Setting	Rate at which the user programmed the pre-filter replacement fluid flow	1.1.6.30	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL/min	0..1	C8
MDCX_HDIALY_RF_PRE_FILTER_FLOW_RATE_AVG		Replacement Fluid	Average Pre-Filter Replacement Fluid Flow Rate	Average rate of pre-filter replacement fluid flow during course of treatment.	1.1.6.31	Intradialytic	Status	Periodic	Numeric	XXXX	mL/min	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_TEMP		Replacement Fluid	Actual Pre-Filter Replacement Fluid Temperature	Temperature of pre-filter replacement fluid	1.1.6.32	Intradialytic	Status	Periodic	Numeric	XX.X	Cel	1..1	C7
MDCX_HDIALY_RF_PRE_FILTER_TEMP_SETTING		Replacement Fluid	Pre-Filter Replacement Fluid Temperature Setting	Temperature at which pre-filter replacement fluid is to be delivered to patient	1.1.6.33	Intradialytic	Parameter	Episodic	Numeric	XX.X	Cel	1..1	C7
MDCX_HDIALY_RF_PRE_FILTER_VOL		Replacement Fluid	Total Pre-Filter Replacement Fluid Volume	Total amount of pre-filter replacement fluid used during a treatment	1.1.6.34	Intradialytic	Status	Periodic	Numeric	XXX.XX	L	0..1	O
MDCX_HDIALY_RF_PRE_FILTER_VOL_SETTING		Replacement Fluid	Pre-Filter Replacement Fluid Volume Setting	Programmed total amount of pre -filter replacement fluid volume	1.1.6.35	Intradialytic	Parameter	Episodic	Numeric	XXX.XX	L	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_RF_CONV_CLEARANCE		Replacement Fluid	Convective Clearance - Urea	The product of the sieving coefficient and the total water flux. KX = QUF x SX where KX is Convective Clearance for solute "X", QUF is Ultrafiltration flow across membrane and SX is Membrane sieving coefficient for solute "X"	1.1.6.36	Intradialytic	Status	Periodic	Numeric	XXX.X	mL/min	0..1	O
MDCX_EVT_HDIALY_RF_EXCESS_DELIV	tech-high	Replacement Fluid	Excessive Replacement Fluid Bolus Delivered	Notification that more replacement fluid has been delivered than intended	1.1.6.37	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_EVT_HDIALY_RF_INSUFF_DELIV	tech-low	Replacement Fluid	Insufficient Replacement Fluid Bolus Delivered	Notification that less replacement fluid has been delivered than intended	1.1.6.38	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_RF_PREPOST_FLOW_RATIO_SETTING		Replacement Fluid	Pre-Post Replacement Fluid Flow Rate Ratio Setting	The ratio desired between the pre and post Replacement Fluid Flow Rate Settings	1.1.6.39	Intradialytic	Parameter	Episodic	Numeric	XX.X	N/A	0..1	O
MDCX_EVT_HDIALY_RF_BOLUS_RATE	tech high low	Replacement Fluid	Replacement Fluid Bolus Rate	User defined replacement fluid bolus rate.	1.1.6.40	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL/min	0..1	C9
MDCX_HDIALY_RF_BOLUS_VOL_SETTING		Replacement Fluid	Replacement Fluid Bolus Volume Setting	User defined amount of replacement fluid to be delivered in this single large dose	1.1.6.41	Intradialytic	Parameter	Episodic	Numeric	XXX	mL	0..1	C9
MDCX_HDIALY_RF_BOLUS_VOL_DELIVERED	tech high low	Replacement Fluid	Replacement Fluid Bolus Volume Delivered	Amount of replacement fluid delivered in a single large dose.	1.1.6.42	Intradialytic	Parameter	Episodic	Numeric	XXX	mL	0..1	C9
MDCX_HDIALY_RF_FLOW_MODE		Replacement Fluid	Replacement Fluid Flow Mode	The mode of controlling replacement fluid flow rate	1.1.6.43	All	Parameter	Episodic	_TBL_07	Alphanumeric	N/A	1..1	C8
MDCX_HDIALY_RF_SOURCE		Replacement Fluid	Replacement Fluid Source	The source of the replacement fluid.	1.1.6.44	Intradialytic	Status	Episodic	_TBL_09	Alphanumeric	N/A	0..1	O
MDCX_DEV_HDIALY_SAFETY_SYSTEMS_CHAN					1.1.7								
MDCX_EVT_HDIALY_SAFETY_ART_AIR_DETECT	tech	Safety Systems	Arterial Air Detector	Notification that air is sensed in the arterial blood lines	1.1.7.1	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_EVT_HDIALY_SAFETY_DIALYSATE_AIR_DETECT	tech	Safety Systems	Dialysate Air Detector	Notification that air is sensed in dialysate line	1.1.7.2	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	C6
MDCX_EVT_HDIALY_SAFETY_DIALYSATE_COMPOSITION	tech	Safety Systems	Dialysis Composition	Notification that the dialyzing fluid composition is incorrect	1.1.7.3	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_EVT_HDIALY_SAFETY_SYSTEM_GENERAL	tech	Safety Systems	General System	Notification of any dialysis machine malfunction that are not related to patient treatment	1.1.7.4	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_EVT_HDIALY_SAFETY_SELF_TEST	tech	Safety Systems	Self-Test	Notification that a self-test failed	1.1.7.5	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_EVT_HDIALY_SAFETY_VEN_ACCESS	tech	Safety Systems	Venous Access	Notification that the venous access may have become disconnected	1.1.7.6	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_EVT_HDIALY_SAFETY_VEN_AIR_DETECT	tech	Safety Systems	Venous Air Detector	Notification that air is sensed in the venous blood lines	1.1.7.7	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ALERT	tech	Safety Systems	Wetness Detector	Notification that the wetness detector has sensed fluid	1.1.7.8	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ERROR	tech	Safety Systems	Wetness Detector Error	Notification of a detected error on the wetness detector device	1.1.7.9	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_DEV_HDIALY_THERAPY_OUTCOMES_CHAN					1.1.8								
MDCX_HDIALY_THERAPY_MASS_TRF_AREA_COEFF		Therapy Outcomes	KoA	Overall mass transfer coefficient multiplied by surface area of a dialyzer	1.1.8.1	Intradialytic	Status	Episodic	Numeric	XXXX	N/A	0..1	O
MDCX_HDIALY_THERAPY_ACCESS_FLOW		Therapy Outcomes	Access Flow	Flow through a fistula or other conduit connecting an artery to a vein	1.1.8.2	Intradialytic	Status	Periodic	Numeric	±XXXX	mL/min	0..1	O
MDCX_HDIALY_THERAPY_RATIO_EKT_OVER_V_DELIVERED		Therapy Outcomes	Delivered Equilibrated Kt/V	Fractional solute clearance for an intermittent treatment corrected to take account of solute redistribution following the treatment	1.1.8.4	Intradialytic	Status	Episodic	Numeric	X.XX	%	0..1	O
MDCX_HDIALY_THERAPY_KT_DELIVERED		Therapy Outcomes	Delivered Kt	Volume of fluid cleared of urea by dialysis treatment	1.1.8.5	Intradialytic	Status	Episodic	Numeric	X.XX	L	0..1	O
MDCX_HDIALY_THERAPY_RATIO_SPKT_OVER_V_DELIVERED		Therapy Outcomes	Delivered Single Pool Kt/V	Fractional solute clearance, calculated on the basis of solute concentration before and after an intermittent treatment, and on the basis of effluent dialysate volume and anthropometric estimate of volume of urea distribution for peritoneal dialysis	1.1.8.6	Intradialytic	Status	Episodic	Numeric	X.XX	%	0..1	O
MDCX_EVT_HDIALY_THERAPY_PAT_TX_GENERAL	phys	Therapy Outcomes	General Patient Treatment	Notification that the machine has detected an alarm condition not covered by other alarms that is related to patient treatment.	1.1.8.7	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_HDIALY_THERAPY_RATIO_KT_OVER_V_GOAL		Therapy Outcomes	Kt/V Goal	Target dialysis treatment adequacy measurement	1.1.8.10	Intradialytic	Parameter	Episodic	Numeric	X.XX	%	0..1	O
MDCX_HDIALY_THERAPY_MEAN_UREA_CLEARANCE		Therapy Outcomes	Mean Urea Clearance	Time weighted average of the individual effective conductivity clearance measurements	1.1.8.11	Intradialytic	Status	Episodic	Numeric	XXX	N/A	0..1	O
MDCX_HDIALY_THERAPY_BODY_START_WT		Therapy Outcomes	Patient Start Weight	Weight of patient at beginning of treatment	1.1.8.12	All	Parameter	Episodic	Numeric	XXX.X	Kg	0..1	O
MDCX_HDIALY_THERAPY_PCT_RECIRC		Therapy Outcomes	Percent Recirculation	Percentage of the blood flow in the access line that is coming directly from the venous line without going through the patient's body.	1.1.8.13	Intradialytic	Status	Periodic	Numeric	±XXX.X	%	0..1	O
MDCX_HDIALY_THERAPY_PLASMA_NA_CONC		Therapy Outcomes	Plasma Sodium Concentration	Sodium ion concentration measured in patient plasma.	1.1.8.14	Intradialytic	Status	Episodic	Numeric	XXX	mmol/L	0..1	O
MDCX_HDIALY_THERAPY_RATIO_SPKT_OVER_V_PROJECTED		Therapy Outcomes	Projected Single Pool Kt/V	Fractional solute clearance, calculated on the basis of solute concentration before and after an intermittent treatment, and on the basis of effluent dialysate volume and anthropometric estimate of volume of urea distribution for peritoneal dialysis, projected on the basis of data available before the treatment occurs or before the treatment is completed	1.1.8.16	Intradialytic	Status	Episodic	Numeric	X.XX	%	0..1	O
MDCX_HDIALY_THERAPY_COMPLETE_METHOD		Therapy Outcomes	Treatment Completion Method	Parameter used to determine when treatment has been complete	1.1.8.18	Intradialytic	Status	Episodic	_TBL_12	Alphanumeric	N/A	1..1	M
MDCX_EVT_HDIALY_THERAPY_TX_END_TIME	tech	Therapy Outcomes	Treatment Time Ended	Notification that target time of treatment has been met	1.1.8.19	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_THERAPY_UREA_CLEARANCE		Therapy Outcomes	Urea Clearance	Calculated clearance based on the change in conductivity of the pre-dialyzer vs post dialyzer dialysate.	1.1.8.20	Intradialytic	Status	Episodic	Numeric	XXX	N/A	0..1	O
MDCX_HDIALY_THERAPY_BODY_END_WT		Therapy Outcomes	Patient End Weight	Weight of patient at end of treatment	1.1.8.21	All	Parameter	Episodic	Numeric	XXX.X	Kg	0..1	O
MDCX_DEV_HDIALY_UF_CHAN					1.1.9								
MDCX_EVT_HDIALY_UF_LO	tech-low	Ultrafiltration System	Insufficient Ultrafiltration	Notification that the ultrafiltration system is not removing as much fluid as intended	1.1.9.1	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDCX_EVT_HDIALY_UF_NEG	tech-low	Ultrafiltration System	Negative Ultrafiltration	Notification that the ultrafiltration system has resulted in a negative ultrafiltration rate	1.1.9.2	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_UF_TIME_TO_TARGET		Ultrafiltration System	Remaining Ultrafiltration Time	Time remaining until target volume to remove will have been achieved.	1.1.9.3	Intradialytic	Status	Periodic	Numeric	XXXX	min	0..1	O
MDCX_HDIALY_UF_TARGET_VOL_TO_REMOVE		Ultrafiltration System	Target Volume To Remove	Planned fluid removal	1.1.9.4	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL	1..1	M
MDCX_HDIALY_UF_ACTUAL_REMOVED_VOL		Ultrafiltration System	Total Fluid Volume Removed	Measured fluid removed from patient	1.1.9.5	Intradialytic	Status	Periodic	Numeric	XXXX	mL	1..1	M
MDCX_HDIALY_UF_RATE_MAX_NON_ALARM_SETTING		Ultrafiltration System	Ultrafiltration Flow Rate Maximum Limit	Maximum ultrafiltration rate allowed without an alarm generation	1.1.9.6	Intradialytic	Parameter	Episodic	Numeric	XXXX	mL/h	0..1	O
MDCX_EVT_HDIALY_UF_GOAL_MET	tech	Ultrafiltration System	Ultrafiltration Goal Met	Notification that target ultrafiltration volume has been met	1.1.9.7	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_HDIALY_UF_MODE		Ultrafiltration System	Ultrafiltration Mode	Process used to remove fluid volume from the patient	1.1.9.8	Intradialytic	Parameter	Episodic	_TBL_13	Alphanumeric	N/A	1..1	M
MDCX_HDIALY_UF_RATE		Ultrafiltration System	Ultrafiltration Rate	Rate at which fluid is removed from the patient	1.1.9.9	Intradialytic	Status	Periodic	Numeric	XXXX	mL/h	1..1	M
MDCX_EVT_HDIALY_UF_RATE_RANGE	tech-range	Ultrafiltration System	Ultrafiltration Rate Out of Range	Notification that the ultrafiltration rate has deviated from the intended ultrafiltration rate more than is safe	1.1.9.10	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	M
MDC_DEV_PRESS_BLD_NONINVASIVE_VMD				Non-Invasive Blood Pressure Monitor	1.2								
MDC_DEV_PRESS_BLD_NONINVASIVE_CHAN					1.2.1								
MDC_ATTR_TIME_PD_MS_MT		Blood Pressure Monitor	Blood Pressure Measurement Interval Setting	Time allowed to elapse between automatic blood pressure measurements. A value of zero indicates that automatic measurements are not being taken.	1.2.1.1	All	Parameter	Episodic	Numeric	XXX	min	0..1	O
MDC_ATTR_BODY_POSN		Blood Pressure Monitor	Blood Pressure Measurement Position	Patient position at time of blood pressure measurement	1.2.1.2	All	Blood Pressure	Episodic	_TBL_04	Alphanumeric	N/A	0..1	O
MDC_PRESS_BLD_NONINVASIVE_DIA	phys high low thr	Blood Pressure Monitor	Diastolic Pressure	Minimum arterial pressure during relaxation and dilatation of the ventricles of the heart when the ventricles fill with blood	1.2.1.3	All	Blood Pressure	Episodic	Numeric	XXX	mm[Hg]	0..1	O

Dialysis Machine Implementation Guide

Table 1 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Observation Sub ID	Phase	Message	Temporal	Data Type	Format	UCUM UOM	CARD	USE
MDC_PULS_RATE_NON_INV	phys high low thr	Blood Pressure Monitor	Heart Rate	Number of heart beats per minute	1.2.1.4	All	Blood Pressure	Episodic	Numeric	XXX	{beats}/min	0..1	O
MDC_PRESS_BLD_NONINV_MEAN		Blood Pressure Monitor	Mean Arterial Pressure	Average pressure in a patient's arteries during one cardiac cycle	1.2.1.5	All	Blood Pressure	Episodic	Numeric	XXX	mm[Hg]	0..1	O
MDC_PRESS_BLD_NONINV_SYS	phys high low thr	Blood Pressure Monitor	Systolic Pressure	Maximum arterial pressure during contraction of the left ventricle of the heart	1.2.1.6	All	Blood Pressure	Episodic	Numeric	XXX	mm[Hg]	0..1	O
MDC_DEV_ANALY_SAT_O2_VMD					1.3								
MDC_DEV_ANALY_SAT_O2_CHAN					1.3.1								
MDC_PULS_OXIM_SAT_O2	phys high low thr	Pulse Oximeter Monitor	SpO2 Oxygen Saturation	Hemoglobin oxygen saturation in patient blood measured by pulse oximetry	1.3.1.1	Intradialytic	SpO2	Periodic	Numeric	XXX	%	0..1	O
MDC_PULS_OXIM_PULS_RATE	phys high low thr	Pulse Oximeter Monitor	SpO2 Pulse Rate	Number of heart beats per minute measured by pulse oximetry	1.3.1.2	Intradialytic	SpO2 Pulse Rate	Periodic	Numeric	XXX	{beats}/min	0..1	O
MDCX_EVT_PULS_OXIM_ERROR	tech	Safety Systems	Pulse Oximeter Error	Notification of a detected error on the pulse oximeter device	1.3.1.3	Intradialytic	Alert	Episodic	Bool	T / F	N/A	0..1	O
MDCX_DEV_BLOOD_CHEM_VMD					1.4								
MDCX_DEV_BLOOD_CHEM_CHAN				Blood Chemistry Channel	1.4.1								
MDCX_HDIALY_PLASMA_VOLUME_MARKER		Therapy Outcomes	Plasma Volume Marker	An indicator applied to the blood volume percent reading to designate an event	1.4.1.1	Intradialytic	Status	Periodic	_TBL_10	Alphanumeric	N/A	0..1	O
MDCX_HDIALY_PLASMA_VOLUME_PROFILE		Therapy Outcomes	Plasma Volume Profile	The rate of change in relative blood volume.	1.4.1.2	Intradialytic	Status	Periodic	_TBL_11	Alphanumeric	N/A	0..1	O
MDC_CONC_HCT	tech low	Therapy Outcomes	Hematocrit	100*(sample red blood cell volume/blood sample volume)	1.4.1.3	Intradialytic	Status	Periodic	Numeric	XX.X	{vol}	0..1	O
MDC_CONC_HB		Therapy Outcomes	Hemoglobin	Concentration of hemoglobin in a sample expressed as mass/volume	1.4.1.4	Intradialytic	Status	Periodic	Numeric	±XXX.X	g/dL	0..1	O
MDCX_HDIALY_REL_BLOOD_VOLUME	tech- low	Therapy Outcomes	Relative Blood Volume	The percent change in blood volume since the start of the treatment. The BV change value is calculated from the following equation: $\Delta BV\% = [(H0/H1) - 1] \times 100$ Where H0 = initial Hct or Hb H1 = current Hct or Hb	1.4.1.5	Intradialytic	Alert	Episodic	Numeric	XXX.X	%	0..1	O
MDC_SAT_O2		Therapy Outcomes	SpO2 Oxygen Saturation	Hemoglobin oxygen saturation in patient blood	1.4.1.6	Intradialytic	Status	Periodic	Numeric	XXX	%	0..1	O

7.3 Alarms/Alerts

Table 2 – Dialysis Machine Alarms/Alerts

Source Identifier	Event Identifier	Alert Type	Dialysis Common Name	Definition	Use
MDCX_DEV_HDIALY_MACHINE_MDS					
MDCX_DEV_HDIALY_VMD					
MDCX_DEV_HDIALY_MACH_CONFIG_CHAN					
MDCX_DEV_HDIALY_ANTICOAG_PUMP_CHAN					
MDCX_HDIALY_ANTICOAG_PUMP_CHAN	MDCX_EVT_HDIALY_ANTICOAG_STOP	tech	Anticoagulant Pump Stop	Notification that the anticoagulant pump has stopped	O
MDCX_HDIALY_ANTICOAG_PUMP_CHAN	MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_EMPTY	tech	Syringe Empty	Notification that the syringe plunger has reached the end of its functional range	O
MDCX_HDIALY_ANTICOAG_PUMP_CHAN	MDCX_EVT_HDIALY_ANTICOAG_SYRINGE_SIZE	tech	Wrong Syringe Size	Notification that the size of the syringe is not the size defined for anticoagulant administration	O
MDCX_DEV_HDIALY_BLOOD_PUMP_CHAN					
MDCX_HDIALY_BLD_PRESS_ART	MDC_EVT_HI	tech high thr	Arterial Pressure High	Notification that the arterial pressure has exceeded maximum setting	M
MDCX_HDIALY_BLD_PRESS_ART	MDC_EVT_LO	tech low thr	Arterial Pressure Low	Notification that the arterial pressure has fallen below the minimum setting	M
MDCX_HDIALY_BLOOD_PUMP_CHAN	MDCX_EVT_HDIALY_BLD_PUMP_STOP	tech	Blood Pump Stop	Notification that the blood pump has stopped	M
MDCX_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP	MDC_EVT_HI	tech high thr	Post Pump Arterial Pressure High	Notification that the post-pump arterial pressure is above limit.	O
MDCX_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP	MDC_EVT_LO	tech low thr	Post Pump Arterial Pressure Low	Notification that the post-pump arterial pressure is below limit.	O
MDCX_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS	MDC_EVT_HI	tech high thr	Needle Pressure High	Notification that the needle pressure is above limit.	O
MDCX_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS	MDC_EVT_LO	tech low thr	Needle Pressure Low	Notification that the needle pressure is below limit.	O
MDCX_HDIALY_BLD_PUMP_PRESS_VEN	MDC_EVT_HI	tech high thr	Venous Pressure High	Notification that the venous pressure has exceeded maximum setting	M
MDCX_HDIALY_BLD_PUMP_PRESS_VEN	MDC_EVT_LO	tech low thr	Venous Pressure Low	Notification that the venous pressure has fallen below the minimum setting	M
MDCX_DEV_HDIALY_FLUID_CHAN					
MDCX_HDIALY_DIALYSATE_COND	MDC_EVT_HI	tech high	Conductivity High	Notification that the conductivity of the dialysate has exceeded the maximum setting	O
MDCX_HDIALY_DIALYSATE_COND	MDC_EVT_LO	tech low	Conductivity Low	Notification that the conductivity of the dialysate has fallen below the minimum setting	O
MDCX_HDIALY_DIALYSATE_FLOW_RATE	MDC_EVT_LO	tech low	Dialysate Flow Low	Notification that the dialysate flow is below the dialysate flow rate minimum setting	O
MDCX_HDIALY_DIALYSATE_TEMP	MDC_EVT_LO	tech low	Dialysate Temperature Low	Notification that the dialysate temperate has fallen below the minimum temperate setting	C (if mode is HD or HDF)
MDCX_HDIALY_DIALYSATE_TEMP	MDC_EVT_HI	tech high	Dialysate Temperature High	Notification that the dialysate temperate exceeds the maximum temperate setting	C (if mode is HD or HDF)

Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Alarms/Alerts

Source Identifier	Event Identifier	Alert Type	Dialysis Common Name	Definition	Use
MDCX_HDIALY_FLUID_CHAN	MDCX_EVT_HDIALY_BLOOD_LEAK	tech	Blood Leak	Notification that blood has been detected in the dialysate fluid	M
MDCX_DEV_HDIALY_FILTER_CHAN					
MDCX_HDIALY_FILTER_TRANSMEMBRANE_PRESS	MDC_EVT_LO	tech low	Transmembrane Pressure Low	Notification that transmembrane pressure has fallen below the minimum setting	M
MDCX_HDIALY_FILTER_TRANSMEMBRANE_PRESS	MDC_EVT_HI	tech high	Transmembrane Pressure High	Notification that transmembrane pressure has exceeded the maximum setting	M
MDCX_DEV_HDIALY_CONVECTIVE_CHAN					
MDCX_HDIALY_CONVECTIVE_CHAN	MDCX_EVT_HDIALY_RF_EXCESS_DELIV	tech-high	Excessive Replacement Fluid Bolus Delivered	Notification that more replacement fluid has been delivered than intended	O
MDCX_HDIALY_CONVECTIVE_CHAN	MDCX_EVT_HDIALY_REPLACE_FLUID_IN_SUFF_DELIV	tech-low	Insufficient Replacement Fluid Bolus Delivered	Notification that less replacement fluid has been delivered than intended	O
MDCX_HDIALY_RF_BOLUS_RATE	MDC_EVT_HI	tech high	Replacement Fluid Bolus Rate High	Notification that the replacement fluid bolus is delivering more fluid than intended	O
MDCX_HDIALY_RF_BOLUS_RATE	MDC_EVT_LO	tech low	Replacement Fluid Bolus Rate Low	Notification that the replacement fluid bolus is delivering less than the intended.	O
MDCX_HDIALY_RF_BOLUS_VOL_DELIVERED	MDC_EVT_HI	tech high	Replacement Fluid Bolus Volume High	Notification that the replacement fluid bolus delivered more fluid than intended	O
MDCX_HDIALY_RF_BOLUS_VOL_DELIVERED	MDC_EVT_LO	tech low	Replacement Fluid Bolus Volume Low	Notification that the replacement fluid bolus delivered less fluid than intended.	O
MDCX_DEV_HDIALY_SAFETY_SYSTEMS_CHAN					
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_ART_AIR_DETECT	tech	Arterial Air Detector	Notification that air is sensed in the arterial blood lines	M
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_DIALYSATE_AIR_DETECT	tech	Dialysate Air Detector	Notification that air is sensed in dialysate line	C (if mode is HD or HDF)
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_DIALYSATE_COMPOSITION	tech	Dialysis Composition	Notification that the dialyzing fluid composition is incorrect	O
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_SYSTEM_GENERAL	tech	General System	Notification of any dialysis machine malfunction that are not related to patient treatment	M
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_SELF_TEST	tech	Self-Test	Notification that a self-test failed	M
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_VEN_ACCESS	tech	Venous Access	Notification that the venous access may have become disconnected	O
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_VEN_AIR_DETECT	tech	Venous Air Detector	Notification that air is sensed in the venous blood lines	M
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ALERT	tech	Wetness Detector	Notification that the wetness detector has sensed fluid	O
MDCX_HDIALY_SAFETY_SYSTEMS_CHAN	MDCX_EVT_HDIALY_SAFETY_WETNESS_DETECT_ERROR	tech	Wetness Detector Error	Notification of a detected error on the wetness detector device	O
MDCX_DEV_HDIALY_THERAPY_OUTCOMES_CHAN					
MDCX_HDIALY_THERAPY_OUTCOMES_CHAN	MDCX_EVT_HDIALY_THERAPY_PAT_TX_GENERAL	phys	General Patient Treatment	Notification that the machine has detected an alarm condition not covered by other alarms that is related to patient treatment.	O
MDCX_HDIALY_THERAPY_OUTCOMES_CHAN	MDCX_EVT_HDIALY_THERAPY_TX_END_TIME	tech	Treatment Time Ended	Notification that target time of treatment has been met	O

Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Alarms/Alerts

Source Identifier	Event Identifier	Alert Type	Dialysis Common Name	Definition	Use
MDCX_DEV_HDIALY_UF_CHAN					
MDCX_HDIALY_UF_CHAN	MDCX_EVT_HDIALY_UF_LO	tech low	Insufficient Ultrafiltration	Notification that the ultrafiltration system is not removing as much fluid as intended	O
MDCX_HDIALY_UF_CHAN	MDCX_EVT_HDIALY_UF_NEG	tech low	Negative Ultrafiltration	Notification that the ultrafiltration system has resulted in a negative ultrafiltration rate	O
MDCX_HDIALY_UF_CHAN	MDCX_EVT_HDIALY_UF_GOAL_MET	tech	Ultrafiltration Goal Met	Notification that target ultrafiltration volume has been met	O
MDCX_HDIALY_UF_CHAN	MDCX_EVT_HDIALY_UF_RATE_RANGE	tech range	Ultrafiltration Rate Out of Range	Notification that the ultrafiltration rate has deviated from the intended ultrafiltration rate more than is safe	M
MDC_DEV_PRESS_BLD_NONINV_VMD					
MDC_DEV_PRESS_BLD_NONINV_CHAN					
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_HI	phys high thr	Diastolic High	Notification that diastolic pressure exceeded the diastolic maximum setting	O
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_LOW	phys low thr	Diastolic Low	Notification that diastolic pressure fell below the diastolic minimum setting	O
MDC_PULS_RATE_NON_INV	MDC_EVT_HI	phys high thr	Heart Rate High	Heart rate exceeded heart rate maximum setting	O
MDC_PULS_RATE_NON_INV	MDC_EVT_LOW	phys low thr	Heart Rate Low	Heart rate fell below heart rate minimum setting	O
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_HI	phys high thr	Systolic High	Notification that systolic pressure exceeded the systolic maximum setting	O
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_LOW	phys low thr	Systolic Low	Notification that systolic pressure fell below the systolic minimum setting	O
MDC_DEV_ANALY_SAT_O2_VMD					
MDC_DEV_ANALY_SAT_O2_CHAN					
MDC_PULS_OXIM_SAT_O2	MDC_EVT_HI	phys high thr	Oximeter Oxygen Saturation High	Notification that Oxygen Saturation is above limit.	O
MDC_PULS_OXIM_SAT_O2	MDC_EVT_LO	phys low thr	Oximeter Oxygen Saturation Low	Notification that Oxygen Saturation is below limit.	
MDC_PULS_OXIM_PULS_RATE	MDC_EVT_HI	phys high thr	Oximeter Pulse High	Notification that the pulse rate measure by the pulse oximeter is above limit.	O
MDC_PULS_OXIM_PULS_RATE	MDC_EVT_LO	phys low thr	Oximeter Pulse Low	Notification that the pulse rate measure by the pulse oximeter is below limit.	
MDC_DEV_ANALY_SAT_O2_CHAN	MDCX_EVT_PULS_OXIM_ERROR	tech	Pulse Oximeter Error	Notification of a detected error on the pulse oximeter device	O
MDCX_DEV_BLOOD_CHEM_VMD					
MDCX_DEV_BLOOD_CHEM_CHAN					
MDC_CONC_HCT	MDC_EVT_HI	phys high thr	Hematocrit High	Notification that patient hematocrit has risen above limit.	O
MDC_CONC_HCT	MDC_EVT_LO	phys low thr	Hematocrit Low	Notification that patient hematocrit has fallen below limit.	O
MDC_CONC_HB	MDC_EVT_HI	phys high thr	Hemoglobin High	Notification that patient hemoglobin has risen above limit.	O
MDC_CONC_HB	MDC_EVT_LO	phys low thr	Hemoglobin Low	Notification that patient hemoglobin has fallen below limit.	O
MDCX_HDIALY_REL_BLOOD_VOL	MDC_EVT_LO	phys low thr	Blood Volume	Notification that the relative blood volume has dropped below the defined setting	O
MDC_SAT_O2	MDC_EVT_HI	phys high thr	Oxygen Saturation High	Notification that Patient Oxygen Saturation has risen above limit.	O
MDC_SAT_O2	MDC_EVT_LO	phys low thr	Oxygen Saturation Low	Notification that Patient Oxygen Saturation has fallen below limit.	O