

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12

# Dialysis Machine HL7 Implementation Guidelines

## 13 Copyright Licenses – HL7®

14 HL7 licenses its standards and select IP free of charge. If you did not acquire a free license from HL7 for this  
15 document, you are not authorized to access or make any use of it. To obtain a free license, please visit  
16 <http://www.HL7.org/implement/standards/index.cfm>.

17 If you are the individual that obtained the license for this HL7 Standard, specification or other freely licensed work  
18 (in each and every instance "Specified Material"), the following describes the permitted uses of the Material.

19 **A. HL7 INDIVIDUAL, STUDENT AND HEALTH PROFESSIONAL MEMBERS**, who register and agree to the terms  
20 of HL7's license, are authorized, without additional charge, to read, and to use Specified Material to develop and  
21 sell products and services that implement, but do not directly incorporate, the Specified Material in whole or in part  
22 without paying license fees to HL7.

23 **INDIVIDUAL, STUDENT AND HEALTH PROFESSIONAL MEMBERS** wishing to incorporate additional items of  
24 Special Material in whole or part, into products and services, or to enjoy additional authorizations granted to HL7  
25 ORGANIZATIONAL MEMBERS as noted below, must become ORGANIZATIONAL MEMBERS of HL7.

26 **B. HL7 ORGANIZATION MEMBERS**, who register and agree to the terms of HL7's License, are authorized,  
27 without additional charge, on a perpetual (except as provided for in the full license terms governing the Material),  
28 non-exclusive and worldwide basis, the right to (a) download, copy (for internal purposes only) and share this  
29 Material with your employees and consultants for study purposes, and (b) utilize the Material for the purpose of  
30 developing, making, having made, using, marketing, importing, offering to sell or license, and selling or licensing,  
31 and to otherwise distribute, Compliant Products, in all cases subject to the conditions set forth in this Agreement  
32 and any relevant patent and other intellectual property rights of third parties (which may include members of HL7).  
33 No other license, sublicense, or other rights of any kind are granted under this Agreement.

34 **C. NON-MEMBERS**, who register and agree to the terms of HL7's IP policy for Specified Material, are authorized,  
35 without additional charge, to read and use the Specified Material for evaluating whether to implement, or in  
36 implementing, the Specified Material, and to use Specified Material to develop and sell products and services that  
37 implement, but do not directly incorporate, the Specified Material in whole or in part.

38 **NON-MEMBERS** wishing to incorporate additional items of Specified Material in whole or part, into products  
39 and services, or to enjoy the additional authorizations granted to HL7 ORGANIZATIONAL MEMBERS, as noted  
40 above, must become ORGANIZATIONAL MEMBERS of HL7.

41 Please see <http://www.HL7.org/legal/ippolicy.cfm> for the full license terms governing the Material.  
42

## 43 Copyright Licenses - IHE

44 IHE International hereby grants to each Member Organization, and to any other user of these documents, an  
45 irrevocable, worldwide, perpetual, royalty-free, nontransferable, nonexclusive, non-sublicensable license under its  
46 copyrights in any IHE profiles and Technical Framework documents, as well as any additional copyrighted materials  
47 that will be owned by IHE International and will be made available for use by Member Organizations, to reproduce  
48 and distribute (in any and all print, electronic or other means of reproduction, storage or transmission) such IHE  
49 Technical Documents.

50 The licenses covered by this Copyright License are only to those copyrights owned or controlled by IHE International  
51 itself. If parts of the Technical Framework are included in products that also include materials owned or controlled  
52 by other parties, licenses to use those products are beyond the scope of this IHE document and would have to be  
53 obtained from that other party.

## 55 Copyright Licenses - IEEE

56 IEEE holds copyright in the standards referenced. IEEE, as part of its support of the RTMMS database and on-going,  
57 royalty-free agreement with the NIST, makes these terms available for the development of IEEE11073 compliant  
58 products and supporting material (e.g. in user documentation, collateral, etc.). Any use of IEEE terms beyond  
59 compliant products and support material may require prior approval from IEEE. Please notify IEEE of any request to  
60 use, modify, or reproduce these terms in any manner beyond the permitted use described above. To request

61 permission, please submit your request to [stds-ipr@ieee.org](mailto:stds-ipr@ieee.org).

62 The following information is provided free of charge to all NIST RTMMS users via the IEEE-SA and NIST  
63 Royalty Free Agreement:

64

- 65 • Reference ID ('REFID' in RTMMS)
- 66 • Terminology Code ('CODE10', 'CF\_CODE10', 'UCODE10', 'CF\_UCODE10', 'ECODE10', and  
67 'CF\_ECODE10' in RTMMS)
- 68 • Description ('Term Description' in RTMMS)
- 69 • Systematic Name ('Systematic Name' in RTMMS)
- 70 • Common Term ('Common term' in RTMMS)

71

## 72 Copyright of Base Standards

73 IHE technical documents refer to and make use of a number of standards developed and published by several  
74 standards development organizations. All rights for their respective base standards are reserved by these  
75 organizations. This agreement does not supersede any copyright provisions applicable to such base standards.

76 Health Level Seven<sup>®</sup> has granted permission to IHE to reproduce tables from the HL7 standard. The HL7 tables in  
77 this document are copyrighted by Health Level Seven. All rights reserved. Material drawn from these documents is  
78 credited where used.

79

## 80 Trademark

81 IHE<sup>®</sup> and the IHE logo are trademarks of the Healthcare Information Management Systems Society in the United  
82 States and trademarks of IHE Europe in the European Community. They may only be used with the written consent  
83 of the IHE International Board Operations Committee, which may be given to a Member Organization in broad  
84 terms for any use that is consistent with the IHE mission and operating principles.

85

## 86 Disclaimer Regarding Patent Rights

87 Attention is called to the possibility that implementation of the specifications in this document may require use of  
88 subject matter covered by patent rights. By publication of this document, no position is taken with respect to the  
89 existence or validity of any patent rights in connection therewith.

90

## Revision History

Version	Change Description	Date
1.0	Initial Release	January 2020
2.0	<p>Added appendix summarizing HL7 Batch Protocol for capturing an entire treatment</p> <p>Added support for multiple filters during a treatment</p> <p>Added description of IEEE private terms for supporting manufacturer specific data.</p> <p>Added support for open-ended treatments that have to specific terminating conditions.</p> <p>Clarified the support for two types of reported alarms/alerts. The guide now supports the IEEE events and the original proposed Boolean values.</p>	October 2020
3.0	<p>Changed "MDCX" to "MDC".</p> <p>Added IEEE Reference IDs</p>	January 2021
3.1	<p>Expanded Patient Identification Section to include IHE PDQ Transaction.</p> <p>Added Prescription Transfer Section.</p> <p>Added Prescription Objects</p> <p>Added Usage Column for all objects.</p>	January 2023
4.0	<p>Added Message Transport section.</p> <p>Added note to alarm limit range report that limits can be fixed or machine generated.</p>	March 2023
5.0	Added support for Peritoneal Dialysis.	January 2025

91

92

93

94

95	<b>Table of Contents</b>		
96	1	Introduction & Scope	7
97	1.1	Purpose	7
98	1.2	Audience	7
99	1.3	Scope	7
100	1.4	Assumptions	7
101	1.5	Conventions	7
102	1.5.1	Other Standards	7
103	1.5.2	Usage/Cardinality	8
104	1.5.3	Object Reference Identifiers	9
105	1.6	Referenced Standards	9
106	1.7	Glossary	9
107	2	Time Synchronization	10
108	3	Message Transport	10
109	4	Patient Identification	10
110	4.1	Overview	10
111	4.2	Message Structure	12
112	4.2.1	Patient Demographics Query	12
113	4.2.2	Patient Demographics Response	12
114	4.3	Sample Messages	13
115	4.3.1	Example 1 – Query using Medical Record Number	13
116	4.3.2	Example 2 – Query Response with No Matches	13
117	4.3.3	Example 3 – Query using First Name and Last Name	13
118	4.3.4	Example 4 – Query Response with Multiple Matches	13
119	4.3.5	Example 5 – Query using Person Number	13
120	5	Prescription Transfer	14
121	5.1	Overview	14
122	5.2	Message Structure	14
123	5.2.1	Dialysis Prescription Query	14
124	5.2.2	Dialysis Prescription Response	15
125	5.3	Profile Representation	15
126	5.3.1	Vendor Profile	17
127	5.3.2	Constant Profile	18
128	5.3.3	Linear Profile	18
129	5.3.4	Exponential Profile	19
130	5.3.5	Step Profile	19
131	5.4	PD Prescription Representation	20
132	5.5	Sample Messages	21
133	5.5.1	Example 1 – HD Prescription Request	21
134	5.5.2	Example 2 – HD Prescription Response	21
135	5.5.3	Example 2 – HD Prescription Response for HF Therapy with a Linear UF Profile	22
136	5.5.4	Example 3 – Prescription Response returning No Prescription	23
137	5.5.5	Example 4 – PD Prescription Request	23
138	5.5.6	Example 5 – PD Prescription Response, Exchange Based	23
139	5.5.7	Example 6 – PD Prescription Response, Treatment Based	24
140	6	Reporting Treatment Information	25
141	6.1	Overview	25
142	6.2	Device Request	26
143	6.2.1	Message Structure	26
144	6.2.2	Observation Results Structure	27

145	6.2.3	Event Status	27
146	6.2.4	Example 1 – Minimal Message when Idle	29
147	6.2.5	Example 2 – Minimal HD Message when Treating	29
148	6.2.6	Example 3 – Full HD Message	30
149	6.2.7	Example 4 – PD Message	35
150	6.3	EMR Response	37
151	6.3.1	Message Structure	37
152	6.3.2	Example 1 – Accepted Response	37
153	6.3.3	Example 2 – Rejected Response	37
154	7	Reporting Alarm Information	37
155	7.1	Overview	37
156	7.2	Device Request	38
157	7.2.1	Message Structure	39
158	7.2.2	Example 1 - Alarm Initiation	46
159	7.2.3	Example 2 - Alarm Cessation	46
160	7.2.4	Example 3 - Alarm Mute	46
161	7.2.5	Example 4 – Vendor Specific Alarm	47
162	7.3	EMR Response	47
163	7.3.1	Message Structure	48
164	7.3.2	Example 1 – Accepted Response	48
165	8	HL7 Data Elements	48
166	8.1	Tables	48
167	8.2	Segments	51
168	8.2.1	MSA – Message Acknowledgement Segment	51
169	8.2.2	MSH – Message Header Segment	52
170	8.2.3	ORC – Order Common	54
171	8.2.4	OBR– Observation Request Segment	55
172	8.2.5	OBX– Observation Results Segment	58
173	8.2.6	PID – Patient Identification Segment	61
174	8.2.7	PV1 – Patient Visit Segment	62
175	8.2.8	QAK – Query Acknowledge	64
176	8.2.9	QPD – Query Parameter Definition	65
177	8.2.10	RCP – Response Control Parameters	65
178	9	Hemodialysis Data Elements	67
179	9.1	Tables	67
180	9.2	Data Objects	73
181	9.3	Events	92
182	10	Peritoneal Dialysis Data Elements	95
183	10.1	Tables	95
184	10.2	Data Object	98
185	10.3	Events	108
186	Appendix A – HL7 Batch Protocol		113
187	A.1	Segments	113
188	A.1.1	BHS – Batch Header Segment	113
189	A.1.2	BTS – Batch Trailer Segment	113
190	A.1.3	FHS – File Header Segment	114
191	A.1.4	FTS – File Trailer Segment	114
192			
193			

# 194 1 Introduction & Scope

## 195 1.1 Purpose

196 The Dialysis Machine Implementation Guide is a messaging specification intended to standardize the electronic  
197 reporting of patient, clinical, and machine data to electronic health record (EHR) systems and/or electronic  
198 medical record (EMR) systems found in most care settings today. The goal of this specification is to provide an  
199 applicable set of data definitions and structures that can be adopted as a worldwide industry standard, thereby  
200 obviating the need for proprietary solutions for dialysis healthcare interoperability.

## 201 1.2 Audience

202 The standards defined in this guide are intended for use by data exchange systems found in Renal Care groups  
203 listed below:

- 204 • Renal Care Providers
- 205 • Dialysis System Manufacturers
- 206 • Other groups or agencies that employ a dialysis machine to EHR/EMR interface

## 207 1.3 Scope

208 This specification covers the reporting of acute and chronic hemodialysis and peritoneal dialysis treatments to  
209 the EMR and the receipt of dialysis prescription information from the EMR.

## 210 1.4 Assumptions

211 This document assumes that the user is familiar with the HL7 v2.6 standard as well as the IHE PCD and IHE ITI  
212 specifications.

## 213 1.5 Conventions

### 214 1.5.1 Other Standards

215 This specification is primarily based on the Integrating the Healthcare Enterprise (IHE), Patient Care Device  
216 (PCD) technical framework and IT Infrastructure (ITI) Technical Framework, which is based on a profile as  
217 defined in the Health Level 7 (HL7)'s v2.6 Normative Standard. A message profile is an unambiguous  
218 specification of an HL7 message type intended for a particular use case, as defined in Section 2.B of the HL7  
219 v2.6 Standard. Additionally, the message structures and terminology are based on the ISO/IEEE 11073  
220 Standard for Medical Device Communication.

221 An HL7 message profile defines both the *dynamic* aspects of information interchange (i.e., the systems that  
222 participate in such interchanges and the real-world events that trigger messaging) as well as the *static*  
223 aspects of messaging (i.e., the structure and contents of the electronic messages that are exchanged).

224 The IHE Patient Care Device Technical Framework uses an information model and a nomenclature from the  
225 IEEE 11073. The information model is defined in ISO/IEEE 11073- 10201 Health Informatics – Point-of-care  
226 medical device communication – Part 10201: Domain Information Model. The nomenclature is defined in  
227 ISO/IEEE 11073-10101 Health Informatics – Point -of-care medical device communication – Part 10101:  
228 Nomenclature. Familiarity with these standards is necessary for implementers of the Device Observation  
229 Reporter and Device Observation Consumer Actors.

230 HL7 V2.6 Chapter 7 Observation Reporting defines the general HL7 syntax and coding requirements related  
231 to observation reporting, used for PCD data communications in the PCD TF. Familiarity with HL7 Chapter 7  
232 is necessary for implementers of the PCD TF transactions.

233 This PCD Technical Framework specifies conventions that are used to represent the information model  
234 hierarchy for medical devices embodied in the IEEE 11073 Domain Information Model within the syntactic

235 and semantic conventions of HL7 v. 2.6

236 The values reported in the observations should be the value as shown to the user on the device. The  
 237 purpose of this is to create a run sheet that could be referenced back to events as seen by the users. As  
 238 such, sampling rates, average methods, or other mathematical manipulations of signals is beyond the scope  
 239 of this standard. Likewise, logic for initiating or terminating events is not defined in this standard. Since  
 240 these items may vary between manufacturers or machines, if this information is relevant to understanding  
 241 the reported data, it should be disclosed in the manufacturer's technical documentation.

242 The IHE IT Integration Technical Framework identifies a subset of the functional components of the  
 243 healthcare enterprise, called IHE actors, and specifies their interactions in terms of a set of coordinated,  
 244 standards-based transactions. The Framework offers a common language that healthcare professionals and  
 245 vendors can use to discuss integration needs of healthcare enterprises and the integration capabilities of  
 246 information systems in precise terms. Integration Profiles specify implementations of standards that are  
 247 designed to meet identified clinical needs.

248 This implementation guide uses the IHE Patient Demographics Query (PDQ) for retrieving information  
 249 about the dialysis patient. The information includes first name, last name, date of birth, gender, and  
 250 medical record number. This information is used when requesting prescription information from the EMR.

251 This implementation guide creates a new query based upon the HL7 Query By Parameter/Segment Pattern  
 252 Response model. The dialysis machine sends patient identifying information to the EMR which in turn  
 253 responds with parameters needed to perform the dialysis treatment specific to the patient.

### 254 1.5.2 Usage/Cardinality

255 Optionality (aka. Usage) is defined as one of the values in the following table.

Usage Values

Value	Description
B	Left for backwards compatibility.
C	Conditional
CE	Conditional but maybe empty.
O	Optional
R	Required
RE	Required but maybe Empty. The field or data type component description must stipulate when the field or data type component may be empty.
W	Withdrawn
X	Not used

256

257 The following rules are used when converting usage to cardinality.

258 1. Any object listed as Mandatory will have a cardinality of "1..X". Where X is the maximum number  
 259 possible or "N" if the number is unlimited/unknown.

260 2. Any object listed as Optional will have a cardinality of "0..X".

261 3. Any object listed as Conditional will have a cardinality of "0..X". If the condition is TRUE then the  
 262 effective cardinality is "1..X". If the condition is FALSE then the effective cardinality is "0..0".

263 4. If a channel is conditional and the condition is FALSE, then none of the channel's metrics will be  
 264 present regardless of their MOC/cardinality. This is another way to say cardinality scope applies  
 265 within the parent container.

266 5. Metric should not be listed as conditional unless their condition is different from the channel's  
 267 condition.



268

269 

### 1.5.3 Object Reference Identifiers

270 This standard defines objects to support a wide range of dialysis machines and therapies. Each object is  
 271 assigned an IEEE 11073 Reference Identifier or Ref ID for short. Although the number of objects is large it is  
 272 not exhaustive. If there is a need to capture an object that is not defined, 11073 Private Terms may be  
 273 used. 11073 Terms consists of a Partition and Term Code. The Term is equal to  $\text{Partition} * 2^{16} + \text{Term Code}$ .  
 274 Private metrics are placed in partition 2 and private events are placed in partition 3. Private Terms use Term  
 275 Codes that are in the range of 0xF000 to 0xFFFF. For Dialysis Machine Private Metrics, these values are in  
 276 the range of 2:0xF000 and 2:0xFFFF. For example, a manufacturer that needed to define a term for a Blood  
 277 Warmer Settings could use the following:

278 192548^MDCACME\_BLOOD\_TEMP\_SETTING^MDC

279 The value 129548 is partition 2 and item 0xF024.

280 Text based REFID should be descriptive and unique enough to allow the manufacturer to be identified. For  
 281 example, "MDCACME".

282 It is the responsibility of the manufacturer to disclose the meaning of any private terms that they report if  
 283 they are needed by an external entity.

284

285 

## 1.6 Referenced Standards

Ref	Title
[1]	IHE Patient Care Device, Technical Framework Volume 1 - Profiles, Rev 9.0 – Final Text, December 12, 2019
[2]	IHE Patient Care Device, Technical Framework Volume 2 - Transactions, Rev 9.0 – Final Text, December 12, 2019
[3]	IHE Patient Care Device, Technical Framework, Volume 3 – Sematic Content, Rev 9.0 – Final Text, December 12, 2019
[4]	IHE IT Infrastructure (ITI) Technical Framework Volume 1 - Integration Profiles, Rev 14.0, July 21, 2017
[5]	HL7 v2.6 Normative Standard, October 2007
[6]	ISO/IEEE 11073-10101R Health informatics – Point-of-care, Medical Device Communication – Part 10101: Nomenclature, 2020
[7]	ISO/IEEE 11073-10201 Health informatics – Point-of-care, Medical Device Communication – Part 10201: Domain information model, First Edition, 2004-12-15
[8]	ISO/IEEE 11073-10101R Health informatics – Point-of-care, Medical Device Communication – Part 10101b: Nomenclature, 2023
[9]	ISO/IEEE 11073-10301 Health informatics – Point-of-care, Medical Device Communication – Part 10201: Domain information model, First Edition, 2004-12-15
[10]	ISO/IEEE 11073-20101 Health informatics – Point-of-care, Medical Device Communication – Part 20101: Application Profiles – Base Standard, First Edition, 2004-12-15

286

287 

## 1.7 Glossary

288 **DEV:** An IHE domain formed to address the integration of medical devices into the healthcare  
 289 enterprise.

290 **DOC:** Device Observation Consumer

291 **DOR:** Device Observation Reporter

292 **HD:** Hemodialysis

293 **IHE:** Integrating the Healthcare Enterprise – An initiative by healthcare professionals and

294 industry to improve the way computer systems in healthcare share information.  
295 **ITI:** Information Technology Infrastructure  
296 **MDS:** Medical Device System  
297 **MDC:** Medical Device Communications  
298 **MLLP:** Minimal Lower Layer Protocol – is used for transferring HL7 messages over Ethernet. It  
299 defines delimiters which identify the beginning and ends of the HL7 message.  
300 **MRN:** Medical Record Number  
301 **PCD:** Patient Care Device – An IHE DEV Program formed to address the integration of medical  
302 devices into the healthcare enterprise.  
303 **PD:** Peritoneal Dialysis  
304 **VMD:** Virtual Medical Device  
305 **UOM:** Unit of Measure  
306 **CARD:** Cardinality

307

## 308 2 Time Synchronization

309 The data reported by a dialysis machine must be reconcilable with data from other parts of the health care  
310 enterprise. For this to occur, time must be synchronized between the machines involved in the transactions. To this  
311 end, dialysis machines compliant with this standard will use the IHE Consistent Time (CT) Protocol. The Consistent  
312 Time Profile specifies the use of the Network Time Protocol (NTP) defined in RFC1305.

## 313 3 Message Transport

314 It is not the intent of this standard to specify the manner in which messages are transported.

315 That said, the default and minimal transport for messages is the HL7 defined “Minimum Lower Layer Protocol”  
316 (MLLP) for sending unencrypted HL7 messages over a TCP/IP connection without additional security provisions.  
317 MLLP is widely used at IHE Connectathons and other test venues. It is frequently used for “behind the hospital  
318 firewall” Local Area Network (LAN) installations or when using Virtual Private Networks (VPNs) or other security  
319 solutions provided by the vendor, customer or third party solution over public networks.

320 There are other means that could be used for transferring messages. Regardless of the transport method, it is  
321 critical that security considerations be taken into account.


322

## 323 4 Patient Identification

### 324 4.1 Overview

325 Patient Identification is perhaps the most essential infrastructural component of any interoperability and  
326 communication process, particularly when prescriptions are downloaded from the EMR. It is the key element in  
327 medical device, communication, data analysis, reporting and record keeping. Automation of the entry of patient  
328 identification to patient care device has the potential for improving throughput, reducing errors, increasing safety  
329 and device and drug effectiveness, and efficiency. It is strongly recommended that implementations use IHE Patient  
330 Demographics Query (PDQ) Integration Profile compliant transactions for acquisition of a unique Patient  
331 Identification. Other mechanisms such as bar code or RFID are also valid alternatives or complements. Examples of  
332 Patient Identifiers include Medical Record Number (MR), Person Number (PN), and Social Security Number (SS).

333

**Figure 1 – Patient Demographics Query [ITI-21]**

334

335

336 There are 6 use cases that are specifically identified.

337 Case #1 – Patient Identifier Scanned by Dialysis Machine from Wrist Band

338 In the case where the Dialysis Machine can scan the Patient Identifier from a wrist band then no further  
 339 verification activities are needed. This is based on the precedent set by infusion pump manufacturers.

340 Case #2 – Patient Identifier Scanned from Something other than Wrist Band

341 There may be cases where some other means of obtaining the Patient Identifier is performed. An example  
 342 would be an ID card containing a photo of the patient and a barcode containing the Patient Identifier. In  
 343 this case, the caregiver can verify the link between the artifact used to identify the patient. If this link can  
 344 be established, then the scanned Patient Identifier could be used with no further verification.

345 Case #3 – Patient Demographics Loaded from External Device

346 In this case the patient demographics are loaded from some external device such as a memory card or USB  
 347 stick. Once the demographics are loaded, this use case is identical to the one in which demographics are  
 348 manually entered by the caregiver.

349 Case #4 – Patient Identifier Manually Entered by Caregiver

350 In the case where the care giver enters the Patient Identifier into the machine manually, the system must  
 351 verify that no mistake was made. In this case the dialysis machine uses a PDQ transaction to get the  
 352 demographics information, first name, last name, gender, and date of birth, from the EMR. The dialysis  
 353 machine displays this information to the user so that they can confirm that the entered Patient Identifier  
 354 matches the patient.

355 Case #5 – Patient Demographics Manually Entered by Care Giver

356 In the case where no Patient Identifier is known, the care giver manually enters the first name, last name,  
 357 gender, and date of birth. This information is sent to the EMR in a PDQ transaction. The EMR returns a list  
 358 of matching patients. The user must choose from this list of patients to get the correct Patient Identifier.

359 Case #6 – No Patient Demographics or Patient Identifier is Available

360 If no information is available to identify the patient, then the dialysis machine will report the Machine  
 361 Model and Serial Number as the Patient Identifier. If this is done, then the patient prescription cannot be  
 362 loaded from the EMR.

363 4.2 Message Structure

364 4.2.1 Patient Demographics Query

365 In the Patient Demographics Query, the Dialysis Machine provides the search criteria for the desired  
 366 patient.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
QPD	Query Parameter Definition	R	[1..1]
RCP	Response Control Parameters	R	[1..1]
[DSC]	Continuation Pointer	X	[0..0]

367 The dialysis machine use of this message does not deviate from the IHE Standard.

368 For simplicity, this implementation guide specifies use of a Single Domain Patient Demographics Supplier  
 369 (See IHE ITI Technical Framework Vol 2x, Appendix M). Therefore, Receiving Application (MSH-5) and  
 370 Receiving Facility (MSH-6) are left blank.

371 4.2.2 Patient Demographics Response

372 The EMR replies to the Patient Demographics Request with a Patient Demographics Response. This  
 373 message will contain 0 or more PIDs segments containing information about the patients matching the  
 374 criteria specified in the query.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
MSA	Message Acknowledgment	R	[1..1]
[ {ERR} ]	Error		
QAK	Query Acknowledge	R	[1..1]
QPD	Query Parameter Definition	R	[1..1]
[ {		O	[0..1]
PID	Patient Identifier	O	[0..1]
[PD1]			
[QRI]			
} ]			
[DSC]	Continuation Pointer	X	[0..0]

375 The dialysis machine use of this message does not deviate from the IHE Standard.

376 Upon receiving the response, the dialysis machine shall verify the following items:

- 377 • The Control ID (MSA-2) in the response matches the Control ID (MSH-10) in the request.
- 378 • The Query Tag (QAK-1) in the response matches the Query Tag (QPD-2) in the request.

- 379
- The Query Name (QAK-3) in the response matches the Query Name (QPD-1) in the request.
- 380
- If the dialysis machine does not support continuation, then Hits Remaining (QAK-6) in the response
- 381 must be 0.

## 382 4.3 Sample Messages

### 383 4.3.1 Example 1 – Query using Medical Record Number

384 The query is for the patient with a medical record number of 555444222111.

```
385 MSH|^~\&|ACME Dialysis Machine^00059AFFFE3C7A00^EUI-
386 64|||202204120831230000||QBP^Q22^QBP_Q21|20220412083123173|P|2.6||AL|NE|||
387 QPD|IHE PDQ Query|20220412083123174|@PID.3^555444222111^^^MR
388 RCP|I||R|
```

### 389 4.3.2 Example 2 – Query Response with No Matches

390 This response to the query for medical record number 555444222111 returns no data.

```
391 MSH|^~\&||ACME Dialysis Machine^00059AFFFE3C7A00^EUI-
392 64|||202204120831230000||RSP^K22^RSP_K21|20220412083123174|P|2.6||NE|NE|||
393 MSA|AA|20220412083123173
394 QAK|20220412083123174|NF|IHE PDQ Query|0|0|0
395 QPD|IHE PDQ Query|20220412083123174|@PID.3^555444222111^^^MR
```

### 396 4.3.3 Example 3 – Query using First Name and Last Name

397 The query is for patient John Smith.

```
398 MSH|^~\&|ACME Dialysis Machine^00059AFFFE3C7A00^EUI-
399 64|||202204120831230000||QBP^Q22^QBP_Q21|20220412083123138|P|2.6||AL|NE|||
400 QPD|IHE PDQ Query|20220412083123153|@PID.5.1^Smith~@PID.5.2^John
401 RCP|I||R|
```

### 402 4.3.4 Example 4 – Query Response with Multiple Matches

403 This request for all patients named John Smith returns two patients, one born in 1964 and the other in

404 2000.

```
405 MSH|^~\&||ACME Dialysis Machine^00059AFFFE3C7A00^EUI-
406 64|||202204120831230000||RSP^K22^RSP_K21|20220412083123170|P|2.6||NE|NE|||
407 MSA|AA|20220412083123138
408 QAK|20220412083123153|OK|IHE PDQ Query|2|2|0
409 QPD|IHE PDQ Query|20220412083123153|@PID.5.1^Smith~@PID.5.2^John
410 PID|||555444222111^^^MR||Smith^John^^^^U||19640306
411 PID|||555444999999^^^MR||Smith^John^^^^U||02000921
```

### 412 4.3.5 Example 5 – Query using Person Number

413 The query is for the patient with a person number of 010199-000H.

```
414 MSH|^~\&|ACME Dialysis Machine^00059AFFFE3C7A00^EUI-
415 64|||202204120831230000||QBP^Q22^QBP_Q21|20220412083123173|P|2.6||AL|NE|||
416 QPD|IHE PDQ Query|20220412083123174|@PID.3^010199-000H^^^PN
417 RCP|I||R|
```

## 418 5 Prescription Transfer

### 419 5.1 Overview

420 Prescription information is retrieved from the EMR using a custom query that aligns with the HL7 Query By  
 421 Parameter with a Segment Pattern Response (QBP/RSP). The dialysis machine provides the patient's MRN and  
 422 the EMR returns a series of OBX records containing the parameters that define the therapy. The parameters in  
 423 the response are grouped hierarchically following the containment models defined in Section 9.2 and 10.2.  
 424 Which parameters may be returned are defined by the Rx Use column in Table 2 and Table 4.

425 Although it is unlikely, it is possible that a patient may have more than one active prescription in the EMR. This  
 426 could occur when a patient that normally performs dialysis in a home setting is admitted to the hospital for  
 427 some condition. In this case, the hospital prescription could be different than the home prescription. The dialysis  
 428 machine expects the EMR to send back one and only-one prescription. It is up to the EMR to decide what is the  
 429 appropriate prescription to send to the dialysis machine. If the EMR is unable to select a single prescription it  
 430 should send back an error code in the Dialysis Prescription Response's ERR Segment.

431 Any setting value that is sent from the EMR to the Dialysis Machine will be sent back in the PCD-01 messages. If  
 432 the value remains unchanged from the value sent from the EMR, the setting will be identified as a Remote  
 433 Setting (RSET) in the corresponding OBX-17 field. If a value is changed by the user, this field will change to  
 434 Manual Setting (MSET). If the setting changes because of internal algorithms within the Dialysis Machine, the  
 435 OBX-17 field will be Automatic Setting (ASET). Once a remote setting changes to a manual or automatic setting,  
 436 it will not return to a remote setting even if the user or system returns it to its original value.

### 437 5.2 Message Structure

#### 438 5.2.1 Dialysis Prescription Query

439 The dialysis prescription query is used by the dialysis machine to request a prescription for a particular MRN  
 440 from the EMR.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
QPD	Query Parameter Definition	R	[1..1]
RCP	Response Control Parameters	R	[1..1]
[DSC]	Continuation Pointer	X	[0..0]

441

#### 442 QPD-1 Message Query Name

443 The name of the query shall be the CWE "MDC\_HDIALY\_RX\_QUERY^Hemodialysis Prescription  
 444 Query^MDC" for a Hemodialysis Prescription and "MDC\_PDIALY\_RX\_QUERY^Peritoneal Dialysis  
 445 Prescription Query^MDC" for a Peritoneal Dialysis Prescription. The EMR system copies this value into the  
 446 QAK-3 field in the response.

#### 447 QPD-2 Query Tag

448 The dialysis machine will populate the Query Tag field with an identifier that uniquely identifies this query  
 449 for the dialysis machine. It does not need to be globally unique. A simple method would be to use a  
 450 timestamp in the format of "YYYYMMDDHHMMSSZZZ" where "ZZZ" is milliseconds. Assuming no more than  
 451 one query is sent every millisecond this format will result in a unique value. The EMR system copies this  
 452 value into the QAK-1 field in the response.

#### 453 QPD-3 User Parameters

454 The Dialysis Prescription Query specifies parameters in the same manner as the IHE PDQ message. Initially,

455 only the MRN field will be included but using this format allows for future expansion.

## 456 5.2.2 Dialysis Prescription Response

457 The EMR responds to a Dialysis Prescription Query with a Dialysis Prescription Response. The response  
458 contains the current prescription for the patient with the MRN specified in the query.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
MSA	Message Acknowledgment	R	[1..1]
[ {ERR} ]	Error		
QAK	Query Acknowledge	R	[1..1]
QPD	Query Parameter Definition	R	[1..1]
[ {		O	[0..1]
ORC	Common Order Segment	O	[0..1]
OBX	Observation Results	O	[0..N]
} ]			
[DSC]	Continuation Pointer	X	[0..0]

459 Upon receiving the response, the dialysis machine shall verify the following items:

- 460 • The Control ID (MSA-2) in the response matches the Control ID (MSH-10) in the request.
- 461 • The Query Tag (QAK-1) in the response matches the Query Tag (QPD-2) in the request.
- 462 • The Query Name (QAK-3) in the response matches the Query Name (QPD-1) in the request.
- 463 • The MRN (QPD-3) in the response matches the MRN (QPD-3) in the request.

464

465 In cases where the dialysis machine receives a valid Dialysis Prescription Response, but the machine is  
466 unable to meet all of the parameters contained in the response, the dialysis machine has several options.

467 Option 1 – The dialysis machine may discard the prescription in its entirety. Either the treatment is  
468 not performed, or the caregiver manually enters the prescription. This is no different than  
469 performing a treatment without requesting the prescription. In this case, the Filler Order Number  
470 (OBR-3) reported in the PCD-01 Status Messages will be empty.

471 Option 2 – The caregiver can use the Ordering Provider (ORC-12) or Call back Phone Number (ORC-  
472 14) to request an updated prescription. In this case, the treatment is delayed. The Filler Order  
473 Number (OBR-3) reported and setting parameters in the PCD-01 Status message will match the Filler  
474 Order Number (ORC-3) and parameters specified in the Dialysis Prescription Response.

475 Option 3 – The caregiver can accept the parameters from the prescription that can be performed by  
476 the dialysis machine. They can manually override parameters that do not align with the specific  
477 dialysis machine. In this case, the Filler Order Number (OBR-3) in the PCD-01 Status message will  
478 match Filler Order Number (ORC-3) in the Dialysis Prescription Response but the setting parameters  
479 may differ.


480

## 481 5.3 Profile Representation

482 Dialysis machines support two modes of pump operation, constant speed and profiled speed. For constant  
483 speed, the pump always runs at the same speed when it is running. For profiled speed, the speed of the pump  
484 changes depending on phase of the treatment. The mechanism used to transfer the prescription from the EMR  
485 to the dialysis machines needs to support both modes of operation.

486

Figure 2 –Profile Flow Chart



487

488

489

The following table identifies the possible profiles in a prescription.

Pump	Mode Parameter	Constant Speed Parameter
Blood	N/A because blood pump does not support profiled speed.	MDC_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_SETTING
Dialysate	MDC_HDIALY_DIALYSATE_FLOW_MODE	MDC_HDIALY_DIALYSATE_FLOW_RATE_SETTING
Post-Filter Substitution (RF) Fluid	MDC_HDIALY_RF_FLOW_MODE	MDC_HDIALY_RF_POST_FILTER_FLOW_RATE_SETTING
Pre-Filter Substitution (RF) Fluid	MDC_HDIALY_RF_FLOW_MODE	MDC_HDIALY_RF_PRE_FILTER_FLOW_RATE_SETTING
Ultra Filtration	MDC_HDIALY_UF_MODE	MDC_HDIALY_UF_RATE_SETTING
Anticoagulant	MDC_HDIALY_ANTICOAG_MODE	MDC_HDIALY_ANTICOAG_INFUS_RATE_SETTING
Sodium	MDC_HDIALY_DIALYSATE_CONC_NA_MODE	MDC_HDIALY_DIALYSATE_CONC_NA_SETTING

490

491

The dialysis machine can represent two types of profiles. The first is a manufacturer specific profile. In this type



492 of profile, the details of the profile are specific to a particular manufacturer. The second type of profile is a  
 493 generic profile. The Implementation Guide supports 4 type of generic profiles, Constant, Linear, Exponential, and  
 494 Step. For a generic profile the parameters of the profile are specified as Facet Level 5 objects of the Mode  
 495 Parameter.

496

ID	Object Name	Description
1	MDCX_HDIALY_PROFILE_TYPE	Contains the type of the profile. This will be VENDOR, CONSTANT, LINEAR, EXPONENTIAL, or STEP.
2	MDCX_HDIALY_PROFILE_VALUE	Contains the control points of the profile. The control points are a numeric array. For a constant profile there is only one point. For a linear or exponential profile there are two points, one for the starting value and one for the ending value. For step profile there is one value for each step.
3	MDCX_HDIALY_PROFILE_TIME	This parameter is optional. If it is not present, the control points are assumed to be equally space across the total treatment time. If it is present, it is an array of numerical values that are the treatment times at which a control point takes effect.
4	MDCX_HDIALY_PROFILE_EXP_HALF_TIME	This parameter is optional and only applies to the Exponential Profile. It specifies the time at which 50% of the change has been complete.
5	MDCX_HDIALY_PROFILE_NAME	This parameter is optional and only applies to the Vendor Profile. It uniquely identifies the profile.

497

### 498 5.3.1 Vendor Profile

499 In Vendor Profile, the details of the profile are specific to a particular manufacturer. The profile is identified  
 500 using the MDCX\_HDIALY\_PROFILE\_NAME object and consists of three components, Manufacturer, Model,  
 501 and Name. The Manufacturer is a string that uniquely identifies the maker of the dialysis machine. Model is  
 502 the model of the dialysis machine. If the profile is universal across a manufacturer's product line, then this  
 503 field may be left blank. For example, "FMC^^Profile 1" or "NxStage^SystemOne^Gentle". It is the  
 504 manufacturer's responsibility to make sure that these profile identifiers are compatible or unique across  
 505 device models and software versions. When a manufacturer specific profile is used, it is imperative that  
 506 this profile be clearly defined. The best way of accomplishing this is by defining the profile in the device's  
 507 instructions for use.

508 The segments below represent a B-Braun Dialog+ UF Profile 2. Note, to profile name is provided for  
 509 example only and may not represent an actual profile.

510 OBX|66|ST|158619^MDC\_HDIALY\_UF\_MODE^MDC|1.1.4.27|PRO-WOT|||||F

511 OBX|67|ST|0^MDCX\_HDIALY\_PROFILE\_TYPE^MDC|1.1.4.27.1|VENDOR|||||F

512 OBX|68|NA|0^MDCX\_HDIALY\_PROFILE\_NAME^MDC|1.1.4.27.2|BBRAUN^DIALOGPLUS^UF1|||||F

513

514  
515  
516  
517  
518

### 5.3.2 Constant Profile

A constant profile starts at one level and stays at this level for the course of the treatment. A constant profile could also be represented as a linear or exponential profile where the starting and ending values are the same or as step profile with only a single step.

**Figure 3 – Constant Profile**



519  
520  
521  
522  
523  
524

The segments below represent a constant UF profile of 10 ml/hour.


```
OBX|66|ST|158619^MDC_HDIALY_UF_MODE^MDC|1.1.4.27|PRO-WOT|||||F
OBX|67|ST|0^MDCX_HDIALY_PROFILE_TYPE^MDC|1.1.4.27.1|CONSTANT|||||F
OBX|68|NA|0^MDCX_HDIALY_PROFILE_VALUE^MDC|1.1.4.27.2|10|ml/h^milliliter per|||||F
```

525  
526  
527

### 5.3.3 Linear Profile

A linear profile starts at one level and linearly changes to another level.

**Figure 4 –Linear Profile**



528  
529  
530  
531  
532

The segments below define a linear profile that starts at 10 ml/hour and ends at 60 ml/hour.

```
OBX|66|ST|158619^MDC_HDIALY_UF_MODE^MDC|1.1.4.27|PRO-WOT|||||F
OBX|67|ST|0^MDCX_HDIALY_PROFILE_TYPE^MDC|1.1.4.27.1|LINEAR|||||F
OBX|68|NA|0^MDCX_HDIALY_PROFILE_VALUE^MDC|1.1.4.27.2|10^30|ml/h^milliliter per hour^UCUM|||||F
```

533 **5.3.4 Exponential Profile**


534 An exponential profile starts at one level and exponentially transitions to another level. The general  
535 equation for the exponential curve is

$$536 \quad y = (A - B) e^{-kt} + B$$

537 The standard has three different ways of specifying the time constant, k.

- 538 1) The time constant can be specified by specifying the Half Time,  $T_{\text{half}}$ , of the change using the  
539 MDCX\_HDIALY\_PROFILE\_EXP\_HALF\_TIME. This value is the time at which 50% of the change has  
540 been completed. In this case, k is  $0.693/T_{\text{half}}$ .
- 541 2) If a time is specified using MDCX\_HDIALY\_PROFILE\_TIME, then the standard defines this time as the  
542 point when 99.7% of the change has been completed. In this case, k is  $0.003/T$ .
- 543 3) If neither MDCX\_HDIALY\_PROFILE\_TIME nor MDC\_HDIALY\_EXP\_HALF\_TIME is used, then the curve  
544 will be computed assuming that 99.7% of the change is completed at the end of the treatment. In  
545 this case k is  $0.003/T_{\text{tx}}$ .

546 **Figure 5 – Exponential Profile**



547

548

549 The segments below define a profile that starts at 30 ml/hr and reaches 10 ml/hour after 60 minutes. The  
550 transition between these two points is exponential.

551 `OBX|66|ST|158619^MDC_HDIALY_UF_MODE^MDC|1.1.4.27|PRO-WOT|||||F`

552 `OBX|67|ST|0^MDCX_HDIALY_PROFILE_TYPE^MDC|1.1.4.27.1|EXPONENTIAL|||||F`

553 `OBX|68|NA|0^MDCX_HDIALY_PROFILE_VALUE^MDC|1.1.4.27.2|30^10|ml/h^milliliter per hour^UCUM|||||F`

554 `OBX|69|NA|0^MDCX_HDIALY_PROFILE_TIME^MDC|1.1.4.27.3|0^60.0|min^minutes^UCUM|||||F`


555

556 **5.3.5 Step Profile**

557 A step profile contains one or more discrete values. Each value is run for a specific amount of time before  
558 switching to the next value.

559

Figure 6 – Step Profile



560

561

562

563

The segments below define a Step Profile with three steps. It starts at 10 ml/hr. At Tx Time of 30 minutes, it changes to 20 ml/hr. At Tx Time of 60 minutes, it changes to 30 ml/hr. The rate remains at 30 ml/hr for the remainder of treatment.

564

```
OBX|66|ST|158619^MDC_HDIALY_UF_MODE^MDC|1.1.4.27|PRO-WOT|||||F
```

565

```
OBX|67|ST|0^MDCX_HDIALY_PROFILE_TYPE^MDC|1.1.4.27.1|STEP|||||F
```

566

567

```
OBX|68|NA|0^MDCX_HDIALY_PROFILE_VALUE^MDC|1.1.4.27.2|10^20^30|ml/h^milliliter per  
hour^UCUM|||||F
```

568

```
OBX|69|NA|0^MDCX_HDIALY_PROFILE_TIME^MDC|1.1.4.27.3|0^30.0^60.0|min^minutes^UCUM|||||F
```

569

## 570 5.4 PD Prescription Representation

571

572

There are 3 ways in which a PD treatment can be specified. Which one of the options used is reported in the MDCX\_PDIALY\_PRESCRIPTION\_TYPE object.

573

### Option 1) Exchange Based

574

575

576

577

578

579

580

When specifying an Exchange Based Treatment, all of the information is conveyed in the Exchange Channels. Each channel must have a Fill Volume, Fill Duration, Dwell Duration, Drain Volume, and Drain Duration. The following terms should NOT be present in the prescription:  
MDCX\_PDIALY\_APD\_START\_TIME\_SETTING, MDCX\_PDIALY\_APD\_END\_TIME\_SETTING,  
MDCX\_PDIALY\_APD\_TOTAL\_VOLUME\_SETTING, MDCX\_PDIALY\_APD\_CYCLES\_SETTING,  
MDCX\_PDIALY\_CAPD\_START\_TIME\_SETTING, MDCX\_PDIALY\_CAPD\_END\_TIME\_SETTING,  
MDCX\_PDIALY\_CAPD\_TOTAL\_VOLUME\_SETTING, and MDCX\_PDIALY\_CAPD\_CYCLES\_SETTING.

581

### Option 2a) Treatment Based – End Time Specified

582

583

584

585

586

587

588

589

When specifying a Treatment Based Treatment, exchange channels are not present in the prescription. The exchanges are computed by the machine based upon the values on  
MDCX\_PDIALY\_APD\_END\_TIME\_SETTING, MDCX\_PDIALY\_APD\_TOTAL\_VOLUME\_SETTING,  
MDCX\_PDIALY\_APD\_CYCLES\_SETTING, MDCX\_PDIALY\_CAPD\_END\_TIME\_SETTING,  
MDCX\_PDIALY\_CAPD\_TOTAL\_VOLUME\_SETTING, and MDCX\_PDIALY\_CAPD\_CYCLES\_SETTING. The machine computes the duration based upon the actual start time and the specified end time. The dwell for each exchange is computed to fit within the computed duration. If a sufficient dwell time cannot be achieved the cycler will display an alert and have the user correct the prescription.

590

### Option 2b) Treatment Based – Duration Specified

591

When specifying a Treatment Based Treatment, exchange channels are not present in the prescription. The

592 exchanges are computed by the machine based up on the values on  
 593 MDCX\_PDIALY\_APD\_START\_TIME\_SETTING, MDCX\_PDIALY\_APD\_END\_TIME\_SETTING,  
 594 MDCX\_PDIALY\_APD\_TOTAL\_VOLUME\_SETTING, MDCX\_PDIALY\_APD\_CYCLES\_SETTING,  
 595 MDCX\_PDIALY\_CAPD\_START\_TIME\_SETTING, MDCX\_PDIALY\_CAPD\_END\_TIME\_SETTING,  
 596 MDCX\_PDIALY\_CAPD\_TOTAL\_VOLUME\_SETTING, and MDCX\_PDIALY\_CAPD\_CYCLES\_SETTING. The machine  
 597 computes the duration as the amount of time between specified end time and specified start time. The  
 598 dwell for each exchange is computed to fit within the computed duration. The actual treatment end time  
 599 will be the actual start time plus the duration of each exchange. Said another way, the actual end time will  
 600 change based upon the actual start time.

601

	1	2a	2b	
	Exchange Based	Treatment Based - Fixed End Time	Treatment Based – Fixed Duration	Condition
MDCX_PDIALY_PRESCRIPTION_TYPE	1..1	1..1	1..1	
MDCX_PDIALY_APD_START_TIME_SETTING	0..0	0..0	0..1	C5
MDCX_PDIALY_APD_END_TIME_SETTING	0..0	0..1	0..1	C5
MDCX_PDIALY_APD_TOTAL_VOLUME_SETTING	0..0	0..1	0..1	C5
MDCX_PDIALY_APD_CYCLES_SETTING	0..0	0..1	0..1	C5
MDCX_PDIALY_CAPD_START_TIME_SETTING	0..0	0..0	0..1	C6
MDCX_PDIALY_CAPD_END_TIME_SETTING	0..0	0..1	0..1	C6
MDCX_PDIALY_CAPD_TOTAL_VOLUME_SETTING	0..0	0..1	0..1	C6
MDCX_PDIALY_CAPD_CYCLES_SETTING	0..0	0..1	0..1	C6
MDCX_PDIALY_FILL_DURATION_SETTING	1..1	0..0	0..0	
MDCX_PDIALY_FILL_VOLUME_SETTING	1..1	0..0	0..0	
MDCX_PDIALY_DWELL_DURATION_SETTING	1..1	0..0	0..0	
MDCX_PDIALY_DRAIN_DURATION_SETTING	1..1	0..0	0..0	
MDCX_PDIALY_DRAIN_VOLUME_SETTING	1..1	0..0	0..0	

602

## 603 5.5 Sample Messages

### 604 5.5.1 Example 1 – HD Prescription Request

605 The following message requests the prescription for the patient with medical record number  
 606 555444222111.

```
607 MSH|^~\&|ACME_Dialysis_Machine^00059AFFFE3C7A00^EUI-
608 64|||20220330125317+0000||QBP^D01^QBP_D01|PQ20211216144700|P|2.6||AL|NE||||
609 QPD|0^MDC_HDIALY_RX_QUERY^MDC|Q001|@PID.3^555444222111^^^MR
610 RCP|I||R|
```

### 611 5.5.2 Example 2 – HD Prescription Response

612 The following message is the prescription for patient 555444222111. The therapy is HD using RFP-204  
 613 bagged dialysate. Blood flow rate is 250 ml/min, Dialysate Flow Rate is 120 ml/min, UF Rate is 400 ml/hr,  
 614 and the weight to remove is 1 L.

```
615 MSH|^~\&|ACME_Dialysis_Machine^00059AFFFE3C7A00^EUI-
616 64|||20220330125317+0000||RSP^K22^RSP_K21|20220330125317627|P|2.6||AL|NE||||
617 MSA|AA|PQ20211216144700
618 QAK|Q001|OK|0^MDC_HDIALY_RX_QUERY^MDC|1|1|0
```

619 QPD|0^MDC\_HDIALY\_RX\_QUERY^MDC|Q001|@PID.3^555444222111^^^MR  
 620 OBC|NW|A226677^PC||||N||||444-44-4444^HIPPOCRATES^HAROLD^^^MD  
 621 OBX|1|ST|70929^MDC\_DEV\_HDIALY\_MACHINE\_MDS^MDC|1|||||F  
 622 OBX|2|ST|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.1|||||F  
 623 OBX|3|ST|70939^MDC\_DEV\_HDIALY\_MACH\_CONFIG\_CHAN^MDC|1.1.1|||||F  
 624 OBX|4|ST|158598^MDC\_HDIALY\_MACH\_TX\_MODALITY^MDC|1.1.1.1|HD|||||F  
 625 OBX|5|ST|70967^MDC\_DEV\_HDIALY\_THERAPY\_OUTCOMES\_CHAN^MDC|1.1.2|||||F  
 626 OBX|6|ST|158618^MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD^MDC|1.1.2.1|UF|||||F  
 627 OBX|7|ST|70947^MDC\_DEV\_HDIALY\_BLOOD\_PUMP\_CHAN^MDC|1.1.3|||||F  
 628 OBX|8|NM|16935956^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE\_SETTING^MDC|1.1.3.1|250|ml/min^ml/min^UCUM||||  
 629 |F  
 630 OBX|9|ST|158604^MDC\_HDIALY\_BLD\_PUMP\_MODE^MDC|1.1.3.2|2N|||||F  
 631 OBX|10|ST|70951^MDC\_DEV\_HDIALY\_FLUID\_CHAN^MDC|1.1.4|||||F  
 632 OBX|11|ST|158606^MDC\_HDIALY\_DIALYSATE\_FLOW\_MODE^MDC|1.1.4.1|CONST|||||F  
 633 OBX|12|NM|16936008^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE\_SETTING^MDC|1.1.4.2|120|ml/min^ml/min^UCUM|||||F  
 634 OBX|13|NM|0^MDC\_HDIALY\_DIALYSATE\_VOL\_SETTING^MDC|1.1.4.3|25|L^L^UCUM|||||F  
 635 OBX|14|ST|158608^MDC\_HDIALY\_DIALYSATE\_NAME^MDC|1.1.4.4|RFP-204|||||F  
 636 OBX|15|ST|70971^MDC\_DEV\_HDIALY\_UF\_CHAN^MDC|1.1.5|||||F  
 637 OBX|16|ST|158619^MDC\_HDIALY\_UF\_MODE^MDC|1.1.5.1|CONST-WT|||||F  
 638 OBX|17|NM|16936252^MDC\_HDIALY\_UF\_RATE\_SETTING^MDC|1.1.5.2|400|ml/h^ml/h^UCUM|||||F  
 639 OBX|18|NM|159028^MDC\_HDIALY\_UF\_TARGET\_VOL\_TO\_REMOVE^MDC|1.1.5.3|1000|ml^ml^UCUM|||||F

### 5.5.3 Example 2 – HD Prescription Response for HF Therapy with a Linear UF Profile

The following message is the prescription for patient 555444222111. The therapy is Pre-Filter HF using RFP-400 bagged dialysate. Blood flow rate is 250 ml/min, RF Flow Rate is 120 ml/min, UF Rate starts at 300 and decreases linearly to 100 ml/hr at end of treatment, and the weight to remove is 1 L.

644 MSH|^~\&|ACME\_Dialysis\_Machine^00059AFFFE3C7A00^EUI-  
 645 64||||20220330132820+0000||RSP^K22^RSP\_K21|20220330132820292|P|2.6|||AL|NE||||  
 646 MSA|AA|PQ20211216144700  
 647 QAK|Q001|OK|0^MDC\_HDIALY\_RX\_QUERY^MDC|1|1|0  
 648 QPD|0^MDC\_HDIALY\_RX\_QUERY^MDC|Q001|@PID.3^555444222111^^^MR  
 649 OBC|NW|A226677^PC||||N||||444-44-4444^HIPPOCRATES^HAROLD^^^MD  
 650 OBX|1|ST|70929^MDC\_DEV\_HDIALY\_MACHINE\_MDS^MDC|1|||||F  
 651 OBX|2|ST|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.1|||||F  
 652 OBX|3|ST|70939^MDC\_DEV\_HDIALY\_MACH\_CONFIG\_CHAN^MDC|1.1.1|||||F  
 653 OBX|4|ST|158598^MDC\_HDIALY\_MACH\_TX\_MODALITY^MDC|1.1.1.1|HF|||||F  
 654 OBX|5|ST|70967^MDC\_DEV\_HDIALY\_THERAPY\_OUTCOMES\_CHAN^MDC|1.1.2|||||F  
 655 OBX|6|ST|158618^MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD^MDC|1.1.2.1|UF|||||F  
 656 OBX|7|ST|70947^MDC\_DEV\_HDIALY\_BLOOD\_PUMP\_CHAN^MDC|1.1.3|||||F  
 657 OBX|8|NM|16935956^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE\_SETTING^MDC|1.1.3.1|250|ml/min^ml/min^UCUM||||  
 658 |F  
 659 OBX|9|ST|158604^MDC\_HDIALY\_BLD\_PUMP\_MODE^MDC|1.1.3.2|2N|||||F  
 660 OBX|10|ST|70959^MDC\_DEV\_HDIALY\_CONVECTIVE\_CHAN^MDC|1.1.4|||||F  
 661 OBX|11|ST|158613^MDC\_HDIALY\_RF\_DILUTION\_LOCATION^MDC|1.1.4.1|PREF|||||F  
 662 OBX|12|NM|16936156^MDC\_HDIALY\_RF\_PRE\_FILTER\_FLOW\_RATE\_SETTING^MDC|1.1.4.2|120|ml/min^ml/min^UCUM||||  
 663 |F  
 664 OBX|13|NM|16936164^MDC\_HDIALY\_RF\_PRE\_FILTER\_VOL\_SETTING^MDC|1.1.4.3|25|L^L^UCUM|||||F  
 665 OBX|14|ST|158615^MDC\_HDIALY\_RF\_PRE\_FILTER\_NAME^MDC|1.1.4.4|RFP-400|||||F

666 OBX|15|ST|70971^MDC\_DEV\_HDIALY\_UF\_CHAN^MDC|1.1.5|||||F  
 667 OBX|16|ST|158619^MDC\_DEV\_HDIALY\_UF\_MODE^MDC|1.1.5.1|PRO-WT|||||F  
 668 OBX|17|ST|0^MDCX\_HDIALY\_PROFILE\_TYPE^MDC|1.1.5.1.1|LINEAR|||||F  
 669 OBX|18|ST|1^MDCX\_HDIALY\_PROFILE\_VALUE^MDC|1.1.5.1.2|300^100|ml/h^ml/h^UCUM|||||F  
 670 OBX|19|NM|159028^MDC\_DEV\_HDIALY\_UF\_TARGET\_VOL\_TO\_REMOVE^MDC|1.1.5.2|1000|ml^ml^UCUM|||||F

671

672 **5.5.4 Example 3 – Prescription Response returning No Prescription**

673 The following message is returned by the EMR is there is no prescription associated with the requested  
 674 MRN.

675 MSH|^~\&|ACME\_Dialysis\_Machine^00059AFFFE3C7A00^EUI-  
 676 64|||||20220330125317+0000||RSP^K22^RSP\_K21|20220330125317737|P|2.6||AL|NE|||||  
 677 MSA|AA|PQ20211216144700  
 678 QAK|Q001|NF|0^MDC\_DEV\_HDIALY\_RX\_QUERY^MDC|0|0|0  
 679 QPD|0^MDC\_DEV\_HDIALY\_RX\_QUERY^MDC|Q001|555444222111

680

681 **5.5.5 Example 4 – PD Prescription Request**

682 The following message requests the prescription for the patient with medical record number  
 683 555444222111.

684 MSH|^~\&|ACME\_Dialysis\_Machine^00059AFFFE3C7A00^EUI-  
 685 64|||||20220330125317+0000||QBP^D01^QBP\_D01|PQ20211216144700|P|2.6||AL|NE|||||  
 686 QPD|0^MDC\_PDIALY\_RX\_QUERY^MDC|Q001|@PID.3^555444222111^^^MR  
 687 RCP|I||R|

688

689 **5.5.6 Example 5 – PD Prescription Response, Exchange Based**

690 The following message is a prescription for PD Therapy with 3 fluids and 3 exchanges. Exchange #1 uses  
 691 Fluid 1, Exchange #2 uses Fluids 1 & 2, and Exchange #3 uses Fluid 3

692 MSH|^~\&|ACME\_Dialysis\_Machine^00059AFFFE3C7A00^EUI-  
 693 64|||||20220330132820+0000||RSP^K22^RSP\_K21|20220330132820292|P|2.6||AL|NE|||||  
 694 MSA|AA|PQ20211216144700  
 695 QAK|Q001|OK|0^MDC\_PDIALY\_RX\_QUERY^MDC|1|1|0  
 696 QPD|0^MDC\_PDIALY\_RX\_QUERY^MDC|Q001|@PID.3^555444222111^^^MR  
 697 OBC|NW|A226677^PC||||N||||444-44-4444^HIPPOCRATES^HAROLD^^^MD  
 698 OBX|1|ST|27^MDCX\_DEV\_PDIALY\_MACHINE\_MDS^MDC|1|||||F  
 699 OBX|2|ST|0^MDCX\_DEV\_PDIALY\_VMD^MDC|1.1|||||F  
 700 OBX|3|ST|4^MDCX\_DEV\_PDIALY\_TREATMENT\_CHAN^MDC|1.1.2|||||F  
 701 OBX|4|ST|5^MDCX\_PDIALY\_TREAT\_TYPE^MDC|1.1.2.1|CAPD|||||F  
 702 OBX|5|ST|6^MDCX\_PDIALY\_TREAT\_LOCATION^MDC|1.1.2.2|HOME|||||F  
 703 OBX|6|ST|7^MDCX\_PDIALY\_PRESCRIPTION\_TYPE^MDC|1.1.2.3|EXCH|||||F  
 704 OBX|7|ST|51^MDCX\_DEV\_PDIALY\_FLUID\_1\_CHAN^MDC|1.1.3|||||F  
 705 OBX|8|ST|10^MDCX\_PDIALY\_FLUID\_SOURCE^MDC|1.1.3.1|BAG|||||F  
 706 OBX|9|NM|9^MDCX\_PDIALY\_FLUID\_BAG\_VOLUME^MDC|1.1.3.2|5.000|L^liter^UCUM|||||F  
 707 OBX|10|ST|52^MDCX\_DEV\_PDIALY\_FLUID\_2\_CHAN^MDC|1.1.4|||||F  
 708 OBX|11|ST|10^MDCX\_PDIALY\_FLUID\_SOURCE^MDC|1.1.4.1|BAG|||||F  
 709 OBX|12|NM|9^MDCX\_PDIALY\_FLUID\_BAG\_VOLUME^MDC|1.1.4.2|5.000|L^liter^UCUM|||||F

710 OBX|13|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.5|||||F  
 711 OBX|14|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.5.1|1|||||F  
 712 OBX|15|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.5.2|PENDING|||||F  
 713 OBX|16|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.5.3|13.3|min^minutes^UCUM|||||F  
 714 OBX|17|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.5.4|2.667^0.000|L^liter^UCUM|||||F  
 715 OBX|18|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.5.5|153.3|min^minutes^UCUM|||||F  
 716 OBX|19|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.5.6|13.3|min^minutes^UCUM|||||F  
 717 OBX|20|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.5.7|2.667|L^liter^UCUM|||||F  
 718 OBX|21|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.6|||||F  
 719 OBX|22|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.6.1|2|||||F  
 720 OBX|23|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.6.2|PENDING|||||F  
 721 OBX|24|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.6.3|13.3|min^minutes^UCUM|||||F  
 722 OBX|25|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.6.4|2.333^0.333|L^liter^UCUM|||||F  
 723 OBX|26|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.6.5|153.3|min^minutes^UCUM|||||F  
 724 OBX|27|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.6.6|13.3|min^minutes^UCUM|||||F  
 725 OBX|28|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.6.7|2.667|L^liter^UCUM|||||F  
 726 OBX|29|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.7|||||F  
 727 OBX|30|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.7.1|3|||||F  
 728 OBX|31|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.7.2|PENDING|||||F  
 729 OBX|32|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.7.3|13.3|min^minutes^UCUM|||||F  
 730 OBX|33|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.7.4|0.000^2.667|L^liter^UCUM|||||F  
 731 OBX|34|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.7.5|153.3|min^minutes^UCUM|||||F  
 732 OBX|35|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.7.6|13.3|min^minutes^UCUM|||||F  
 733 OBX|36|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.7.7|2.667|L^liter^UCUM|||||F

734

### 5.5.7 Example 6 – PD Prescription Response, Treatment Based

The following message is a prescription for PD Therapy with 2 fluids and 3 exchanges lasting 9 hours.

735  
 736  
 737 MSH|^~\&|ACME\_Dialysis\_Machine^00059AFFFE3C7A00^EUI-  
 738 64||||20220330132820+0000||RSP^K22^RSP\_K21|20220330132820292|P|2.6|||AL|NE|||||  
 739 MSA|AA|PQ20211216144700  
 740 QAK|Q001|OK|0^MDC\_PDIALY\_RX\_QUERY^MDC|1|1|0  
 741 QPD|0^MDC\_PDIALY\_RX\_QUERY^MDC|Q001|@PID.3^555444222111^^^MR  
 742 OBC|NW|A226677^PC||||N||||444-44-4444^HIPPOCRATES^HAROLD^^^MD  
 743 OBX|1|ST|27^MDCX\_DEV\_PDIALY\_MACHINE\_MDS^MDC|1|||||F  
 744 OBX|2|ST|0^MDCX\_DEV\_PDIALY\_VMD^MDC|1.1|||||F  
 745 OBX|3|ST|4^MDCX\_DEV\_PDIALY\_TREATMENT\_CHAN^MDC|1.1.2|||||F  
 746 OBX|4|ST|5^MDCX\_PDIALY\_TREAT\_TYPE^MDC|1.1.2.1|CAPD|||||F  
 747 OBX|5|ST|6^MDCX\_PDIALY\_TREAT\_LOCATION^MDC|1.1.2.2|HOME|||||F  
 748 OBX|6|ST|7^MDCX\_PDIALY\_PRESCRIPTION\_TYPE^MDC|1.1.2.3|TXDUR|||||F  
 749 OBX|7|TM|25^MDCX\_PDIALY\_CAPD\_START\_TIME\_SETTING^MDC|1.1.2.4|22:00:00.000|||||F  
 750 OBX|8|TM|26^MDCX\_PDIALY\_CAPD\_END\_TIME\_SETTING^MDC|1.1.2.5|07:00:00.000|||||F  
 751 OBX|9|NM|23^MDCX\_PDIALY\_CAPD\_TOTAL\_VOLUME\_SETTING^MDC|1.1.2.6|8.000|L^liter^UCUM|||||F  
 752 OBX|10|ST|24^MDCX\_PDIALY\_CAPD\_CYCLES\_SETTING^MDC|1.1.2.7|3|||||F  
 753 OBX|11|ST|51^MDCX\_DEV\_PDIALY\_FLUID\_1\_CHAN^MDC|1.1.3|||||F  
 754 OBX|12|ST|10^MDCX\_PDIALY\_FLUID\_SOURCE^MDC|1.1.3.1|BAG|||||F



```


755     OBX|13|NM|9^MDCX_PDIALY_FLUID_BAG_VOLUME^MDC|1.1.3.2|5.000|L^liter^UCUM|||||F
756     OBX|14|ST|52^MDCX_DEV_PDIALY_FLUID_2_CHAN^MDC|1.1.4|||||||F
757     OBX|15|ST|10^MDCX_PDIALY_FLUID_SOURCE^MDC|1.1.4.1|BAG|||||||F
758     OBX|16|NM|9^MDCX_PDIALY_FLUID_BAG_VOLUME^MDC|1.1.4.2|5.000|L^liter^UCUM|||||F
    
```

## 6 Reporting Treatment Information

### 6.1 Overview

The sending of Treatment Information to the EMR System is accomplished using the Patient Care Device (PCD) domain in the Device Enterprise Communication (DEC) IHE profile. In doing so, the Dialysis Machine acts as the Device Observation Report (DOR) and the EMR is the Device Observation Consumer (DOC).

**Figure 7 - Communicate PCD Data Interaction Diagram**



In general, the Dialysis Machine sends periodic reports at an interval of between several times per minute (high acuity) and a maximum interval of 24 hours (chronic, home health). The minimum and maximum intervals are configured at implementation. The Dialysis Machine may also send aperiodic reports for "event type" information.

This specification includes different observation types based upon differing real world triggers. The objects in each observation are indicated in the Phase and Temporal columns of Table 2 and Table 4. The phase column specifies the machine state under which the object is sent. The possible values are:

- All – The object is always sent.
- Intradialytic – The object is only sent when the cyclor is delivering therapy

- 776 The Temporal column specifies what machine condition causes the object to be sent. The possible values are:
- 777
- All – The object is sent in every PCD-01 Message.
  - Episodic – The object is sent when the associated alert is active or changing state.
- 778

779

780 Since PCD-01 messages are sent periodically, at a minimum, objects list as “All” will be sent at the update rate.

781 These will also be sent when a change to an Episodic object causes a PCD-01 Message to be sent.

782 A Status Report will contain all mandatory elements, the applicable conditional elements, and the optional

783 elements supported by the dialysis machine. In addition, the dialysis machine may report other data objects not

784 identified by this standard. This allows the device manufacturer to report information that is not common across

785 all dialysis machines and is therefore not part of this standard. The hierarchy must be maintained but ordering

786 within a section can change. The ordering in Table 2 and Table 4 do not prescribe a particular order in the PCD-

787 01 message.

788 The nature of the PCD-01 message is that it contains a snapshot in time of the treatment. The capturing of an

789 entire treatment into a "run sheet" can be accomplished using the HL7 Batch Protocol described in Appendix A.

## 790 6.2 Device Request

### 791 6.2.1 Message Structure

792 Segments shown in gray are not used in this standard.

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
[[SFT]]	Software Segment	X	[0..0]
[UAC]	User Authentication Credential	X	[0..0]
{	--- PATIENT_RESULT begin		
[	--- PATIENT begin		
PID	Patient Identification	R	[1..1]
[PD1]	Additional Demographics	X	[0..0]
[[PRT]]	Participation	X	[0..0]
[[NTE]]	Notes and Comments	X	[0..0]
[[NK1]]	Next of Kin/Associated Parties	O	[0..3]
[	--- VISIT begin		
PV1	Patient Visit	R	[1..1]
[PV2]	Patient Visit – Additional Info	X	[0..0]
[[PRT]]	Participation	X	[0..0]
]	--- VISIT end		
]	--- PATIENT end		
{	---ORDER_OBSERVATION begin		
[ORC]	Order Common	X	[0..0]
OBR	Observation Request	R	[1..1]
[[NTE]]	Notes and Comments	O	[0..1]
[[PRT]]	Participation	X	[0..0]
{	--- TIMING_QTY begin		
TQ1	Timing/Quantity	X	[0..0]
[[TQ2]]	Timing/Quantity Order Sequence	X	[0..0]
}}	--- TIMING_QTY end		

Segment	Meaning	Use	Card
[CTD]	Contact Data	X	[0..0]
{	--- OBSERVATION begin		
OBX	Observation Result	R	[1..N]
{{PRT}}	Participation	X	[0..0]
{{NTE}}	Notes and comments	X	[0..0]
}	--- OBSERVATION end		
{{FT1}}	Financial Transaction	X	[0..0]
{{CTI}}	Clinical Trial Identification	X	[0..0]
{	--- SPECIMEN begin		
SPM	Specimen	X	[0..0]
{{OBX}}	Observation related to Specimen	X	[0..0]
}	--- SPECIMEN end		
}	--- ORDER_OBSERVATION end		
}	--- PATIENT_RESULT end		
[DSC]	Continuation Pointer	X	[0..0]

793

794 **6.2.2 Observation Results Structure**

795 Observation results are grouped in the hierarchy of Medical Device System (MDS), Virtual Medical Device  
 796 (VMD), Channel, and Metric/Observation. When a VMD is not present in the dialysis machine, the VMD  
 797 information is not reported. For example, when a dialysis machine does not have a blood pressure cuff, the  
 798 MDC\_DEV\_PRESS\_BLD\_NONINV\_VMD is not reported. If a Channel does not contain any  
 799 Metrics/Observations due to dialysis machine therapy or operating mode, the channel is not reported. For  
 800 example, when a dialysis machine is performing a Hemofiltration Therapy, the  
 801 MDC\_DEV\_HDIALY\_FLUID\_CHAN will not be reported. Because of this dynamic nature, sub-observation IDs  
 802 will change from treatment to treatment and should not be used for parsing messages from the dialysis  
 803 machine.

804 **6.2.3 Event Status**

805 This standard allows the use of inline events. These non-actionable events are captured in one of two ways.  
 806 For events associated with a MDS, VMD, or Channel, the event is shown as a level 4 metric observation. For  
 807 events associated with a metric, the event is reported as a level 5 facet of the metric.

808 Example, an event associated with a channel.

```
809 OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|F|||||F
```

810 Example, an event associated with a metric.

```
811 OBX|4|NM|158776^MDC_HDIALY_BLD_PUMP_PRESS_VEN^MDC|1.1.3.15|15|mm[Hg]^Millimeters of  
812 Mercury^UCUM |20-400|||F
```

```
813 OBX|4|CWE|196670^MDC_EVT_LO^MDC|1.1.3.15.1|158776^MDC_HDIALY_BLD_PUMP_PRESS_VEN^MDC|||PH~SP||  
814 |F
```

815

816 This standard allows event information to be transmitted in one of two methods.

- 817 • True/False
- 818 • Start/Continue/End

819

820 When True/False is used, the dialysis machine communicates the state of each event in every PCD-01  
 821 message. The event is either active (i.e. True) or Inactive (i.e. False). Since the dialysis machine is sending  
 822 PCD-01 messages every time status changes instead of just periodically, the time associated with an event  
 823 state transition is the time of the observation as specified by time point OBX-14 or the default time point  
 824 specified by OBR-7.

825

826 Example, Blood Leak not detected

827 `OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|F|||||F`

828 Example, Blood Leak detected

829 `OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|T|||||F`

830

831 When Start/Continue/End is used the dialysis machine communicates the starting and ending time of the  
 832 event. The time associated with these events is expressly indicated in the OBX-14 portion of the event. In  
 833 PCD-01 messages where the event did not start or end but is still active the dialysis machine reports  
 834 continue. In the case where a dialysis machine does not know the starting time of an event, it will transmit  
 835 continue as the first state for the event.

836

837 Example, Blood Leak detected

838 `OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|start|||||||20191003092005+0000`

839 Example, Blood Leak ongoing

840 `OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|continue|||||||20191003092015+0000`

841 Example, Blood Leak done


842 `OBX|20|ST|198244^MDC_EVT_HDIALY_BLOOD_LEAK^MDC|1.1.4.15|end|||||||20191003092032+0000`

843

844 The following figure shows how the dialysis machine would report an event for each of the two methods.

845

**Figure 8 – Event Status in PCD-01 Messages**



Using Start/Continue/End

- \* PCD01 Msg at Time 9:20:00 would not mention the event.
- \* PCD01 Msg at Time 9:20:05 would contain the Start with a Time of 9:20:05.
- \* PCD01 Msg at Time 9:20:15 and 9:20:25 would contain the Continue.
- \* PCD01 Msg at Time 9:20:32 would contain the End with a timestamp of 9:20:32.
- \* PCD01 msg at Time 9:20:42 and beyond would not mention the event

Using True/False without Timestamp in OBX

- \* PCD01 Msg at Time 9:20:00 would have an event value of False.
- \* PCD01 Msg at Time 9:20:05, 9:20:15, & 9:20:25 would have an event value of True.
- \* PCD01 Msg at Time 9:20:32 and beyond would have an event value of false.

846

## 847

### 848 6.2.4 Example 1 – Minimal Message when Idle

849 This message is a periodic report of the dialysis machine status. It is reported when the machine is Idle or in  
850 Service.

```
851 MSH|^~\&|ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
852 64|||20191003092006+0000||ORU^R01^ORU_R01|20191003092005|P|2.6||AL|NE||||IHE_PCD_001^IH
853 EPCD^1.3.6.1.4.1.19376.1.6.1.1.1^ISO
854 PID|||Scrubber 2000/SC678932^^^"^^U||^"^^^U
855 OBR|1||080019FFFE3ED02D20110602045842^ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
856 64|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|||20191003092005+0000
857 OBX|1|ST|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|1.0.0|||||F
858 OBX|2|ST|531970^MDC_ID_MODEL_MANUFACTURER^MDC|1.0.0.1|NxStage|||||F
859 OBX|3|ST|531969^MDC_ID_MODEL_NUMBER^MDC|1.0.0.2|System One|||||F
860 OBX|4|ST|531972^MDC_ID_PROD_SPEC_SERIAL^MDC|1.0.0.3|1000478|||||F
861 OBX|5|ST|531975^MDC_ID_PROD_SPEC_SW^MDC|1.0.0.4|1.2.3.4|||||F
862 OBX|6|ST|70934^MDC_DEV_HDIALY_VMD^MDC|1.1|||||F
863 OBX|7|ST|70939^MDC_DEV_HDIALY_MACH_CONFIG_CHAN^MDC|1.1.1|||||F
864 OBX|8|DTM|67975^MDC_ATTR_TIME_ABS^MDC|1.1.1.1|20191003092005+0000|||||F
865 OBX|9|ST|158594^MDC_HDIALY_MACH_MODE_OF_OPERATION^MDC|1.1.1.3|IDL|||||F
```

## 866

### 867 6.2.5 Example 2 – Minimal HD Message when Treating

868 This message is a periodic report of a HD treatment status. It contains the minimum number of reported  
869 items.

870 Note, this standard reports Boolean event values in the PCD-01 transaction. This is done so that a complete  
871 record of the machine state is captured in the PCD-01. The state of the events is not considered clinically  
872 actionable. Clinically actionable alerts and events are communicated in the PCD-04 transaction.

```
873 MSH|^~\&|ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
874 64|||20191003092006+0000||ORU^R01^ORU_R01|20191003092005|P|2.6||AL|NE||||IHE_PCD_001^IH
875 EPCD^1.3.6.1.4.1.19376.1.6.1.1.1^ISO
876 PID|||Scrubber 2000/SC678932^^^"^^U||^"^^^U
877 OBR|1||080019FFFE3ED02D20110602045842^ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
878 64|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|||20191003092005+0000
879 OBX|1|ST|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|1.0.0|||||F
880 OBX|2|ST|531970^MDC_ID_MODEL_MANUFACTURER^MDC|1.0.0.1|NxStage|||||F
881 OBX|3|ST|531969^MDC_ID_MODEL_NUMBER^MDC|1.0.0.2|System One|||||F
882 OBX|4|ST|531972^MDC_ID_PROD_SPEC_SERIAL^MDC|1.0.0.3|1000478|||||F
883 OBX|5|ST|531975^MDC_ID_PROD_SPEC_SW^MDC|1.0.0.4|1.2.3.4|||||F
884 OBX|6|ST|70934^MDC_DEV_HDIALY_VMD^MDC|1.1|||||F
885 OBX|7|ST|70939^MDC_DEV_HDIALY_MACH_CONFIG_CHAN^MDC|1.1.1|||||F
886 OBX|8|DTM|67975^MDC_ATTR_TIME_ABS^MDC|1.1.1.1|20191003092005+0000|||||F
887 OBX|9|ST|158594^MDC_HDIALY_MACH_MODE_OF_OPERATION^MDC|1.1.1.3|TX|||||F
888 OBX|10|ST|158596^MDC_HDIALY_MACH_BLD_PUMP_ON^MDC|1.1.1.7|T|||||F
889 OBX|11|ST|158597^MDC_HDIALY_MACH_TX_FLUID_BYPASS^MDC|1.1.1.8|F|||||F
890 OBX|12|ST|158598^MDC_HDIALY_MACH_TX_MODALITY^MDC|1.1.1.9|HD|||||F
891 OBX|13|NM|158720^MDC_HDIALY_MACH_THERAPY_TIME^MDC|1.1.1.10|180|min^min^UCUM|||||F
892 OBX|14|NM|158724^MDC_HDIALY_MACH_TIME_REMAIN^MDC|1.1.1.11|600|min^min^UCUM|||||F
```

893 OBX|15|ST|70947^MDC\_DEV\_HDIALY\_BLOOD\_PUMP\_CHAN^MDC|1.1.3|||||F

894 OBX|16|NM|16935956^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE\_SETTING^MDC|1.1.3.2|250|ml/min^ml/min^

895 UCUM|||||F

896 OBX|17|NM|158744^MDC\_HDIALY\_BLD\_PRESS\_ART^MDC|1.1.3.4|-75|mm[Hg]^mm[Hg]^UCUM|< -200|||||F

897 OBX|18|ST|158604^MDC\_HDIALY\_BLD\_PUMP\_MODE^MDC|1.1.3.5|2N|||||F

898 OBX|19|ST|198242^MDC\_EVT\_HDIALY\_BLD\_PUMP\_STOP^MDC|1.1.3.6|F|||||F

899 OBX|20|NM|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|1.1.3.15|200|mm[Hg]^mm[Hg]^UCUM|20-

900 400|||||F

901 OBX|21|ST|70951^MDC\_DEV\_HDIALY\_FLUID\_CHAN^MDC|1.1.4|||||F

902 OBX|22|ST|158606^MDC\_HDIALY\_DIALYSATE\_FLOW\_MODE^MDC|1.1.4.1|CONST|||||F

903 OBX|23|NM|16936008^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE\_SETTING^MDC|1.1.4.2|100|ml/min^ml/min^UCUM|

904 |||||F

905 OBX|24|NM|158792^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE^MDC|1.1.4.3|99|ml/min^ml/min^UCUM|||||F

906 OBX|25|NM|158788^MDC\_HDIALY\_DIALYSATE\_COND^MDC|1.1.4.4|13.81|mS/cm^mS/cm^UCUM|||||F

907 OBX|26|ST|198244^MDC\_EVT\_HDIALY\_BLOOD\_LEAK^MDC|1.1.4.5|F|||||F

908 OBX|27|ST|70955^MDC\_DEV\_HDIALY\_FILTER\_CHAN^MDC|1.1.5|||||F

909 OBX|28|NM|158852^MDC\_HDIALY\_FILTER\_TRANSMEMBRANE\_PRESS^MDC|1.1.5.2|35|mm[Hg]^mm[Hg]^UCUM|||||

910 F

911 OBX|29|ST|70963^MDC\_DEV\_HDIALY\_SAFETY\_SYSTEMS\_CHAN^MDC|1.1.7|||||F

912 OBX|30|ST|198252^MDC\_EVT\_HDIALY\_SAFETY\_ART\_AIR\_DETECT^MDC|1.1.7.1|F|||||F

913 OBX|31|ST|198254^MDC\_EVT\_HDIALY\_SAFETY\_DIALYSATE\_AIR\_DETECT^MDC|1.1.7.2|F|||||F

914 OBX|32|ST|198258^MDC\_EVT\_HDIALY\_SAFETY\_SYSTEM\_GENERAL^MDC|1.1.7.3|F|||||F

915 OBX|33|ST|198216^MDC\_EVT\_SELFTEST\_FAILURE^MDC|1.1.7.4|F|||||F

916 OBX|34|ST|198262^MDC\_EVT\_HDIALY\_SAFETY\_VEN\_AIR\_DETECT^MDC|1.1.7.5|F|||||F

917 OBX|35|ST|70967^MDC\_DEV\_HDIALY\_THERAPY\_OUTCOMES\_CHAN^MDC|1.1.8|||||F

918 OBX|36|ST|158618^MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD^MDC|1.1.8.18|UF|||||F

919 OBX|37|ST|70971^MDC\_DEV\_HDIALY\_UF\_CHAN^MDC|1.1.9|||||F

920 OBX|38|NM|159028^MDC\_HDIALY\_UF\_TARGET\_VOL\_TO\_REMOVE^MDC|1.1.9.1|2000|ml^ml^UCUM|||||F

921 OBX|39|NM|159032^MDC\_HDIALY\_UF\_ACTUAL\_REMOVED\_VOL^MDC|1.1.9.2|555|ml^ml^UCUM|||||F

922 OBX|40|ST|158619^MDC\_HDIALY\_UF\_MODE^MDC|1.1.9.3|CONST-WT|||||F

923 OBX|41|NM|159036^MDC\_HDIALY\_UF\_RATE^MDC|1.1.9.4|100|ml/h^ml/h^UCUM|||||F

924 OBX|42|NM|16936252^MDC\_HDIALY\_UF\_RATE\_SETTING^MDC|1.1.9.5|100|ml/h^ml/h^UCUM|||||F

925 OBX|43|ST|198276^MDC\_EVT\_HDIALY\_UF\_RATE\_RANGE^MDC|1.1.9.6|F|||||F

926

### 6.2.6 Example 3 – Full HD Message

This message reports all mandatory and optional fields for a pre/post dilution HDF therapy.

929 MSH|^~\&|ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-

930 64|||||20191003092025+0000||ORU^R01^ORU\_R01|20191003092024|P|2.6|||AL|NE|||||IHE\_PCD\_001^IH

931 EPCD^1.3.6.1.4.1.19376.1.6.1.1.1^ISO

932 PID|||Scrubber 2000/SC678932^^^"^^U||^^^^^^U

933 OBX|1||080019FFFE3ED02D20110602045842^ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-

934 64|70929^MDC\_DEV\_HDIALY\_MACHINE\_MDS^MDC|||20191003092024+0000

935 OBX|1|ST|70929^MDC\_DEV\_HDIALY\_MACHINE\_MDS^MDC|1.0.0|||||F

936 OBX|2|ST|531970^MDC\_ID\_MODEL\_MANUFACTURER^MDC|1.0.0.1|NxStage|||||F

937 OBX|3|ST|531969^MDC\_ID\_MODEL\_NUMBER^MDC|1.0.0.2|System One|||||F

938 OBX|4|ST|531972^MDC\_ID\_PROD\_SPEC\_SERIAL^MDC|1.0.0.3|1000478|||||F

939 OBX|5|ST|531975^MDC\_ID\_PROD\_SPEC\_SW^MDC|1.0.0.4|1.2.3.4|||||F

940 OBX|6|ST|67916^MDC\_ATTR\_ID\_UDI^MDC|1.0.0.6|+M535NX10003A0/\$\$+735241/16D20180305J|||||F

941 OBX|7|ST|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.1|||||F

942 OBX|8|ST|70939^MDC\_DEV\_HDIALY\_MACH\_CONFIG\_CHAN^MDC|1.1.1|||||F

943 OBX|9|DTM|67975^MDC\_ATTR\_TIME\_ABS^MDC|1.1.1.1|20191003092024+0000|||||F

944 OBX|10|ST|158593^MDC\_HDIALY\_MACH\_MODE\_DESCRIPTION^MDC|1.1.1.2|NxStage System One|||||F

945 OBX|11|ST|158594^MDC\_HDIALY\_MACH\_MODE\_OF\_OPERATION^MDC|1.1.1.3|TX|||||F

946 OBX|12|NM|184195^MDC\_TIME\_PD\_MAINTENANCE\_TO\_NEXT\_SERVICE^MDC|1.1.1.4|12|h^h^UCUM|||||F

947 OBX|13|DTM|184199^MDC\_MAINTENANCE\_NEXT\_SERVICE\_DATE^MDC|1.1.1.5|20191003|||||F

948 OBX|14|NM|158595^MDC\_HDIALY\_MACH\_MAINT\_TX\_REMAIN^MDC|1.1.1.6|45|||||F

949 OBX|15|ST|158596^MDC\_HDIALY\_MACH\_BLD\_PUMP\_ON^MDC|1.1.1.7|T|||||F

950 OBX|16|ST|158597^MDC\_HDIALY\_MACH\_TX\_FLUID\_BYPASS^MDC|1.1.1.8|F|||||F

951 OBX|17|ST|158598^MDC\_HDIALY\_MACH\_TX\_MODALITY^MDC|1.1.1.9|HDF|||||F

952 OBX|18|NM|158720^MDC\_HDIALY\_MACH\_THERAPY\_TIME^MDC|1.1.1.10|180|min^min^UCUM|||||F

953 OBX|19|NM|158724^MDC\_HDIALY\_MACH\_TIME\_REMAIN^MDC|1.1.1.11|600|min^min^UCUM|||||F

954 OBX|20|NM|188508^MDC\_TEMP\_ROOM^MDC|1.1.1.12|20|Cel^Cel^UCUM|||||F

955 OBX|21|ST|70943^MDC\_DEV\_HDIALY\_ANTICOAG\_PUMP\_CHAN^MDC|1.1.2|||||F

956 OBX|22|ST|198236^MDC\_EVT\_HDIALY\_ANTICOAG\_STOP^MDC|1.1.2.4|F|||||F

957 OBX|23|ST|158599^MDC\_HDIALY\_ANTICOAG\_NAME^MDC|1.1.2.5|heparin|||||F

958 OBX|24|ST|158600^MDC\_HDIALY\_ANTICOAG\_MODE^MDC|1.1.2.6|CON|||||F

959 OBX|25|NM|0^MDC\_HDIALY\_ANTICOAG\_INFUS\_RATE\_SETTING^MDC|1.1.2.7|50.0|ml/hr^ml/hr^UCUM|||||F

960 OBX|26|NM|158736^MDC\_HDIALY\_ANTICOAG\_INFUS\_RATE^MDC|1.1.2.8|50.0|ml/hr^ml/hr^UCUM|||||F

961 OBX|27|NM|158728^MDC\_HDIALY\_ANTICOAG\_ACCUM\_DELIV^MDC|1.1.2.9|60.0|ml^ml^UCUM|||||F

962 OBX|28|NM|68142^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.2.10|1|||||F

963 OBX|29|ST|198238^MDC\_EVT\_HDIALY\_ANTICOAG\_SYRINGE\_EMPTY^MDC|1.1.2.11|F|||||F

964 OBX|30|ST|158602^MDC\_HDIALY\_ANTICOAG\_SYRINGE\_BRAND^MDC|1.1.2.12|Fishman|||||F

965 OBX|31|NM|158603^MDC\_HDIALY\_ANTICOAG\_SYRINGE\_VOL^MDC|1.1.2.13|60|ml^ml^UCUM|||||F

966 OBX|32|ST|198240^MDC\_EVT\_HDIALY\_ANTICOAG\_SYRINGE\_SIZE^MDC|1.1.2.14|F|||||F

967 OBX|33|ST|70947^MDC\_DEV\_HDIALY\_BLOOD\_PUMP\_CHAN^MDC|1.1.3|||||F

968 OBX|34|NM|158740^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE^MDC|1.1.3.1|250|ml/min^ml/min^UCUM|||||F

969 OBX|35|NM|16935956^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE\_SETTING^MDC|1.1.3.2|250|ml/min^ml/min^UCUM|||||F

970

971 OBX|36|NM|158743^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_FLOW\_RATE\_MEAN^MDC|1.1.3.3|250|ml/min^ml/min^UCUM|

972 |||F

973 OBX|37|NM|158744^MDC\_HDIALY\_BLD\_PRESS\_ART^MDC|1.1.3.4|-75|mm[Hg]^mm[Hg]^UCUM|< -200|||F

974 OBX|38|ST|158604^MDC\_HDIALY\_BLD\_PUMP\_MODE^MDC|1.1.3.5|2N|||||F

975 OBX|39|ST|198242^MDC\_EVT\_HDIALY\_BLD\_PUMP\_STOP^MDC|1.1.3.6|F|||||F

976 OBX|40|ST|158605^MDC\_HDIALY\_BLD\_PUMP\_TUBING\_SIZE^MDC|1.1.3.7|8 mm|||||F

977 OBX|41|NM|158748^MDC\_HDIALY\_BLOOD\_TEMP\_ART^MDC|1.1.3.8|39.1|Cel^Cel^UCUM|||||F

978 OBX|42|NM|158752^MDC\_HDIALY\_BLD\_PUMP\_CHANGE\_IN\_ENERGY^MDC|1.1.3.9|1|kJ/h^kJ/h^UCUM|||||F

979 OBX|43|NM|158756^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_ART\_POST\_PUMP^MDC|1.1.3.10|100|mm[Hg]^mm[Hg]^UCUM|

980 |||F

981 OBX|44|NM|158760^MDC\_HDIALY\_BLD\_PUMP\_PRIMING\_VOL^MDC|1.1.3.11|191|ml^ml^UCUM|||||F

982 OBX|45|NM|158764^MDC\_HDIALY\_BLD\_PUMP\_SINGLE\_NEEDLE\_PRESS^MDC|1.1.3.12|200|mm[Hg]^ml[Hg]^UCUM|

983 |||F

984 OBX|46|NM|158772^MDC\_HDIALY\_BLD\_PUMP\_BLOOD\_PROCESSED\_TOTAL^MDC|1.1.3.14|120|L^L^UCUM|||||F

985 OBX|47|NM|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|1.1.3.15|200|mm[Hg]^mm[Hg]^UCUM|20-

986 400|||F

987 OBX|48|NM|158780^MDC\_HDIALY\_BLOOD\_TEMP\_VEN^MDC|1.1.3.16|39.1|Cel^Cel^UCUM|||||F

988 OBX|49|NM|16935996^MDC\_HDIALY\_BLOOD\_TEMP\_VEN\_SETTING^MDC|1.1.3.17|39.0|Cel^Cel^UCUM|||||F

989 OBX|50|ST|70951^MDC\_DEV\_HDIALY\_FLUID\_CHAN^MDC|1.1.4|||||F

990 OBX|51|NM|158784^MDC\_HDIALY\_BICARB\_COND^MDC|1.1.4.1|13.81|mS/cm^mS/cm^UCUM||||F

991 OBX|52|NM|16936000^MDC\_HDIALY\_BICARB\_COND\_SETTING^MDC|1.1.4.2|13.81|mS/cm^mS/cm^UCUM||||F

992 OBX|53|NM|158788^MDC\_HDIALY\_DIALYSATE\_COND^MDC|1.1.4.3|13.81|mS/cm^mS/cm^UCUM||||F

993 OBX|54|NM|16936004^MDC\_HDIALY\_DIALYSATE\_COND\_SETTING^MDC|1.1.4.4|13.81|mS/cm^mS/cm^UCUM||||F

994 OBX|55|NM|158792^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE^MDC|1.1.4.5|99|ml/min^ml/min^UCUM||||F

995 OBX|56|NM|16936008^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE\_SETTING^MDC|1.1.4.6|100|ml/min^ml/min^UCUM|  
996 |||F

997 OBX|57|ST|158606^MDC\_HDIALY\_DIALYSATE\_FLOW\_MODE^MDC|1.1.4.7|CONST|||||F

998 OBX|58|NM|158800^MDC\_HDIALY\_DIALYSATE\_AMMONIA^MDC|1.1.4.10|0.00|[ppm]^ [ppm]^UCUM||||F

999 OBX|59|NM|158795^MDC\_HDIALY\_DIALYSATE\_FLOW\_RATE\_MEAN^MDC|1.1.4.11|100|ml/min^ml/min^UCUM||||  
1000 F

1001 OBX|60|NM|16936020^MDC\_HDIALY\_CONC\_HCO3\_SETTING^MDC|1.1.4.13|32.0|mmol/L^mmol/L^UCUM||||F

1002 OBX|61|ST|158607^MDC\_HDIALY\_CONC\_HCO3\_MODE^MDC|1.1.4.14|PRO|||||F

1003 OBX|62|ST|198244^MDC\_EVT\_HDIALY\_BLOOD\_LEAK^MDC|1.1.4.15|F|||||F

1004 OBX|63|ST|158608^MDC\_HDIALY\_DIALYSATE\_NAME^MDC|1.1.4.16|FMC smartbag 111.5||Bbraun Duosol  
1005 35||||F

1006 OBX|64|NM|158808^MDC\_HDIALY\_DIALYSATE\_CONC\_ACETATE^MDC|1.1.4.17|3.00|mmol/L^mmol/L^UCUM||||F

1007 OBX|65|NM|158812^MDC\_HDIALY\_DIALYSATE\_CONC\_HCO3^MDC|1.1.4.18|35.0|mmol/L^mmol/L^UCUM||||F

1008 OBX|66|NM|158816^MDC\_HDIALY\_DIALYSATE\_CONC\_CHLORIDE^MDC|1.1.4.19|108.00|mmol/L^mmol/L^UCUM|||  
1009 ||F|||20191003085024+0000

1010 OBX|67|NM|158820^MDC\_HDIALY\_DIALYSATE\_CONC\_MG^MDC|1.1.4.20|0.5|mmol/L^mmol/L^UCUM||||F|||201  
1011 91003085024+0000

1012 OBX|68|NM|158824^MDC\_HDIALY\_DIALYSATE\_CONC\_CA^MDC|1.1.4.21|1.50|mmol/L^mmol/L^UCUM||||F|||20  
1013 191003085024+0000

1014 OBX|69|NM|158828^MDC\_HDIALY\_DIALYSATE\_CONC\_CITRATE^MDC|1.1.4.22|0.0|mmol/L^mmol/L^UCUM||||F

1015 OBX|70|NM|158832^MDC\_HDIALY\_DIALYSATE\_CONC\_GLU^MDC|1.1.4.23|1.00|mmol/L^mmol/L^UCUM||||F

1016 OBX|71|NM|158836^MDC\_HDIALY\_DIALYSATE\_CONC\_K^MDC|1.1.4.24|1.00|mmol/L^mmol/L^UCUM||||F|||201  
1017 91003085024+0000

1018 OBX|72|NM|158840^MDC\_HDIALY\_DIALYSATE\_CONC\_NA^MDC|1.1.4.25|138|mmol/L^mmol/L^UCUM||||F|||201  
1019 91003085024+0000

1020 OBX|73|NM|16936056^MDC\_HDIALY\_DIALYSATE\_CONC\_NA\_SETTING^MDC|1.1.4.26|140|mmol/L^mmol/L^UCUM||  
1021 |||F

1022 OBX|74|ST|158609^MDC\_HDIALY\_DIALYSATE\_CONC\_NA\_MODE^MDC|1.1.4.27|CONST|||||F

1023 OBX|75|NM|158844^MDC\_HDIALY\_DIALYSATE\_CONC\_PH^MDC|1.1.4.28|7.0|[pH]^ [pH]^UCUM||||F

1024 OBX|76|NM|158848^MDC\_HDIALY\_DIALYSATE\_VOL\_DELIV^MDC|1.1.4.29|24.34|L^L^UCUM|||||F

1025 OBX|77|ST|70955^MDC\_DEV\_HDIALY\_FILTER\_CHAN^MDC|1.1.5|||||F

1026 OBX|78|ST|158610^MDC\_HDIALY\_FILTER\_NAME^MDC|1.1.5.1|NxStage CAR125|||||F

1027 OBX|79|NM|158852^MDC\_HDIALY\_FILTER\_TRANSMEMBRANE\_PRESS^MDC|1.1.5.2|35|mm[Hg]^mm[Hg]^UCUM||||  
1028 F

1029 OBX|80|NM|158611^MDC\_HDIALY\_FILTER\_NUM^MDC|1.1.5.3|1||| |F

1030 OBX|81|ST|158612^MDC\_HDIALY\_FILTER\_UDI^MDC|1.1.5.4|  
1031 (01)00842289101845(17)201200(10)LOT00606|| |F

1032 OBX|82|ST|70959^MDC\_DEV\_HDIALY\_CONVECTIVE\_CHAN^MDC|1.1.6|||||F

1033 OBX|83|ST|158613^MDC\_HDIALY\_RF\_DILUTION\_LOCATION^MDC|1.1.6.1|PREF-POSTF|||||F

1034 OBX|84|ST|158614^MDC\_HDIALY\_RF\_POST\_FILTER\_NAME^MDC|1.1.6.2|multiBIC|| |F

1035 OBX|85|NM|158856^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_ACETATE^MDC|1.1.6.3|0|mmol/L^mmol/L^UCUM||||  
1036 F

1037 OBX|86|NM|158860^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_HCO3^MDC|1.1.6.4|35.0|mmol/L^mmol/L^UCUM||||  
1038 F

1039 OBX|87|NM|158864^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_CHLORIDE^MDC|1.1.6.5|111.00|mmol/L^mmol/L^UCU



1040 M||||F

1041 OBX|88|NM|158868^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_MG^MDC|1.1.6.6|0.5|mmol/L^mmol/L^UCUM||||F

1042 OBX|89|NM|158824^MDC\_HDIALY\_DIALYSATE\_CONC\_CA^MDC|1.1.6.7|1.50|mmol/L^mmol/L^UCUM||||F

1043 OBX|90|NM|158876^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_CITRATE^MDC|1.1.6.8|0.0|mmol/L^mmol/L^UCUM|||

1044 |F

1045 OBX|91|NM|158880^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_GLU^MDC|1.1.6.9|5.55|mmol/L^mmol/L^UCUM||||F

1046 OBX|92|NM|158884^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_K^MDC|1.1.6.10|2.00|mmol/L^mmol/L^UCUM||||F

1047 OBX|93|NM|158888^MDC\_HDIALY\_RF\_POST\_FILTER\_CONC\_NA^MDC|1.1.6.11|140|mmol/L^mmol/L^UCUM||||F

1048 OBX|94|NM|158892^MDC\_HDIALY\_RF\_POST\_FILTER\_FLOW\_RATE^MDC|1.1.6.12|60|ml/min^ml/min^UCUM||||F

1049 OBX|95|NM|16936108^MDC\_HDIALY\_RF\_POST\_FILTER\_FLOW\_RATE\_SETTING^MDC|1.1.6.13|60|ml/min^ml/min^

1050 UCUM||||F

1051 OBX|96|NM|158895^MDC\_HDIALY\_RF\_POST\_FILTER\_FLOW\_RATE\_MEAN^MDC|1.1.6.14|60|ml/min^ml/min^UCUM|

1052 |||F

1053 OBX|97|NM|158896^MDC\_HDIALY\_RF\_POST\_FILTER\_TEMP^MDC|1.1.6.15|39.9|Cel^Cel^UCUM||||F

1054 OBX|98|NM|16936112^MDC\_HDIALY\_RF\_POST\_FILTER\_TEMP\_SETTING^MDC|1.1.6.16|41.0|Cel^Cel^UCUM||||

1055 F

1056 OBX|99|NM|158900^MDC\_HDIALY\_RF\_POST\_FILTER\_VOL^MDC|1.1.6.17|6.00|L^L^UCUM||||F

1057 OBX|100|NM|16936116^MDC\_HDIALY\_RF\_POST\_FILTER\_VOL\_SETTING^MDC|1.1.6.18|12.00|L^L^UCUM||||F

1058 OBX|101|ST|158615^MDC\_HDIALY\_RF\_PRE\_FILTER\_NAME^MDC|1.1.6.19|multiBIC|||||F

1059 OBX|102|NM|158904^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_ACETATE^MDC|1.1.6.20|0|mmol/L^mmol/L^UCUM|||

1060 |F

1061 OBX|103|NM|158908^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_HCO3^MDC|1.1.6.21|35.0|mmol/L^mmol/L^UCUM|/L|

1062 |||F

1063 OBX|104|NM|158912^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_CHLORIDE^MDC|1.1.6.22|111.00|mmol/L^mmol/L^UC

1064 UM||||F

1065 OBX|105|NM|158916^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_MG^MDC|1.1.6.23|0.5|mmol/L^mmol/L^UCUM||||F

1066 OBX|106|NM|158920^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_CA^MDC|1.1.6.24|1.50|mmol/L^mmol/L^UCUM||||F

1067 OBX|107|NM|158924^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_CITRATE^MDC|1.1.6.25|0.0|mmol/L^mmol/L^UCUM||

1068 |||F

1069 OBX|108|NM|158928^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_GLU^MDC|1.1.6.26|5.55|mmol/L^mmol/L^UCUM|||

1070 F

1071 OBX|109|NM|158932^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_K^MDC|1.1.6.27|2.00|mmol/L^mmol/L^UCUM||||F

1072 OBX|110|NM|158936^MDC\_HDIALY\_RF\_PRE\_FILTER\_CONC\_NA^MDC|1.1.6.28|140|mmol/L^mmol/L^UCUM||||F

1073 OBX|111|NM|158940^MDC\_HDIALY\_RF\_PRE\_FILTER\_FLOW\_RATE^MDC|1.1.6.29|60|ml/min^ml/min^UCUM||||F

1074 OBX|112|NM|16936156^MDC\_HDIALY\_RF\_PRE\_FILTER\_FLOW\_RATE\_SETTING^MDC|1.1.6.30|60|ml/min^ml/min^

1075 UCUM||||F

1076 OBX|113|NM|158943^MDC\_HDIALY\_RF\_PRE\_FILTER\_FLOW\_RATE\_MEAN^MDC|1.1.6.31|60|ml/min^ml/min^UCUM|

1077 |||F

1078 OBX|114|NM|158944^MDC\_HDIALY\_RF\_PRE\_FILTER\_TEMP^MDC|1.1.6.32|39.9|Cel^Cel^UCUM||||F

1079 OBX|115|NM|16936160^MDC\_HDIALY\_RF\_PRE\_FILTER\_TEMP\_SETTING^MDC|1.1.6.33|41.0|Cel^Cel^UCUM||||

1080 F

1081 OBX|116|NM|158948^MDC\_HDIALY\_RF\_PRE\_FILTER\_VOL^MDC|1.1.6.34|1.23|L^L^UCUM||||F

1082 OBX|117|NM|16936164^MDC\_HDIALY\_RF\_PRE\_FILTER\_VOL\_SETTING^MDC|1.1.6.35|3.00|L^L^UCUM||||F

1083 OBX|118|NM|158952^MDC\_HDIALY\_RF\_CONV\_CLEARANCE^MDC|1.1.6.36|31.6|ml/min^ml/min^UCUM||||F

1084 OBX|119|ST|198246^MDC\_EVT\_HDIALY\_RF\_EXCESS\_DELIV^MDC|1.1.6.37|F|||||F

1085 OBX|120|ST|198248^MDC\_EVT\_HDIALY\_RF\_INSUFF\_DELIV^MDC|1.1.6.38|F|||||F

1086 OBX|121|NM|16936172^MDC\_HDIALY\_RF\_PRE\_POST\_FLOW\_RATIO\_SETTING^MDC|1.1.6.39|3.00|||||F

1087 OBX|122|NM|158960^MDC\_HDIALY\_RF\_BOLUS\_RATE^MDC|1.1.6.40|0|ml/min^ml/min^UCUM||||F

1088 OBX|123|NM|16936180^MDC\_HDIALY\_RF\_BOLUS\_VOL\_SETTING^MDC|1.1.6.42|3.00|mL^mL^UCUM||||F

1089 OBX|124|NM|158968^MDC\_HDIALY\_RF\_BOLUS\_VOL\_DELIVERED^MDC|1.1.6.41|0.00|mL^mL^UCUM||||F

1090 OBX|125|ST|158616^MDC\_HDIALY\_RF\_FLOW\_MODE^MDC|1.1.6.43|CONST|||||F

1091 OBX|126|ST|158617^MDC\_HDIALY\_RF\_SOURCE^MDC|1.1.6.44|BAG||||||F

1092 OBX|127|ST|70963^MDC\_DEV\_HDIALY\_SAFETY\_SYSTEMS\_CHAN^MDC|1.1.7|||||F

1093 OBX|128|ST|198252^MDC\_EVT\_HDIALY\_SAFETY\_ART\_AIR\_DETECT^MDC|1.1.7.1|F|||||F

1094 OBX|129|ST|198254^MDC\_EVT\_HDIALY\_SAFETY\_DIALYSATE\_AIR\_DETECT^MDC|1.1.7.2|F|||||F

1095 OBX|130|ST|198256^MDC\_EVT\_HDIALY\_SAFETY\_DIALYSATE\_COMPOSITION^MDC|1.1.7.3|F|||||F

1096 OBX|131|ST|198258^MDC\_EVT\_HDIALY\_SAFETY\_SYSTEM\_GENERAL^MDC|1.1.7.4|F|||||F

1097 OBX|132|ST|198216^MDC\_EVT\_SELFTEST\_FAILURE^MDC|1.1.7.5|F|||||F

1098 OBX|133|ST|198260^MDC\_EVT\_HDIALY\_SAFETY\_VEN\_ACCESS^MDC|1.1.7.6|F|||||F

1099 OBX|134|ST|198262^MDC\_EVT\_HDIALY\_SAFETY\_VEN\_AIR\_DETECT^MDC|1.1.7.7|F|||||F

1100 OBX|135|ST|198264^MDC\_EVT\_HDIALY\_SAFETY\_WETNESS\_DETECT\_ALERT^MDC|1.1.7.8|F|||||F

1101 OBX|136|ST|198266^MDC\_EVT\_HDIALY\_SAFETY\_WETNESS\_DETECT\_ERROR^MDC|1.1.7.9|F|||||F

1102 OBX|137|ST|68489^MDC\_ATTR\_ALERT\_ID\_NUM^MDC|1.1.7.10|0|||||F

1103 OBX|138|ST|68546^MDC\_ATTR\_ALERT\_TEXT^MDC|1.1.7.11|||||F

1104 OBX|139|ST|70967^MDC\_DEV\_HDIALY\_THERAPY\_OUTCOMES\_CHAN^MDC|1.1.8|||||F

1105 OBX|140|NM|158972^MDC\_HDIALY\_THERAPY\_MASS\_TRF\_AREA\_COEFF^MDC|1.1.8.1|700|||||F

1106 OBX|141|NM|158976^MDC\_HDIALY\_THERAPY\_ACCESS\_FLOW^MDC|1.1.8.2|250|ml/min^ml/min^UCUM|||||F

1107 OBX|142|NM|158980^MDC\_HDIALY\_THERAPY\_RATIO\_EKT\_OVER\_V\_DELIVERED^MDC|1.1.8.4|1.1|%^%^UCUM

1108 |||||F

1109 OBX|143|NM|158984^MDC\_HDIALY\_THERAPY\_KT\_DELIVERED^MDC|1.1.8.5|42.0|L^L^UCUM|||||F

1110 OBX|144|NM|158988^MDC\_HDIALY\_THERAPY\_RATIO\_SPKT\_OVER\_V\_DELIVERED^MDC|1.1.8.6|1.1|%^%^UCUM

1111 |||||F

1112 OBX|145|ST|198268^MDC\_EVT\_HDIALY\_THERAPY\_PAT\_TX\_GENERAL^MDC|1.1.8.7|F|||||F

1113 OBX|146|NM|158992^MDC\_HDIALY\_THERAPY\_RATIO\_KT\_OVER\_V\_GOAL^MDC|1.1.8.10|1.21|%^%^UCUM |||||F

1114 OBX|147|NM|159019^MDC\_HDIALY\_THERAPY\_UREA\_CLEARANCE\_MEAN^MDC|1.1.8.11|70|||||F

1115 OBX|148|NM|159000^MDC\_HDIALY\_THERAPY\_BODY\_START\_WT^MDC|1.1.8.12|75.9|kg^kg^UCUM|||||F

1116 OBX|149|NM|159004^MDC\_HDIALY\_THERAPY\_PCT\_RECIRC^MDC|1.1.8.13|20|%^%^UCUM|||||F

1117 OBX|150|NM|159008^MDC\_HDIALY\_THERAPY\_PLASMA\_NA\_CONC^MDC|1.1.8.14|140|mmol/L^mmol/L^UCUM|||||F

1118 OBX|151|NM|159012^MDC\_HDIALY\_THERAPY\_RATIO\_SPKT\_OVER\_V\_PROJECTED^MDC|1.1.8.16|1.1|%^%^UCUM

1119 |||||F

1120 OBX|152|ST|158618^MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD^MDC|1.1.8.18|UF|||||F

1121 OBX|153|ST|198270^MDC\_EVT\_HDIALY\_THERAPY\_TX\_END\_TIME^MDC|1.1.8.19|F|||||F

1122 OBX|154|NM|159016^MDC\_HDIALY\_THERAPY\_UREA\_CLEARANCE^MDC|1.1.8.20|196|||||F

1123 OBX|155|NM|159020^MDC\_HDIALY\_THERAPY\_BODY\_END\_WT^MDC|1.1.8.21|75.9|kg^kg^UCUM|||||F

1124 OBX|156|ST|70971^MDC\_DEV\_HDIALY\_UF\_CHAN^MDC|1.1.9|||||F

1125 OBX|157|ST|198272^MDC\_EVT\_HDIALY\_UF\_LO^MDC|1.1.9.1|F|||||F

1126 OBX|158|ST|198274^MDC\_EVT\_HDIALY\_UF\_NEG^MDC|1.1.9.2|F|||||F

1127 OBX|159|NM|159024^MDC\_HDIALY\_UF\_TIME\_TO\_TARGET^MDC|1.1.9.3|55|min^min^UCUM|||||F

1128 OBX|160|NM|159028^MDC\_HDIALY\_UF\_TARGET\_VOL\_TO\_REMOVE^MDC|1.1.9.4|2000|ml^ml^UCUM|||||F

1129 OBX|161|NM|159032^MDC\_HDIALY\_UF\_ACTUAL\_REMOVED\_VOL^MDC|1.1.9.5|555|ml^ml^UCUM|||||F

1130 OBX|162|ST|198276^MDC\_EVT\_HDIALY\_UF\_RATE\_RANGE^MDC|1.1.9.6|F|||||F

1131 OBX|163|ST|198278^MDC\_EVT\_HDIALY\_UF\_GOAL\_MET^MDC|1.1.9.7|F|||||F

1132 OBX|164|ST|158619^MDC\_HDIALY\_UF\_MODE^MDC|1.1.9.8|CONST-WT|||||F

1133 OBX|165|NM|159036^MDC\_HDIALY\_UF\_RATE^MDC|1.1.9.9|100|ml/h^ml/h^UCUM|||||F

1134 OBX|166|NM|16936252^MDC\_HDIALY\_UF\_RATE\_SETTING^MDC|1.1.9.9|100|ml/h^ml/h^UCUM|||||F

1135 OBX|167|ST|198276^MDC\_EVT\_HDIALY\_UF\_RATE\_RANGE^MDC|1.1.9.10|F|||||F

1136 OBX|168|NM|16936257^MDC\_HDIALY\_UF\_RATE\_LIMIT\_HIGH\_SETTING^MDC|1.1.9.11|150|ml/h^ml/h^UCUM|||

1137 |F

1138 OBX|169|NM|16936259^MDC\_HDIALY\_UF\_RATE\_LIMIT\_LOW\_SETTING^MDC|1.1.9.12|90|ml/h^ml/h^UCUM|||||F  
 1139 OBX|170|ST|70686^MDC\_DEV\_PRESS\_BLD\_NONINV\_VMD^MDC|1.2|||||F  
 1140 OBX|171|ST|70687^MDC\_DEV\_PRESS\_BLD\_NONINV\_CHAN^MDC|1.2.1|||||F  
 1141 OBX|172|NM|67979^MDC\_ATTR\_TIME\_PD\_MSMT^MDC|1.2.1.1|30|min^min^UCUM|||||F||20191003085024+000  
 1142 0  
 1143 OBX|173|ST|68135^MDC\_ATTR\_PT\_BODY\_POSN^MDC|1.2.1.2|SUPINE|||||F||20191003085024+0000  
 1144 OBX|174|NM|150022^MDC\_PRESS\_BLD\_NONINV\_DIA^MDC|1.2.1.3|80|mm[Hg]^mm[Hg]^UCUM|||||F||20191003  
 1145 085024+0000  
 1146 OBX|175|NM|149546^MDC\_PULS\_RATE\_NON\_INV^MDC|1.2.1.4|70|{beats}/min^{beats}/min^UCUM|||||F||2  
 1147 0191003085024+0000  
 1148 OBX|176|NM|150023^MDC\_PRESS\_BLD\_NONINV\_MEAN^MDC|1.2.1.5|100|mm[Hg]^mm[Hg]^UCUM|||||F||201910  
 1149 03085024+0000  
 1150 OBX|177|NM|150021^MDC\_PRESS\_BLD\_NONINV\_SYS^MDC|1.2.1.6|120|mm[Hg]^mm[Hg]^UCUM|||||F||2019100  
 1151 3085024+0000  
 1152 OBX|178|ST|69642^MDC\_DEV\_ANALY\_SAT\_O2\_VMD^MDC|1.3|||||F  
 1153 OBX|179|ST|69643^MDC\_DEV\_ANALY\_SAT\_O2\_CHAN^MDC|1.3.1|||||F  
 1154 OBX|180|NM|150456^MDC\_PULS\_OXIM\_SAT\_O2^MDC|1.3.1.1|98|%%^UCUM|||||F  
 1155 OBX|181|NM|149530^MDC\_PULS\_OXIM\_PULS\_RATE^MDC|1.3.1.2|67|{beats}/min^{beats}/min^UCUM|||||F  
 1156 OBX|182|ST|196638^MDC\_EVT\_ERR^MDC|1.3.1.3|F|||||F  
 1157 OBX|183|ST|70974^MDC\_DEV\_BLOOD\_CHEM\_VMD^MDC|1.4|||||F  
 1158 OBX|184|ST|70975^MDC\_DEV\_BLOOD\_CHEM\_CHAN^MDC|1.4.1|||||F  
 1159 OBX|185|ST|158620^MDC\_HDIALY\_PLASMA\_VOL\_MARKER^MDC|1.4.1.1|NONE|||||F  
 1160 OBX|186|ST|158621^MDC\_HDIALY\_PLASMA\_VOL\_PROFILE^MDC|1.4.1.2|FLAT|||||F  
 1161 OBX|187|NM|160132^MDC\_CONC\_HCT\_GEN^MDC|1.4.1.3|0.45|%{vol}^{vol}^UCUM|||||F  
 1162 OBX|188|NM|160120^MDC\_CONC\_HB\_GEN^MDC|1.4.1.4|13.6|g/dL^g/dL^UCUM|||||F  
 1163 OBX|189|NM|159044^MDC\_HDIALY\_REL\_BLOOD\_VOL^MDC|1.4.1.5|35|%%^UCUM|||||F  
 1164 OBX|190|NM|150316^MDC\_SAT\_O2^MDC|1.4.1.6|98|%%^UCUM|||||F

1165 **6.2.7 Example 4 – PD Message**

1166 In this message the PD treatment is dwelling in the second of three exchanges.

1167 MSH|^~\&|ACME\_Dialysis\_Machine^025041FFFE000001^EUI-  
 1168 64|||20241217125311+0000||ORU^R01^ORU\_R01|20241217125311437|P|2.6||AL|NE||||IHE\_PCD\_001^IHEPCD^1.  
 1169 3.6.1.4.12559.11.1.1.129^ISO  
 1170 PID||Shifter 100/19640306^^^U||^^^U||  
 1171 OBR|1||025041FFFE00000120241217125311^ACME\_Dialysis\_Machine^025041FFFE000001^EUI-  
 1172 64|27^MDCX\_DEV\_PDIALY\_MACHINE\_MDS^MDC||20241217125311+0000  
 1173 OBX|1|ST|27^MDCX\_DEV\_PDIALY\_MACHINE\_MDS^MDC|1|||||F  
 1174 OBX|2|ST|531969^MDC\_ID\_MODEL\_NUMBER^MDC|1.0.0.1|Shifter 100|||||F  
 1175 OBX|3|ST|531970^MDC\_ID\_MODEL\_MANUFACTURER^MDC|1.0.0.2|Acme|||||F  
 1176 OBX|4|ST|531972^MDC\_ID\_PROD\_SPEC\_SERIAL^MDC|1.0.0.3|19640306|||||F  
 1177 OBX|5|ST|531975^MDC\_ID\_PROD\_SPEC\_SW^MDC|1.0.0.4|1.2.3.4|||||F  
 1178 OBX|6|ST|0^MDCX\_DEV\_PDIALY\_VMD^MDC|1.1|||||F  
 1179 OBX|7|ST|1^MDCX\_DEV\_PDIALY\_MACH\_CONFIG\_CHAN^MDC|1.1.1|||||F  
 1180 OBX|8|DTM|67975^MDC\_ATTR\_TIME\_ABS^MDC|1.1.1.1|20241217125311+0000|||||F  
 1181 OBX|9|ST|4^MDCX\_DEV\_PDIALY\_TREATMENT\_CHAN^MDC|1.1.2|||||F  
 1182 OBX|10|ST|5^MDCX\_DEV\_PDIALY\_TREAT\_TYPE^MDC|1.1.2.1|CAPD|||||F  
 1183 OBX|11|ST|8^MDCX\_DEV\_PDIALY\_CURRENT\_PHASE^MDC|1.1.2.2|DWELL|||||F  
 1184 OBX|12|ST|6^MDCX\_DEV\_PDIALY\_TREAT\_LOCATION^MDC|1.1.2.3|HOME|||||F  
 1185 OBX|13|ST|7^MDCX\_DEV\_PDIALY\_PRESCRIPTION\_TYPE^MDC|1.1.2.4|EXCH|||||F

1186 OBX|14|ST|51^MDCX\_DEV\_PDIALY\_FLUID\_1\_CHAN^MDC|1.1.3|||||F

1187 OBX|15|ST|29^MDCX\_PDIALY\_FLUID\_NAME^MDC|1.1.3.1|Dialneal Low Calcium 2.5%|||||F

1188 OBX|16|ST|10^MDCX\_PDIALY\_FLUID\_SOURCE^MDC|1.1.3.2|BAG|||||F

1189 OBX|17|NM|9^MDCX\_PDIALY\_FLUID\_BAG\_VOLUME^MDC|1.1.3.3|5.000|L^liter^UCUM|||||F

1190 OBX|18|ST|52^MDCX\_DEV\_PDIALY\_FLUID\_2\_CHAN^MDC|1.1.4|||||F

1191 OBX|19|ST|29^MDCX\_PDIALY\_FLUID\_NAME^MDC|1.1.4.1|Dialneal Low Calcium 2.5%|||||F

1192 OBX|20|ST|10^MDCX\_PDIALY\_FLUID\_SOURCE^MDC|1.1.4.2|BAG|||||F

1193 OBX|21|NM|9^MDCX\_PDIALY\_FLUID\_BAG\_VOLUME^MDC|1.1.4.3|5.000|L^liter^UCUM|||||F

1194 OBX|22|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.5|||||F

1195 OBX|23|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.5.1|1|||||F

1196 OBX|24|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.5.2|COMPLETE|||||F

1197 OBX|25|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.5.3|13.3|min^minutes^UCUM|||||F

1198 OBX|26|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.5.4|2.667^0.000|L^liter^UCUM|||||F

1199 OBX|27|DTM|37^MDCX\_PDIALY\_FILL\_START\_TIME^MDC|1.1.5.5|20230913220000+0000|||||F

1200 OBX|28|DTM|38^MDCX\_PDIALY\_FILL\_END\_TIME^MDC|1.1.5.6|20230913221320+0000|||||F

1201 OBX|29|NM|11^MDCX\_PDIALY\_FILL\_DURATION^MDC|1.1.5.7|13.3|min^minutes^UCUM|||||F

1202 OBX|30|NA|13^MDCX\_PDIALY\_FILL\_VOLUME^MDC|1.1.5.8|2.667^0.000|L^liter^UCUM|||||F

1203 OBX|31|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.5.9|153.3|min^minutes^UCUM|||||F

1204 OBX|32|DTM|39^MDCX\_PDIALY\_DWELL\_START\_TIME^MDC|1.1.5.10|20230913221320+0000|||||F

1205 OBX|33|DTM|40^MDCX\_PDIALY\_DWELL\_END\_TIME^MDC|1.1.5.11|20230914004640+0000|||||F

1206 OBX|34|NM|15^MDCX\_PDIALY\_DWELL\_DURATION^MDC|1.1.5.12|153.3|min^minutes^UCUM|||||F

1207 OBX|35|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.5.13|13.3|min^minutes^UCUM|||||F

1208 OBX|36|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.5.14|2.667|L^liter^UCUM|||||F

1209 OBX|37|DTM|41^MDCX\_PDIALY\_DRAIN\_START\_TIME^MDC|1.1.5.15|20230914004640+0000|||||F

1210 OBX|38|DTM|42^MDCX\_PDIALY\_DRAIN\_END\_TIME^MDC|1.1.5.16|20230914010000+0000|||||F

1211 OBX|39|NM|17^MDCX\_PDIALY\_DRAIN\_DURATION^MDC|1.1.5.17|13.3|min^minutes^UCUM|||||F

1212 OBX|40|NM|19^MDCX\_PDIALY\_DRAIN\_VOLUME^MDC|1.1.5.18|2.667|L^liter^UCUM|||||F

1213 OBX|41|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.6|||||F

1214 OBX|42|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.6.1|2|||||F

1215 OBX|43|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.6.2|DWELL|||||F

1216 OBX|44|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.6.3|13.3|min^minutes^UCUM|||||F

1217 OBX|45|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.6.4|2.333^0.333|L^liter^UCUM|||||F

1218 OBX|46|DTM|37^MDCX\_PDIALY\_FILL\_START\_TIME^MDC|1.1.6.5|20230914010000+0000|||||F

1219 OBX|47|DTM|38^MDCX\_PDIALY\_FILL\_END\_TIME^MDC|1.1.6.6|20230914011320+0000|||||F

1220 OBX|48|NM|11^MDCX\_PDIALY\_FILL\_DURATION^MDC|1.1.6.7|13.3|min^minutes^UCUM|||||F

1221 OBX|49|NA|13^MDCX\_PDIALY\_FILL\_VOLUME^MDC|1.1.6.8|2.333^0.333|L^liter^UCUM|||||F

1222 OBX|50|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.6.9|153.3|min^minutes^UCUM|||||F

1223 OBX|51|DTM|39^MDCX\_PDIALY\_DWELL\_START\_TIME^MDC|1.1.6.10|20230914011320+0000|||||F

1224 OBX|52|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.6.11|13.3|min^minutes^UCUM|||||F

1225 OBX|53|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.6.12|2.667|L^liter^UCUM|||||F

1226 OBX|54|ST|3^MDCX\_DEV\_PDIALY\_EXCHANGE\_CHAN^MDC|1.1.7|||||F

1227 OBX|55|ST|100^MDC\_ATTR\_CHAN\_NUM\_LOGICAL^MDC|1.1.7.1|3|||||F

1228 OBX|56|ST|8^MDCX\_PDIALY\_CURRENT\_PHASE^MDC|1.1.7.2|PENDING|||||F

1229 OBX|57|NM|12^MDCX\_PDIALY\_FILL\_DURATION\_SETTING^MDC|1.1.7.3|13.3|min^minutes^UCUM|||||F

1230 OBX|58|NA|14^MDCX\_PDIALY\_FILL\_VOLUME\_SETTING^MDC|1.1.7.4|0.000^2.667|L^liter^UCUM|||||F

1231 OBX|59|NM|16^MDCX\_PDIALY\_DWELL\_DURATION\_SETTING^MDC|1.1.7.5|153.3|min^minutes^UCUM|||||F

OBX|60|NM|18^MDCX\_PDIALY\_DRAIN\_DURATION\_SETTING^MDC|1.1.7.6|13.3|min^minutes^UCUM|||||F

OBX|61|NM|20^MDCX\_PDIALY\_DRAIN\_VOLUME\_SETTING^MDC|1.1.7.7|2.667|L^liter^UCUM|||||F

## 6.3 EMR Response

The EMR responds with a simple Acknowledgement Message (ACK\_R01\_ACK). Although the response to a missing or rejected response is machine specific, this standard recommends the following actions:

- 1) If the EMR does not respond, then the dialysis machine should retry the message. If no response is received for the retry, the dialysis machine should generate an alert locally.
- 2) If the EMR responds but rejects the PCD-01 message, the dialysis machine should generate an alert locally.

### 6.3.1 Message Structure

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
{{SFT}}	Software Segment	X	[0..0]
[UAC]	User Authentication Credential	X	[0..0]
MSA	Message Acknowledgement	R	[1..1]
ERR	Error	RE	[0..100]

### 6.3.2 Example 1 – Accepted Response

The following message is an example of an accepted PCD-01 message.

```
MSH|^~\&|EMR|||20191003092025+0000||ACK^R01^ACK|XX3657|P|2.6|||NE|NE
MSA|AA|20191003092024
```

### 6.3.3 Example 2 – Rejected Response

The following message is an example of a rejected PCD-01 message.


```
MSH|^~\&|EMR|||20191003092025+0000||ACK^R01^ACK|XX3657|P|2.6|||NE|NE
MSA|AR|20191003092024
ERR||PID^1^11^^9|103^Table Value Not Found|E
```

## 7 Reporting Alarm Information

### 7.1 Overview

The sending of Alarm Information to the EMR System is accomplished using the PCD-04 Transaction of the Patient Care Device (PCD) domain in the Device Enterprise Communication (DEC) IHE profile. In doing so, the Dialysis Machine acts as the Alert Reporter (AR) and the EMR is the Alert Manager (AM) and/or the Alert Consumer (ACON).

1262

**Figure 9 - AR sends Report Alert to AM and/or ACON as an HL7 ORU message.**

1263

1264

1265 The Dialysis Machine will report the alarm when the alarm state changes, when the alarm activity state changes,  
 1266 and periodically while it is active. The period of the "keep-alive" PCD-04 messages should be disclosed in the  
 1267 manufacturer's documents. It is recommended that this be between 10 and 30 seconds. The dialysis machine  
 1268 will send the appropriate parameters and measurements with each alarm. For example, when an Arterial  
 1269 Pressure High Alarm is reported, the message will contain the Arterial pressure measurement and upper alarm  
 1270 limit.

1271 The ORU^R40 message represents a unitary alert, which is to be acknowledged, in whole, by an ORA message.  
 1272 Multiple alerts requiring separate acknowledgement will be sent as individual messages.

## 1273 7.2 Device Request

1274 The PCD-04 Report Alert message is used to communicate ACM data from an Alert Reporter (AR) to Alert  
 1275 Manager (AM) and/or Alert Consumer (ACON). In addition, the Dialysis Machine sends a PCD-04 Message to  
 1276 the Device Observation Consumer (DOC) so that the information can be logged with the PCD-01 Message.

1277 Common HL7 segments are defined in HL7 Data Elements. There are sections discussing considerations  
 1278 specific to PCD-04 where applicable.

1279 While there can be multiple OBR segments per PCD-04 transaction (in support of inclusion of alert common  
 1280 containment and evidentiary data) there is at most one alert per PCD-04 transaction.

1281  
1282

### 7.2.1 Message Structure

**Message Type - ORU^R40^ORU\_R40**

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
[[SFT]]	Software Segment	X	[0..0]
{	--- ALERT_begin		
[	--- PATIENT begin		
PID	Patient Identification	R	[1..1]
[	--- LOCATION begin		
PV1	Alert Location	O	[0..1]
]	--- LOCATION end		
]	--- PATIENT end		
{	--- ALERT_IDENTIFICATION begin		
[ORC]	Alert Order Common	X	[0..0]
{OBR}	Alert Identification	R	[1..*]
[{	--- ALERT_OBSERVATION begin		
{OBX}	Alert observations relative to OBR	R	[1..*]
[[NTE]]	Notes and Comments	X	[0..0]
}}	--- ALERT_OBSERVATION end		
}	--- ALERT_IDENTIFICATION end		
}	--- ALERT end		

1283  
1284  
1285  
1286  
1287  
1288  
1289  
1290  
1291  
1292

A single Report Alert [PCD-04] transaction contains at most one alert for a given patient. The PCD-04 message contains a single OBR segment that contains five OBX segments in a specific order.

- OBX 1 identifies the type of alarm.
- OBX 2 identifies the event source. For numeric values this includes the object value and limits. This is the Source Column of Table 3 and 5.
- OBX 3 identifies the event phase.
- OBX 4 identifies the alarm state.
- OBX 5 identifies the alarm activity state.

The first OBX (ID=1) identifies the type of alert.

**HL7 Attribute Table - OBX Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

1293

**OBX-1 Set ID**

1294

This field contains the sequence number. For this segment it will always be 1.

1295

**OBX-2 Value Type**

1296

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "CWE".

1297

1298

**OBX-3 Observation Identifier**

1299

This field will be one of the following.

Field Value	Usage
196670^MDC_EVT_LO^MDC	This value is used when the alarm is generated because a measure value is below a numeric limit.
196648^MDC_EVT_HI^MDC	This value is used when the alarm is generated because a measure value is above a numeric limit.
196616^MDC_EVT_ALARM^MDC	This value is used when the alarm is not generated due to a numeric value.

1300

**OBX-4 Observation Sub-ID**

1301

This value will always be "1.0.0.0.1".

1302

**OBX-5 Observation Value**

1303

This field identifies the alarm being generated. For numeric values, this field is the value in the Source Column in Table 3. For non-numeric alarms, this field is the Event Identifier Column in Table 3.

1304

1305

1306

For example, "198240^MDC\_EVT\_HDIALY\_ANTICOAG\_SYRINGE\_SIZE^MDC"

1307

"198242^MDC\_EVT\_HDIALY\_BLD\_PUMP\_STOP^MDC"



1308

**OBX-8 Interpretation Codes**

1309

This field is used to provide the Alarm Priority and Alarm Type, and optionally the Abnormality Type.

1310

Supported values for these items can be found in the definition of [OBX-8](#).

1311

Examples,

1312

High Priority, Technical, Upper Limit Alarm H~ST~PH

1313

Non-Numeric High Priority, Technical Alarm ST~PH

1314

1315

**OBX-11 Observation Result Status**

1316

For this segment the value will always be "F".

1317

1318

The second OBX (ID=2) identifies the parameter and the reference range.

**HL7 Attribute Table - OBX Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

1319

**OBX-1 Set ID**

1320

For this segment the value will always be 2.

1321

**OBX-2 Value Type**

1322

This field defines the data type of OBX-5, Observation Value. For numeric alarms, this value will be "NM".

1323

For non-numeric alarms, this value will be "CWE".

1324

1325

**OBX-3 Observation Identifier**

1326  
1327

For numeric alarms, this field identifies the object being monitored. This object is identified in the Source Column in Table 3.

1328

For non-numeric alarms, this field will always be "68480^MDC\_ATTR\_ALERT\_SOURCE^MDC".

1329

**OBX-4 Observation Sub-ID**

1330  
1331

For numeric alarms, this field identifies the object being measured. This value is the Observation Sub ID of the object identified in the Source Column in Table 3.

1332  
1333

For non-numeric alarms, this field identifies the VMD of the device generating the alarm. This value is the Observation Sub-ID of the object identified in the Source Column in Table 3.

1334

In both cases, the observation sub-ID has ".2" appended to signify facet 2.

1335

For example, the Therapy Outcomes VMD would be reported as "1.1.8.2".

1336

**OBX-5 Observation Value**

1337

For numeric alarms, this field contains the value of the object that resulted in the alarm.

1338  
1339

For non-numeric alarms, this field contains the identifier of the VMD generating the alarm. This is identified in the Source Column in Table 3.

1340

Example, "70967^MDC\_DEV\_HDIALY\_THERAPY\_OUTCOMES\_CHAN^MDC".

1341

**OBX-6 Observation Units**

1342

For numeric alarms, this field contains the units of measure for the object.

1343

For non-numeric alarms, this field is blank.

1344

**OBX-7 Reference Range**

1345  
1346

This field contains the alarm limits for numeric alarms. For non-numeric alarms this field is empty. Ranges can be specified in one of three ways.

	Format	Example
Upper Limit Only	> Upper	> 400
Lower Limit Only	< Lower	< -200
Upper and Lower Limit	Lower – Upper	20 - 400

1347

Note, alarm limits may be fixed or machine generated.

1348

**OBX-11 Observation Result Status**

1349

For this segment the value will always be "F".

1350

1351

The third OBX (ID=3) identifies the event phase.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

**OBX-1 Set ID**

This field contains the sequence number. For this segment it will always be 3.

**OBX-2 Value Type**

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "ST".

**OBX-3 Observation Identifier**

This field will always be "68481^MDC\_ATTR\_EVENT\_PHASE^MDC".

**OBX-4 Observation Sub-ID**

This value will always be "1.0.0.0.3".

**OBX-5 Observation Value**

This value will always be one of the following.

- start - start of an interval event/alert, an end is expected.
- continue - continuation of an ongoing interval event/alert.
- End - end of an interval event/alert.

**OBX-11 Observation Result Status**

For this segment the value will always be "F".

The third OBX (ID=4) identifies the alarm state.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

**OBX-1 Set ID**

This field contains the sequence number. For this segment it will always be 4.

**OBX-2 Value Type**

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "ST".

**OBX-3 Observation Identifier**

This field will always be "68482^MDC\_ATTR\_ALARM\_STATE^MDC".

**OBX-4 Observation Sub-ID**

This value will always be "1.0.0.0.4".

**OBX-5 Observation Value**

This value will one of the following:

- off
- inactive
- active
- latched

**OBX-11 Observation Result Status**

For this segment the value will always be "F".

The third OBX (ID=5) identifies the alarm activity state.

HL7 Attribute Table - OBX Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R		Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R		Units
7	ST	X		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	X		Date/Time of the Observation
15	CWE	X		Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X		Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

1387

**OBX-1 Set ID**

1388

This field contains the sequence number. For this segment it will always be 5.

1389

**OBX-2 Value Type**

1390

This field defines the data type of OBX-5, Observation Value. For this segment the value will always be "ST".

1391

**OBX-3 Observation Identifier**

1392

This field will always be "68483^MDC\_ATTR\_ALARM\_INACTIVATION\_STATE^MDC".

1393

**OBX-4 Observation Sub-ID**

1394

This value will always be "1.0.0.0.5".

1395

**OBX-5 Observation Value**

1396

This value will one of the following:

1397

- enabled

1398

- audio-paused

1399

- audio-off

1400

- alarm-paused

1401

- alarm-off

1402

- alert-acknowledged (IEC 60601-1-8 Amendment 1; CP126)

**OBX-11 Observation Result Status**

For this segment the value will always be "F".

**7.2.2 Example 1 - Alarm Initiation**

The following message is an example of a Venous Pressure Low Alarm.

```

MSH|^~\&|ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1408 64|||20191003092025+0000||ORU^R40^ORU_R40|20191003092024|P|2.6||AL|NE||||IHE_PCD_001^IH
1409 EPCD^1.3.6.1.4.1.19376.1.6.1.4.1^ISO
1410
PID|||Scrubber 2000/SC678932^^^"U||^^^^^U
1411
OBR|1||080019FFFE3ED02D20110602045842^ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1412 64|196616^MDC_EVT_ALARM^MDC||20191003092024+0000
1413
OBX|1|ST|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|1.0.0|||||F
1414
OBX|2|ST|70934^MDC_DEV_HDIALY_VMD^MDC|1.1|||||F
1415
OBX|3|CWE|196670^MDC_EVT_LO^MDC|1.0.0.0.1|158776^MDC_HDIALY_BLD_PUMP_PRESS_VEN^MDC||PH~SP~L|
1416 ||F
1417
OBX|4|NM|158776^MDC_HDIALY_BLD_PUMP_PRESS_VEN^MDC|1.1.3.15.2|15|mm[Hg]^Millimeters of
1418 Mercury^UCUM |20-400|||F
1419
OBX|5|ST|68481^MDC_ATTR_EVENT_PHASE^MDC|1.0.0.0.3|start|||||F
1420
OBX|6|ST|68482^MDC_ATTR_ALARM_STATE^MDC|1.0.0.0.4|active|||||F
1421
OBX|7|ST|68483^MDC_ATTR_ALARM_INACTIVATION_STATE^MDC|1.0.0.0.5|enabled|||||F
1422

```

**7.2.3 Example 2 - Alarm Cessation**

The following message shows the cessation of a Blood Pump Stopped alarm.

```

MSH|^~\&|ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1426 64|||20191003092025+0000||ORU^R40^ORU_R40|20191003092024|P|2.6||AL|NE||||IHE_PCD_001^IH
1427 EPCD^1.3.6.1.4.1.19376.1.6.1.4.1^ISO
1428
PID|||Scrubber 2000/SC678932^^^"U||^^^^^U
1429
OBR|1||080019FFFE3ED02D20110602045842^ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1430 64|196616^MDC_EVT_ALARM^MDC||20191003092024+0000
1431
OBX|1|ST|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|1.0.0|||||F
1432
OBX|2|ST|70934^MDC_DEV_HDIALY_VMD^MDC|1.1|||||F
1433
OBX|3|CWE|196616^MDC_EVT_ALARM^MDC|1.0.0.0.1|198242^MDC_EVT_HDIALY_ALARM_BLD_PUMP_STOP^MDC|||
1434 ||F
1435
OBX|4|CWE|68480^MDC_ATTR_ALERT_SOURCE^MDC|1.1.3.0.2|||||F
1436
OBX|5|ST|68481^MDC_ATTR_EVENT_PHASE^MDC|1.0.0.0.3|end|||||F
1437
OBX|6|ST|68482^MDC_ATTR_ALARM_STATE^MDC|1.0.0.0.4|inactive|||||F
1438
OBX|7|ST|68483^MDC_ATTR_ALARM_INACTIVATION_STATE^MDC|1.0.0.0.5|enabled|||||F
1439

```

**7.2.4 Example 3 - Alarm Mute**

The following message shows a VP Low Alarm that has been muted by the user.

```

MSH|^~\&|ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1443 64|||20191003092025+0000||ORU^R40^ORU_R40|20191003092024|P|2.6||AL|NE||||IHE_PCD_001^IH
1444 EPCD^1.3.6.1.4.1.19376.1.6.1.4.1^ISO
1445
PID|||Scrubber 2000/SC678932^^^"U||^^^^^U
1446
OBR|1||080019FFFE3ED02D20110602045842^ACME_Dialysis_Machine^080019FFFE3ED02D^EUI-
1447 64|196616^MDC_EVT_ALARM^MDC||20191003092024+0000
1448
OBX|1|ST|70929^MDC_DEV_HDIALY_MACHINE_MDS^MDC|1.0.0|||||F
1449

```

1450 OBX|2|ST|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.1|||||F  
 1451 OBX|3|CWE|196670^MDC\_EVT\_LO^MDC|1.0.0.0.1|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|||PH~SP|||  
 1452 F  
 1453 OBX|4|NM|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|1.1.3.15.2|15|mm[Hg]^Millimeters of  
 1454 Mercury^UCUM |20-400|||F  
 1455 OBX|5|ST|68481^MDC\_ATTR\_EVENT\_PHASE^MDC|1.0.0.0.3|continue|||||F  
 1456 OBX|6|ST|68482^MDC\_ATTR\_ALARM\_STATE^MDC|1.0.0.0.4|active|||||F  
 1457 OBX|7|ST|68483^MDC\_ATTR\_ALARM\_INACTIVATION\_STATE^MDC|1.0.0.0.5|audio-paused|||||F

1458

1459 The following message shows a VP Low Alarm whose mute has expired.

1460 MSH|^~\&|ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-  
 1461 64|||20191003092025+0000||ORU^R40^ORU\_R40|20191003092024|P|2.6||AL|NE||||IHE\_PCD\_001^IH  
 1462 EPCD^1.3.6.1.4.1.19376.1.6.1.4.1^ISO  
 1463 PID|||Scrubber 2000/SC678932^^^U||^U  
 1464 OBR|1||080019FFFE3ED02D20110602045842^ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-  
 1465 64|196616^MDC\_EVT\_ALARM^MDC|||20191003092024+0000  
 1466 OBX|1|ST|70929^MDC\_DEV\_HDIALY\_MACHINE\_MDS^MDC|1.0.0|||||F  
 1467 OBX|2|ST|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.1|||||F  
 1468 OBX|3|CWE|196670^MDC\_EVT\_LO^MDC|1.0.0.0.1|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|||PH~SP|||  
 1469 F  
 1470 OBX|4|NM|158776^MDC\_HDIALY\_BLD\_PUMP\_PRESS\_VEN^MDC|1.1.3.15.2|15|mm[Hg]^Millimeters of  
 1471 Mercury^UCUM |20-400|||F  
 1472 OBX|5|ST|68481^MDC\_ATTR\_EVENT\_PHASE^MDC|1.0.0.0.3|continue|||||F  
 1473 OBX|6|ST|68482^MDC\_ATTR\_ALARM\_STATE^MDC|1.0.0.0.4|active|||||F  
 1474 OBX|7|ST|68483^MDC\_ATTR\_ALARM\_INACTIVATION\_STATE^MDC|1.0.0.0.5|enabled|||||F

1475

### 7.2.5 Example 4 – Vendor Specific Alarm

1476 The following shows how to report a vendor specific alarm. In this case a Venous Air Alarm. The vendor  
 1477 specific alarm information is reported following the mandatory PCD-04 alarm information. In this case, the  
 1478 Alarm ID is “10” and the Alarm Text is “Venous Air”.

1479 MSH|^~\&|ACME\_Dialysis\_Machine^025041FFFE000001^EUI-  
 1480 64|||20241121154324+0000||ORU^R40^ORU\_R40|20241121154324827|P|2.6||AL|NE||||IHE\_PCD\_ACM\_00  
 1481 1^IHEPCD^1.3.6.1.4.1.19376.1.6.4.4^ISO  
 1482 PID|||Scrubber 2000/SC678932^^^U||^U  
 1483 OBR|1||025041FFFE00000120241121154324^ACME\_Dialysis\_Machine^025041FFFE000001^EUI-  
 1484 64|196616^MDC\_EVT\_ALARM^MDC|||20241121154324+0000  
 1485 OBX|1|CWE|61439^MDC\_EVT\_NOS^MDC|1.0.0.1.1|||||F  
 1486 OBX|2|CWE|70934^MDC\_DEV\_HDIALY\_VMD^MDC|1.0.0.1.2|||||F  
 1487 OBX|3|ST|68481^MDC\_ATTR\_EVENT\_PHASE^MDC|1.0.0.1.3|start|||||F  
 1488 OBX|4|ST|68482^MDC\_ATTR\_ALARM\_STATE^MDC|1.0.0.1.4|active|||||F  
 1489 OBX|5|ST|68483^MDC\_ATTR\_ALARM\_INACTIVATION\_STATE^MDC|1.0.0.1.5|enabled|||||F  
 1490 OBX|6|ST|999999^MDC\_ATTR\_VENDOR\_ALERT\_ID^MDC|1.0.0.2|10|||||F  
 1491 OBX|7|ST|999999^MDC\_ATTR\_VENDOR\_ALERT\_TEXT^MDC|1.0.0.3|Venous Air|||||F

1492

## 7.3 EMR Response

1493 The EMR responds with a Report Alert Acknowledgement (ORA^R41^ORA\_R41).  
 1494  
 1495

### 7.3.1 Message Structure

Message Type - ORA^R41^ORA\_R41

Segment	Meaning	Use	Card
MSH	Message Header	R	[1..1]
{{SFT}}	Software Segment	X	[0..0]
[UAC]	User Authentication Credential	X	[0..0]
MSA	Message Acknowledgement	R	[1..1]
ERR	Error	RE	[0..100]
{{PRT}}	Participation (Acknowledging User)	O	[0..1]

### 7.3.2 Example 1 – Accepted Response

The following message is an example of an accepted PCD-04 message.

```
MSH|^~\&|EMR|||20191003092025+0000||ACK^R40^ACK|XX3657|P|2.6|||NE|NE
MSA|AA|20191003092024
```

## 8 HL7 Data Elements

The information contained in this section is a subset of the full HL7 Specification. Only segments used by this standard are shown. Restrictions that are specific to this specification are identified for each segment.

### 8.1 Tables

The values contained in the tables in this section may be a subset of the values supported by the HL7 specification. Only values supported or used by Dialysis Machines are shown.

HL7 Table 0003 – Event

Value	Description
R01	PCD-01, PCD-01 Ack
R40	PCD-04, PCD-04 Ack
D01	Prescription Query Request
K22	Prescription Query Response, Patient Query Response
Q22	Patient Query Request
	Other values not used in this specification.

HL7 Table 0004 – Patient Class

Value	Description
B	Obstetrics (Not Used)
C	Commercial (Not User)
E	Emergency
I	Inpatient (Used for acute treatments)
N	Not Applicable
O	Outpatient (Used for chronic incenter treatments)
P	Preadmit (Not Used)
R	Recurring Patient (Used for chronic home treatments)
U	Unknown



1515

**HL7 Table 0008 - Acknowledgment Code**

Value	Description
AA	Original mode: Application Accept - Enhanced mode: Application acknowledgment: Accept
AE	Original mode: Application Error - Enhanced mode: Application acknowledgment: Error
AR	Original mode: Application Reject - Enhanced mode: Application acknowledgment: Reject
CA	Enhanced mode: Accept acknowledgment: Commit Accept
CE	Enhanced mode: Accept acknowledgment: Commit Error
CR	Enhanced mode: Accept acknowledgment: Commit Reject

1516

1517

**HL7 Table 0038 - Order status**

Value	Description
A	Some, but not all, results available
CA	Order was canceled
CM	Order is completed
DC	Order was discontinued
ER	Error, order not found
HD	Order is on hold
IP	In process, unspecified
RP	Order has been replaced
SC	In process, scheduled

1518

1519

**HL7 Table 0076 – Message Type**

Value	Description
ORU	PCD-01, PCD-04
ACK	PCD-01 Ack
QBP	Prescription Query Request, Patient Query Request
RSP	Prescription Query Response, Patient Query Response
	Other values not used in this specification.

1520

1521

**HL7 Table 0085 - Observation result status codes interpretation**

Value	Description
F	Final results; Can only be changed with a corrected result.

1522

1523

**HL7 Table 0091 – Query Priority**

Value	Description
D	Deferred
I	Immediate

1524

1525

**HL7 Table 0123 - Result status for OBR segment**

Value	Description
O	Order received; specimen not yet received
I	No results available; specimen received, procedure incomplete
S	No results available; procedure scheduled, but not done
A	Some, but not all, results available
P	Preliminary: A verified early result is available, final results not yet obtained
C	Correction to results

Value	Description
R	Results stored; not yet verified
F	Final results; results stored and verified. Can only be changed with a corrected result.
X	No results available; Order canceled.
Y	No order on record for this test. (Used only on queries)
Z	No record of this patient. (Used only on queries)

1526

1527

**HL7 Table 0125 - Value type**

Value	Description
DT	Date
NM	Numeric
ST	String Data.
TM	Time

1528

1529

**HL7 Table 0155 - Accept/application acknowledgment conditions**

Value	Description
AL	Always

1530

1531

**HL7 Table 0208 – Query Response Status**

Value	Description
OK	Data found, no errors. This is the default value.
NF	No data found, no errors.
AE	Application Error
AR	Application Reject

1532

1533

**HL7 Table 0354 – Message Structure**

Value	Description
ORU_R01	PCD-01
ORU_R40	PCD-04
ACK	PCD-01 Ack, PCD-03 Ack
QBP_DO1	Prescription Request
RSP_K22	Prescription Response
QBP_Q21	Patient Query Request
RSP_K21	Patient Query Response
	Other values not used in the standard

1534

1535

**HL7 table 0357 – Message Error Condition Codes**

Value	Description	Comment
0	Message accepted	Success. Optional, as the AA conveys success. Used for systems that shall always return a status code.
100	Segment sequence error	Error: The message segments were not in the proper order, or required segments are missing.
101	Required field missing	Error: A required field is missing from a segment
102	Data type error	Error: The field contained data of the wrong data type, e.g. an NM field contained "FOO".
103	Table value not found	Error: A field of data type ID or IS was compared against the corresponding table, and no match was found.
200	Unsupported message type	Rejection: The Message Type is not supported.

Value	Description	Comment
201	Unsupported event code	Rejection: The Event Code is not supported.
202	Unsupported processing id	Rejection: The Processing ID is not supported.
203	Unsupported version id	Rejection: The Version ID is not supported.
204	Unknown key identifier	Rejection: The ID of the patient, order, etc., was not found. Used for transactions <i>other than</i> additions, e.g. transfer of a non-existent patient.
205	Duplicate key identifier	Rejection: The ID of the patient, order, etc., already exists. Used in response to addition
206	Application record locked	Rejection: The transaction could not be performed at the application storage level, e.g., database locked.
207	Application internal error	Rejection: A catchall for internal errors not explicitly covered by other codes.

1536

1537

**HL7 Table 0394 – Response Modality**

Value	Description
R	Real Time
T	Bolus (a series of responses sent at the same time without use of batch formatting)
B	Batch

1538

1539

**HL7 Table 0396 – Coding Systems**

Value	Description
MDC	ISO/IEEE
UCUM	Unified Code for Units of Measure (UCUM)

1540

1541

**HL7 table 0516 – Error Severity**

Value	Description	Comment
W	Warning	Transaction successful, but there may issues
I	Information	Transaction was successful but includes information e.g., inform patient
E	Error	Transaction was unsuccessful

1542

1543

## 8.2 Segments

1544

### 8.2.1 MSA – Message Acknowledgement Segment

1545

1546

1547

The MSA segment contains INFORMATION sent while acknowledging another message. This message is sent from the EMR System to the Dialysis Machine. Although other fields may be present, the Dialysis Machine processes MSA-1 Acknowledgement Code and MSA-2 Message Control ID.

1548

**HL7 Attribute Table - MSA Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ID	R	0008	Acknowledgment Code
2	ST	R		Message Control ID
3		W		Text Message
4	NM	O		Expected Sequence Number
5		W		Delayed Acknowledgment Type
6		W		Error Condition
7	NM	O		Message Waiting Number

HL7 Attribute Table - MSA Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
8	ID	O	0520	Message Waiting Priority

1549

1550

**MSA-1 Acknowledgment Code (ID)**

1551

HL7 Definition: This field contains an acknowledgment code, see message processing rules. Refer to HL7 Table 0008 - Acknowledgment Code for valid values.

1552

1553

In this table, the value "CA" corresponds to "Enhanced mode: Accept acknowledgment: Commit Accept."

1554

This value means that the receiving system has committed the message to safe storage in a manner that

1555

releases the sending system from the need to resend the message. "CR" means that a message was

1556

rejected for a message uniquely identified by the information in MSH-9, MSH-12, MSH-11, MSH-21. "CE" is

1557

a message processing error or failure condition.

1558

**This Specification:** The Dialysis Machine accepts both original and enhanced mode responses. The response of the Dialysis Machine to an Error or Rejected message is device and manufacturer specific and is not covered by this guide.

1559

1560

1561

**MSA-2 Message Control ID (ST)**

1562

HL7 Definition: This field contains the message control ID of the message sent by the sending system. It allows the sending system to associate this response with the message for which it is intended.

1563

1564

**This Specification:** Note that the combination of the values in *MSA-2 Message Control ID* and *MSH-6 Receiving Facility* should be unique to the recipient of the acknowledgement message (i.e., the Order Fulfiller).

1565

1566

1567

**Note on Element Length:** The element length for MSA-2 has been extended to 50 characters from the HL7-prescribed length of 20 characters. This element has been extended to accommodate the extended element length of MSH-10 Message Control ID in the ORU message types.

1568

1569

1570

1571

**8.2.2 MSH – Message Header Segment**

1572

The MSH segment defines the intent, source, destination, and some specifics of the syntax of a message.

1573

HL7 Attribute Table - MSH Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	R		Field Separator
2	ST	R		Encoding Characters
3	HD	R	0361	Sending Application
4	HD	X	0362	Sending Facility
5	HD	X	0361	Receiving Application
6	HD	X	0362	Receiving Facility
7	DTM	R		Date/Time of Message
8	ST	X		Security
9	MSG	R		Message Type
10	ST	R		Message Control ID
11	PT	R		Processing ID
12	VID	R		Version ID
13	NM	X		Sequence Number
14	ST	X		Continuation Pointer
15	ID	X	0155	Accept Acknowledgment Type

HL7 Attribute Table - MSH Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
16	ID	X	0155	Application Acknowledgment Type
17	ID	X	0399	Country Code
18	ID	X	0211	Character Set
19	CWE	X		Principal Language Of Message
20	ID	X	0356	Alternate Character Set Handling Scheme
21	EI	R		Message Profile Identifier
22	XON	X		Sending Responsible Organization
23	XON	X		Receiving Responsible Organization
24	HD	X		Sending Network Address
25	HD	X		Receiving Network Address

1574

1575

**MSH-1 Field separator (ST)**

1576

**HL7 Definition:** This field contains the separator between the segment ID and the first real field. As such it serves as the separator and defines the character to be used as a separator for the rest of the message. Recommended value is | (ASCII 124).

1577

1578

1579

**This Specification:** This value will be | (ASCII 124).

1580

**MSH-2 Encoding characters (ST)**

1581

**HL7 Definition:** This field contains four characters in the following order: the component separator, repetition separator, escape character, and subcomponent separator. Recommended values are ^~\& (ASCII 94, 126, 92, and 38, respectively).

1582

1583

1584

**This Specification:** This value will be recommended values are ^~\& (ASCII 94, 126, 92, and 38, respectively).

1585

1586

**MSH-3 Sending Application (HD)**

1587

**HL7 Definition:** This field uniquely identifies the sending application among all other applications within the network enterprise. The network enterprise consists of all those applications that participate in the exchange of HL7 messages within the enterprise.

1588

1589

1590

**This Specification:** This value will be the name of the dialysis machine and the Extended Unique Identified (EUI-64). For example,

1591

1592

ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-64

1593

The EUI-64 is created from the MAC Address of the interface sending the messages. If the MAC address of the interface is 12-34-56-78-AB-CD then the EUI-64 is 123456FFFE78ABCD.

1594

1595

**MSH-7 Date/time of message (TS)**

1596

**HL7 Definition:** This field contains the date/time that the sending system created the message. If the time zone is specified, it will be used throughout the message as the default time zone.

1597

1598

**MSH-9 Message type (MSG)**

1599

**HL7 Definition:** This element contains the message type, trigger event, and the message structure ID for the message.

1600

1601

CM Components: (ID)> <message type (ID)> ^ <trigger event (ID)> ^ <message structure

1602

1603

The allowed components of this element are listed in several tables maintained by HL7. See HL7 Table 0076 – Message Type, HL7 Table 0003 – Event, and HL7 Table 0354 – Message Structure.

1604

1605 **This Specification:** In the MT-ACK-1 message type, this element shall be hard coded to ACK^R01^ACK

1606 **MSH-10 Message control ID (ST)**

1607 **HL7 Definition:** This element contains a number or other identifier that uniquely identifies the message.

1608 **This Specification:** The sending system shall assign an identifier for the message that is unique within the

1609 namespace of the sending facility and/or application.

1610 Note that the value of this element is *not* a reference to the message that is being acknowledged. The

1611 Message control ID of the acknowledged message appears in MSA-2 Message Control ID.

1612 **Note on Element Length:** The element length for MSH-10 has been extended to 50 characters from the

1613 HL7-prescribed length of 20 characters. This extension allows sending systems to use globally unique

1614 identifiers (such as GUIDs) for Message IDs, an increasingly common practice.

1615 **MSH-11 Processing ID (PT)**

1616 **HL7 Definition:** This field is used to decide whether to process the message as defined in HL7 Application

1617 (level 7) Processing rules.

1618 **This Specification:** This value will be "P".

1619 **MSH-12 Version ID (VID)**

1620 **HL7 Definition:** This field identifies the version.

1621 **This Specification:** This value will be "2.6".

1622 **MSH-21 Message Profile Identifier (EI)**

1623 **HL7 Definition:** Sites may use this element to assert adherence to a Conformance Statement published by

1624 HL7 or by a site. Conformance Statements contain detailed explanations of grammar, syntax, and usage for

1625 a particular message or set of messages. This element is treated like all EI data type information, if the

1626 underscores or dashes are used as encoding characters in an interface, those characters shall be escaped to

1627 be sent in the above element.

1628 **This Specification:** In the MT-ACK-1 message type, MSH-21 Conformance statement ID should be:

1629 DIALY\_MT-ACK-1\_R1

1630

1631 8.2.3 ORC – Order Common

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ID	R	0119	Order Control
2	EI	R		Placer Order Number
3	EI	X		Filler Order Number
4	EI	X		Placer Group Number
5	ID	X	0038	Order Status
6	ID	O	0121	Response Flag
7		X		Quantity/Timing
8	EIP	X		Parent
9	DTM	X		Date/Time of Transaction
10	XCN	X		Entered By
11	XCN	X		Verified By
12	XCN	O		Ordering Provider
13	PL	X		Enterer's Location
14	XTN	O		Call Back Phone Number
15	DTM	X		Order Effective Date/Time
16	CWE	X	9999	Order Control Code Reason

SEQ	DT	OPT	TBL#	ELEMENT NAME
17	CWE	X	9999	Entering Organization
18	CWE	X	9999	Entering Device
19	XCN	X		Action By
20	CWE	X	0339	Advanced Beneficiary Notice Code
21	XON	X		Ordering Facility Name
22	XAD	X		Ordering Facility Address
23	XTN	X		Ordering Facility Phone Number
24	XAD	X		Ordering Provider Address
25	CWE	X	9999	Order Status Modifier
26	CWE	X	0552	Advanced Beneficiary Notice Override Reason
27	DTM	X		Filler's Expected Availability Date/Time
28	CWE	X	0177	Confidentiality Code
29	CWE	X	0482	Order Type
30	CNE	X	0483	Enterer Authorization Mode
31	CWE	X		Parent Universal Service Identifier
32	DT	X		Advanced Beneficiary Notice Date
33	CX	X		Alternate Placer Order Number

**ORC-1 Order Control**

Determines the function of the order segment. For the Prescription Query Response, the Order Control will always be "NW" (new order).

**ORC-2 Placer Order Number**

This field is the placer application's order number. The value is echoed back in the OBR-2 of the PCD-1 and PCD-4 messages.

**ORC-6 Response Flag**

This field allows the placer (sending) application to determine the amount of information to be returned from the filler. For the Prescription Query Response, the Response Flag will always be "N" (Only the MSA Segment is returned).

**ORC-12 Ordering Provider**

This field contains the identity of the person who is responsible for creating the request (i.e., ordering physician). If this field is present it may be used by the dialysis machine for prescription confirmation. In other words, the care-giver could ask if Dr. Smith is the patient's nephologist.

**ORC-14 Call Back Phone Number**

This field contains the telephone number to call for clarification of a request or other information regarding the order. If this field is present, a dialysis machine could be present its value to the user in case of a problem with prescription.

**8.2.4 OBR– Observation Request Segment**

The ORU message may include discrete OBX segments for individual observations reported. An OBR Segment will be used for each set of such OBX segments to establish the equipment context for the observations. All observation dates and times reported here should match OBX segments that report the same information.

HL7 Attribute Table - OBR Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBR
2	EI	C		Placer Order Number
3	EI	R		Filler Order Number
4	CWE	R	9999	Universal Service Identifier
5		X		Priority
6		X		Requested Date/Time
7	DTM	R		Observation Date/Time #
8	DTM	X		Observation End Date/Time #
9	CQ	X		Collection Volume *
10	XCN	X		Collector Identifier *
11	ID	X	0065	Specimen Action Code *
12	CWE	X	9999	Danger Code
13	ST	X		Relevant Clinical Information
14		X		Specimen Received Date/Time *
15		X		Specimen Source
16	XCN	X		Ordering Provider
17	XTN	X		Order Callback Phone Number
18	ST	X		Placer Field 1
19	ST	X		Placer Field 2
20	ST	X		Filler Field 1 +
21	ST	X		Filler Field 2 +
22	DTM	X		Results Rpt/Status Chng – Date/Time +
23	MOC	X		Charge to Practice +
24	ID	X	0074	Diagnostic Serv Sect ID
25	ID	C	0123	Result Status +
26	PRL	X		Parent Result +
27		X		Quantity/Timing
28	XCN	X		Result Copies To
29	EIP	X		Parent
30	ID	X	0124	Transportation Mode
31	CWE	X	9999	Reason for Study
32	NDL	X		Principal Result Interpreter +
33	NDL	X		Assistant Result Interpreter +
34	NDL	X		Technician +
35	NDL	X		Transcriptionist +
36	DTM	X		Scheduled Date/Time +
37	NM	X		Number of Sample Containers *
38	CWE	X	9999	Transport Logistics of Collected Sample *
39	CWE	X	9999	Collector's Comment *
40	CWE	X	9999	Transport Arrangement Responsibility
41	ID	X	0224	Transport Arranged
42	ID	X	0225	Escort Required
43	CWE	X	9999	Planned Patient Transport Comment
44	CNE	X	0088	Procedure Code
45	CNE	X	0340	Procedure Code Modifier
46	CWE	X	0411	Placer Supplemental Service Information
47	CWE	X	0411	Filler Supplemental Service Information
48	CWE	X	0476	Medically Necessary Duplicate Procedure Reason
49	CWE	X	0507	Result Handling
50	CWE	X		Parent Universal Service Identifier



HL7 Attribute Table - OBR Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
51	EI	X		Observation Group ID
52	EI	X		Parent Observation Group ID
53	CX	X		Alternate Placer Order Number

1656

1657

**OBR-2 Placer Order Number**

1658

If the dialysis machine downloaded a prescription, then this field will contain the value from OBC-2. Otherwise, the field is blank.

1659

1660

**OBR-3 Filler Order Number**

1661

**HL7 Definition:** This field is the order number associated with the filling application. This is a permanent identifier for an order and its associated observations.

1662

1663

**This specification:** This value is used to convey a unique identifier for the therapy being performed. It is a concatenation of the Dialysis Machine EUI-64 (see MSH-3) and a timestamp. The format is

1664

1665

Therapy\_ID^Machine\_Name^Extended\_Unique\_Identifier^EUI-64

1666

The Machine Name and Extended Unique identified are the same as conveyed in the MSH-3. The

1667

Therapy\_ID is the Extended Unique Identifier concatenated with the timestamp using the following format.

1668

XXXXXXXXXXXXXXXXXXYYYYMMDDhhmmss

1669

Where XXXXXXXXXXXXXXXXXXXX is the Extended Unique Identifier

1670

YYYY is the year

1671

MM is the month with zero padding.

1672

DD is the day of the month with zero padding.

1673

hh is the hours (0-23) with zero padding.

1674

mm is the minutes with zero padding

1675

ss is the seconds with zero padding

1676

Example of a Therapy\_ID, 080019FFFE3ED02D2011060204584

1677

Example of the OBR-3 Field,

1678

080019FFFE3ED02D20110602045842^ACME\_Dialysis\_Machine^080019FFFE3ED02D^EUI-64

1679

1680

**OBR-4 Universal Service ID**

1681

**HL7 Definition:** This field contains the identifier code for the requested observation/test/battery.

1682

**This specification:** Identifier and Text can identify unique OBR segments that partition observations. The values for this field will be taken from the 11073\_10103 MDC\_IDC\_SESS\_TYPE enumerator MDC\_IDC\_ENUM\_SESS\_TYPE.

1683

1684

1685

1686

**OBR-7 Observation Date/Time**

1687

**HL7 Definition:** This field is the clinically relevant date/time of the observation. In the case of observations taken directly from a subject, it is the actual date and time the observation was obtained. In the case of a specimen associated study, this field shall represent the date and time the specimen was collected or obtained. (This is a results-only field except when the placer or a third party has already drawn the

1688

1689

1690

specimen.) This field is conditionally required. When the OBR is transmitted as part of a report message, the field must be filled in. If it is transmitted as part of a request and a sample has been sent along as part of the request, this field must be filled in because this specimen time is the physiologically relevant date/time of the observation.

**This Specification:** This value will be the same as the message time, MSH-7.

**OBR-25 Result Status**

**HL7 Definition:** This field contains the status of results for this order. This conditional field is required whenever the OBR is contained in a report message. See HL7 Table 0085 - Observation result status codes interpretation.

**This Specification:** This value will always be Final ("F").

**8.2.5 OBX– Observation Results Segment**

Discrete OBX segments for individual observations will be encoded into separate OBX segments as individual observations or measurements. These OBX segments will be preceded by an appropriate OBR segment to set the context for observations dealing with the dialysis device, characteristics, or treatment measurements.

**HL7 Attribute Table - OBX Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID – OBX
2	ID	C	0125	Value Type
3	CWE	R	9999	Observation Identifier
4	ST	R		Observation Sub-ID
5	varies	R		Observation Value
6	CWE	R	9999	Units
7	ST	CE		References Range
8	CWE	X		Interpretation Codes
9	NM	X		Probability
10	ID	X	0080	Nature of Abnormal Test
11	ID	R	0085	Observation Result Status
12	DTM	X		Effective Date of Reference Range
13	ST	X		User Defined Access Checks
14	DTM	CE		Date/Time of the Observation
15	CWE	X	9999	Producer's ID
16	XCN	X		Responsible Observer
17	CWE	X	9999	Observation Method
18	EI	X		Equipment Instance Identifier
19	DTM	X		Date/Time of the Analysis
20	CWE	X	0163	Observation Site
21	EI	X		Observation Instance Identifier
22	CNE	X	0725	Mood Code
23	XON	X		Performing Organization Name
24	XAD	X		Performing Organization Address
25	XCN	X		Performing Organization Medical Director
26	ID	X	0909	Patient Results Release Category

1708 **OBX-1 Set ID**

1709 HL7 Specification: This field contains the sequence number.

1710 **OBX-2 Value Type**

1711 HL7 Specification: This field defines the data type of OBX-5, Observation Value.

1712 This specification: Only value identified in the table below will be used. The value is based upon the value  
1713 type as identified in Table 2 – Dialysis Machine Data Objects.

1714 **IEEE to HL7 Data Type Matching**

Applicable IEEE 11073 MDC_IDC types	HL7 v2 data type
String	ST
Enumerated	CWE or CNE
Date Time	DTM
Numeric	NM

1715

1716 **OBX-3 Observation Identifier**

1717 HL7 Specification: This field contains a unique identifier for the observation.

1718 This specification: The format is that of the Coded Element (CWE). Example:  
1719 "158606^MDC\_HDIALY\_DIALYSATE\_FLOW\_MODE^MDC".

1720 **OBX-4 Observation Sub-ID**

1721 HL7 Specification: This field is used to distinguish between multiple OBX segments with the same  
1722 observation ID organized under one OBR.

1723 This specification: This value is expressed in a dotted notation consisting of 4 parts, MDS, VMD, CHAN,  
1724 and Parameter.

1725 **OBX-5 Observation Value**

1726 HL7 Specification: This field contains the value observed by the observation producer.

1727 This specification: This is the actual value is expressed a string value.

1728 **OBX-6 Unit**

1729 HL7 Specification: This field contains the units of measurement for the value.

1730 This specification: Will be coded with the MDC\_IDC Nomenclature (based on UCUM) Unit for associated  
1731 observation. Example, "ml/min^milliliter per minute^UCUM"

1732 **OBX-7 Reference Range**

1733 HL7 Specification: Reference range for the value.

- 1734 a) lower limit-upper limit (when both lower and upper limits are defined, e.g., for Venous Pressure, 20-400)
- 1735 b) > lower limit (if no upper limit, e.g., >10)
- 1736 c) < upper limit (if no lower limit, e.g., <15)

1737 This specification: Dialysis Machines use this field in two ways.

- 1738 1) For numeric measurements, this field may be used to provide the alarm limits.
- 1739 2) For settings, this field may be used to specify legal values. For example, to require blood pump rate  
1740 remain below 200 ml/min, a value of "< 200" would be used. To lock the blood pump rate at 250, a  
1741 value of "250-250" would be used.

1742

**OBX-8 Interpretation Codes**

1743

HL7 Specification: One or more codes specifying a categorical assessment of the observation value (OBX-5), such as "Normal", "Abnormal", "Positive", "Negative", "Resistant", "Susceptible", etc.

1744

1745

This specification: This field is used to provide zero or more codes separated by the tilde ("~"). A dialysis machine provides two mandatory codes, Alarm priority and Alarm Type, and one optional code, Abnormality Type.

1746

1747

1748

The alarm priority is dialysis machine specific and may vary by manufacturer.

Alarm Priority	Code
High	PH
Medium	PM
Low	PL
Informational	PI
Priority Not Indicated	PN
Priority Unknown	PU

1749

The Alarm Type is specific in the Alert Type Column in Table 2 – Dialysis Machine Data Objects.

Alarm Type	Code
Physiological	SP
Technical	ST
Advisory	SA

1750

Abnormality Type	Code
Below low normal	L
Above high normal	H

1751

1752

**OBX-11 Observation Result Status**

1753

HL7 Specification: This field contains the observation result status. See HL7 Table 0085 - Observation result status codes interpretation.

1754

1755

This specification: This field will always be Final ("F").

1756

**OBX-14 Date/Time of Observation**

1757

HL7 Specification: This field is required in two circumstances. The first is when the observations reported beneath one report header (OBR) have different dates/times. This could occur in the case of queries, timed test sequences, or clearance studies where one measurement within a battery may have a different time than another measurement.

1758

1759

1760

1761

This specification: Dialysis Machines use this field to report the time and date of spot check values such as non-invasive blood pressure.

1762

1763

**OBX-17 Observation Method**

1764

HL7 Specification: This optional field can be used to transmit the method or procedure by which an observation was obtained when the sending system wishes to distinguish among one measurement obtained by different methods and the distinction is not implicit in the test ID.

1765

1766

1767

This specification: Dialysis Machines will follow the IHE PCD recommendations for this field.

OBX-17	Description
AMEAS^auto-measurement^MDC	A measurement automatically taken by the machine. Examples are Arterial Pressure and Dialysate Temperature.
MMEAS^manual-measurement^MDC	A measurement taken due to user interaction. An example would be a Blood Pressure measurement taken using the NIBP system.
ASET^auto-setting^MDC	A setting that is being automatically determined by the machine.
MSET^manual-setting^MDC	A setting that was manually set or changed by the user. A prescription setting that was changed by the user would fall into this category.
RSET^remote-setting^MDC	A setting that was remotely set through a prescription downloaded from an EMR.

1768 The field is optional for measurements. It is required for settings.

1769

### 1770 8.2.6 PID – Patient Identification Segment

1771 The PID segment is used by all applications as the primary means of communicating patient identification  
 1772 information. This segment contains permanent patient identifying and demographic information that, for  
 1773 the most part, is not likely to change frequently.

HL7 Attribute Table - PID Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID - PID
2		X		Patient ID
3	CX	R		Patient Identifier List
4		X		Alternate Patient ID - PID
5	XPN	R	0200	Patient Name
6	XPN	X		Mother's Maiden Name
7	DTM	X		Date/Time of Birth
8	CWE	X	0001	Administrative Sex
9		X		Patient Alias
10	CWE	X	0005	Race
11	XAD	X		Patient Address
12		X		County Code
13	XTN	X		Phone Number - Home
14	XTN	X		Phone Number - Business
15	CWE	X	0296	Primary Language
16	CWE	X	0002	Marital Status
17	CWE	X	0006	Religion
18	CX	X	0061	Patient Account Number
19		X		SSN Number - Patient
20		X		Driver's License Number - Patient
21	CX	X	0061	Mother's Identifier
22	CWE	X	0189	Ethnic Group
23	ST	X		Birth Place
24	ID	X	0136	Multiple Birth Indicator
25	NM	X		Birth Order
26	CWE	X	0171	Citizenship
27	CWE	X	0172	Veterans Military Status
28		X		Nationality

HL7 Attribute Table - PID Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
29	DTM	X		Patient Death Date and Time
30	ID	X	0136	Patient Death Indicator
31	ID	X	0136	Identity Unknown Indicator
32	CWE	X	0445	Identity Reliability Code
33	DTM	X		Last Update Date/Time
34	HD	X		Last Update Facility
35	CWE	X	0446	Species Code
36	CWE	X	0447	Breed Code
37	ST	X		Strain
38	CWE	X	0429	Production Class Code
39	CWE	X	0171	Tribal Citizenship
40	XTN	X		Patient Telecommunication Information

1774  
1775  
1776  
1777  
1778  
1779  
1780  
1781

**PID-3.1 Patient Identifier List**

ID Number contains a unique identifier for the patient assigned by the Device or entered by the User. Sample Identifier Type Codes are shown Table 0203 listed below (others can be included as defined in the HL7 Standard). The last identifier will always be the unique model/serial number of the Dialysis Machine with an identifier type code of U (see table following).

HL7 Table 0203

Code	Description	Notes	Use
MR	Medical Record Number		O
PN	Person Number		O
SS	Patient Social Security Number		O
U	Machine Identifier	Model, MDC_ID_MODEL_NUMBER, and Serial number, MDC_ID_PROD_SPEC_SERIAL, will be concatenated together to uniquely identify the machine. The format of the ID will be following: "model/serial"	R

1782  
1783  
1784  
1785  
1786  
1787  
1788  
1789

Example Patient Identifier List with just model and serial number.

`Scrubber 2000/SC678932^^^U`

Example Patient Identifier List with Medical Record Number

`555444222111^^^MR~Scrubber 2000/SC678932^^^U`

Example Patient Identifier List with Person Number

`010199-000H^^^PN~Scrubber 2000/SC678932^^^U`

**8.2.7 PV1 – Patient Visit Segment**

The PV1 segment is used by Registration/Patient Administration applications to communicate information on an account or visit-specific basis. Because it is required in the PCD-01 Message, it is included in this standard. The Dialysis Machine Provides the bare minimum of data to be compliant.

1794

HL7 Attribute Table - PV1 Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	SI	R		Set ID - PV1
2	CWE	R	0004	Patient Class
3	PL	O		Assigned Patient Location
4	CWE	X	0007	Admission Type
5	CX	X		Preadmit Number
6	PL	X		Prior Patient Location
7	XCN	X	0010	Attending Doctor
8	XCN	X	0010	Referring Doctor
9	XCN	X		Consulting Doctor
10	CWE	X	0069	Hospital Service
11	PL	X		Temporary Location
12	CWE	X	0087	Preadmit Test Indicator
13	CWE	X	0092	Re-admission Indicator
14	CWE	X	0023	Admit Source
15	CWE	X	0009	Ambulatory Status
16	CWE	X	0099	VIP Indicator
17	XCN	X	0010	Admitting Doctor
18	CWE	X	0018	Patient Type
19	CX	RE		Visit Number
20	FC	X	0064	Financial Class
21	CWE	X	0032	Charge Price Indicator
22	CWE	X	0045	Courtesy Code
23	CWE	X	0046	Credit Rating
24	CWE	X	0044	Contract Code
25	DT	X		Contract Effective Date
26	NM	X		Contract Amount
27	NM	X		Contract Period
28	CWE	X	0073	Interest Code
29	CWE	X	0110	Transfer to Bad Debt Code
30	DT	X		Transfer to Bad Debt Date
31	CWE	X	0021	Bad Debt Agency Code
32	NM	X		Bad Debt Transfer Amount
33	NM	X		Bad Debt Recovery Amount
34	CWE	X	0111	Delete Account Indicator
35	DT	X		Delete Account Date
36	CWE	X	0112	Discharge Disposition
37	DLD	X	0113	Discharged to Location
38	CWE	X	0114	Diet Type
39	CWE	X	0115	Servicing Facility
40		X		Bed Status
41	CWE	X	0117	Account Status
42	PL	X		Pending Location
43	PL	X		Prior Temporary Location
44	DTM	X		Admit Date/Time
45	DTM	X		Discharge Date/Time
46	NM	X		Current Patient Balance
47	NM	X		Total Charges
48	NM	X		Total Adjustments
49	NM	X		Total Payments
50	CX	X	0203	Alternate Visit ID

HL7 Attribute Table - PV1 Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
51	CWE	X	0326	Visit Indicator
52		X		Other Healthcare Provider
53	ST	X		Service Episode Description
54	CX	X		Service Episode Identifier

1795

**PV1-2 Patient Class (CWE)**

HL7 Definition: This field is used by systems to categorize patients by site. See HL7 Table 0004.

This specification: Dialysis machines will always use the Unknown Patient Class.

1796

1797

1798

1799

**PV1-3 Patient Class (PL)**

HL7 Definition: This field contains the patient's initial assigned location or the location to which the patient is being moved.

This specification: This field contains the patient's current location.

The fields of the Person Location (PL) data type are complex and the reader is referred to the HL7 Chapter on Data Types. Several examples from that chapter are shown below.

1800

1801

1802

1803

1804

Example: Nursing Unit

A nursing unit at Community Hospital: 4 East, room 136, bed B

4E^136^B^CommunityHospital^^N^^^

1806

1807

Example: Clinic

A clinic at University Hospitals: Internal Medicine Clinic located in the Briones building, 3rd floor.

InternalMedicine^^^UniversityHospitals^^C^Briones^3^

1808

1809

1810

Example: Home

The patient was treated at his home.

^^^^^H^^^

1812

1813

**PV1-19 Visit Number (CX)**

HL7 Definition: This field contains the unique number assigned to each patient visit.

This specification: This field will be empty.

1814

1815

1816

1817

**8.2.8 QAK – Query Acknowledge**

HL7 Attribute Table - QAK Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	R		Query Tag
2	ID	R	0208	Query Response Status
3	CWE	O		Message Query Name
4	NM	O		Hit Count Total
5	NM	O		Hits In this Message
6	NM	R		Hits Remaining

**QAK-1 Query Tag**

HL7 Specification: This field may be valued by the initiating system to identify the query and may be used to match response messages to the originating query.

1818

1819

1820



1821  
1822  
1823  
1824  
1825  
1826  
1827  
1828  
1829  
1830  
1831  
1832  
1833  
1834  
1835  
1836  
1837  
1838  
1839  
1840  
1841  
1842  
1843  
1844

**QAK-2 Query Response Status**

HL7 Specification: This field allows the responding system to return a precise response status. It is especially useful in the case where no data is found that matches the query parameters, but where there is also no error.

**QAK-3 Message Query Name**

HL7 Specification: This field contains the name of the query.

**QAK-4 Hit Count Total**

HL7 Specification: This field contains the total number of records found by the EMR that matched the query.

**QAK-5 Hits In This Message**

HL7 Specification: This field contains the total number of matching records that the Server sent in the current response.

**QAK-6 Hits Remaining**

HL7 Specification: This field contains the number of matching records found by the Server that have yet to be sent.

This specification: This value is used by the dialysis machine to determine when it has received all the results of a query.

**8.2.9 QPD – Query Parameter Definition**

HL7 Attribute Table - QPD Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	CWE	R		Message Query Name
2	ST	C		Query Tag
3-N	varies	R		Query Parameters

This segment is used in both the Patient Demographics Query (PDQ) transaction and the Patient Prescription Query (PPQ). The details of the fields are specified in those sections.

**8.2.10RCP – Response Control Parameters**

HL7 Attribute Table - RCP Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ID	R	0091	Query Priority
2	CQ	X	0126	Quantity Limited Request
3	CNE	R	0394	Response Modality
4	DTM	X		Execution and Delivery Time
5	ID	X	0395	Modify Indicators
6	SRT	X		Sort By Field

HL7 Attribute Table - RCP Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
7	ID	X	0391	Segment Group Inclusion

1849

1850

**RCP-1 Query Priority**

1851

HL7 Specification: This field contains the time frame in which the response is expected, See HL7 Table 0091

1852

– Query Priority.

1853

This specification: This field will always be Immediate ("I").

1854

1855

**RCP-3 Response Modality**

1856

HL7 Specification: This field specifies the timing and grouping of the response message(s), see HL7 Table

1857

0394 – Response Modality.

1858

This specification: This field will always be Real Time ("R").

1859

## 9 Hemodialysis Data Elements

### 9.1 Tables

This section defines tables of values used by the message elements.

**HD\_TBL\_01 – Mode of Operation**

<b>Value</b>	<b>Description</b>
PRETX	Pre-Treatment - Preparing for dialysis but no blood in extracorporeal circuit
TX	Treatment - Blood in the extracorporeal circuit
POSTTX	Post-Treatment - After dialysis but no longer any blood in extracorporeal circuit
DIS	Disinfection or rinse of dialysis fluid path.
IDL	Idle
SVC	Service mode

**HD\_TBL\_02 – Treatment Modality**

<b>Value</b>	<b>Description</b>
HD	Hemodialysis – dialysis with diffusive transport of solutes. AKA CVVHD
HDF	Hemodiafiltration – dialysis with both diffusive and convective transport of solutes, CVVHDF
HF	Hemofiltration – dialysis with convective transport of solutes. CVVH
SLED	Sustained low efficiency dialysis.
IUF	Isolated Ultrafiltration - Removal of plasma water without dialysis, SCUF
HP	Hemoperfusion – Removal of solutes from blood by adsorption.

**HD\_TBL\_03 – Anticoagulation Mode**

<b>Value</b>	<b>Description</b>
BOL	Bolus - Administration of a fixed amount of anticoagulant drug during a specific time
CON	Continuous - Anticoagulant drug administered at a constant rate
BOLCON	Combination of Bolus and Continuous mode
PRO	Profile - Anticoagulant drug administered at a variable rate according to a programmed profile
BOLPRO	Combination of Bolus and Profile mode
NONE	None

**HD\_TBL\_04 – Patient Position**

<b>Value</b>	<b>Description</b>
SITTING	Sitting
STANDING	Standing
SUPINE	Supine

**HD\_TBL\_05 – Blood Pump Mode**

<b>Value</b>	<b>Description</b>
2N	Double Needle - Separate needles/catheter lumens for arterial and venous blood flow
1N1P	Single Needle/Single Pump - The same needle/catheter lumen used for both arterial and venous blood flow, blood flow controlled by one pump

**HD\_TBL\_05 – Blood Pump Mode**

<b>Value</b>	<b>Description</b>
1N2P	Single Needle/Double Pump - The same needle/catheter lumen used for both arterial and venous blood flow, blood flow controlled by two pumps

**HD\_TBL\_06 – Dialysis Flow Mode**

<b>Value</b>	<b>Description</b>
CONST	Constant dialysis fluid flow
AUTO	Auto - Dialysis fluid flow is automatically controlled by the dialysis machine
PRO	Profile - Variable dialysis fluid flow according to a programmed profile
STBY	Standby - dialysis fluid flow is reduced when dialysis machine is not providing treatment
NONE	Dialysis fluid does not flow.

**HD\_TBL\_07 – Replacement Fluid Delivery Mode**

<b>Value</b>	<b>Description</b>
CONST	Constant replacement fluid flow
PRO	Profile - Variable replacement fluid flow according to a programmed profile
NONE	None, Not Applicable.
BOLUS	Bolus - administration of a fixed amount of replacement fluid during a specific time
AUTO	Automatic - Replacement fluid flow is automatically controlled by the dialysis machine

**HD\_TBL\_08 – RF Dilution Location**

<b>Value</b>	<b>Description</b>
PREF	Pre-Filter
POSTF	Post-Filter
PREF-POSTF	Both Pre-Filter and Post-Filter
PREP-POSTF	Both Pre-Pump and Post-Filter
MIDF	Mid Filter

**HD\_TBL\_09 – Fluid Source**

<b>Value</b>	<b>Description</b>
BAG	Bag – Replacement fluid
ONLINE	Online - Replacement fluid prepared by dialysis machine

**HD\_TBL\_10 – Plasma Volume Marker**

<b>Value</b>	<b>Description</b>
NONE	None – No selection made at time of reading/reporting
SYMP	Symptomatic – Noting a symptom for continue monitoring.
INTER	Intervention – Action was taken or a machine setting was changed

**HD\_TBL\_11 – Plasma Volume Profile**

<b>Value</b>	<b>Description</b>
FLAT	Change less than or equal to 3% per hour
GRAD	Change greater than 3% and less than or equal to 6.5% per hour

**HD\_TBL\_11 – Plasma Volume Profile**

<b>Value</b>	<b>Description</b>
STEEP	Change greater than 6.5% per hour

**HD\_TBL\_12 – Completion Method**

<b>Value</b>	<b>Description</b>
CT	Clock Time
AT	Actual Treatment Time
UF	UF Removed
KTV	Kt/V
USER	The treatment runs until stopped by the user.

**HD\_TBL\_13 – Ultrafiltration Mode**

<b>Value</b>	<b>Description</b>
NONE	None
CONST-WT	UF removed at a constant rate with a target volume
PRO-WT	UF removed according to a programmed profile with a Target Volume
CONST-WOT	UF removal at a constant rate without a target Volume
PRO-WOT	UF removed according to programmed profile without a target volume.

**HD\_TBL\_14 – Bicarb Delivery Mode**

<b>Value</b>	<b>Description</b>
CONST	Constant – Bicarbonate concentration in dialysis fluid is constant.
PRO	Profile - Bicarbonate concentration in

**HD\_TBL\_14 – Bicarb Delivery Mode**

<b>Value</b>	<b>Description</b>
	dialysis fluid is adjusted according to a programmed profile
NA	Not Applicable

**HD\_TBL\_15 – Sodium Delivery Mode**

<b>Value</b>	<b>Description</b>
CONST	Constant - Sodium concentration in dialysis fluid is constant
PRO	Profile - Sodium concentration in dialysis fluid is adjusted according to a programmed profile

**HD\_TBL\_16 – Anticoagulant Delivery Location**


<b>Value</b>	<b>Description</b>
PREP	Before the blood pump
PREF	After the blood pump but before the filter.
POSTF	After the filter.
PAT	At the patient

**HD\_TBL\_17 – Profile Type**

<b>Value</b>	<b>Description</b>
CONST	Constant rate
LINEAR	Linear profile
EXPONENTIAL	Exponential profile
STEP	Step profile
VENDOR	Vendor Specific Profile



## 9.2 Data Objects



In the table below when usage is conditional (C), a number is given. The number refers to one of the following conditions. If the condition is not met then the cardinality is "0..0".

- 1) If Anticoagulation Mode is not None
- 2) If Anticoagulation Mode includes Bolus
- 3) If Anticoagulation Mode includes Continuous or Profile
- 4) One for each anticoagulant pump
- 5) If Blood Pump Mode is not double needle
- 6) If Mode is HD or HDF
- 7) If Mode is HF or HDF and Dilution is Post or Pre-Post
- 8) If Mode is HDF or HF
- 9) If Replacement Fluid Mode is Bolus
- 10) If Completion Method is not None
- 11) if UF removal includes a target volume
- 12) If Mode is HD or HDF and the dialysis machine knows the conductivity through direct measurement or by knowing what type of bagged fluid is being used.
- 13) If alarm limit is being set.
- 14) If Completion Method (MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD) is Clock Time.
- 15) If Completion Method (MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD) is Treatment Time.
- 16) If Mode is HF or HDF and Dilution is Pre or Pre-Post
- 17) If Completion Method (MDC\_HDIALY\_THERAPY\_COMPLETE\_METHOD) is Kt/V.
- 18) If UF Target (MDC\_HDIALY\_UF\_TARGET\_VOL\_TO\_REMOVE) or Treatment Time (MDC\_HDIALY\_THERAPY\_TX\_TIME\_SETTING) is not present.
- 19) If UF Mode (MDC\_HDIALY\_UF\_MODE) is not constant.
- 20) If Anticoagulant Mode (MDC\_HDIALY\_ANTICOAG\_MODE) is not constant.
- 22) If Dialysate Mode (MDC\_HDIALY\_DIALYSATE\_FLOW\_MODE) is not constant.
- 23) If Sodium Mode (MDC\_HDIALY\_DIALYSATE\_CONC\_NA\_MODE) is not constant.
- 24) If RF Delivery Mode (MDC\_HDIALY\_RF\_FLOW\_MODE) is not constant.
- 91) If Profile Type is not Vendor then 'M' otherwise 'X'.
- 92) If Profile Type is not Vendor then 'O' otherwise 'X'.
- 93) If Profile Type is Exponential then 'O' otherwise 'X'.
- 94) If Profile Type is Vendor then 'M' otherwise 'X'.
- 95) If Profile Type is Vendor then 'O' otherwise 'X'.

In addition to the above conditions, the following general rule applies.

- 1) Any prescription setting sent from the EMR to the dialysis machine, will be reported in the PCD-01 Status message.

Any object below that is a cumulative total is the total since the start of treatment and not since the start of the filter.

Table 1 – Profile Parameters contains the objects that are needed to define a profile. These objects are referenced in Table 2 – Dialysis Machine Data Objects where appropriate.

Table 1 – Profile Parameters

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDCX_HDIALY_PROFILE_TYPE			Profile Type	Contains the anticoagulant profile type.	All	Episodic	Enum	HD_TBL_17	N/A	1..1	M	M
... MDCX_HDIALY_PROFILE_VALUE			Profile Values	Contains the control points of the profile. The control points are a numeric array. For a constant profile there is only one point. For a linear or exponential profile there are two points, one for the starting value and one for the ending	All	Episodic	Array	Numeric	mL/h	0..1	C91	C91
... MDCX_HDIALY_PROFILE_TIME			Profile Times	This parameter is optional. If it is not present, the control points are assumed to be equally space across the total treatment time. If it is present, it is an array of numerical values that are the treatment times at which a control point takes effect.	All	Episodic	Array	Numeric	min	0..1	C92	C92
... MDCX_HDIALY_PROFILE_EXP_HALF_TIME			Exponential Half Time	This parameter is optional and only applies to the Exponential Profile. It specifies the time at which 50% of the change has been complete.	All	Episodic	Numeric	Numeric	min	0..1	C93	C93
... MDCX_HDIALY_PROFILE_NAME			Profile Name	Vendor specific profile name.	All	Episodic	String	Alphanumeric	N/A	0..1	C94	C94
... MDCX_HDIALY_PROFILE_PARAMETERS			Profile Parameters	Vendor specific parameters for a vendor specific profile.	All	Episodic	Array	Numeric	N/A	0..1	C95	C95

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
MDC_DEV_HDIALY_MACHINE_MDS (70929)										1..1	M	M
... MDC_ID_MODEL_MANUFACTURER (531970)		Dialysis Device	Dialysis Device Manufacturer	String containing device manufacturer's name	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ID_MODEL_NUMBER (531969)		Dialysis Device	Dialysis Device Model	String containing device manufacturer's model identifier for the device	All	All	String	Alphanumeric	N/A	1..1	M	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_ID_PROD_SPEC_SERIAL (531972)		Dialysis Device	Dialysis Device Serial Number	String containing the device's serial number	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ID_PROD_SPEC_SW (531975)		Dialysis Device	Device Software Version	String defining the version of the software on the machine	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ATTR_ID_UDI (67916)		Dialysis Device	Dialysis Device UDI	String containing the device's FDA required Universal Device Identifier [Note: MDC_ATTR_ID_UDI contains three components: UdiAuthority, UdiIssuer, UdiLabel to support non-FDA entities, but currently there is no approved way of encoding this using HL7 V2; PRT-10 and PRT-16-20 is used instead; UdiLabel = PRT-10.]	All	All	String	Alphanumeric	N/A	1..1	O	X
. MDC_DEV_HDIALY_VMD (70934)										1..1	M	M
.. MDC_DEV_HDIALY_MACH_CONFIG_CHAN (70939)										1..1	M	M
... MDC_ATTR_TIME_ABS (67975)		Dialysis Device	Dialysis Device Time	Date and time as recorded on the dialysis device's internal clock with offset to UTC (based on location and DST)	All	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	M	X
... MDC_HDIALY_MACH_MODE_DESCRIPTION (158593)		Dialysis Device	Machine Mode Description	Manufacturer-specific description of the machine mode of operation	All	Episodic	String	Alphanumeric	N/A	0..1	O	X
... MDC_HDIALY_MACH_MODE_OF_OPERATION (158594)		Dialysis Device	Machine Mode Of Operation	Process the machine is currently performing	All	Episodic	Enum	HD_TBL_01	N/A	1..1	M	X
... MDC_TIME_PD_MAINTENANCE_TO_NEXT_SERVICE (184195)		Dialysis Device	Preventive Maintenance Time Remaining	Hours of operation remaining until next preventative maintenance is required	All	Episodic	Numeric	XXXX	hours	0..1	O	X
... MDC_MAINTENANCE_NEXT_SERVICE_DATE (184199)		Dialysis Device	Preventative Maintenance Due Date	The date on which preventative maintenance is due.	All	Episodic	Date	YYYYMMDD	N/A	0..1	O	X
... MDC_HDIALY_MACH_MAINT_TX_REMAIN (158595)		Dialysis Device	Preventive Maintenance Treatments Remaining	Number of treatments remaining until next preventative maintenance is required	All	Episodic	Numeric	XXXX	N/A	0..1	O	X
... MDC_HDIALY_MACH_BLD_PUMP_ON (158596)		Dialysis Device	Treatment Blood Pump On	Status indicating if blood pump is running	Intradialytic	Episodic	Bool	T / F	N/A	1..1	M	X
... MDC_HDIALY_MACH_TX_FLUID_BYPASS (158597)		Dialysis Device	Treatment Fluid Bypass	Fluid pathway state to indicate if fluid is NOT flowing through dialyzer	Intradialytic	Episodic	Bool	T / F	N/A	1..1	M	X
... MDC_HDIALY_MACH_TX_MODALITY (158598)		Dialysis Device	Treatment Modality	The way a patient receives dialysis	All	Episodic	Enum	HD_TBL_02	N/A	1..1	M	M

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use	
... MDC_HDIALY_MACH_THERAPY_TIME (158720)		Dialysis Device State	Elapsed Treatment Time	Accumulated time that the patient has been actively treated. This does not include time spend in alarms, recirculating or rinsing back.	Intradialytic		Periodic	Numeric	XXXX	minutes	1..1	M	X
... MDC_HDIALY_MACH_TIME_REMAIN (158724)		Dialysis Device State	Remaining Treatment Time	Number of minutes left for current dialysis session	Intradialytic		Periodic	Numeric	XXXX	minutes	1..1	C10	X
... MDC_TEMP_ROOM (188508)		Dialysis Device State	Room Temperature	Room Temperature in which the Dialysis Machine is operating	All		Periodic	Numeric	XX.X	°C	0..1	O	X
.. MDC_DEV_HDIALY_ANTICOAG_PUMP_CHAN (70943)											1..*	C1	C1
... MDC_ATTR_CHAN_NUM_LOGICAL (68142)		Anticoagulant Pump	Anticoagulant Pump Number	The number associated with the anticoagulant pump in this channel. The value starts at 1 and is incremented for each anticoagulant pump in use.	Intradialytic		Episodic	Numeric	X	N/A	1..1	M	M
... MDC_HDIALY_ANTICOAG_MODE (158600)		Anticoagulant Pump	Anticoagulation Mode	Method of anticoagulant administration during dialysis treatment	Intradialytic		Episodic	Enum	HD_TBL_03	N/A	1..1	C4	M
... Profile Parameters				Anticoagulant Profile								C20	C20
... MDC_HDIALY_ANTICOAG_INFUS_RATE_SETTING (0)		Anticoagulant Pump	Anticoagulant Infusion Rate Setting	User defined amount of anticoagulation drug to be administered per unit time	Intradialytic		Episodic	Numeric	XX.X	mL/h	0..1	C3	M
... MDC_HDIALY_ANTICOAG_INFUS_RATE (158736)		Anticoagulant Pump	Anticoagulant Infusion Rate	User defined amount of anticoagulation drug to be administered per unit time	Intradialytic		Episodic	Numeric	XX.X	mL/h	0..1	C3	X
... MDC_HDIALY_ANTICOAG_ACCUM_DELIV (158728)		Anticoagulant Pump	Accumulated Delivered Anticoagulant	Total amount of anticoagulation drug administered to patient since the start of treatment	Intradialytic		Both	Numeric	XXX.X	mL	0..1	C1	X
... MDC_HDIALY_ANTICOAG_BOLUS_VOL (158732)		Anticoagulant Pump	Anticoagulant Bolus Volume	User defined amount of anticoagulation drug to be administered in this single large dose	Intradialytic		Episodic	Numeric	XX.X	mL	0..1	C2	O
... MDC_HDIALY_INFUSION_TIME_SETTING (0)		Anticoagulant Pump	Anticoagulant Infusion Time Setting	User specified infusion time.	Intradialytic		Episodic	Numeric	XXXX	minutes	0..1	O	O
... MDC_EVT_HDIALY_ANTICOAG_STOP (198236)	tech	Anticoagulant Pump	Anticoagulant Pump Stop	Notification that the anticoagulant pump has stopped	Intradialytic		Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_ANTICOAG_NAME (158599)		Anticoagulant Pump	Anticoagulant Type	Name of the anticoagulant	Intradialytic		Episodic	String	Alphanumeric	N/A	0..1	O	M

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temp-oral	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_EVT_HDIALY_ANTICOAG_SYRINGE_EMPTY (198238)	tech	<b>Anticoagulant Pump</b>	Syringe Empty	Notification that the syringe plunger has reached the end of its functional range	Intradialytic	Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_ANTICOAG_SYRINGE_BRAND (158602)		<b>Anticoagulant Pump</b>	Syringe Name	Name of the syringe	Intradialytic	Episodic	String	Alphanumeric	N/A	0..1	O	O
... MDC_HDIALY_ANTICOAG_SYRINGE_VOL (158603)		<b>Anticoagulant Pump</b>	Syringe Volume	Volume of the syringe with anticoagulant.	Intradialytic	Episodic	Numeric	XXX	mL	0..1	O	O
... MDC_EVT_HDIALY_ANTICOAG_SYRINGE_SIZE (198240)	tech	<b>Anticoagulant Pump</b>	Wrong Syringe Size	Notification that the size of the syringe is not the size defined for anticoagulant administration	Intradialytic	Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_ANTICOAG_LOCATION (0)		<b>Anticoagulant Pump</b>	Anticoagulant Deliver Location	Location for delivery of the anticoagulant.	Intradialytic	Episodic	Enum	<b>HD_TBL_16</b>	N/A	0..1	O	M
<b>.. MDC_DEV_HDIALY_BLOOD_PUMP_CHAN (70947)</b>										<b>1..1</b>	<b>M</b>	<b>M</b>
... MDC_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_SETTING (16935956)		<b>Blood Pump</b>	Blood Flow Rate Setting	The rate at which the user programmed the blood flow	Intradialytic	Episodic	Numeric	XXX	mL/min	1..1	M	M
... MDC_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE (158740)		<b>Blood Pump</b>	Actual Blood Flow Rate	Adjusted blood flow rate based on the blood flow rate setting and pressure drop caused by blood line, needle and vascular access	Intradialytic	Periodic	numeric	XXX	mL/min	0..1	O	X
... MDC_HDIALY_BLD_PUMP_BLOOD_FLOW_RATE_MEAN (158743)		<b>Blood Pump</b>	Average Blood Flow Rate	Average of the actual blood flow rate over the course of treatment.	Intradialytic	Periodic	Numeric	XXX	mL/min	0..1	O	X
... MDC_HDIALY_BLD_PRESS_ART (158744)	phys tech high low thr	<b>Blood Pump</b>	Arterial Pressure	Pressure of arterial access line pre blood pump	Intradialytic	Periodic	Numeric	±XXX	mmHg	1..1	M	X
... MDC_HDIALY_BLD_PUMP_MODE (158604)		<b>Blood Pump</b>	Blood Pump Mode	Therapy method in which blood is retrieved and returned to the patient.	Intradialytic	Episodic	Enum	<b>HD_TBL_05</b>	N/A	1..1	M	M
... MDC_EVT_HDIALY_BLD_PUMP_STOP (198242)	tech	<b>Blood Pump</b>	Blood Pump Stop	Notification that the blood pump has stopped	Intradialytic	Episodic	Bool	T / F	N/A	0..1	M	X
... MDC_HDIALY_BLD_PUMP_TUBING_SIZE (158605)		<b>Blood Pump</b>	Blood Tubing Size	Inner diameter of pump segment of tubing used in extracorporeal circulation	Intradialytic	Episodic	String	alphanumeric	N/A	0..1	O	X
... MDC_HDIALY_BLOOD_TEMP_ART (158748)		<b>Blood Pump</b>	Arterial Blood Temperature	Arterial blood temperature as measured by the device.	Intradialytic	Periodic	Numeric	XX.X	°C	0..1	O	X
... MDC_HDIALY_BLD_PUMP_CHANGE_IN_ENERGY (158752)		<b>Blood Pump</b>	Change In Energy	Measurement of the energy (temperate) flux to and from the patient	Intradialytic	Periodic	Numeric	±XXX.X	kJ/h	0..1	O	X
... MDC_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP (158756)	tech high low thr	<b>Blood Pump</b>	Post Pump Arterial Pressure	Pressure of arterial access line post blood pump	Intradialytic	Periodic	Numeric	±XXX	mmHg	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_BLD_PU MP_PRIMING_VOL (158760)		<b>Blood Pump</b>	Priming Volume	Volume of fluid used to prime the extracorporeal circuit	Intradialytic		Numeric	XXX	mL	0..1	O	X
... MDC_HDIALY_BLD_PU MP_SINGLE_NEEDLE_PRES S (158764)	tech high low thr	<b>Blood Pump</b>	Single Needle Switching Pressure	In single needle mode, the pressure used to switch phases at high pressure	Intradialytic		Numeric	± XXX	mmHg	0..1	O	X
... MDC_HDIALY_BLD_PU MP_STROKE_VOL (158768)		<b>Blood Pump</b>	Stroke Volume	The volume of blood passing through the dialyzer during each single needle cycle.	Intradialytic		Numeric	XXX	mL	0..1	C5	X
... MDC_HDIALY_BLD_PU MP_BLOOD_PROCESSED_T OTAL (158772)		<b>Blood Pump</b>	Total Blood Processed	Accumulated volume of blood circulated by extracorporeal circulation.	Intradialytic		Numeric	XXXX.XX	L	0..1	O	X
... MDC_HDIALY_BLD_PU MP_PRESS_VEN (158776)	phys tech high low thr	<b>Blood Pump</b>	Venous Pressure	Pressure of the venous access line	Intradialytic		Numeric	±XXX	mmHg	1..1	M	X
... MDC_HDIALY_BLOOD_ TEMP_VEN (158780)		<b>Blood Pump</b>	Venous Temperature	Temperature of the blood measured in the venous access line	Intradialytic		Numeric	XX.X	°C	0..1	O	X
... MDC_HDIALY_BLOOD_ TEMP_VEN_SETTING (16935996)		<b>Blood Pump</b>	Venous Temperature Setting	Desired temperature of the blood in the venous access line	Intradialytic		Numeric	XX.X	°C	0..1	O	O
<b>.. MDC_DEV_HDIALY_FLUID D_CHAN (70951)</b>										<b>1..1</b>	<b>C6</b>	<b>C6</b>
... MDC_HDIALY_DIALYSA TE_FLOW_MODE (158606)		<b>Dialysis Fluid</b>	Dialysate Flow Mode	The mode of controlling dialysate flow rate	All		Enum	<b>HD_TBL_06</b>	N/A	1..1	C6	M
... Profile Parameters				Dialysate pump profile parameters							C22	C22
... MDC_HDIALY_DIALYSA TE_FLOW_RATE_SETTING (16936008)		<b>Dialysis Fluid</b>	Dialysate Flow Rate Setting	Rate at which the user programmed the dialysate flow	Intradialytic		Numeric	XXXX	mL/min	0..1	C6	M
... MDC_HDIALY_DIALYSA TE_FLOW_RATE (158792)	tech low	<b>Dialysis Fluid</b>	Actual Dialysate Flow Rate	Rate of dialysate flow at any given point during a treatment	Intradialytic		Numeric	XXXX	mL/min	1..1	C6	X
... MDC_HDIALY_DIALYSA TE_VOL_DELIV (158848)		<b>Dialysis Fluid</b>	Total Dialysate Volume	Total amount of dialysate delivered to patient during a treatment	Intradialytic		Numeric	XXX.XX	L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_VOL_DELIV_SETTING (0)		<b>Dialysis Fluid</b>	Dialysate Volume Setting	User programmed dialysate volume to deliver	Intradialytic		Numeric	XXX.XX	L	0..1	O	O
... MDC_HDIALY_DIALYSA TE_NAME (158608)		<b>Dialysis Fluid</b>	Dialysate Name	String containing a descriptive name of the dialysate	Intradialytic		String	Alphanumeric	N/A	1..1	O	M
... MDC_HDIALY_BICARB_ COND (158784)		<b>Dialysis Fluid</b>	Actual Bicarbonate Conductivity	Measured conductivity of bicarbonate	Intradialytic		Numeric	XX.XX	mS/cm	0..1	O	X
... MDC_HDIALY_BICARB_ COND_SETTING (16936000)		<b>Dialysis Fluid</b>	Bicarbonate Conductivity Setting	Bicarbonate conductivity	Intradialytic		Numeric	XX.XX	mS/cm	0..1	O	X



## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_DIALYSA TE_COND (158788)	tech high low	Dialysis Fluid	Actual Conductivity	Measured conductivity of dialysate	Intradialytic		Numeric	XX.XX	mS/cm	1..1	C12	X
... MDC_HDIALY_DIALYSA TE_COND_SETTING (16936004)		Dialysis Fluid	Dialysate Conductivity Setting	Acid and Bicarbonate conductivity combined	Intradialytic		Numeric	XX.XX	mS/cm	0..1	O	X
... MDC_HDIALY_DIALYSA TE_TEMP (158796)	tech high low	Dialysis Fluid	Actual Dialysate Temperature	Temperature of dialysate	Intradialytic		Numeric	XX.X	°C	1..1	O	X
... MDC_HDIALY_DIALYSA TE_TEMP_SETTING (16936012)		Dialysis Fluid	Dialysate Temperature Setting	Temperature at which dialysate is to be delivered to patient	Intradialytic		Numeric	XX.X	°C	1..1	O	O
... MDC_HDIALY_DIALYSA TE_AMMONIA (158800)		Dialysis Fluid	Ammonia Concentration	Measured ammonia concentration of the dialysate	Intradialytic		Numeric	XX.XX	ppm	0..1	O	X
... MDC_HDIALY_DIALYSA TE_FLOW_RATE_MEAN (158795)		Dialysis Fluid	Average Dialysate Flow Rate	Average rate of dialysate flow during course of treatment	Intradialytic		Numeric	XXXX	mL/min	0..1	O	X
... MDC_HDIALY_CONC_H CO3_SETTING (16936020)		Dialysis Fluid	Bicarbonate Concentration Setting	User programmed Bicarbonate concentration in dialysate.	Intradialytic		Numeric	XXX.X	mmol/L	1..1	O	O
... MDC_HDIALY_CONC_H CO3_MODE (158607)		Dialysis Fluid	Bicarbonate Mode	Process for delivering adjusted bicarbonate concentration in dialysate	Intradialytic		Enum	HD_TBL_14	N/A	1..1	O	O
... MDC_EVT_HDIALY_BLO OD_LEAK (198244)	tech	Dialysis Fluid	Blood Leak	Notification that blood has been detected in the dialysate fluid	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_HDIALY_DIALYSA TE_CONC_ACETATE (158808)		Dialysis Fluid	Dialysate Acetate Concentration	Acetate concentration in dialysate	Intradialytic		Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_HCO3 (158812)		Dialysis Fluid	Dialysate Bicarbonate Concentration	Bicarbonate concentration in dialysate	Intradialytic		Numeric	XXX.X	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_CHLORIDE (158816)		Dialysis Fluid	Dialysate Chloride Concentration	Chloride concentration in dialysate	Intradialytic		Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_MG (158820)		Dialysis Fluid	Dialysate Magnesium Concentration	Magnesium concentration in dialysate	Intradialytic		Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_CA (158824)		Dialysis Fluid	Dialysate Calcium Concentration	Calcium concentration in dialysate	Intradialytic		Numeric	X.XX	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_CITRATE (158828)		Dialysis Fluid	Dialysate Citrate Concentration	Citrate concentration in dialysate	Intradialytic		Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_GLU (158832)		Dialysis Fluid	Dialysate Glucose Concentration	Glucose concentration in dialysate	Intradialytic		Numeric	XXX	mmol/L	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_DIALYSA TE_CONC_K (158836)		Dialysis Fluid	Dialysate Potassium Concentration	Potassium concentration in dialysate	Intradialytic	Episodic	Numeric	X.X	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_NA (158840)		Dialysis Fluid	Dialysate Sodium Concentration	Sodium concentration in dialysate	Intradialytic	Episodic	Numeric	XXXX	mmol/L	0..1	O	X
... MDC_HDIALY_DIALYSA TE_CONC_NA_SETTING (16936056)		Dialysis Fluid	Dialysate Sodium Concentration Setting	User programmed sodium concentration in dialysate	Intradialytic	Episodic	Numeric	XXXX	mmol/L	1..1	O	O
... MDC_HDIALY_DIALYSA TE_CONC_NA_MODE (158609)		Dialysis Fluid	Dialysate Sodium Mode	Process for delivering adjusted sodium concentration in dialysate	Intradialytic	Episodic	Enum	HD_TBL_15	N/A	1..1	O	O
... Profile Parameters				Sodium Mode							C23	C23
... MDC_HDIALY_DIALYSA TE_CONC_PH (158844)		Dialysis Fluid	pH	Measured pH of the dialysate	Intradialytic	Periodic	Numeric	XX.X	pH	0..1	O	X
.. MDC_DEV_HDIALY_FILTER_CHAN (70955)										1..1	M	M
... MDC_HDIALY_FILTER_NAME (158610)		Filter (Dialyzer)	Dialyzer Name	Dialyzer manufacturer and model	Intradialytic	Episodic	String	Alphanumeric	N/A	0..1	O	M
... MDC_HDIALY_FILTER_TRANSMEMBRANE_PRESS (158852)	tech high low	Filter (Dialyzer)	Actual Transmembrane Pressure	Difference in pressure between blood compartment and dialysate compartment of dialyzer	Intradialytic	Periodic	Numeric	±XXX	mmHg	1..1	M	X
... MDC_HDIALY_FILTER_NUMBER (158611)		Filter (Dialyzer)	Filter Number	The filter number since the start of treatment. This value starts at 1 and increments with each new filter used during the current treatment.	Intradialytic	Periodic	Numeric	X	N/A	0..1	O	X
... MDC_HDIALY_FILTER_UDI (158612)		Filter (Dialyzer)	Filter UDI	The UDI for the current filter.	Intradialytic	Periodic	String	N/A	N/A	0..1	O	X
.. MDC_DEV_HDIALY_CONVECTIVE_CHAN (70959)										1..1	C8	C8
... MDC_HDIALY_RF_DILUTION_LOCATION (158613)		Replacement Fluid	Replacement Fluid Dilution Location	Delivery location of the replacement fluid	Intradialytic	Episodic	Enum	HD_TBL_08	N/A	1..1	C8	M
... MDC_HDIALY_RF_FLOW_MODE (158616)		Replacement Fluid	Replacement Fluid Flow Mode	The mode of controlling replacement fluid flow rate	All	Episodic	Enum	HD_TBL_07	N/A	1..1	C8	C8
... Profile Parameters				RF Delivery profile parameters							C24	C24
... MDC_HDIALY_RF_POST_FILTER_FLUID_NAME (158614)		Replacement Fluid	Post Filter Replacement Fluid Name	String containing a descriptive name of the post-filter replacement fluid	Intradialytic	Episodic	String	Alphanumeric	N/A	1..1	O	C7
... MDC_HDIALY_RF_POST_FILTER_CONC_ACETATE (158856)		Replacement Fluid	Post-Filter RF Acetate Concentration	Acetate concentration in the post-filter replacement fluid	Intradialytic	Episodic	Numeric	XX.X	mmol/L	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_RF_POST_FILTER_CONC_HCO3 (158860)		Replacement Fluid	Post-Filter RF Bicarbonate Concentration	Bicarbonate concentration in the post-filter replacement fluid	Intradialytic		Numeric	XXX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_CHLORIDE (158864)		Replacement Fluid	Post-Filter RF Chloride Concentration	Chloride concentration in the post-filter replacement fluid	Intradialytic		Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_MG (158868)		Replacement Fluid	Post-Filter RF Magnesium Concentration	Magnesium concentration in the post-filter replacement fluid	Intradialytic		Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_CA (158872)		Replacement Fluid	Post-Filter RF Calcium Concentration	Calcium concentration in the post-filter replacement fluid	Intradialytic		Numeric	X.XX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_CITRATE (158876)		Replacement Fluid	Post-Filter RF Citrate Concentration	Citrate concentration in the post-filter replacement fluid	Intradialytic		Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_GLU (158880)		Replacement Fluid	Post-Filter RF Glucose Concentration	Glucose concentration in the post-filter replacement fluid	Intradialytic		Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_K (158884)		Replacement Fluid	Post-Filter RF Potassium Concentration	Potassium concentration in the post-filter replacement fluid	Intradialytic		Numeric	X.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_CONC_NA (158888)		Replacement Fluid	Post-Filter RF Sodium Concentration	Sodium concentration in the post-filter replacement fluid	Intradialytic		Numeric	XXXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_FLOW_RATE_SETTING (16936108)		Replacement Fluid	Post-Filter Replacement Fluid Flow Rate Setting	Rate at which the user programmed the post-filter replacement fluid flow	Intradialytic		Numeric	XXXX	mL/min	0..1	C7	C7
... MDC_HDIALY_RF_POST_FILTER_FLOW_RATE (158892)		Replacement Fluid	Actual Post-Filter Replacement Fluid Flow Rate	Rate of post-filter replacement fluid flow at any given point during a treatment	Intradialytic		Numeric	XXXX	mL/min	1..1	C7	X
... MDC_HDIALY_RF_POST_FILTER_FLOW_RATE_MEAN (158895)		Replacement Fluid	Average Post-Filter Replacement Fluid Flow Rate	Average rate of post-filter replacement fluid flow during course of treatment.	Intradialytic		Numeric	XXXX	mL/min	0..1	O	X
... MDC_HDIALY_RF_POST_FILTER_VOL_SETTING (16936116)		Replacement Fluid	Post-Filter Replacement Fluid Volume setting	Programmed total amount of post-filter replacement fluid volume	Intradialytic		Numeric	XXX.XX	L	0..1	O	C7
... MDC_HDIALY_RF_POST_FILTER_VOL (158900)		Replacement Fluid	Total Post-Filter Replacement Fluid Volume	Total amount of post-filter replacement fluid used during a treatment	Intradialytic		Numeric	XXX.XX	L	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_RF_POST_FILTER_TEMP_SETTING (16936112)		Replacement Fluid	Post-Filter Replacement Fluid Temperature Setting	Temperature at which post-filter replacement fluid is to be delivered to patient	Intradialytic		Numeric	XX.X	°C	1..1	C7	O
... MDC_HDIALY_RF_POST_FILTER_TEMP (158896)		Replacement Fluid	Actual Post-Filter Replacement Fluid Temperature	Temperature of post-filter replacement fluid	Intradialytic	Periodic	Numeric	XX.X	°C	1..1	C7	X
... MDC_HDIALY_RF_PRE_FILTER_FLUID_NAME (158615)		Replacement Fluid	Pre-Filter RF Name	String containing a descriptive name of the pre-filter replacement fluid	Intradialytic	Episodic	String	Alphanumeric	N/A	1..1	O	C16
... MDC_HDIALY_RF_PRE_FILTER_CONC_ACETATE (158904)		Replacement Fluid	Pre-Filter RF Acetate Concentration	Acetate concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_HCO3 (158908)		Replacement Fluid	Pre-Filter RF Bicarbonate Concentration	Bicarbonate concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XXX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_CHLORIDE (158912)		Replacement Fluid	Pre-Filter RF Chloride Concentration	Chloride concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_MG (158916)		Replacement Fluid	Pre-Filter RF Magnesium Concentration	Magnesium concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_CA (158920)		Replacement Fluid	Pre-Filter RF Calcium Concentration	Calcium concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	X.XX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_CITRATE (158924)		Replacement Fluid	Pre-Filter RF Citrate Concentration	Citrate concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XX.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_GLU (158928)		Replacement Fluid	Pre-Filter RF Glucose Concentration	Glucose concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_K (158932)		Replacement Fluid	Pre-Filter RF Potassium Concentration	Potassium concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	X.X	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_CONC_NA (158936)		Replacement Fluid	Pre-Filter RF Sodium Concentration	Sodium concentration in the pre-filter replacement fluid	Intradialytic	Episodic	Numeric	XXXX	mmol/L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_FLOW_RATE_SETTING (16936156)		Replacement Fluid	Pre-Filter Replacement Fluid Flow Rate Setting	Rate at which the user programmed the pre-filter replacement fluid flow	Intradialytic	Episodic	Numeric	XXXX	mL/min	0..1	C8	C16
... MDC_HDIALY_RF_PRE_FILTER_FLOW_RATE (158940)		Replacement Fluid	Actual Pre-Filter Replacement Fluid Flow Rate	Rate of pre-filter replacement fluid flow at any given point during a treatment	Intradialytic	Periodic	Numeric	XXXX	mL/min	1..1	C7	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_RF_PRE_FILTER_FLOW_RATE_MEAN (158943)		Replacement Fluid	Average Pre-Filter Replacement Fluid Flow Rate	Average rate of pre-filter replacement fluid flow during course of treatment.	Intradialytic	Periodic	Numeric	XXXX	mL/min	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_VOL_SETTING (16936164)		Replacement Fluid	Pre-Filter Replacement Fluid Volume Setting	Programmed total amount of pre -filter replacement fluid volume	Intradialytic	Episodic	Numeric	XXX.XX	L	0..1	O	C16
... MDC_HDIALY_RF_PRE_FILTER_VOL (158948)		Replacement Fluid	Total Pre-Filter Replacement Fluid Volume	Total amount of pre-filter replacement fluid used during a treatment	Intradialytic	Periodic	Numeric	XXX.XX	L	0..1	O	X
... MDC_HDIALY_RF_PRE_FILTER_TEMP_SETTING (16936160)		Replacement Fluid	Pre-Filter Replacement Fluid Temperature Setting	Temperature at which pre-filter replacement fluid is to be delivered to patient	Intradialytic	Episodic	Numeric	XX.X	°C	1..1	C7	O
... MDC_HDIALY_RF_PRE_FILTER_TEMP (158944)		Replacement Fluid	Actual Pre-Filter Replacement Fluid Temperature	Temperature of pre-filter replacement fluid	Intradialytic	Periodic	Numeric	XX.X	°C	1..1	C7	X
... MDC_HDIALY_RF_CONV_CLEARANCE (158952)		Replacement Fluid	Convective Clearance - Urea	The product of the sieving coefficient and the total water flux. KX = QUF x SX where KX is Convective Clearance for solute "X", QUF is Ultrafiltration flow across membrane and SX is Membrane sieving coefficient for solute "X"	Intradialytic	Periodic	Numeric	XXX.X	mL/min	0..1	O	X
... MDC_EVT_HDIALY_RF_EXCESS_DELIV (198246)	tech-high	Replacement Fluid	Excessive Replacement Fluid Bolus Delivered	Notification that more replacement fluid has been delivered than intended	Intradialytic	Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_EVT_HDIALY_RF_INSUFF_DELIV (198248)	tech-low	Replacement Fluid	Insufficient Replacement Fluid Bolus Delivered	Notification that less replacement fluid has been delivered than intended	Intradialytic	Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_RF_POST_FLOW_RATIO_SETTING (16936172)		Replacement Fluid	Pre-Post Replacement Fluid Flow Rate Ratio Setting	The ratio desired between the pre and post Replacement Fluid Flow Rate Settings	Intradialytic	Episodic	Numeric	XX.X	N/A	0..1	O	O
... MDC_HDIALY_RF_BOLUS_RATE (158960)	tech-high low	Replacement Fluid	Replacement Fluid Bolus Rate	User defined replacement fluid bolus rate.	Intradialytic	Episodic	Numeric	XXXX	mL/min	0..1	C9	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_RF_BOLUS_VOL_SETTING (16936180)		<b>Replacement Fluid</b>	Replacement Fluid Bolus Volume Setting	User defined amount of replacement fluid to be delivered in this single large dose	Intradialytic		Numeric	XXX	mL	0..1	C9	O
... MDC_HDIALY_RF_BOLUS_VOL_DELIVERED (158968)	tech high low	<b>Replacement Fluid</b>	Replacement Fluid Bolus Volume Delivered	Amount of replacement fluid delivered in a single large dose.	Intradialytic		Numeric	XXX	mL	0..1	C9	X
... MDC_HDIALY_RF_SOURCE (158617)		<b>Replacement Fluid</b>	Replacement Fluid Source	The source of the replacement fluid.	Intradialytic		Enum	<b>HD_TBL_09</b>	N/A	0..1	O	X
... MDC_EVT_HDIALY_RF_BOLUS (198250)		<b>Replacement Fluid</b>	Replacement Fluid Bolus Event	True when an RF Bolus is occurring.	Intradialytic		Bool	T / F	N/A	0..1	C9	X
... MDC_DEV_HDIALY_SAFETY_SYSTEMS_CHAN (70963)										<b>1..1</b>	<b>M</b>	<b>X</b>
... MDC_EVT_HDIALY_SAFETY_ART_AIR_DETECT (198252)	tech	<b>Safety Systems</b>	Arterial Air Detector	Notification that air is sensed in the arterial blood lines	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_EVT_HDIALY_SAFETY_DIALYSATE_AIR_DETECT (198254)	tech	<b>Safety Systems</b>	Dialysate Air Detector	Notification that air is sensed in dialysate line	Intradialytic		Bool	T / F	N/A	0..1	C6	X
... MDC_EVT_HDIALY_SAFETY_DIALYSATE_COMPOSITION (198256)	tech	<b>Safety Systems</b>	Dialysis Composition	Notification that the dialyzing fluid composition is incorrect	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_EVT_HDIALY_SAFETY_SYSTEM_GENERAL (198258)	tech	<b>Safety Systems</b>	General System	Notification of any dialysis machine malfunction that are not related to patient treatment	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_EVT_SELFTEST_FAILURE (198216)	tech	<b>Safety Systems</b>	Self-Test	Notification that a self-test failed	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_EVT_HDIALY_SAFETY_VEN_ACCESS (198260)	tech	<b>Safety Systems</b>	Venous Access	Notification that the venous access may have become disconnected	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_EVT_HDIALY_SAFETY_VEN_AIR_DETECT (198262)	tech	<b>Safety Systems</b>	Venous Air Detector	Notification that air is sensed in the venous blood lines	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_EVT_HDIALY_SAFETY_WETNESS_DETECT_ALERT (198264)	tech	<b>Safety Systems</b>	Wetness Detector	Notification that the wetness detector has sensed fluid	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_EVT_HDIALY_SAFETY_WETNESS_DETECT_ERROR (198266)	tech	<b>Safety Systems</b>	Wetness Detector Error	Notification of a detected error on the wetness detector device	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_ATTR_ALERT_ID_NUM (68489)		<b>Safety Systems</b>	Alert ID	The number associated with the current alert. If no alert is active then this object is not reported.	Intradialytic		String	XXXX	N/A	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_ATTR_ALERT_TEXT (68546)		<b>Safety Systems</b>	Alert Text	The text name of the current alert. If no alert is active then this object is not reported. IN most case this object will contain the text shown to the user.	Intradialytic		String	N/A	N/A	0..1	O	X
<b>.. MDC_DEV_HDIALY_THERAPY_OUTCOMES_CHANNEL (70967)</b>										<b>1..1</b>	<b>M</b>	<b>M</b>
... MDC_HDIALY_THERAPY_MASS_TRF_AREA_COEFF (158972)		<b>Therapy Outcomes</b>	KoA	Overall mass transfer coefficient multiplied by surface area of a dialyzer	Intradialytic		Numeric	XXXX	N/A	0..1	O	X
... MDC_HDIALY_THERAPY_ACCESS_FLOW (158976)		<b>Therapy Outcomes</b>	Access Flow	Flow through a fistula or other conduit connecting an artery to a vein	Intradialytic		Numeric	±XXXX	mL/min	0..1	O	X
... MDC_HDIALY_THERAPY_RATIO_KT_OVER_V_GOAL (158992)		<b>Therapy Outcomes</b>	Kt/V Goal	Target dialysis treatment adequacy measurement	Intradialytic		Numeric	X.XX	%	0..1	O	C17
... MDC_HDIALY_THERAPY_UREA_CLEARANCE_MEAN (159019)		<b>Therapy Outcomes</b>	Mean Urea Clearance	Time weighted average of the individual effective conductivity clearance measurements	Intradialytic		Numeric	XXX	N/A	0..1	O	X
... MDC_HDIALY_THERAPY_UREA_CLEARANCE (159016)		<b>Therapy Outcomes</b>	Urea Clearance	Calculated clearance based on the change in conductivity of the pre-dialyzer vs post dialyzer dialysate.	Intradialytic		Numeric	XXX	N/A	0..1	O	X
... MDC_HDIALY_THERAPY_UREA_DISTRIBUTION_VOL_SETTING (0)		<b>Therapy Outcomes</b>	Urea Distribution Volume Setting	User setting for urea distribution volume used as V in the equation Kt/V.	All		Numeric	X.XX	L	0..1	O	O
... MDC_HDIALY_THERAPY_KT_DELIVERED (158984)		<b>Therapy Outcomes</b>	Delivered Kt	Volume of fluid cleared of urea by dialysis treatment	Intradialytic		Numeric	X.XX	L	0..1	O	X
... MDC_HDIALY_THERAPY_RATIO_EKT_OVER_V_DELIVERED (158980)		<b>Therapy Outcomes</b>	Delivered Equilibrated Kt/V	Fractional solute clearance for an intermittent treatment corrected to take account of solute redistribution following the treatment	Intradialytic		Numeric	X.XX	%	0..1	O	X
... MDC_HDIALY_THERAPY_RATIO_SPKT_OVER_V_DELIVERED (158988)		<b>Therapy Outcomes</b>	Delivered Single Pool Kt/V	Fractional solute clearance, calculated on the basis of solute concentration before and after an intermittent treatment, and on the basis of effluent dialysate volume and anthropometric estimate of volume of urea distribution for peritoneal dialysis	Intradialytic		Numeric	X.XX	%	0..1	O	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temp-oral	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_HDIALY_THERAPY_RATIO_SPKT_OVER_V_PROJECTED (159012)		<b>Therapy Outcomes</b>	Projected Single Pool Kt/V	Fractional solute clearance, calculated on the basis of solute concentration before and after an intermittent treatment, and on the basis of effluent dialysate volume and anthropometric estimate of volume of urea distribution for peritoneal dialysis, projected on the basis of data available before the treatment occurs or before the treatment is completed	Intradialytic		Numeric	X.XX	%	0..1	O	X
... MDC_EVT_HDIALY_THERAPY_PAT_TX_GENERAL (198268)	phys	<b>Therapy Outcomes</b>	General Patient Treatment	Notification that the machine has detected an alarm condition not covered by other alarms that is related to patient treatment.	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_THERAPY_PCT_RECIRC (159004)		<b>Therapy Outcomes</b>	Percent Recirculation	Percentage of the blood flow in the access line that is coming directly from the venous line without going through the patient's body.	Intradialytic		Numeric	±XXX.X	%	0..1	O	X
... MDC_HDIALY_THERAPY_PLASMA_NA_CONC (159008)		<b>Therapy Outcomes</b>	Plasma Sodium Concentration	Sodium ion concentration measured in patient plasma.	Intradialytic		Numeric	XXX	mmol/L	0..1	O	X
... MDC_HDIALY_THERAPY_COMPLETE_METHOD (158618)		<b>Therapy Outcomes</b>	Treatment Completion Method	Parameter used to determine when treatment has been complete	Intradialytic		Enum	<b>HD_TBL_12</b>	N/A	1..1	M	M
... MDC_EVT_HDIALY_THERAPY_TX_END_TIME (198270)	tech	<b>Therapy Outcomes</b>	Treatment Time Ended	Notification that target time of treatment has been met	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_THERAPY_BODY_START_WT (159000)		<b>Therapy Outcomes</b>	Patient Start Weight	Weight of patient at beginning of treatment	All		Numeric	XXX.X	Kg	0..1	O	X
... MDC_HDIALY_THERAPY_BODY_END_WT (159020)		<b>Therapy Outcomes</b>	Patient End Weight	Weight of patient at end of treatment	All		Numeric	XXX.X	Kg	0..1	O	X
... MDC_HDIALY_THERAPY_BODY_END_WT_SETTING (0)		<b>Therapy Outcomes</b>	Patient End Weight Target	Desired weight of patient at end of treatment. This is present as a means to specify UF when patient weight is changing. Requires patient's pre-treatment weight which could come from patient, scale, or EMR.	All		Numeric	XXX.X	Kg	0..1	O	O



## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use	
... MDC_HDIALY_THERAPY_END_TIME_SETTING (0)		<b>Therapy Outcomes</b>	Treatment End Time	The clock time that the therapy ends when completion method is Clock Time,	All		Episodic	Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	0..1	O	C14
... MDC_HDIALY_THERAPY_TX_TIME_SETTING (0)		<b>Therapy Outcomes</b>	Treatment Time	The duration of the therapy when completion method is Tx Time.	All		Episodic	Numeric	XXXX	minutes	0..1	O	C15
<b>.. MDC_DEV_HDIALY_UF_CHAN (70971)</b>										<b>1..1</b>	<b>M</b>	<b>M</b>	
... MDC_HDIALY_UF_MODE (158619)		<b>Ultrafiltration System</b>	Ultrafiltration Mode	Process used to remove fluid volume from the patient	Intradialytic		Episodic	Enum	<b>HD_TBL_13</b>	N/A	1..1	M	M
... Profile Parameters				UF Profile parameters							C19	C19	
... MDC_HDIALY_UF_RATE_SETTING (16936252)		<b>Ultrafiltration System</b>	Ultrafiltration Rate	Rate at which fluid is removed from the patient For systems that compute this value using Treatment Time and UF Target, this value would be reported as an Automatic Setting (ASET).	Intradialytic		Periodic	Numeric	XXXX	mL/h	1..1	M	C18
... MDC_HDIALY_UF_RATE (159036)		<b>Ultrafiltration System</b>	Ultrafiltration Rate	Current rate at which fluid is removed from the patient	Intradialytic		Periodic	Numeric	XXXX	mL/h	1..1	M	X
... MDC_EVT_HDIALY_UF_LO (198272)	tech-low	<b>Ultrafiltration System</b>	Insufficient Ultrafiltration	Notification that the ultrafiltration system is not removing as much fluid as intended	Intradialytic		Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_EVT_HDIALY_UF_NEG (198274)	tech-low	<b>Ultrafiltration System</b>	Negative Ultrafiltration	Notification that the ultrafiltration system has resulted in a negative ultrafiltration rate	Intradialytic		Episodic	Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_UF_TIME_TO_TARGET (159024)		<b>Ultrafiltration System</b>	Remaining Ultrafiltration Time	Time remaining until target volume to remove will have been achieved.	Intradialytic		Periodic	Numeric	XXXX	minutes	0..1	O	X
... MDC_HDIALY_UF_TARGET_VOL_TO_REMOVE (159028)		<b>Ultrafiltration System</b>	Target Volume To Remove	Planned fluid removal	Intradialytic		Episodic	Numeric	XXXX	mL	1..1	C11	C11
... MDC_HDIALY_UF_ACTUAL_REMOVED_VOL (159032)		<b>Ultrafiltration System</b>	Total Fluid Volume Removed	Measured fluid removed from patient	Intradialytic		Periodic	Numeric	XXXX	mL	1..1	M	X

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
... MDC_EVT_HDIALY_UF_RATE_RANGE (198276)	tech-range	<b>Ultrafiltration System</b>	Ultrafiltration Rate Out of Range	Notification that the ultrafiltration rate has deviated from the intended ultrafiltration rate more than is safe. Specified another way, MDC_HDIALY_UF_RATE is not between MDC_HDIALY_UF_RATE_LIMIT_LOW_SETTING and MDC_HDIALY_UF_RATE_LIMIT_HIGH_SETTING	Intradialytic		Bool	T / F	N/A	0..1	M	X
... MDC_EVT_HDIALY_UF_GOAL_MET (198278)	tech	<b>Ultrafiltration System</b>	Ultrafiltration Goal Met	Notification that target ultrafiltration volume has been met	Intradialytic		Bool	T / F	N/A	0..1	O	X
... MDC_HDIALY_UF_RATE_LIMIT_HIGH_SETTING (16936257)		<b>Ultrafiltration System</b>	UF Rate Upper Limit	Maximum ultrafiltration rate allowed without an alarm generation.	Intradialytic		Numeric	XXXX	mL/h	0..1	O	X
... MDC_HDIALY_UF_RATE_LIMIT_LOW_SETTING (16936259)		<b>Ultrafiltration System</b>	UF Rate Lower Limit	Maximum ultrafiltration rate allowed without an alarm generation.	Intradialytic		Numeric	XXXX	mL/h	0..1	O	X
. MDC_DEV_PRESS_BLD_NONINV_VMD (70686)				<b>Non-Invasive Blood Pressure Monitor</b>						0..1	O	O
.. MDC_DEV_PRESS_BLD_NONINV_CHAN (70687)										0..1	O	O
... MDC_ATTR_TIME_PD_MSMT (67979)		<b>Blood Pressure Monitor</b>	Blood Pressure Measurement Interval Setting	Time allowed to elapse between automatic blood pressure measurements. A value of zero indicates that automatic measurements are not being taken.	All		Numeric	XXX	minutes	0..1	O	X
... MDC_ATTR_PT_BODY_POSN (68135)		<b>Blood Pressure Monitor</b>	Blood Pressure Measurement Position	Patient position at time of blood pressure measurement	All		Enum	HD_TBL_04	N/A	0..1	O	X
... MDC_PRESS_BLD_NONINV_DIA (150022)	phys high low thr	<b>Blood Pressure Monitor</b>	Diastolic Pressure	Minimum arterial pressure during relaxation and dilatation of the ventricles of the heart when the ventricles fill with blood	All		Numeric	XXX	mmHg	0..1	O	C13
... MDC_PULS_RATE_NONINV (149546)	phys high low thr	<b>Blood Pressure Monitor</b>	Heart Rate	Number of heart beats per minute	All		Numeric	XXX	beats/min	0..1	O	C13
... MDC_PRESS_BLD_NONINV_MEAN (150023)		<b>Blood Pressure Monitor</b>	Mean Arterial Pressure	Average pressure in a patient's arteries during one cardiac cycle	All		Numeric	XXX	mmHg	0..1	O	C13
... MDC_PRESS_BLD_NONINV_SYS (150021)	phys high low thr	<b>Blood Pressure Monitor</b>	Systolic Pressure	Maximum arterial pressure during contraction of the left ventricle of the heart	All		Numeric	XXX	mmHg	0..1	O	C13

## Dialysis Machine Implementation Guide

Table 2 – Dialysis Machine Data Objects

REFID	Alert Type	Dialysis Device VMD	Dialysis Common Name	Definition	Phase	Temporal	Data Type	Format	UOM	CARD	USE	Rx Use
. MDC_DEV_ANALY_SAT_O2_VMD (69642)										0..1	O	O
.. MDC_DEV_ANALY_SAT_O2_CHAN (69643)										0..1	O	O
... MDC_PULS_OXIM_SAT_O2 (150456)	phys high low thr	<b>Pulse Oximeter Monitor</b>	SpO2 Oxygen Saturation	Hemoglobin oxygen saturation in patient blood measured by pulse oximetry	Intradialytic	Periodic	Numeric	XXX	%	0..1	O	C13
... MDC_PULS_OXIM_PULS_RATE (149530)	phys high low thr	<b>Pulse Oximeter Monitor</b>	SpO2 Pulse Rate	Number of heart beats per minute measured by pulse oximetry	Intradialytic	Periodic	Numeric	XXX	beats/min	0..1	O	C13
... MDC_EVT_ERR (196638)	tech	<b>Safety Systems</b>	Pulse Oximeter Error	Notification of a detected error on the pulse oximeter device	Intradialytic	Episodic	Bool	T / F	N/A	0..1	O	X
. MDC_DEV_BLOOD_CHEM_VMD (70974)										0..1	O	X
.. MDC_DEV_BLOOD_CHEM_CHAN (70975)				<b>Blood Chemistry Channel</b>						0..1	O	X
... MDC_HDIALY_PLASMA_VOL_MARKER (158620)		<b>Therapy Outcomes</b>	Plasma Volume Marker	An indicator applied to the blood volume percent reading to designate an event	Intradialytic	Periodic	Enum	HD_TBL_10	N/A	0..1	O	X
... MDC_HDIALY_PLASMA_VOL_PROFILE (158621)		<b>Therapy Outcomes</b>	Plasma Volume Profile	The rate of change in relative blood volume.	Intradialytic	Periodic	Enum	HD_TBL_11	N/A	0..1	O	X
... MDC_CONC_HCT_GEN (160132)	phys low high thr	<b>Therapy Outcomes</b>	Hematocrit	100*(sample red blood cell volume/blood sample volume)	Intradialytic	Periodic	Numeric	XX.X	%{vol}	0..1	O	X
... MDC_CONC_HB_GEN (160120)	phys low high thr	<b>Therapy Outcomes</b>	Hemoglobin	Concentration of hemoglobin in a sample expressed as mass/volume	Intradialytic	Periodic	Numeric	±XXX.X	g/dL	0..1	O	X
... MDC_HDIALY_REL_BLOOD_VOL (159044)	phys low thr	<b>Therapy Outcomes</b>	Relative Blood Volume	The percent change in blood volume since the start of the treatment. The BV change value is calculated from the following equation: $\Delta BV\% = [(H0/H1) - 1] \times 100$ Where H0 = initial Hct or Hb H1 = current Hct or Hb	Intradialytic	Episodic	Numeric	XXX.X	%	0..1	O	X
... MDC_SAT_O2 (150316)	phys low high thr	<b>Therapy Outcomes</b>	SpO2 Oxygen Saturation	Hemoglobin oxygen saturation in patient blood	Intradialytic	Periodic	Numeric	XXX	%	0..1	O	X

## 9.3 Events

Table 3 – Hemodialysis Machine Alarms/Alerts

Source Identifier	Event Identifier <sup>1</sup>	Alert Type	Dialysis Common Name	Definition	Use
MDC_DEV_HDIALY_MACHINE_MDS					
MDC_DEV_HDIALY_VMD					
MDC_DEV_HDIALY_MACH_CONFIG_CHAN					
MDC_DEV_HDIALY_ANTICOAG_PUMP_CHAN					
MDC_HDIALY_ANTICOAG_PUMP_CHAN	MDC_EVT_HDIALY_ANTICOAG_STOP	tech	Anticoagulant Pump Stop	Notification that the anticoagulant pump has stopped	O
MDC_HDIALY_ANTICOAG_PUMP_CHAN	MDC_EVT_HDIALY_ANTICOAG_SYRINGE_EMPTY	tech	Syringe Empty	Notification that the syringe plunger has reached the end of its functional range	O
MDC_HDIALY_ANTICOAG_PUMP_CHAN	MDC_EVT_HDIALY_ANTICOAG_SYRINGE_SIZE	tech	Wrong Syringe Size	Notification that the size of the syringe is not the size defined for anticoagulant administration	O
MDC_DEV_HDIALY_BLOOD_PUMP_CHAN					
MDC_HDIALY_BLD_PRESS_ART	MDC_EVT_HI	tech high thr	Arterial Pressure High	Notification that the arterial pressure has exceeded maximum setting	M
MDC_HDIALY_BLD_PRESS_ART	MDC_EVT_LO	tech low thr	Arterial Pressure Low	Notification that the arterial pressure has fallen below the minimum setting	M
MDC_HDIALY_BLOOD_PUMP_CHAN	MDC_EVT_HDIALY_BLD_PUMP_STOP	tech	Blood Pump Stop	Notification that the blood pump has stopped	M
MDC_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP	MDC_EVT_HI	tech high thr	Post Pump Arterial Pressure High	Notification that the post-pump arterial pressure is above limit.	O
MDC_HDIALY_BLD_PUMP_PRESS_ART_POST_PUMP	MDC_EVT_LO	tech low thr	Post Pump Arterial Pressure Low	Notification that the post-pump arterial pressure is below limit.	O
MDC_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS	MDC_EVT_HI	tech high thr	Needle Pressure High	Notification that the needle pressure is above limit.	O
MDC_HDIALY_BLD_PUMP_SINGLE_NEEDLE_PRESS	MDC_EVT_LO	tech low thr	Needle Pressure Low	Notification that the needle pressure is below limit.	O
MDC_HDIALY_BLD_PUMP_PRESS_VEN	MDC_EVT_HI	tech high thr	Venous Pressure High	Notification that the venous pressure has exceeded maximum setting	M
MDC_HDIALY_BLD_PUMP_PRESS_VEN	MDC_EVT_LO	tech low thr	Venous Pressure Low	Notification that the venous pressure has fallen below the minimum setting	M
MDC_DEV_HDIALY_FLUID_CHAN					
MDC_HDIALY_DIALYSATE_COND	MDC_EVT_HI	tech high	Conductivity High	Notification that the conductivity of the dialysate has exceeded the maximum setting	O
MDC_HDIALY_DIALYSATE_COND	MDC_EVT_LO	tech low	Conductivity Low	Notification that the conductivity of the dialysate has fallen below the minimum setting	O
MDC_HDIALY_DIALYSATE_FLOW_RATE	MDC_EVT_LO	tech low	Dialysate Flow Low	Notification that the dialysate flow is below the dialysate flow rate minimum setting	O
MDC_HDIALY_DIALYSATE_TEMP	MDC_EVT_LO	tech low	Dialysate Temperature Low	Notification that the dialysate temperate has fallen below the minimum temperate setting	C (if mode is HD or HDF)
MDC_HDIALY_DIALYSATE_TEMP	MDC_EVT_HI	tech high	Dialysate Temperature High	Notification that the dialysate temperate exceeds the maximum temperate setting	C (if mode is HD or HDF)
MDC_HDIALY_FLUID_CHAN	MDC_EVT_HDIALY_BLOOD_LEAK	tech	Blood Leak	Notification that blood has been detected in the dialysate fluid	M
MDC_DEV_HDIALY_FILTER_CHAN					
MDC_HDIALY_FILTER_TRANSMEMBRANE_PRESS	MDC_EVT_LO	tech low	Transmembrane Pressure Low	Notification that transmembrane pressure has fallen below the minimum setting	M

## Dialysis Machine Implementation Guide

**Table 3 – Hemodialysis Machine Alarms/Alerts**

Source Identifier	Event Identifier <sup>1</sup>	Alert Type	Dialysis Common Name	Definition	Use
MDC_HDIALY_FILTER_TRANSMEMBRANE_PRESS	MDC_EVT_HI	tech high	Transmembrane Pressure High	Notification that transmembrane pressure has exceeded the maximum setting	M
<b>MDC_DEV_HDIALY_CONVECTIVE_CHAN</b>					
MDC_HDIALY_CONVECTIVE_CHAN	MDC_EVT_HDIALY_RF_EXCESS_DELIV	tech-high	Excessive Replacement Fluid Bolus Delivered	Notification that more replacement fluid has been delivered than intended	O
MDC_HDIALY_CONVECTIVE_CHAN	MDC_EVT_HDIALY_REPLACE_FLUID_INSUFF_DELIV	tech-low	Insufficient Replacement Fluid Bolus Delivered	Notification that less replacement fluid has been delivered than intended	O
MDC_HDIALY_RF_BOLUS_RATE	MDC_EVT_HI	tech high	Replacement Fluid Bolus Rate High	Notification that the replacement fluid bolus is delivering more fluid than intended	O
MDC_HDIALY_RF_BOLUS_RATE	MDC_EVT_LO	tech low	Replacement Fluid Bolus Rate Low	Notification that the replacement fluid bolus is delivering less than the intended.	O
MDC_HDIALY_RF_BOLUS_VOL_DELIVERED	MDC_EVT_HI	tech high	Replacement Fluid Bolus Volume High	Notification that the replacement fluid bolus delivered more fluid than intended	O
MDC_HDIALY_RF_BOLUS_VOL_DELIVERED	MDC_EVT_LO	tech low	Replacement Fluid Bolus Volume Low	Notification that the replacement fluid bolus delivered less fluid than intended.	O
<b>MDC_DEV_HDIALY_SAFETY_SYSTEMS_CHAN</b>					
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_ART_AIR_DETECT	tech	Arterial Air Detector	Notification that air is sensed in the arterial blood lines	M
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_DIALYSATE_AIR_DETECT	tech	Dialysate Air Detector	Notification that air is sensed in dialysate line	C (if mode is HD or HDF)
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_DIALYSATE_COMPOSITION	tech	Dialysis Composition	Notification that the dialyzing fluid composition is incorrect	O
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_SYSTEM_GENERAL	tech	General System	Notification of any dialysis machine malfunction that are not related to patient treatment	M
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_SELFTEST_FAILURE	tech	Self-Test	Notification that a self-test failed	M
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_VEN_ACCESS	tech	Venous Access	Notification that the venous access may have become disconnected	O
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_VEN_AIR_DETECT	tech	Venous Air Detector	Notification that air is sensed in the venous blood lines	M
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_WETNESS_DETECT_ALERT	tech	Wetness Detector	Notification that the wetness detector has sensed fluid	O
MDC_HDIALY_SAFETY_SYSTEMS_CHAN	MDC_EVT_HDIALY_SAFETY_WETNESS_DETECT_ERROR	tech	Wetness Detector Error	Notification of a detected error on the wetness detector device	O
<b>MDC_DEV_HDIALY_THERAPY_OUTCOMES_CHAN</b>					
MDC_HDIALY_THERAPY_OUTCOMES_CHAN	MDC_EVT_HDIALY_THERAPY_PAT_TREATMENT_GENERAL	phys	General Patient Treatment	Notification that the machine has detected an alarm condition not covered by other alarms that is related to patient treatment.	O
MDC_HDIALY_THERAPY_OUTCOMES_CHAN	MDC_EVT_HDIALY_THERAPY_TREATMENT_TIME_ENDED	tech	Treatment Time Ended	Notification that target time of treatment has been met	O
<b>MDC_DEV_HDIALY_UF_CHAN</b>					
MDC_HDIALY_UF_CHAN	MDC_EVT_HDIALY_UF_LO	tech low	Insufficient Ultrafiltration	Notification that the ultrafiltration system is not removing as much fluid as intended	O
MDC_HDIALY_UF_CHAN	MDC_EVT_HDIALY_UF_NEG	tech low	Negative Ultrafiltration	Notification that the ultrafiltration system has resulted in a negative ultrafiltration rate	O
MDC_HDIALY_UF_CHAN	MDC_EVT_HDIALY_UF_GOAL_MET	tech	Ultrafiltration Goal Met	Notification that target ultrafiltration volume has been met	O

## Dialysis Machine Implementation Guide

**Table 3 – Hemodialysis Machine Alarms/Alerts**

Source Identifier	Event Identifier <sup>1</sup>	Alert Type	Dialysis Common Name	Definition	Use
MDC_HDIALY_UF_CHAN	MDC_EVT_HDIALY_UF_RATE_RANGE	tech range	Ultrafiltration Rate Out of Range	Notification that the ultrafiltration rate has deviated from the intended ultrafiltration rate more than is safe	M
MDC_DEV_PRESS_BLD_NONINV_VMD					
MDC_DEV_PRESS_BLD_NONINV_CHAN					
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_HI	phys high thr	Diastolic High	Notification that diastolic pressure exceeded the diastolic maximum setting	O
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_LOW	phys low thr	Diastolic Low	Notification that diastolic pressure fell below the diastolic minimum setting	O
MDC_PULS_RATE_NON_INV	MDC_EVT_HI	phys high thr	Heart Rate High	Heart rate exceeded heart rate maximum setting	O
MDC_PULS_RATE_NON_INV	MDC_EVT_LOW	phys low thr	Heart Rate Low	Heart rate fell below heart rate minimum setting	O
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_HI	phys high thr	Systolic High	Notification that systolic pressure exceeded the systolic maximum setting	O
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_LOW	phys low thr	Systolic Low	Notification that systolic pressure fell below the systolic minimum setting	O
MDC_DEV_ANALY_SAT_O2_VMD					
MDC_DEV_ANALY_SAT_O2_CHAN					
MDC_PULS_OXIM_SAT_O2	MDC_EVT_HI	phys high thr	Oximeter Oxygen Saturation High	Notification that Oxygen Saturation is above limit.	O
MDC_PULS_OXIM_SAT_O2	MDC_EVT_LO	phys low thr	Oximeter Oxygen Saturation Low	Notification that Oxygen Saturation is below limit.	
MDC_PULS_OXIM_PULS_RATE	MDC_EVT_HI	phys high thr	Oximeter Pulse High	Notification that the pulse rate measure by the pulse oximeter is above limit.	O
MDC_PULS_OXIM_PULS_RATE	MDC_EVT_LO	phys low thr	Oximeter Pulse Low	Notification that the pulse rate measure by the pulse oximeter is below limit.	
MDC_DEV_ANALY_SAT_O2_CHAN	MDC_EVT_ERR	tech	Pulse Oximeter Error	Notification of a detected error on the pulse oximeter device	O
MDC_DEV_BLOOD_CHEM_VMD					
MDC_DEV_BLOOD_CHEM_CHAN					
MDC_CONC_HCT	MDC_EVT_HI	phys high thr	Hematocrit High	Notification that patient hematocrit has risen above limit.	O
MDC_CONC_HCT	MDC_EVT_LO	phys low thr	Hematocrit Low	Notification that patient hematocrit has fallen below limit.	O
MDC_CONC_HB	MDC_EVT_HI	phys high thr	Hemoglobin High	Notification that patient hemoglobin has risen above limit.	O
MDC_CONC_HB	MDC_EVT_LO	phys low thr	Hemoglobin Low	Notification that patient hemoglobin has fallen below limit.	O
MDC_HDIALY_REL_BLOOD_VOL	MDC_EVT_LO	phys low thr	Blood Volume	Notification that the relative blood volume has dropped below the defined setting	O
MDC_SAT_O2	MDC_EVT_HI	phys high thr	Oxygen Saturation High	Notification that Patient Oxygen Saturation has risen above limit.	O
MDC_SAT_O2	MDC_EVT_LO	phys low thr	Oxygen Saturation Low	Notification that Patient Oxygen Saturation has fallen below limit.	O

### Notes

- 1) In the case where the Event Identified is MDC\_EVT\_LO or MDC\_EVT\_HI, the dialysis machine will report the value of the parameter that is out of bounds in OBX-2 in the second observation. The parameter that is reported is shown in the Source Identifier column. The limit/range that is violated is reported in OBX-7 in the second observation.

## 10 Peritoneal Dialysis Data Elements

### 10.1 Tables

**PD\_TBL\_01 – Treatment Type**

<b>Value</b>	<b>Description</b>
CAPD	Continuous Ambulatory Peritoneal Dialysis
APD	Automated Peritoneal Dialysis
CCPD	Continuous Cycling Peritoneal Dialysis
NIPD	Nocturnal Intermittent Peritoneal Dialysis
PDplus	PD Plus
TPD	Tidal Peritoneal Dialysis
aAPD	Adapted Automated Peritoneal Dialysis
IPD	Intermittent Peritoneal Dialysis
Other	A treatment that does not fit into any of the other defined types.

**PD\_TBL\_02 – Treatment Location**

<b>Value</b>	<b>Description</b>
Home	Home
Clinic	Clinic
Other	Other

**PD\_TBL\_03 – Treatment Phase**

<b>Value</b>	<b>Description</b>
PENDING	Pending
FILL	Fill
DWELL	Dwell
DRAIN	Drain
COMPLETE	Complete

**PD\_TBL\_04 – Osmotic Agent**

<b>Value</b>	<b>Description</b>
GLU	Glucose
DEX	Dextrose
ICO	Icodextrin
OTHER	Other

**PD\_TBL\_05 – Fluid Source**

<b>Value</b>	<b>Description</b>
BAG	Bags
ONLINE	Online
UNKNOWN	Unknown

**PD\_TBL\_06 – Fill Mode**

<b>Value</b>	<b>Description</b>
AUTO	Cycler
MANUAL	Manual

**PD\_TBL\_07 – Drain Mode**

<b>Value</b>	<b>Description</b>
TOBAG	Cycler to Bag
TODRAIN	Cycler to Drain
UNKNOWN	Manual

**PD\_TBL\_08 – Transport Type**

<b>Value</b>	<b>Description</b>
LOW	Low
MEDLOW	Medium-Low
MEDHIGH	Medium-High
HIGH	High

**PD\_TBL\_09 – Catheter Type**

<b>Value</b>	<b>Description</b>
STRAIGHT-TENCK	Straight Tenckhoff
OREO-ZELL	Oreopoulos-Zellerman
TWH	Toronto Western Hospital
COILED-TENCK	Coiled Tenckhoff
STRAIGHT-SWAN-TENCK	Straight Swan-neck Tenckhoff
COILED-SWAN-TENCK	Coiled Swan-neck Tenckhoff
STRAIGHT-MO-SWAN	Straight Missouri Swan-neck
COLIED-MO-SWAN	Coiled Missouri Swan-neck
MON-POP	Moncrief-Popovich
VINCENZA	Vincenza
DO-PAOLA	Do Paolo
VALLI	Valli (Balloon)
ASH-ADV	Ash Advantage (T-Fluted)
CRUZ	Cruz (Pail-Handle)
OTHER	Other/Unknown



**PD\_TBL\_10 – Measurement Phase**

<b>Value</b>	<b>Description</b>
NONE	A measurement is not taken.
PRE	A measurement is taken prior to the start of the treatment.
POST	A measurement is taken after the treatment is complete.
BOTH	A measurement is taken prior to the start of treatment and again when the treatment is complete.

**PD\_TBL\_11 – Air Sensor State**

<b>Value</b>	<b>Description</b>
AIR	The sensor is reading air.
FLUID	The sensor is reading fluid.
UNKNOWN	The state of the sensor is unknown.


**PD\_TBL\_12 – Last Drain Mode**

<b>Value</b>	<b>Description</b>
NEVER	It will never take place
ALWAYS	It will always take place
UFLIMIT	It will take place if the total estimated UF for the therapy is below the Extra Drain UF Limit. The Extra Drain UF Limit is expressed as a percentage of total estimated UF for the therapy.

**PD\_TBL\_13 – Prescription Type**

<b>Value</b>	<b>Description</b>
EXCH	Exchange Based
TXDUR	Treatment Based with Fixed Duration
TXEND	Treatment Based with Fixed End Time

## 10.2 Data Object



### Fluid Channels

A PD PCD-01 Message (Prescription or Treatment Status) will contain one Fluid Channel for each fluid source. A fluid source will be a single bag or online delivery system. The standard supports the use of up to 6 fluid sources. The machine reports each source in a separate channel. These channels are numbered from 1 to 6. The channels' name is abbreviated as MDCX\_DEV\_PDIALY\_FLUID\_x\_CHAN. In the PCD-01 message, "x" is replaced with the source number. For example, a treatment that has 3 bags will report fluid channels MDCX\_DEV\_PDIALY\_FLUID\_1\_CHAN, MDCX\_DEV\_PDIALY\_FLUID\_2\_CHAN, and MDCX\_DEV\_PDIALY\_FLUID\_3\_CHAN.

Exchanges

A PD exchange uses data from one or more fluid sources. The cyclor reports how much volume from each source is used in the MDCX\_PDIALY\_FILL\_VOLUME object. This object is a numeric array. The number of items in the array will be the same as the number of fluid channels. If a particular source is not used during the exchange, then the volume will be 0. For example, given a treatment with three sources, the cyclor would report [0.000^0.000^2.500] for an exchange of 2.5 L from the last bag.

Supply Lines

The supply lines are part of the Disposable Set. The Disposable Set varies from cyclor to cyclor. The supply lines are physically connected to the fluid sources and convey fluid from the bags to the cyclor. The Supply Line Channel represents all lines that are connected to a Fluid Source. If the cyclor can associate an event (low flow, no flow, etc.) with one of the Fluid Sources, then the cyclor will link the event to the appropriate Fluid Source Channel. If the cyclor cannot relate the event to a specific fluid source, the cyclor will link the event to the generic Supply Line Channel.

In the table below when usage is conditional (C), a number is given. The number refers to one of the following conditions. If the condition is not met then the cardinality is "0..0".

- 1) The Exchange Phase is Fill or greater (MDCX\_PDIALY\_CURRENT\_PHASE = FILL or DWELL or DRAIN or COMPLETE).
- 2) The Exchange Phase is Dwell or greater (MDCX\_PDIALY\_CURRENT\_PHASE = DWELL or DRAIN or COMPLETE).
- 3) The Exchange Phase is Drain or greater (MDCX\_PDIALY\_CURRENT\_PHASE = DRAIN or COMPLETE).
- 4) The Fluid Source is bags (MDCX\_PDIALY\_FLUID\_SOURCE = BAG).
- 5) There are APD exchanges (MDCX\_PDIALY\_TREAT\_TYPE = APD).
- 6) There are CAPD exchanges (MDCX\_PDIALY\_TREAT\_TYPE = CAPD).
- 7) The prescription type is exchange based (MDCX\_PDIALY\_PRESCRIPTION\_TYPE = EXCH).

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
MDCX_DEV_PDIALY_MACHINE_MDS		Hemodialysis machine	A device or system used to treat renal failure delivered by removing toxins from the blood through an artificial kidney, called a dialyzer.	All	All				1..1	M	M
... MDC_ID_MODEL_MANUFACTURER (531950)		Dialysis Device Manufacturer	String containing device manufacturer's name	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ID_MODEL_NUMBER (531969)		Dialysis Device Model	String containing device manufacturer's model identifier for the device	All	All	String	Alphanumeric	N/A	1..1	M	X

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDC_ID_PROD_SPEC_SERIAL		Dialysis Device Serial Number	String containing the device's serial number	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ID_PROD_SPEC_SW		Device Software Version	String defining the version of the software on the machine	All	All	String	Alphanumeric	N/A	1..1	M	X
... MDC_ATTR_ID_UDI		Dialysis Device UDI	String containing the device's FDA required Universal Device Identifier. [Note: MDC_ATTR_ID_UDI contains three components: UdiAuthority, UdiIssuer, UdiLabel to support non-FDA entities, but currently there is no approved way of encoding this using HL7 V2; PRT-10 and PRT-16-20 is used instead; UdiLabel = PRT-10.]	All	All	String	Alphanumeric	N/A	1..1	O	X
. MDCX_DEV_PDIALY_VMD		Hemodialysis subsystem	A device or subsystem used to treat renal failure delivered by removing toxins from the blood through an artificial kidney, called a dialyzer.	All	All				1..1	M	M
.. MDCX_DEV_PDIALY_MACH_CONFIG_CHAN		Hemodialysis configuration	The principal dialysis machine operational modes, treatment modality, settings and other information about the dialysis machine.	All	All				1..1	M	X
... MDC_ATTR_TIME_ABS (67975)		Dialysis Device Time	Date and time as recorded on the dialysis device's internal clock with offset to UTC (based on location and DST)	All	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	M	X
.. MDCX_DEV_PDIALY_TREATMENT_CHANNEL		<b>Treatment</b>		<b>Intradialytic</b>	<b>All</b>				1..1	M	M
... MDCX_PDIALY_TREAT_TYPE		Treatment Type	Treatment Type	Intradialytic	All	Enumeration	PD_TBL_01	N/A	1..1	M	M
... MDCX_PDIALY_TREAT_LOCATION		Treatment Location	Treatment Location	Intradialytic	All	Enumeration	PD_TBL_02	N/A	0..1	O	O
... MDCX_PDIALY_CURRENT_PHASE		Current Phase	The current treatment phase.	Intradialytic	All	Enumeration	PD_TBL_03	N/A	1..1	M	X
... MDCX_PDIALY_PRESCRIPTION_TYPE		Prescription Type	Prescription Type	Intradialytic	All	Enumeration	PD_TBL_13	N/A	1..1	M	M
... MDCX_PDIALY_APD_START_TIME_SETTING		APD Start Time	APD Start Time	Intradialytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C5	C5
... MDCX_PDIALY_APD_END_TIME_SETTING		APD End Time	APD End Time	Intradialytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C5	C5
... MDCX_PDIALY_APD_TOTAL_VOLUME_SETTING		APD Total Volume	Total volume in APD exchanges.	Intradialytic	All	Numeric	X.XXX	L	1..1	C5	C5
... MDCX_PDIALY_APD_CYCLES_SETTING		APD Cycles	APD Cycles	Intradialytic	All	Numeric	XX	N/A	1..1	C5	C5

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_PDIALY_CAPD_START_TIME_SETTING		CAPD Start Time	CAPD Start Time	Intradialytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C6	C6
... MDCX_PDIALY_CAPD_END_TIME_SETTING		CAPD End Time	CAPD End Time	Intradialytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C6	C6
... MDCX_PDIALY_CAPD_TOTAL_VOLUME_SETTING		CAPD Exchanges	CAPD Exchanges	Intradialytic	All	Numeric		N/A	1..1	C6	C6
... MDCX_PDIALY_CAPD_CYCLES_SETTING		CAPD Total Volume	Total volume in CAPD exchanges	Intradialytic	All	Numeric	XX	N/A	1..1	C6	C6
... MDCX_PDIALY_TIDAL_PERCENTAGE		Tidal Percentage	The percentage of Night Fill Volume to be drained during a tidal cycle.	Intradialytic	All	Numeric	XXX.X	%	0..1	O	O
... MDCX_PDIALY_FULL_DRAIN_FREQ		Full Drain Frequency	The number of cycles between full drains when performing tidal exchanges. The value is expressed as the number of drains between full drains.	Intradialytic	All	Numeric	X	N/A	0..1	O	O
... MDCX_PDIALY_PERITONEUM_VOL_MAX		Maximum Peritoneum Volume		Intradialytic	All	Numeric	X.XXX	L	0..1	O	O
... MDCX_PDIALY_INIT_MIN_DRAIN_VOL_SETTING		Minimum Initial Drain Volume	The minimum volume to remove during the initial drain.	Intradialytic	All	Numeric	X.XXX	L	0..1	O	O
... MDCX_PDIALY_INIT_MAX_DRAIN_VOL_SETTING		Maximum Initial Drain Time	The maximum volume to remove during the initial drain.	Intradialytic	All	Numeric	X.XXX	L	0..1	O	O
... MDCX_PDIALY_SAMPLE_REMINDER		Sample Reminder	This value is true when the user should be reminded to collect an effluent sample.	Intradialytic	All	Bool	T / F	N/A	0..1	O	O
... MDCX_PDIALY_MEDICATION_SURVEY		Medications Survey	This value is true when the user should be surveyed for medications taken.	Intradialytic	All	Bool	T / F	N/A	0..1	O	O
... MDCX_PDIALY_EXCHANGE_SURVEY		Exchange Survey	This value is true when the user should be surveyed for number of day exchanges performed.	Intradialytic	All	Bool	T / F	N/A	0..1	O	O
... MDCX_PDIALY_DAY_DRAIN_TIME_MIN		Minimum Day Drain Time	Minimum duration for a day time drain.	Intradialytic	All	Numeric	XXX.X	min	0..1	O	O
... MDCX_PDIALY_DAY_DRAIN_VOL_PCT_MIN		Minimum Day Drain Volume Percent	The min amount of fluid to drain expressed as a percentage of the fill volume.	Intradialytic	All	Numeric	XXX.X	%	0..1	O	O
... MDCX_PDIALY_NIGHT_DRAIN_TIME_MIN		Minimum Night Drain Time	Minimum duration for a night time drain.	Intradialytic	All	Numeric	XXX.X	min	0..1	O	O
... MDCX_PDIALY_LAST_DRAIN_MODE		Extra Last Drain Mode	This value specifies how the last drain phase is performed.	Intradialytic	All	Enumeration	<b>PD_TBL_12</b>	N/A	0..1	O	O

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_PDIALY_LAST_DRAIN_UF_LIMIT		Extra Last Drain UF Limit	This limits which determines if the extra drain is performed. If the total estimated UF for the therapy is below the Extra Drain UF Limit the extra drain is performed. The Extra Drain UF Limit is expressed as a percentage of total estimated UF for the therapy.	Intradialytic	All	Numeric	XXX.X	%	0..1	O	O
... MDCX_PDIALY_NIGHT_UF_ESTIMATE		Estimated Night UF	The estimated UF volume for the night therapy. This value is used to keep track of the fluid balance for the patient, so that the fluid that was filled plus the fluid generated as UF is drained from the patient.	Intradialytic	All	Numeric	X.XXX	L	0..1	O	O
... MDCX_PDIALY_TEMP_DIALYSATE_SETTING		Dialysate Temperature Setting	Desired dialysate temperature during treatment	Intradialytic	All	Numeric	XX.X	°C	0..1	O	X
.. MDCX_DEV_PDIALY_FLUID_1_CHAN .. MDCX_DEV_PDIALY_FLUID_2_CHAN .. MDCX_DEV_PDIALY_FLUID_3_CHAN .. MDCX_DEV_PDIALY_FLUID_4_CHAN .. MDCX_DEV_PDIALY_FLUID_5_CHAN .. MDCX_DEV_PDIALY_FLUID_6_CHAN		<b>Fluid Definition</b>		<b>Intradialytic</b>	<b>All</b>				<b>1..1</b>	<b>M</b>	<b>M</b>
... MDCX_PDIALY_FLUID_NAME		Dialysate Name	Dialysate Name	Intradialytic	All	String	Alphanumeric	N/A	0..1	O	O
... MDCX_DIALY_FLUID_GLUCOSE_CONC		Glucose	Glucose	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_ICODEXTRIM_CONC		Icodextrin	Icodextrin	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_AMINO_ACID_CONC		Amino acid	Amino acid	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_OSMOTIC_AGENT		Other osmotic agent type	Other osmotic agent type	Intradialytic	All	Enumeration	<b>PD_TBL_04</b>	N/A	0..1	O	O
... MDCX_DIALY_FLUID_OSMOTIC_AGENT_CONC		Other osmotic agent concentration	Other osmotic agent concentration	Intradialytic	All	Numeric	XX.XX	%	0..1	O	O
... MDCX_DIALY_FLUID_BICARB_CONC		Bicarbonate (HCO3)	Bicarbonate (HCO3)	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_CHLORIDE_CONC		Chloride (CL)	Chloride (CL)	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_LACTATE_CONC		Lactate	Lactate	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_CALCIIUM_CONC		Calcium (CA)	Calcium (CA)	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_CITRATE_CONC		Citrate	Citrate	Intradialytic	All	Numeric	XXX.X	mmol /L	0..1	O	O

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_DIALY_FLUID_POTASSIUM_CONC		Potassium (K)	Potassium (K)	Intradial ytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_DIALY_FLUID_MAGNESIUM_CONC		Magnesium	Magnesium	Intradial ytic	All	Numeric	XXX.X	mmol /L	0..1	O	O
... MDCX_PDIALY_FLUID_BAG_VOLUME		Bag Volume	Bag Volume	Intradial ytic	All	Numeric	X.XXX	L	1..1	C4	C4
... MDCX_PDIALY_FLUID_SOURCE		Fluid Source	Source of the fluid.	Intradial ytic	All	Enumeration	<b>PD_TBL_05</b>	N/A	1..1	M	M
... MDCX_DEV_PDIALY_EXCHANGE_CHAN		Peritoneal Exchange	The values for the current and previous exchanges.	Intradial ytic	All				<b>1..1</b>	<b>M</b>	<b>C7</b>
... MDC_ATTR_CHAN_NUM_LOGICAL		Exchange Number	The number of the exchange in the current treatment. The first exchange is number 1.	Intradial ytic	All	Numeric	X	N/A	1..1	M	C7
... MDCX_PDIALY_CURRENT_PHASE		Current Phase	The current exchange phase.	Intradial ytic	All	Enumeration	<b>PD_TBL_03</b>	N/A	1..1	M	X
... MDCX_PDIALY_FILL_START_TIME		Fill Start time	The time that the fill phase started.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C1	X
... MDCX_PDIALY_FILL_END_TIME		Fill End time	The time that the fill phase ended.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C1	X
... MDCX_PDIALY_FILL_MODE		Fill mode	Fill mode	Intradial ytic	All	Enumeration	<b>PD_TBL_06</b>	N/A	1..1	C1	C7
... MDCX_PDIALY_FILL_DURATION_SETTING		Prescribed Fill duration	Prescribed fill duration	Intradial ytic	All	Numeric	XXX.X	min	1..1	C1	C7
... MDCX_PDIALY_FILL_DURATION		Actual fill duration	Actual fill duration	Intradial ytic	All	Numeric	XXX.X	min	1..1	C1	X
... MDCX_PDIALY_FILL_VOLUME_SETTING		Prescribed Fill volume	Prescribed Fill volume	Intradial ytic	All	Numeric Array	X.XXX	L	1..1	C1	C7
... MDCX_PDIALY_FILL_VOLUME		Actual Fill volume	Actual Fill volume	Intradial ytic	All	Numeric Array	X.XXX	L	1..1	C1	X
... MDCX_PDIALY_FLUID_TEMP_SETTING		Fluid temperature setting	Fluid temperature	Intradial ytic	All	Numeric	XX.X	degC	0..1	O	O
... MDCX_PDIALY_FLUID_TEMP		Fluid temperature	Fluid temperature	Intradial ytic	All	Numeric	XX.X	degC	0..1	O	X
... MDCX_PDIALY_DWELL_START_TIME		Dwell Start time	The time that the dwell phase started.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C2	X
... MDCX_PDIALY_DWELL_END_TIME		Dwell End time	The time that the dwell phase ended.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C2	X
... MDCX_PDIALY_DWELL_DURATION_SETTING		Prescribed Dwell Duration	Prescribed dwell duration	Intradial ytic	All	Numeric	XXX.X	min	1..1	M	C7
... MDCX_PDIALY_DWELL_DURATION		Actual Dwell Duration	Actual dwell duration	Intradial ytic	All	Numeric	XXX.X	min	1..1	M	X
... MDCX_PDIALY_DRAIN_START_TIME		Drain Start time	The time that the drain phase started.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C3	X
... MDCX_PDIALY_DRAIN_END_TIME		Drain End time	The time that the drain phase ended.	Intradial ytic	All	Date Time	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	1..1	C3	X

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_PDIALY_DRAIN_MODE		Drain mode	Drain mode	Intradialytic	All	Enumeration	PD_TBL_07	N/A	1..1	C3	C7
... MDCX_PDIALY_DRAIN_DURATION_SETTING		Prescribed Drain Duration	Prescribed drain duration	Intradialytic	All	Numeric	XXX.X	min	1..1	C3	C7
... MDCX_PDIALY_DRAIN_DURATION		Actual Drain Duration	Actual drain duration	Intradialytic	All	Numeric	XXX.X	min	1..1	C3	X
... MDCX_PDIALY_DRAIN_VOLUME_SETTING		Prescribed Drain volume	Prescribed Drain volume	Intradialytic	All	Numeric	X.XXX	L	1..1	C3	C7
... MDCX_PDIALY_DRAIN_VOLUME		Actual Drain volume	Actual Drain volume	Intradialytic	All	Numeric	X.XXX	L	1..1	C3	X
.. MDCX_DEV_PDIALY_SUPPLY_LINE_CHANNEL				All	All				0..1	O	X
... MDCX_PDIALY_PRESSURE_DIALYSATE		Dialysate Supply Pressure Sensor Value	Dialysate pressure in the supply line(s).	Intradialytic	All	Numeric	XXXX	mmHg	0..1	O	X
.. MDCX_DEV_PDIALY_DRAIN_LINE_CHANNEL				All	All				0..1	O	X
... MDCX_PDIALY_PRESSURE_DRAIN		Drain Pressure Sensor Value	Dialysate pressure in the drain line.	Intradialytic	All	Numeric	XXXX	mmHg	0..1	O	X
.. MDCX_DEV_PDIALY_PATIENT_LINE_CHANNEL				All	All				0..1	O	X
... MDCX_PDIALY_PRESSURE_PATIENT		Patient Pressure Sensor Value	Dialysate pressure in the patient line.	Intradialytic	All	Numeric	XXXX	mmHg	0..1	O	X
.. MDCX_DEV_PDIALY_HEATER_LINE_CHANNEL				All	All				0..1	O	X
.. MDCX_DEV_PDIALY_MACH_STATE_CHANNEL				All	All				0..1	O	X
... MDC_VENDOR_ALERT_NUM		Alert ID	The number associated with the current alert. If no alert is active, then this object is not reported.	All	All	Numeric	XXXX	N/A	0..1	O	X
... MDC_VENDOR_ALERT_TEXT		Alert Text	The text name of the current alert. If no alert is active, then this object is not reported. In most case this object will contain the text shown to the user.	All	All	String	N/A	N/A	0..1	O	X
... MDCX_PDIALY_AIR_SENSOR_STATE		Air sensor state	The state of the air sensor.	All	All	Enumeration	PD_TBL_11	N/A	0..1	O	X
... MDCX_PDIALY_PUMP_SPEED		Pump Speed	The speed of the pump.	All	All	Numeric	XXX.X	ml/min	0..1	O	X
... MDCX_PDIALY_PUMP_PRESSURE		Pump Pressure	The pressure at the pump outlet.	All	All	Numeric	XXX.X	mmHg	0..1	O	X



## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_PDIALY_TEMP_DIALYSATE_IN	tech high low thr	Dialysate Temperature Sensor Input	Dialysate temperature at the inlet.	All	All	Numeric	XX.X	°C	0..1	O	X
... MDCX_PDIALY_TEMP_DIALYSATE_OUT	tech high low thr	Dialysate Temperature Sensor Output	Dialysate temperature at the outlet.	All	All	Numeric	XX.X	°C	0..1	O	X
... MDCX_PDIALY_TEMP_HEATER_BAG		Dialysate Heater Bag Temperature Sensor	Dialysate temperature in the heater bag.	All	All	Numeric	XX.X	°C	0..1	O	X
... MDC_TEMP_DEVICE_INTERNAL		Device Internal Temperature Monitoring	Temperature inside the cyclor.	All	All	Numeric	XX.X	°C	0..1	O	X
... MDCX_PDIALY_AIR_SENSOR_IN		Inlet Air Sensor Value	Measurement of the air going to the patient.	All	All	Numeric	XXX	mL	0..1	O	X
... MDCX_PDIALY_AIR_SENSOR_OUT		Outlet Air Sensor Value	Measurement of the air coming from the patient.	All	All	Numeric	XXX	mL	0..1	O	X
... MDCX_PDIALY_LEAK_SENSOR		Leak Sensor(s)	Detects fluid leaks in the cyclor.	All	All	Bool			0..1	O	X
... MDCX_ATTR_AC_MAINS_RMS		Mains Power Monitoring	Monitors AC line in.	All	All	Numeric	XXX.X	V	0..1	O	X
... MDCX_PDIALY_VOLTAGE_MOTOR		Motor Voltage Monitoring	Motor Voltage Monitoring	All	All	Numeric	XXX.X	V	0..1	O	X
... MDCX_PDIALY_VOLTAGE_HEATER		Heater Input Voltage	The heater input voltage level is monitored.	All	All	Numeric	XXX.X	V	0..1	O	X
... MDC_ATTR_BATT_VOLTAGE		Battery Charge Level Monitoring	The battery voltage level is monitored.	All	All	Numeric	XX.X	V	0..1	O	X
... MDCX_ATTR_POWER_SUPPLY_VOLTAGE		Voltage Monitoring	The cyclor monitors power supply voltages.	All	All	Numeric	XX.X	V	0..1	O	X
... MDCX_PDIALY_PRESSURE_COMPRESSED_AIR		Pneumatic: Compressed air pressure	Air pressure is monitored. Compressed air is positive pressure.	All	All	Numeric	XXXX	mmHg	0..1	O	X
... MDCX_PDIALY_PRESSURE_VACUUM		Pneumatic: Vacuum pressure	Vacuum level is monitored. Vacuum is negative pressure.	All	All	Numeric	XXXX	mmHg	0..1	O	X
... MDCX_PDIALY_PRESSURE_PNEUMATIC		Hydraulic Fluid Pressure	Hydraulic pump fluid pressure is monitored.	All	All	Numeric	XXXX	mmHg	0..1	O	X
... MDCX_PDIALY_VOLUME_SUPPLY_BAG		Dialysate Supply Bag Volume Status Monitoring	The cyclor monitors the volume of solution pumped out of the supply bag(s).	All	All	Numeric	X.XXX	L	0..1	O	X
... MDCX_PDIALY_VOLUME_HEATER_BAG		Dialysate Heater Bag Volume Status Monitoring	The cyclor monitors the volume of solution pumped into and out of the heater bag.	All	All	Numeric	X.XXX	L	0..1	O	X
... MDCX_PDIALY_WEIGHT_HEATER_BAG		Heater Bag Weight	Weight of the heater bag	All	All	Numeric	XXXX.X	g	0..1	O	X
. MDCX_DEV_PDIALY_PATIENT_VMD		Patient "Device"		Intradialytic	All				0..1	O	O
.. MDCX_DEV_PATIENT_INFO		Patient channel	A channel to hold information about the patient.	Intradialytic	All				0..1	O	O

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDCX_PDIALY_PATIENT_TRANSPORT_TYPE		Transport Type	The rate at which small molecules cross the peritoneal membrane.	Intradialytic	All	Enumeration	PD_TBL_08	N/A	0..1	O	O
... MDCX_PDIALY_PATIENT_CATHETER_TYPE		Catheter Type	Catheter Type	Intradialytic	All	Enumeration	PD_TBL_09	N/A	0..1	O	O
... MDCX_PDIALY_PATIENT_CATHETER_DATE		Catheter Placement date	Catheter Placement date	Intradialytic	All	Date	YYYYMMDDHHMMSS[.SSS]+/-ZZZZ	N/A	0..1	O	O
... MDCX_PDIALY_PAT_TARGET_WT		Target Weight	The patient's target weight. This is captured for information purposes only and is not used by the dialysis machine when determining treatment settings.	Intradialytic	All	Numeric	XXX.X	Kg	0..1	O	O
. MDC_DEV_PRESS_BLD_NONINV_VMD		Non-invasive blood pressure	Instrument for the non-invasive measurement of blood pressure.	Intradialytic	All				0..1	O	O
.. MDC_DEV_PRESS_BLD_NONINV_CHAN				Intradialytic	All				0..1	O	O
... MDCX_DIALY_BLD_MEAS_PHASE		Blood Pressure Measurement Phase	The phase of the treatment when the patient should perform a blood pressure measurement.	Intradialytic	All	String	PD_TBL_10	N/A	0..1	O	O
... MDCX_DIALY_PUL_MEAS_PHASE		Pulse Measurement Phase	The phase of the treatment when the patient should perform a pulse measurement.	Intradialytic	All	String	PD_TBL_10	N/A	0..1	O	O
... MDC_ATTR_TIME_PD_MSMT		Blood Pressure Measurement Interval Setting	Time allowed to elapse between automatic blood pressure measurements. A value of zero indicates that automatic measurements are not being taken.	Intradialytic	All	Numeric	XXX	minutes	0..1	O	O
... MDC_ATTR_PT_BODY_POSN		Patient body position	Patient body position at the time of the {blood pressure} measurement.	Intradialytic	All	Enumeration	HD_TBL_04	N/A	0..1	O	O
... MDC_PRESS_BLD_NONINV_DIA	phys high low thr	Diastolic Pressure	Minimum arterial pressure during relaxation and dilatation of the ventricles of the heart when the ventricles fill with blood	Intradialytic	All	Numeric	XXX	mmHg	0..1	O	X
... MDC_ATTR_VAL_RANGE		Diastolic Pressure Entry Limits	The minimum and maximum value that the user can enter for a diastolic pressure.	Intradialytic	All	Numeric Range	XXX^XXX	mmHg	0..1	O	O
... MDC_PULS_RATE_NON_INV	phys high low thr	Heart Rate	Number of heart beats per minute	Intradialytic	All	Numeric	XXX	BPM	0..1	O	X
... MDC_ATTR_VAL_RANGE		Pulse Entry Limits	The minimum and maximum value that the user can enter for a pulse measurement.	Intradialytic	All	Numeric Range	XXX^XXX	BPM	0..1	O	O
... MDC_PRESS_BLD_NONINV_MEAN		Mean Arterial Pressure	Average pressure in a patient's arteries during one cardiac cycle	Intradialytic	All	Numeric	XXX	mmHg	0..1	O	X
... MDC_ATTR_VAL_RANGE		Mean Arterial Pressure Entry Limits	The minimum and maximum value that the user can enter for a mean arterial pressure.	Intradialytic	All	Numeric Range	XXX^XXX	mmHg	0..1	O	O

## Dialysis Machine Implementation Guide

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDC_PRESS_BLD_NONINV_SYS	phys high thr	Systolic Pressure	Maximum arterial pressure during contraction of the left ventricle of the heart	Intradial ytic	All	Numeric	XXX	mmHg	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Systolic Pressure Entry Limits	The minimum and maximum value that the user can enter for a systolic pressure.	Intradial ytic	All	Numeric Range	XXX^XXX	mmHg	0..1	O	O
. MDC_DEV_SPEC_PROFILE_GLUCOSE		<b>Glucose "Device"</b>		Intradial ytic	<b>All</b>				0..1	<b>O</b>	<b>O</b>
.. MDC_DEV_CHAN		Glucose Channel	A channel to hold information about the patient's glucose.	<b>Intradial ytic</b>	All				0..1	<b>O</b>	<b>O</b>
... MDCX_DIALY_PAT_GLU_MEAS_PHASE		Glucose Measurement Phase	The phase of the treatment when the patient should perform a glucose measurement.	Intradial ytic	All	String	<b>PD_TBL_10</b>	N/A	0..1	O	O
... MDC_CONC_GLU_CAPILLARY_PLASMA		Blood Glucose	Blood glucose at the start of treatment	Intradial ytic	All	Numeric	XX.X	mmol /L	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Glucose Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XX.X^XX.X	mmol /L	0..1	O	O
. MDC_DEV_SPEC_PROFILE_SCALE		<b>Scale "Device"</b>		Intradial ytic	<b>All</b>				0..1	<b>O</b>	<b>O</b>
.. MDC_DEV_CHAN		Weight Channel	A channel to hold information about the patient's weight.	<b>Intradial ytic</b>	All				0..1	<b>O</b>	<b>O</b>
... MDCX_DIALY_PAT_WT_MEAS_PHASE		Weight Measurement Phase	The phase of the treatment when the patient should perform a weight measurement.	Intradial ytic	All	String	<b>PD_TBL_10</b>	N/A	0..1	O	O
... MDC_MASS_BODY_ACTUAL		Patient Weight	Weight of patient at beginning or end of treatment	Intradial ytic	All	Numeric	XXX.X	Kg	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Weight Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XXX.X^XXX.X	Kg	0..1	O	O
. MDC_DEV_SPEC_PROFILE_TEMP		<b>Thermometer "Device"</b>		Intradial ytic	<b>All</b>				0..1	<b>O</b>	<b>O</b>
.. MDC_DEV_CHAN		Temperature Channel	A channel to hold information about the patient's temperature.	<b>Intradial ytic</b>	All				0..1	<b>O</b>	<b>O</b>
... MDCX_DIALY_PAT_TEMP_MEAS_PHASE		Temperature Measurement Phase	The phase of the treatment when the patient should perform a temperature measurement.	Intradial ytic	All	String	<b>PD_TBL_10</b>	N/A	0..1	O	O
... MDC_TEMP_TYMP	phys high thr	Patient Temperature	Patient's temperature when measured at the ear drum.	Intradial ytic	All	Numeric	XX.X	°C	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Temperature Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XX.X^XX.X	°C	0..1	O	O
... MDC_TEMP_ORAL	phys high thr	Patient Temperature	Patient's temperature when measured orally.	Intradial ytic	All	Numeric	XX.X	°C	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Temperature Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XX.X^XX.X	°C	0..1	O	O
... MDC_TEMP_BODY	phys high thr	Patient Temperature	Patient's temperature when measured at an unknown location	Intradial ytic	All	Numeric	XX.X	°C	0..1	O	X
.... MDC_ATTR_VAL_RANGE		Temperature Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XX.X^XX.X	°C	0..1	O	O

Table 4 – Peritoneal Dialysis Data Objects

REFID	Alert Type	Common Term	Description	Phase	Temporal	Data Type	Format	UOM	Card	Use	Rx Use
... MDC_TEMP_AXILLA	phys high thr	Patient Temperature	Patient's temperature when measured at the arm pit.	Intradial ytic	All	Numeric	XX.X	°C	0..1	O	X
... MDC_ATTR_VAL_RANGE		Temperature Entry Limits	The minimum and maximum values for entry by the user.	Intradial ytic	All	Numeric Range	XX.X^XX.X	°C	0..1	O	O

## 10.3 Events

Table 5 – Peritoneal Dialysis Alarms/Alerts

Source	Event	Alert Type	Dialysis Common Name	Definition	Usage
MDC_MASS_BODY_ACTUAL	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check the patient weight.	O
MDCX_DEV_PDIALY_DRAIN_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED	tech		A line is blocked; checking is required.	O
MDCX_DEV_PDIALY_DRAIN_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_CHK_REQD	tech		Slow flow in drain line should be checked by user.	O
MDCX_DEV_PDIALY_DRAIN_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_OUTFLOW	tech		Blocked flow in drain line.	O
MDCX_DEV_PDIALY_DRAIN_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW	tech		Slow flow in drain line.	O
MDCX_DEV_PDIALY_DRAIN_LINE_CHAN	MDCX_EVT_PDIALY_LINE_CHK	tech		Check drain line.	O
MDCX_DEV_PDIALY_HEATER_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED	tech		A line is blocked	O
MDCX_DEV_PDIALY_HEATER_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW	tech		Slow flow in heater line.	O
MDCX_DEV_PDIALY_HEATER_LINE_CHAN	MDCX_EVT_PDIALY_HEATER_BAG_POSN_ERR_OR_LINE_CLAMPED_CHK_REQD	tech		Verify that the heater bag is on the heater tray and the heater line is not clamped	O
MDCX_DEV_PDIALY_HEATER_LINE_CHAN	MDCX_EVT_PDIALY_LINE_CHK	tech		Check heater line.	O
MDCX_DEV_PDIALY_MACH_STATE_CHAN	MDCX_EVT_PROGRAM_SYNC_ERROR	tech		Error during monitoring of program execution synchronization	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_AC_POWER_FAIL	tech		AC Power fail (or short circuit in one of the devices connected to an auxiliary outlet)	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_BATT_LO	tech	Low Battery	{device or sensor} Battery low	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_HID_MALF	tech		HID failure	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_INTERNAL_SYSTEM_FAULT	tech		Unspecified system fault.	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_MSMT_OUT_OF_RANGE	tech		Unspecified measurement out of range	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_SELFTEST_FAILURE	tech		Unspecified selftest failure	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_SPEAKER_FAIL	tech		The speaker power and signals are monitored to ensure the speaker is active.	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_TOUCHSCREEN_FAIL	tech		One or more elements touchscreen have failed.	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDC_EVT_UNRECOV_ERR	tech		Unrecoverable state machine error in communication.	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_AC_POWER_RESTORED	tech		AC power restored	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_DEVICE_TILTED	tech		Device tilted	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_EXTERNAL_MEMORY_ERROR	tech		External memory device error (e.g. external memory not available, cannot communicate, cannot be read, cannot be written)	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_DISPLAY_FAIL	tech		One or more elements of the display have failed.	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_PATIENT_ID_INVALID	tech		Invalid Patient ID	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_TREATMENT_ERR	tech		An error occurred during treatment	O
MDCX_DEV_PDIALY_MACHINE_MDS	MDCX_EVT_UNRECOV_TREATMENT_ERR	tech		An unrecoverable error occurred during treatment	O

## Dialysis Machine Implementation Guide

Table 5 – Peritoneal Dialysis Alarms/Alerts

Source	Event	Alert Type	Dialysis Common Name	Definition	Usage
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_DRAIN_COMPLICATION_CHK_REQD	tech		Drain complication encountered: check patient and drain lines	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FILL_COMPLICATION_CHK_REQD	tech		Fill complication encountered: check heater, patient and drain lines	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED	tech		Tubing system blocked during inflow to patient.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_CHK_REQD	tech		Slow flow in patient line should be checked by user.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_INFLOW	tech		Blocked flow in patient inlet line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_OUTFLOW	tech		Blocked flow in patient outlet line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW	tech		Slow flow in patient line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW_INFLOW	tech		Slow flow in patient inlet line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW_OUTFLOW	tech		Slow flow in patient outlet line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_LINE_ABSENT	tech		Missing or disconnected patient line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_PDIALY_LINE_CHK	tech		Check patient line.	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_STAT_PDIALY_PATIENT_DRAIN_PHASE	tech		Patient drain phase	0
MDCX_DEV_PDIALY_PATIENT_LINE_CHAN	MDCX_EVT_STAT_PDIALY_PATIENT_FILL_PHASE	tech		Patient fill or refill phase	0
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED	tech		Tubing system blocked during outflow from patient.	0
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_BLOCKED_CHK_REQD	tech		Slow flow in supply line should be checked by user.	0
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN	MDCX_EVT_PDIALY_FLOW_SLOW	tech		Slow flow in supply line.	0
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN	MDCX_EVT_PDIALY_LINE_ABSENT	tech		Missing or disconnected dialysate line.	0
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN	MDCX_EVT_PDIALY_LINE_CHK	tech		Check dialysate/supply line.	0

## Dialysis Machine Implementation Guide

Table 5 – Peritoneal Dialysis Alarms/Alerts

Source	Event	Alert Type	Dialysis Common Name	Definition	Usage
MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN					
MDCX_DEV_PDIALY_FLUID_1_CHAN MDCX_DEV_PDIALY_FLUID_2_CHAN MDCX_DEV_PDIALY_FLUID_3_CHAN MDCX_DEV_PDIALY_FLUID_4_CHAN MDCX_DEV_PDIALY_FLUID_5_CHAN MDCX_DEV_PDIALY_FLUID_6_CHAN MDCX_DEV_PDIALY_SUPPLY_LINE_CHAN	MDCX_EVT_PDIALY_SOLUTION_NOT_REQD	tech		Dialysate solution detected on a line where it is not required by the prescription.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check the therapy time setting.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check the total UF setting.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check the total volume setting.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check an unspecified treatment setting.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_LO	tech		Total UF volume is too low.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDC_EVT_LO	tech		Initial drain volume is too low.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDCX_EVT_PDIALY_PRESCRIPTION_ERROR	tech		Invalid prescription received.	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDCX_EVT_PDIALY_ULTRAFILTRATION_NEG	tech		Negative Ultrafiltration	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDCX_EVT_PDIALY_ULTRAFILTRATION_POS	tech		Positive Ultrafiltration	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDCX_EVT_TREATMENT_DATA_INVALID	tech		Invalid treatment data (e.g. possibly after a power failure)	0
MDCX_DEV_PDIALY_TREATMENT_CHAN	MDCX_EVT_VOLUME_DEVIATION	tech		Patient volume monitoring has identified a deviation	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_ADVIS_PDIALY_USER_MESSAGE	tech		Machine specific advisory	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_COMM_ERR_EXTERNAL	tech		Internal communication error.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_COMM_ERR_INTERNAL	tech		Internal communication error.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_DISPLAYED_VALUES_INCONSISTENT	tech		Inconsistent safety-relevant information displayed	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_DISPOSABLE_SET_LEAK	tech		Disposable set leak	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_LOAD_NEW_SET	tech		Load a new set	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_LOAD_NEW_SET_AND_BAGS	tech		Load new set and bags	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_ADVIS_TREATMENT_MAY_BE_INSUFFICIENT	tech		Potentially insufficient treatment detected	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_AIR_DETECTED	tech		Maximum allowable air detected in cassette or lines.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_CASSETTE_ABSENT	tech		Cassette not detected	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_CASSETTE_AIR_DETECTED	tech		Maximum allowable air detected in cassette.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_CASSETTE_FAILURE	tech		Cassette failure	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_COMM_FAILURE	tech		The cyclor monitors the information flow to the cloud to determine that the information has been sent.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_INITIAL_DRAIN_VOLUME_UNEXPECTED	tech		Unexpected volume (of fluid) during initial drain	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_PRIMING_ERROR	tech		Priming error	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_PUMP_FAILURE	tech		The cyclor has determined that the pump is no longer functional.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_RELOAD_SET	tech		Reload the disposable set	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_SET_POSITION_WRONG	tech	Disposable Set Position Sensor	The disposable set status is monitored.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_STRAIN_GAUGE_FAILURE	tech	Strain Gauge Sensor Status	The strain gauge state is monitored.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_TREATMENT_PAUSED_TOO_LONG	tech		Treatment delay	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PDIALY_TREATMENT_TIME_EXCEEDED	tech		Treatment time exceeded	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_PT_SECURITY_KEY_ABSENT	tech		Patient security key or other credential is missing.	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_STAT_DOOR_NOT_CLOSED_OR_LOCKED_CORRECTLY	tech		Door is not closed or locked correctly.	0

## Dialysis Machine Implementation Guide

Table 5 – Peritoneal Dialysis Alarms/Alerts

Source	Event	Alert Type	Dialysis Common Name	Definition	Usage
MDCX_DEV_PDIALY_VMD	MDCX_EVT_STAT_PDIALY_HEATER_BAG_FILLING	tech		Filling heater bag	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_STAT_PDIALY_PT_DISCONN	tech		Patient disconnect timeout (intentional)	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_STAT_PDIALY_TREATMENT_PAUSED	tech		Treatment paused (e.g. stop button pressed)	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_SW_ERR	tech		Software error	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_TIMER_ERR	tech		Error during the plausibility check of the time counter	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_TREATMENT_RESULTS_INCONSISTENT	tech		Inconsistent treatment results	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_UNRECOV_TREATMENT_DATA_ERR	tech		An unrecoverable data error occurred during treatment	0
MDCX_DEV_PDIALY_VMD	MDCX_EVT_VALVE_MALF	tech		Valve malfunction	0
MDCX_PDIALY_DRAIN_VOLUME	MDC_EVT_HI	tech		Drain volume too high	0
MDCX_PDIALY_DRAIN_VOLUME	MDC_EVT_LO	tech		Drain volume too low	0
MDCX_PDIALY_DRAIN_VOLUME	MDCX_EVT_PDIALY_DRAIN_PHASE_NOT_COMPLETED	tech		Drain phase not completed	0
MDCX_PDIALY_FILL_VOLUME_SETTING	MDC_EVT_ADVIS_SETTING_CHK	tech		User should check the day fill volume setting.	0
MDCX_PDIALY_FILL_VOLUME_SETTING	MDC_EVT_ADVIS_SETTING_CHK	tech		User should check the fill volume setting.	0
MDCX_PDIALY_FILL_VOLUME_SETTING	MDC_EVT_ADVIS_SETTING_CHK	tech		User should check the last fill volume setting.	0
MDCX_PDIALY_FILL_VOLUME_SETTING	MDC_EVT_ADVIS_SETTING_CHK	tech		User should check the nigh fill volume setting.	0
MDCX_PDIALY_FILL_VOLUME	MDCX_EVT_PDIALY_FILL_PHASE_NOT_COMPLETED	tech		Fill or refill phase not completed	0
MDCX_PDIALY_PRESSURE_DRAIN	MDC_EVT_MSMT_OUT_OF_RANGE	tech		Drain pressure out of range	0
MDCX_PDIALY_PRESSURE_DRAIN	MDC_EVT_SELFTEST_FAILURE	tech		Drain pressure selftest failure	0
MDCX_PDIALY_PRESSURE_PATIENT	MDC_EVT_HI	tech		Patient pressure too high	0
MDCX_PDIALY_PRESSURE_PATIENT	MDC_EVT_LO	tech		Patient Pressure too low	0
MDCX_PDIALY_PRESSURE_PATIENT	MDC_EVT_MSMT_ERR	tech		Unspecified error in patient pressure	0
MDCX_PDIALY_PRESSURE_PATIENT	MDC_EVT_SELFTEST_FAILURE	tech		Patient pressure selftest failure	0
MDCX_PDIALY_PRESSURE_PATIENT	MDCX_EVT_HDIALY_PT_LINE_PRESS_HI_LOW_POSN_CHK_REQD	tech		Patient line pressure too high or low; check patient position	0
MDCX_PDIALY_PRESSURE_PNEUMATIC	MDC_EVT_HI	tech		Pneumatic pressure too high	0
MDCX_PDIALY_PRESSURE_PNEUMATIC	MDC_EVT_LEAK	tech		Leakage in a gas or fluid filled system detected.	0
MDCX_PDIALY_PRESSURE_PNEUMATIC	MDC_EVT_LO	tech		Pneumatic pressure too low	0
MDCX_PDIALY_PRESSURE_PNEUMATIC	MDC_EVT_SELFTEST_FAILURE	tech		Pneumatic pressure selftest failure	0
MDCX_PDIALY_PRESSURE_PNEUMATIC	MDCX_EVT_PDIALY_VALVE_LEAK	tech		Valve or pressure leak	0
MDCX_PDIALY_PRESSURE_VACUUM	MDC_EVT_HI	tech		Vacuum pressure too high	0
MDCX_PDIALY_PRESSURE_VACUUM	MDC_EVT_LEAK	tech		Leakage in a gas or fluid filled system detected.	0
MDCX_PDIALY_PRESSURE_VACUUM	MDC_EVT_LO	tech		Vacuum pressure too low	0
MDCX_PDIALY_PRESSURE_VACUUM	MDC_EVT_MSMT_OUT_OF_RANGE	tech		Vacuum pressure out of range	0
MDCX_PDIALY_PRESSURE_VACUUM	MDC_EVT_SELFTEST_FAILURE	tech		Vacuum pressure selftest failure	0
MDCX_PDIALY_TEMP_DIALYSATE_IN	MDC_EVT_HI	tech		Dialysate inlet temperature too high	0
MDCX_PDIALY_TEMP_DIALYSATE_IN	MDC_EVT_LO	tech		Dialysate inlet temperature too low	0
MDCX_PDIALY_TEMP_DIALYSATE_OUT	MDC_EVT_HI	tech		Dialysate temperature too high	0
MDCX_PDIALY_TEMP_DIALYSATE_OUT	MDC_EVT_LO	tech		Dialysate temperature too low	0
MDCX_PDIALY_TEMP_DIALYSATE_OUT	MDC_EVT_MSMT_ERR	tech		Unspecified error in dialysate temperature	0
MDCX_PDIALY_TEMP_DIALYSATE_OUT	MDC_EVT_MSMT_OUT_OF_RANGE	tech		Dialysate temperature out of range	0
MDCX_PDIALY_TEMP_DIALYSATE_OUT	MDCX_EVT_STAT_SOLUTION_WARMING	tech		Dialysate solution is warming	0
MDCX_PDIALY_TEMP_HEATER_BAG	MDC_EVT_HI	tech		Heater bag temperature too high	0
MDCX_PDIALY_TEMP_HEATER_BAG	MDC_EVT_LO	tech		Heater bag temperature too low	0
MDCX_PDIALY_TEMP_HEATER_BAG	MDC_EVT_MSMT_OUT_OF_RANGE	tech		Heater Bag Temperature out of range	0
MDCX_PDIALY_TEMP_HEATER_BAG	MDC_EVT_MSMT_OUT_OF_RANGE_HI	tech		Heater bag temperature too high	0
MDCX_PDIALY_TEMP_HEATER_BAG	MDC_EVT_SELFTEST_FAILURE	tech		Heater bag temperature selftest failure	0

## Dialysis Machine Implementation Guide

Table 5 – Peritoneal Dialysis Alarms/Alerts

Source	Event	Alert Type	Dialysis Common Name	Definition	Usage
MDCX_PDIALY_TEMP_HEATER_BAG	MDCX_EVT_TEMP_ERR	tech		Temperature sensor error	0
MDCX_PDIALY_TIDAL_PERCENTAGE	MDC_EVT_ADVIS_SETTING_CHK	tech		The user should check the tidal setting.	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDC_EVT_HI_GT_LIM	tech		A metric exceeds a given threshold	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDC_EVT_MSMT_ERR	tech		Unspecified error in heater bag weight	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDC_EVT_MSMT_FAIL	tech		Failure in heater bag weight	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDCX_EVT_PDIALY_ABRUPT_WEIGHT_CHANGE	tech		Abrupt weight change; possible scale interference.	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDCX_EVT_PDIALY_FLUID_BAG_WEIGHT_VOLUME_ERROR	tech		Fluid bag weight scale (volume) measurement error.	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDCX_EVT_PDIALY_HEATER_BAG_ON_HEATER_TRAY_UNEXPECTED	tech		Heater bag detected on the heater tray when not expected.	0
MDCX_PDIALY_WEIGHT_HEATER_BAG	MDCX_EVT_WEIGHT_INVALID_CHK_REQD	tech		The dialysis machine has detected an invalid weight change and the user must verify/confirm the weight.	0
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_HI	phys	Diastolic High	Diastolic pressure is too high.	0
MDC_PRESS_BLD_NONINV_DIA	MDC_EVT_LO	phys	Diastolic Low	Diastolic pressure is too low.	0
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_HI	phys	Systolic High	Systolic pressure is too high.	0
MDC_PRESS_BLD_NONINV_SYS	MDC_EVT_LO	phys	Systolic Low	Systolic pressure is too low.	0
MDC_PULS_RATE_NON_INV	MDC_EVT_HI	phys	Pulse High	Pulse is too high.	0
MDC_PULS_RATE_NON_INV	MDC_EVT_LO	phys	Pulse Low	Pulse is too low.	0
MDC_TEMP_TYMP	MDC_EVT_HI	phys	Temperature High	Temperature is too high	0
MDC_TEMP_ORAL	MDC_EVT_HI	phys	Temperature High	Temperature is too high	0
MDC_TEMP_BODY	MDC_EVT_HI	phys	Temperature High	Temperature is too high	0
MDC_TEMP_AXILLA	MDC_EVT_HI	phys	Temperature High	Temperature is too high	0



## Appendix A – HL7 Batch Protocol

To capture an entire dialysis machine treatment, the HL7 Batch Protocol should be used. This protocol combines a series of intact MSH messages into a single message. The Batch Protocol preserves the content of the original PCD-01 DEC MSH messages. Although they are not required the manufacturer may include the ACK messages in the batch.

Segment	Meaning	Use	Card
[FHS]	File Header Segment	R	[1..1]
{	Batch Begin		
[BHS]	Batch Header Segment	R	[1..1]
{ [	Message begin		
MSH	One or more HL7 Messages	R	[1..N]
...			
...			
] }	MESSAGE end	R	[1..1]
[BTS]	Batch Trailer Segment	R	[1..1]
[FTS]	File Trailer Segment	R	[1..1]

### A.1 Segments

#### A.1.1 BHS – Batch Header Segment

The BHS segment defines the start of a batch.

HL7 Attribute Table - BHS Segment

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	R		Batch Field Separator
2	ST	R		Batch Encoding Characters
3	HD	R		Batch Sending Application
4	HD	O		Batch Sending Facility
5	HD	O		Batch Receiving Application
6	HD	O		Batch Receiving Facility
7	DTM	O		Batch Date/Time of Message
8	ST	O		Batch Security
9	ST	O		Batch Name/ID/Type
10	ST	O		Batch Comment
11	ST	O		Batch Control ID
12	ST	O		Reference Batch Control ID
13	HD	O		Batch Sending Network Address
14	HD	O		Batch Receiving Network Address

This standard does not change any field definitions from the HL7 Standard.

BHS-3 Batch Sending Application

This Standard: This value is the same as the MSH-3 of the observations.

#### A.1.2 BTS – Batch Trailer Segment

The BTS segment defines the end of a batch.

**HL7 Attribute Table - BTS Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	O		Batch Message Count
2	ST	O		Batch Comments
3	NM	X		Batch Totals

**BTS-3 Batch Totals (NM)**

HL7 Definitions: We encourage new users of this field to use the HL7 Version 2.3 data type of NM and to define it as "repeating." This field contains the batch total. If more than a single batch total exists, this field may be repeated.

This Standard: The dialysis machines do not use the batch totals as this information is already conveyed by observation data.

**A.1.3 FHS – File Header Segment**

The FHS segment is used to head a file (group of batches).

**HL7 Attribute Table - FHS Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	R		File Field Separator
2	ST	R		File Encoding Characters
3	HD	R		File Sending Application
4	HD	O		File Sending Facility
5	HD	O		File Receiving Application
6	HD	O		File Receiving Facility
7	DTM	O		File Creation Date/Time
8	ST	O		File Security
9	ST	O		File Name/ID
10	ST	O		File Comment
11	ST	O		File Control ID
12	ST	O		Reference File Control ID
13	HD	O		File Sending Network Address
14	HD	O		File Receiving Network Address

This standard does not change any field definitions from the HL7 Standard.

**FHS-3 File Sending Application**

This Standard: This value is the same as the MSH-3 of the observations.

**A.1.4 FTS – File Trailer Segment**

The FTS segment defines the end of a file.

**HL7 Attribute Table - FHS Segment**

SEQ	DT	OPT	TBL#	ELEMENT NAME
1	ST	O		File Batch Count
2	ST	O		File Trailer Comments

This standard does not change any field definitions from the HL7 Standard.